

Visual Otoscope



Please read this instructions carefully
before using and follow the instructions

Intended Use

This otoscope can be used for daily inspection and cleaning of ear, nose, mouth or throat. You can take photos and videos when you use the otoscope, and send them to healthcare professional for a quick and easy consultation. This ear scope is not only convenient for yourself, but also for your friends, family, even pets.

Safety Warning

- The product can not be disassembled by the user, and can only be repaired by authorized professionals.
- Do not wash with water or completely soak the product in water.
- Recharge the battery regularly when not in use for a long time.
- Do not apply the strong corrosive chemical solvents to clean the product.
- Do not change or transform the use of the product at will.
- Consult health care provider to answer your personal questions before use.
- Do not aggressively shove or stick the camera into the ear canal or scrape the skin of the canal as it may cause bleeding, discomfort, or damage to the eardrum.
- If using on children, hold child's head firmly to avoid jumping or movement.
- The adult supervision is required for use of this device.
- The ear otoscoped is unsuitable for children under 3 year old.
- Discontinue use and contact your health care provider if you feel pain.
- Discontinue use and immediately turn device off if it becomes hot.
- Strictly for use in the outer ear canal (from the eardrum outward).
- Not for use in mouth or throat.
- Do not use near water or bathtubs.

Kit Includes

- Wireless Ear Otoscoped x1
- Child safety cone tip x1
- Adult safety cone tip x1
- Soft ear spoon x2
- USB charging cable x1
- User Manual x1

Directions For Use

- 1.Charge your ear camera by plugging in the micro-USB power cord to the device and a charging socket.
- 2.Scan QR code to download the app.
- 3.To power on the device, remove the cap. A blue light will appear to let you know your device is on.
- 4.Connect to the device's Wi-Fi by exiting the app. Go to your Wi-Fi settings and select youbetterWi-Fi.
- 5.Return to the Eosera app and press start to begin viewing live footage. Select "left ear" to switch between your right and left ear.
- 6.Press the photo button to take a picture. Press the video button to start recording a video. Press the video button again to stop recording.
- 7.Select the album button on the bottom of your screen to view photos and videos. Photos and videos will be saved to your phone or tablet automatically.
- 8.Place the cap back onto the device to turn off.

Maintenance

- 1.The product cannot be disassembled by the user and can only be repaired by authorized professionals.
- 2.If do not use this product for a long time, please full of electricity and charge it timely before using again.
- 3.Wipe the product with a dry towel. Do not rinse the product directly with water or immerse it in water.

Specification

Product Model: OS806
Power Supply: Lithium battery(Built-in)
Product Size:135× ϕ 13.6mm
Input Voltage: DC 5V
Charging Time: 1.5 hours
Operation condition: +5 ~ +40°C, 35 - 70% relative maximum humidity
Storage conditions: -5 ~ +50°C, 30 - 80% relative maximum humidity
This device complies with the requirements of the Medical Device.
Directive EU 2017/745 regulation.
Technical alterations reserved.

Guarantee

This device is covered by a 1 year guarantee from the date of purchase. During this guarantee period, at our discretion, we will repair or replace the defective product free of charge. Opening or altering the device invalidates the guarantee. The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance.

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

Storage conditions and methods

The Products shall be stored in a cool and dry place away from water source and heating.
Keep away from sunlight, and there shall be no acid corrosive gas or other harmful gas in the surrounding air.
The storage time should not exceed five years from the production date.
Recommended storage temperature is -5 ~ +50°C, relative humidity is between 30%~80%.



Scan QR code to download the app.

Trouble Shooting

Problem	Solution
No response after startup	Please follow the correct startup method and try for several times Please use the matching charger to charge, the charging time is about 1.5 hours
The camera is off	Power off and then restartup Please confirm whether the product has electricity
APP exit instantaneously	Uninstall current APP and scan the QR code again for reinstallation
Failed to display the image after successfully connecting to APP	Disconnect the phone from WIFI, restart device and reconnect to WIFI
WiFi connections are erratic, sometimes disconnecting within few seconds	Restart device and reconnect to APP
If the above operation fails to solve the problem, please contact customer service, we will provide you with practical solutions	



Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.



Humiss International B.V.
Joop Geesinkweg 701, 1114 AB
Amsterdam-Duivendrecht, the Netherlands



Ningbo Albert Novosino Co.,Ltd.
No.1 Xinheng 3 Road Cicheng Town, Jiangbei
District, 315036 Ningbo, Zhejiang , China.
www.albertnovosino.com



Made in China

120 mm

60 mm

115 mm

Visual Otoscope



Please read this instructions carefully
before using and follow the instructions

Intended Use

This otoscope can be used for daily inspection and cleaning of ear, nose, mouth or throat. You can take photos and videos when you use the otoscope, and send them to healthcare professional for a quick and easy consultation. This ear scope is not only convenient for yourself, but also for your friends, family, even pets.

Safety Warning

- The product can not be disassembled by the user, and can only be repaired by authorized professionals.
- Do not wash with water or completely soak the product in water.
- Recharge the battery regularly when not in use for a long time.
- Do not apply the strong corrosive chemical solvents to clean the product.
- Do not change or transform the use of the product at will.
- Consult health care provider to answer your personal questions before use.
- Do not aggressively shove or stick the camera into the ear canal or scrape the skin of the canal as it may cause bleeding, discomfort, or damage to the eardrum.
- If using on children, hold child's head firmly to avoid jumping or movement.
- The adult supervision is required for use of this device.
- The ear otoscoped is unsuitable for children under 3 year old.
- Discontinue use and contact your health care provider if you feel pain.
- Discontinue use and immediately turn device off if it becomes hot.
- Strictly for use in the outer ear canal (from the eardrum outward).
- Not for use in mouth or throat.
- Do not use near water or bathtubs.

Kit Includes

Wireless Ear Otoscoped x1
USB charging cable x1
User Manual x1

Directions For Use

- 1.Charge your ear camera by plugging in the micro-USB power cord to the device and a charging socket.
- 2.Scan QR code to download the app.
- 3.Open the lid, press the button, the device lights up and you can start using it.
- 4.Connect to the device's Wi-Fi by exiting the app. Go to your Wi-Fi settings and select youbetterWi-Fi.
- 5.Return to the youbetter app and press start to begin viewing live footage. Select "left ear" to switch between your right and left ear.
- 6.Press the photo button to take a picture. Press the video button to start recording a video. Press the video button again to stop recording.
- 7.Select the album button on the bottom of your screen to view photos and videos. Photos and videos will be saved to your phone or tablet automatically.
- 8.Please press the button to close.

Maintenance

- 1.The product cannot be disassembled by the user and can only be repaired by authorized professionals.
- 2.If do not use this product for a long time, please full of electricity and charge it timely before using again.
- 3.Wipe the product with a dry towel. Do not rinse the product directly with water or immerse it in water.

Specification

Product Model: OS811
Product Size:φ14×13.2×141mm
Input Voltage: DC 5V
Charging Time: 1.5 hours
Operation condition: +5 ~ +40°C, 35 - 70% relative maximum humidity
Storage conditions: -5 ~ +50°C, 30 - 80% relative maximum humidity
This device complies with the requirements of the Medical Device.
Directive EU 2017/745 regulation.
Technical alterations reserved.

Guarantee

This device is covered by a 1 year guarantee from the date of purchase. During this guarantee period, at our discretion, we will repair or replace the defective product free of charge. Opening or altering the device invalidates the guarantee. The following items are excluded from the guarantee:
• Transport costs and risks of transport.
•Damage caused by incorrect application or non-compliance with the instructions for use.
•Damage caused by leaking batteries.
•Damage caused by accident or misuse.
•Packaging/storage material and instructions for use.
•Regular checks and maintenance.
Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

Storage conditions and methods

The Products shall be stored in a cool and dry place away from water source and heating.
Keep away from sunlight, and there shall be no acid corrosive gas or other harmful gas in the surrounding air.
The storage time should not exceed five years from the production date.
Recommended storage temperature is -5 ~ +50°C, relative humidity is between 30%~80%.



Scan QR code to download the app.

Trouble Shooting

Problem	Solution
No response after startup	Please follow the correct startup method and try for several times. Please use the matching charger to charge, the charging time is about 1.5 hours
The camera is off	Power off and then restartup Please confirm whether the product has electricity
APP exit instantaneously	Uninstall current APP and scan the QR code again for reinstallation
Failed to display the image after successfully connecting to APP	Disconnect the phone from WIFI, restart device and reconnect to WIFI
WIFI connections are erratic, sometimes disconnecting within few seconds	Restart device and reconnect to APP
If the above operation fails to solve the problem, please contact customer service, we will provide you with practical solutions	



Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.



Humiss International B.V.
Joop Geesinkweg 701, 1114 AB
Amsterdam-Duivendrecht, the Netherlands



Ningbo Albert Novosino Co.,Ltd.
No.1 Xinheng 3 Road Cicheng Town, Jiangbei
District, 315036 Ningbo, Zhejiang , China.
www.albertnovosino.com



Made in China

57 mm

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction