
TracWatch

User Manual

V1.0

1.SYMBOLS

2.Install and login


- (1) Download "TracWatch" APP on phone and install.
- (2) Turn on the Bluetooth , switch on the "TracWatch" , and enter the APP home page.

3.Connection

(1)Bluetooth Scanning


Enter the Bluetooth connection page, showing "No device

connected" , click "Bluetooth Scan" to search for the device.


(2) Bluetooth Connection

Click Bluetooth name in the bluetooth list to start the connection.


If disconnection is required, click "Disconnect" to

disconnect the APP from the device.

4. Measurement

Wearing the watch correctly, view the measurement data in real time on APP, including date and time, alcohol value, voltage, and contact value. Data is updated every 30 minutes and a data is recorded.


5. Battery and Charging

Battery : The device is powered by 3.7V 520mAh rechargeable Li-ion battery.

Charging :

- (1) Mate the metal contacts of the device with the metal probe on the charging station and connect the USB connector to the power supply to charge the device.
- (2) During charging, the red light on charging station is on.
- (3) When charging finished, the yellow light is on.

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF Exposure Information

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

ISED Statement

English: This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

(1) This device may not cause interference.

(2) This device must accept any interference, including interference that may cause undesired operation of the device.

The digital apparatus complies with Canadian CAN ICES-3 (B)/NMB-3(B).

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This device meets the exemption from the routine evaluation limits in section 2.5 of RSS 102 and compliance with RSS 102 RF exposure, users can obtain Canadian information on RF exposure and compliance.

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This equipment complies with Canada radiation exposure limits set forth for an uncontrolled environment.

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