

9. Do not use the Product in high-temperature, inflammable or damp place or place with electromagnetic radiation.

10. Do not dismount the Product without permission; ensure the Product free of crashing and water, for fear of malfunction.

>>Product scrap processing

Please dispose of this product according to the relevant urban environmental protection regulations, otherwise it may cause environmental pollution.

>>Product Specifications

Name: Eye Massager
Model: EYE 8
Power supply: 5 V === 1000mA
Lithium battery:3.7V === 700mAh
Input power: 5W
Size of packaging case: 200x91x85mm
Net weight: Around 190g

>> Scientific Knowledge: Prevent “Computer Vision Syndrome”

People will feel uncomfortable if using computer for a long time. Ophthalmologist and optometrist call it as “Computer Vision Syndrome (CVS)”. Details on CVS are as follows:

1.Eye fatigue will be aggravated if not blinking for a long time upon watching fluorescent screen.

This circumstance is often when reading lots of files, for computer user does not blink or blink as less as possible to avoid missing content. Blinking is critical for eyes protection. It offers position change opportunity to tense crystalline lens and iris muscle, so as to lubricate eyes and prevent eyes from generating itching and burning sensations and other symptoms.

2.The abnormal reaction “McCullough effect” will be generated after seeing green words in black fluorescent screen during working before fluorescent screen. Though harmless to body, this symptom can last for several days. It can make white words get pink slightly: A piece of paper pasted on white wall will look like the paper with pink edge.

3.This is related to the operator wearing bifocals. For them, the best position to see the target should be at or below the chest, with a distance of 12-18 inches.If height of fluorescent screen is not appropriate or the distance is too near, computer user will uncomfortably bend or have his or her head keeping a quite unnatural angle, which can cause dizziness, neck sore, shoulder sour or general malaise. According to expert, computer operating is a kind of work required for straining eyes, and it can reduce secretion of lubricant and enzyme in eyes. General speaking, there will be eyes dryness and fatigue, ghosting, blurred vision, head and neck soar and other symptoms if times of blinking is less than 5 and interval is relatively long.

>> Scientific Knowledge: Prevent “Computer Vision Syndrome”

With computer popularized, computer operating is an indispensable part in modern life almost. Many people must watch fluorescent screen of computer all day long for work need. Operating like this for a long time will cause visual fatigue easily and even induce some diseases, such as glaucoma. Therefore, those operating computer often should pay attention to eye care. Items below should be highlighted generally:

1.Avoid operating computer for a long time and pay attention to rest

Have a 5-10 min rest after operating computer for 1 h. Computer user can focus on a distant object or do eye exercises.

2.Keep excellent work posture

Keep the most appropriate work posture to make eyes look at fluorescent screen horizontally or look downwards the screen slightly. The optimal position is center of fluorescent screen is the part 20°below horizontal sight line. That position can relax neck muscles and minimize the area that eyes being expose to air.

3.Keep appropriate work distance
Eyes should be 60 cm at least far way from fluorescent screen of computer.

4.Create and keep excellent work condition
Surrounding light should be gentle; fluorescent screen of computer should have appropriate lightness and excellent clearness. Do not operate in dark environment, because black and white contrast can damage eyes. Heights of desk and chair should be matched with that of computer.

5.Go to hospital to see ophthalmologist if eye care message and rest does not work in case of dry or red eyes, burning or foreign body sensation, heavy eyelids, blurred vision and even eyeballs swelling pain or headache.

Warning:
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.
If this equipment does cause harmful interference to radio or television reception which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
--Reorient or relocate the receiving antenna.
-- Increase the separation between the equipment and receiver
--Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
--Consult the dealer or an experienced radio/TV technician for help
NOTE: This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

RF exposure:
The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

Shenzhen Pango Medical Electronics Co., Ltd
Main Site: No.25 1st Industry Zone, Fenghuang Road, Xikeng Village, Henggang Town, Longgang District, Shenzhen, Guangdong China
Tel: +86-0755-33825988
Fax: +86-0755-33825989
Web: www.pan-go.com



Eye Massager
EYE 8



Please read these instructions carefully prior to use!

>>>Product introduction:

Designed according to ergonomic and physical massage techniques, the EYE 8 touch type eye massager (hereinafter referred to as "the Product") massages multiple eye points and areas around the eye through the unique vibration amplitude and frequency of 16 vibration motors. Four massage modes can be repeated to relieve eye fatigue and relax the body and mind.

This product can be connected to the mobile phone through Bluetooth to listen to music, while massaging, making the body relax naturally.

The product adopts glasses shape, built-in rechargeable battery, integral folding design, stylish and beautiful appearance, comfortable to wear, easy to carry, anytime and anywhere to use, urban gold collar, white-collar workers, women and students stay up late eye care nurse.

>>> Features

- Massage Yuyao point, Cuanzhu point and other ocular acupuncture points high frequency with multiple motors;
- Built-in rechargeable battery: Can be used anytime and anywhere;
- Integrated and collapsible design: Able to contact ocular region more closely, convenient to carry;
- One-key operating: Simple;
- LED indicator;
- Listen to music through Bluetooth.


>>> Parts

1. One integrated mainframe (including the Eye Massager)
2. One USB cable
3. One operation introduction

>>>Product Picture



Front Elevation of Mainframe



1. Tightness adjusting buckle
2. Elastic band
3. On-off/mode key
4. TYPE-C USB charging interface
5. LED indicator

>>> Normal working conditions, transportation and storage conditions of products

1. The normal working conditions of the product
Ambient temperature: +5℃~+40℃;
Relative humidity: ≤80%;
Atmospheric pressure: 86kPa~106kPa.
2. Product transportation and storage conditions
Ambient temperature range: -20℃~+55℃;
Relative humidity range: ≤93%;
Atmospheric pressure range: 50 kPa~106 kPa.

>>> Operating Method

Preparation prior to using

1. Built-in rechargeable battery
Please use the USB cable to connect the USB power supply and the Eye Massager to charge the built-in battery. It is recommended to use after filling.
2. Using upon wearing the Eye Massager
First clean the inside of the Eye Massager on the host, then put the Eye Massager on your eyes, and adjust the elastic band to comfortable position.

Operating Introduction

- In the shutdown state, long press the power button to turn on the product.
- In the power-on state, long press the power button to turn off the product.

1. long press the power button to turn on the product, the three indicator light will all be on. Voice prompt: "Welcome to use the Eye Massager", "integrated mode"; Bluetooth on; Enable Bluetooth function; Remaining time display is 15 min.

2. After the product is powered on, the power button is the mode button. Press it briefly once to change a mode (all mode conversions have voice prompts):

- 1» Mode 1 (Default after power on) : comprehensive Mode (Multifrequency vibration) #1 Light is on, #2 Light is on, #3 Light is on
- 2» Mode 2: Eye Caring Mode (ring sensitive vibration) + music, #1 Light is on
- 3» Mode 3: Eye Protection Mode (rhythm vibration) + music, #2 Light is on
- 4» Mode 4: Sleep Mode (gentle soft vibration) + music, #3 Light is on

| |
|---|
| 1» Mode 1 (Default after power on): Integrated mode (Multifrequency vibration) #1 Light is on, #2 Light is on, #3 Light is on |
| 2» Mode 2: Love eye mode (ring sensitive vibration) + music, #1 Light is on |
| 3» Mode 3: Eye care mode (rhythm vibration) + music, #2 Light is on |
| 4» Mode 4: Sleep Mode (gentle soft vibration) + music, #3 Light is on |

3. The Bluetooth is turned on by default when it is turned on. After connecting with the mobile phone (Bluetooth name is "EYE 8"), the voice prompt "Bluetooth is connected".
4. After the product is powered on, double click the mode button to pause or start the music.

>>> Cleaning and Maintenance


1. For cleaning, wipe the Product with slightly wet soft cloth or towel. No wiping with volatile preparation.
2. Keep the Product in dry and airy place.
3. Do not expose it to direct sunlight or high-temperature or damp place.

Using Method for Rechargeable Battery

1. Low batter of in-built lithium battery if you hear buzzer ticks for three consecutive times while using your eye massager. At this time, you are required to use USB cable to charge battery. When charging, the charging indicator flashes. When fully charged, The charging indicator is always on, and it takes about 2-3 hours to fully charge.
2. The fully charged battery can be used for 4-5 days (15 min for each day).
3. Normally, service life of the chargeable battery is ≥500 times of charge and discharge, around 3 years.
4. The service life will be shortened if the chargeable battery is not be used for a long time. It is suggested to charge the battery every month at least.

Removal method of rechargeable battery

1. The built-in lithium battery of this product contains substances that may pollute the environment. Before discarding this product, the built-in lithium battery should be taken out.
2. Before disassembling the product and removing the lithium battery, be sure to turn off the power or cut off the power.
3. The detailed steps for removing the lithium battery are shown in the figure below. Disassemble the product in the order of ①—② and take out the lithium battery.
4. The taken-out lithium battery should be safely disposed of in a place specified by the government or a recycling unit.



>>> Precautions

1. Please read and understand the Introductions entirely before operating the Product.
2. Please clean eyeshade of the Product prior to using.
3. Children and others without expression ability are forbidden to use the Product.
4. Those who have sensitive skin or insensitive to heat are forbidden to use the Product.
5. Pregnant women, patients suffering heart disease, dysarteriotomy, malignant tumor, cranial vascular disease and other acute diseases and those being treated by doctor are forbidden to use the Product.
6. Stop using the Product immediately if you feel uncomfortable or your skin is abnormal during using and consult doctor.
7. No using the Product or using the Product under guidance of doctor if you are in eye trouble. Those ever provided with eye operation and patients suffering glaucoma, cataract, retinal detachment and other eye iseases are forbidden to use the Product.
8. Those wearing contact lenses cannot use the Product until taking out contact lenses.