



Please read this user manual thoroughly before using the product.

**1.Medical purpose**-The hearing aid intends to compensate for impaired hearing.



3.Side effects - For patients with sensitive skin, the silicone material may irritate your ear canals .
- First use may cause slight pain to the ear canal due to maladjestment.







# 

In-Ear tips	<b>→</b>		<b>→</b>	
N TO WEAD VO	LIDHEN	DINC AID		

Please choose the appropriate size earplugs, and assemble the ear tips as below picture.

-Hold the hearing aid with your index finger and thumb, with the little pore facing up

	Touch Panel LED Indicator Light In-Ear Tips Microphone Port	
Left Hearing Aid		Right Hearing Aid

Number	Picture	Function	Tone Notice	Operation
1		Turn on	Music cue	Open the changing case , wear the hearing aids on your ears, automatically enter hearing aid mode after 10 second
2		Turn off	1	Put the hearing aids back to charging case, close the case cover, then Hearing aids automatically turn off
3		Volume adjustment	Music cue	Press once a total of seven volume levels, which will be increased and cyclically adjusted by the volume buttons
4		Minimum volume	Di	* Press once
5		Maximum volume	DiDiDi	* Press once
6	•	Middle volume	DiDi	Press once the beeps of the middle five gears are the same
7		Manual shut down	Music cue	Press for 5 seconds
8		Manual power on	Music cue	Press for 2 seconds
9	Touch Planel LSD Indicator Light Microphiese Plan	Answer the phone	1	* Press once
10		Hang up the phone	1	Press for 2 seconds
11		Reject the call	1	Press for 2 seconds
12		Low battery in charging compartment		Blinking blue LED light
13	7	Hearing aid fully charged		Red LED light is off
14	Charging Indicator Light	Factory reset the charging case		Press and hold the reset button of the charging compartment for 5 seconds

Symptom	Possible Causes	Possible Causes
	Blocked tubing/earbud	Clean or replace ear wax guard as needed.
Not Loud Enough	Hearing change	Contact your hearing professional.
	Too much dust	Clean both microphone and receiver with brush .
Inconsistent Performance	Restart required	Put the hearing aid into charging case, close the top cover to turn off .Then take out the hearing aids ,press and hold each side for 2 seconds to restart.
	Blocked tubing/earbud	Clean or replace ear wax guard as needed.
Unclear, Distorted	Restart required	Put the hearing aid into charging case, close the top cover to turn off. Then take out the hearing aids, press and hold each side for 2 seconds to restart.
	Blocked tubing/ earbud	Clean or replace ear wax guard as needed.
	Defective hearing aids	Contact your hearing professional.
Dead Hearing aids	Charge required	Put the hearing aid into charging case, close the top cover and wait for charging .
	Blocked earbud/ tubing	Clean or replace ear wax guard as needed.
	Ear canal fulfilled with oil wax	Contact your hearing professional.
No LED indicator when hearing aids are inserted into the charging port	Low battery	Connect the Type-c cable to your charging case and plug it into the wall. The on-board battery LEDs will flash for a few seconds to indicate the connection to the power source.
Blinking blue LED light	Low battery	Insert your type-c cord to charge the case
Hearing aids whistling in the charging case	Dead battery	Put the hearing aid into charging case .close the top cover, and insert your type -c cord to charge the case.

	Project	Possible Causes	Solution		
		The button time is too short to start	Press and hold touch panel for 3 seconds to restart		
	Hearing aid has no sound	Low battery	Put the hearing aid into the charging case to charging		
		Blocked up the earplugs' sound port or microphone port	Clean the sound port or microphone port		
		Low battery	Put the hearing aid into the charging case to charging		
	The voice is very weak and unclear	The volume is too low	Double click to adjust the volume level		
		Earplugs loose	Re-wearing the hearing aid		
	Hearing aid failure self-checking table				

## TECHNICAL SPECIFICATION

Model	H006A
Output.OSPL90 Max.	≤105dB
FOG50 MAX Gain	≤33dB
EQ Input Noise	≤33dB
Frequency Range(Not narrower than)	200Hz-4500Hz W
Working Current	≤8mA/
Total Harmonic Distortion	<5%
Working Voltage	D.C 3.7V
Play Time	8hours

Notice: The above parameter index testing is on basis of FONIX8000Hearing Aid Test System; Inspection standard: IEC 60118-7 2005, 2CC coupling cavity HA-1

### HOW TO DOWNLOAD THE APP

## IOS app download methods:

1.Website: https://apps.apple.com/cn/app/widersonic/ id1577782600

2. Searching "Widersonic" on app store for IOS. The icon is as below:



3. Scan the QR code below:



### HOW TO DOWNLOAD THE APP

## Android app download methods:

### 1. Website:

https://play.google.com/store/apps/details? id=com.usopp.widersonic 2. Searching "Amplifier" in google play. The icon is as below:



3.Scan the QR code below:



# HOW TO DOWNLOAD THE APP

If you can't download the app in google play, please scan the QR code below use the brower app:




If you meet any problem, please contact us any time. Online Service or Customer Service Email: service@mirasing.com

### HOW TO USE THE APP


1 Connect the device. Open the app, press the search device, start scanning, power on you device with the battery, then the device can be scanned, press the device's name, and select the left ear or right ear, you will enter in the main page.





## HOW TO USE THE APP

2 Adjust the volume. Move the volume slider, you will get different volume, select comfortable volume.



## HOW TO USE THE APP

best sound quality

3 Adjust the program.

changing the gain for different frequence, select the comfortable sound. ■ Homepage 
 ■ Voice cor

4 Adjust the Equalizer band. Move the slider for



## **∆**1.CONTRAINDICATION

Press the square for changing the program. Select the mode close to the audiometric curve of yourself according to the icon or Select the comfortable program for you. If you utic value
itial atresia of external ear or deformity of external auditory meatus, a large number of foreignbodies
hal acoustic meatus;
ss is unstable, and the degree of hearing loss often fluctuates;
int with headache, dizziness, earache, linnitus, and other symptoms that ere not suitable for adaptat

### ▲ 3.NOTICE

in section are source anymerer, judy shocking and dropping the product. Do not keep or use it in high humidity environments ( such as the sauna room), not throw the scrapped battery as standard rubbish, not throw the scrapped battery as standard rubbish, not be the production of the standard rubbish.

ine MANOPACTORER. I setting up, using, or maintaining the hearing aid: To report unexpected operations or events.

use the product in a potential explosive or dangerous environment .

hare your hearing aids with others, Your hearing aids setting is personal and individual, share with someone sault in others' hearing damage.

with any nuasuls skir resertings. The careful with any unusual skin reactions. In some instances, wearing a hearing aid can cause an allergic reaction sading to it. Chiness or rashes. If you notice any uncommon skin condition, stop wearing the hearing aids, and check ith a medical specialist.

with a medical specialist.

Do not place the hearing aids near the magnetic objects.

Be sure to take off the hearing aids when getting a diagnosis or treatment.

Never pry the inside of the hearing aids with a needle or similar object; otherwise, it wil damage the device.

Do not put the hearing aid near heating appliances or in direct sunlight, and you should also avoid high moisture areas such as bathrooms or near a humidifier.

In a closed car in summer, the temperature can rise to very high levels. Do not leave your hearing aids in a high-temperature car.

Do not use the hearing aids close to a cell phone, otherwise it may increase the noise and affect the sound quality.

Please contact the MANUFACTURER:

Ny questions is earlier un-visited and affect the sound quality.

Not only are you likely to violate any applicable warranties or insurance, you may also cause further damage. 2. If your hearing aids and charger fail or perform poorly, check the guide on the next page for possible solutions.

guidance of a doctor or an experienced fitting technician of hearing aid.

2. Do not immerse the hearing aid in any liquid. Remove the hearing aid before swimming an After swimming or washing, wait for the ear canal to dry before wearing the hearing aid.

## charging time interval of hearing aid is 28 days;

Do not place hearing aids near high-temperature heating areas such as stoves.

Do not leave the hearing aid in a humid place for a long time.

Cleaning the earplug
 Remove ear wax plugged in the earpiece by using a stick and brush.
 Cleaning hook

Using a soft cloth or tissue, wipe of any moisture, sweat, or debris from the surface of the hook

3. Please use the function and volume control switch directly, open and close as lightly as possible to avoid

When you do not use hearing aids, please charge the product regularly to avoid battery damage; the re-

# 2. Care and maintenance

# 3.Cleaning

-Fully charge hearing aid batteries every night.

# Please put it in where children could not reach, dry and no exposed in the sun

-Do not take apart your hearing aids or insert the cleaning tools inside them.
(When not in use, place your hearing aids in the charger and store them.) - In a dry, safe place . Away from direct sunlight or heat to avoid extreme temperatures.

-Where you can easily find them.
-Safely out of reach from children and pets. Syrergy Charger Care Keep your charger clean. Heat, moisture, and breigh materials can result in poor per -Using the provided cleaning brush, keep the chrging ports clear of debris. -Do not use water, solvents, or cleaning fluids to clean the charging ports -Store your charger in a clean and dry location

1)We do not authorize any institution or individual to maintain and repair the product. If you suspect that the

1)If you Do not use it, turn off the hearing aid, and then put it into the storage case
2)Do not put the device under the sunshine, hign temperature, and moist environment get in touch with fire or is vulneable to vibration. 3)Keep the device storage in following environmental ranges:-20~55℃,95%R700-1060hPa



Please dispose of the device under the directive 2002/96/EC- WEEE(Waste Electrical and Electronic Equipment).and Don't dispose with domestic waste. If you have any queries, please refer to the local

Keep your hearing aids clean. Heat, moisture, and foreign substances can result in poor performance.

-Use a cleaning brush or soft cloth to clean debris from around the user control and microphone; inspect the the receiver, earbud, or ear tip and wax guard regularly - Never use water, solvents, cleaning fluids, or oil to cean your hearing aids. on additional maintenance procedures for your hearing aids, if needed

## that may affect hearing are identified and treated before the hearing aid is purchased.

l observation, or review of any other available information concerning the prospective user, that the prospect

(I) Visible congenital or traumatic deformity of the ear.

(ii)History of active drainage from the ear within the previous 90 days.

(iii)History of sudden or rapidly progressive hearing loss within the previous 90 days.

(iv)Auction ethornic diszinses.

(iv)Unliateral hearing loss of sudden or recent onset within the previous 90 days.

(iv)Audiometric air- bone gap eaual to or greater than 15 decibels at 500 hertz(Hz),1.000 Hz, and 2.000 Hz

(iv)Visible evidence of significant cerumen accumulation or a foreign body in the ear canal. (ivi) Pain or discomfor the ear Special care should be exercised in selecting and fitting hearing aid whose maximum sound pressure exceeds 132 decibels because there may be risk of impairina the remainina hearina of the hearing aid user.

(This provision is required only for those hearing aids with a maximum sound pressure capability greater hearing aids with a maximum sound pressure capability greater hearing aids with a maximum sound pressure capability greater hearing aids with a maximum sound pressure capability greater.

[]i Consult instructions for use B type applied part Manufacturer

## Date of manufacture SN Serial number

Marning/Danger: Improper use might cause serious injuries.

FCC STATEMENT:

Operation is subject to the following two conditions: (1) This device may not cause harmful interference,

(2) This device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: Reorient or

relocate the receiving antenna. Increase the separation between the equipment and receiver. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.