

Peally 3.8 User manual



LED color indicator

When the device is in the docking station:



Yellow

The device is searching for a WIFI connection



Red

The device encountered a connection failure. Disconnect the device from the docking station for 5 seconds. Try to find a new connection.

While charging and data upload:



Blue

Uploading data and charging in progress and are not completed yet.



Green

Uploading data and charging are complete. The device is ready to use. Be aware that the uploading process can take approx' 60 minutes

When device is attached to baby's shirt:



Yellow

To determine if the device is on, tap the device twice. A single yellow indicator led light will appear



Yellow

If you see a blinking yellow indicator led light, the device's battery is low, and the device must be charged



Red

When abnormal movement is detected

Product Manual

We wrote the User Manual in a way that allows you a simple and immediate use of our solution




<https://LittleOne.Care>

Peally 3.8 User manual




What's in your box?


1. Device




2. Ring + spare ring



3. Docking station



4. Cable & wall charger



Start here: Application Setup

1. Scan the QR code to download LittleOne.Care's App



2. Click on the App store icon and start downloading LittleOne.Care's App
3. Now, please follow the App's instructions

FCC compliance statement



FCC Compliance Statements

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

US Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End user must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. The portable device is designed to meet the requirements for exposure to radio waves established by the Federal Communications Commission (USA). These requirements set a SAR limit of 1.6 W/kg averaged over one gram of tissue. The highest SAR value is measured at a distance of 0mm between the radiator and the body.