


# I-SOL SMART INSOLE Product Manual

I-SOL Standard

## Product Configuration



## Preparation and confirmation

### 1. Prepare an APP

Scan the QR code below or search for 'I-SOL Care' in the store and install the app.

\* IOS app will be distributed sequentially



Supports Android 5.0 or higher


### 2. Product specifications

- Item name : I-SOL Standard
- Weight : Less than 45g
- Sensor : Pressure sensor, 3-axis accelerometer
- Communication method : Bluetooth LE
- Operating temperature : -10 to 50°C
- Battery : Replaceable (CR2032, 220mA)


## Product connection

### 1. How to install a battery


As shown in the figure below, install the battery with the + side facing up on the lower part of the smart insole, and then attach the battery cover by turning it clockwise.




- ① Make sure to insert the battery under the smart insole with the + part facing up.




- ② How the battery is mounted



- ③ Align the arrow on the battery cover with the direction of the 11 o'clock groove, and then press the entire frame to install it.



- ④ Insert a coin into the central groove and turn it clockwise to fix it.




- ⑤ The battery is successfully mounted and the cover is fixed

§15.19(a)(3) : All other devices shall bear the following statement in a conspicuous location on the device: This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions : (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

§15.21 : Information to user. The users manual or instruction manual for an intentional or unintentional radiator shall caution the user that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. In cases where the manual is provided only in a form other than paper, such as on a computer disk or over the Internet, the information required by this section may be included in the manual in that alternative form, provided the user can reasonably be expected to have the capability to access information in that form.

## 2. Connecting with a smartphone


- ① Turn on the location of the smartphone and Bluetooth.
- ② After logging in to the I-SOL app, register your member information.
- ③ Register the smart insole in Device Registration.
- ④ After placing the smart insole and smartphone close together, touch the shoe image in the center.
- ⑤ Shake the smart insole you want to register up and down or left and right according to the guide to recognize it.



- ⑥ Save the device after entering the device name.
- ⑦ If the connection is not recognized, move to an environment where there is no other Bluetooth nearby and try again.

## 3. Inserting a smart insole

- ① Check the size of the back of the smart insole to match the size of the shoe you want to wear, and cut it with scissors along the line.
- ② After removing the existing insole, put the smart insole connected to the smartphone to fit the shoe.



Now, walk towards the correct path.

## Precautions for use

- Do not arbitrarily modify the product.
- Do not disassemble the product.
- Do not place the product near fire.
- Parents need to pay attention that children do not swallow the battery and nor touch the product alone.
- For stable operation, keep the product away from water or moisture.
- It is only a light water-proof product; be careful not to submerge the product.
- Do not soak in water for washing, and remove foreign substances from the upper part with a wet towel.

## Warranty

- Product warranty period :  
1 year from the date of purchase
- Any faults caused by customer negligence are not warrantable.
  - If disassembled intentionally by the consumer
  - If damaged due to external impact
  - Inserting the battery or using the product in an abnormal way
  - In case of excessive use that exceeds the durability of the product
- A/S contact information (customer service center)
  - Address: 3rd floor, 12, Yanghyeon-ro 405beon-gil, Jungwon-gu, Seongnam-si, Gyeonggi-do
  - TEL : 031-602-7600

## Other notes

### Details of health responsibilities

Do not use the information provided by this product for diagnosis or treatment purposes.

### Certification

This device is registered for electromagnetic compatibility for home use and can be used in all areas.  
KC Certification No. MSIP-CRM-405-TLI-103N