

Bluetooth device settings

For most mobile phones

1. Enter the settings, find the Bluetooth menu, search for the device starting with "Flow-EZ...", pair this device until the pairing is successful.
2. Pixel 3a with Android OS version 12
 - Search directly after logging into your phone

Register and log in to the system

1. Turn on the measuring device
2. Click Flow-EZ to open the APP. After entering the system, you need to register an account for the first time, click "Register"

4:59



Above Care Inc.

Sign Up

Sign In

By using this application, I accept the Terms and
Condition and Privacy Policy.

3. Create a new account

- Enter a new account email and password
- Press Sign Up

5:00   



Sign In Flow-Ez

Please enter your email and password.





[Forget password](#)

New User? [Sign Up](#)

Flow-EZ does not diagnose, treat, mitigate or prevent any disease, symptom or abnormal physical state.

4. Login to the system

- Enter a new account email and password
- Press Sign In

5:00



Sign In Flow-Ez

Please enter your email and password.



andychen1981@gmail.com



.....

[Forget password](#)


Sign In

New User? [Sign Up](#)

Flow-EZ does not diagnose, treat, mitigate or prevent
any disease, symptom or abnormal physical state.

5. After entering the system, it will automatically find the measuring equipment

After finding the Bluetooth device, click to connect it



Click if the Bluetooth device
is not found the first time

Now start using the system

1. Fill in "Personal Information"

Flow-Ez

Usage ☒ Home ☐ Center

First Name
三

Last Name
張

Age < 30 Race Asian

Device
FB1E0D68-BF29-C2B7-8E3D-E3FF84A52397

Latest Stenosis < 50%

Shunt ☒ AV ☐ Neck

Position LeftUp


Gender ☒ Female ☐ Male ☐ Other

Weight ☐ pounds ☒ KG 65.0


Height ☐ feet inches ☒ meter 1.72

Home Recent profile Share

2. After entering the measurement screen, press the record



3. The result will be displayed after recording



Analysis
Outcome

Play the recording using
a Bluetooth earphone

Remark

- Align vessels with grooves
- Press harder for a clearer signal
- Hold your finger and device steady until the signal appears on the graph

LED Display

White - An indication of Bluetooth connected state

Blue - An indication of Bluetooth transmitting

Green - An indication of battery power

Warning

"For low-power radio-frequency devices that have obtained certification, companies, firms, or users are not allowed to change the frequency, increase power, or change the characteristics and functions of the original design without approval. The use of low-power radio-frequency devices must not affect flight safety and Interfering with legal communications; if any interference is found, it should be stopped immediately, and it can only be used when it is improved to no interference.

The aforementioned legal communication refers to radio communication operated in accordance with the provisions of the Telecommunications Management Act. Low-power radio frequency equipment must endure the interference of legal communication or industrial, scientific and medical radio wave radiation electrical equipment. "

Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- . Reorient or relocate the receiving antenna.
- . Increase the separation between the equipment and receiver.
- . Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- . Consult the dealer or an experienced radio/TV technician for help.

FCC Caution: To assure continued compliance, any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. (Example - use only shielded interface cables when connecting to computer or peripheral devices).

FCC Radiation Exposure Statement

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 0.5 centimeters between the radiator and your body.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

The antennas used for this transmitter must be installed to provide a separation distance of at least 0.5 cm from all persons and must not be co-located or operating in conjunction with any other antenna or transmitter.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.