


Thank you for purchasing Dolittle.
Please read the warnings before using this product.
*The content of this manual can change
due to further development of the product.




* The battery charge level of the product can be low when you first open the package. Please use it after charge.
※ Band, diaphragm, and diaphragm rim are consumables.

Specification	
Name	Dolittle
Material	PC(Polycarbonate)
Size(WxDxH)	47mm x 33mm x 14mm
Weight	13g
Bluetooth	Bluetooth 5.0 (min 4.0)
Operation	Room Temperature
Battery	Lithium 3.7v
Recommended OS	Android 7.0 / iOS 13.0

zentry

Scheduling and Sync



① Touch the 'Monitor' icon.
② Touch the 'clock' icon.

③ Set the measurement start time.
④ Set the measurement end time.
⑤ Reservation setting completed.

⑥ Save measurement data
after selecting a pet.

⑦ Check the data by touching
the 'Record' icon.

Product name




Device

Customized band


**① Band ring
(connected with the device)**

1. Charging



When charging starts, the LED light turns on.
When charging is complete, the LED light turns off.
(Approximately 30 minutes to full charge)

3. Connect and fasten the band



① Insert the band loop into the device slot.
② The attached band wraps around the chest of the pet.
③ Thread the end of the band into the buckle.
④ Adjust the length of the band and fasten the buckle.

4. On/Off



On: If you press the power button for more than 3 seconds, the device turns on while the LED light blinks.


Off: If the power button is pressed for more than 3 seconds while device is on, the device turns off.

* Turn off the device when not in use for a long time.

Solution

Can't measure the heart rate/respiration

1. Make sure the product is ON.
- If it is off, turn it on or charge it to use it.
2. Make sure Bluetooth is connected.
- Turn the Bluetooth on and connect with the app.
3. Check if the diaphragm is broken or deformed.
- If it is damaged or deformed, purchase it from our online mall and replace it.
4. Make sure the device near the heart of the pet.



- > Adjust the band to position the device near the animal's heart.
- 5. Make sure the device is in close contact with the pet's heart.
- Adjust the length of the band so that the device fits close to your pet's heart.


Can't connect to Bluetooth

1. Check the Bluetooth activation.
2. Make sure the product is ON.
3. Check if the product connected to another smart device.
- Briefly press the power button once.

Purpose of the product

Dolittle monitors your pet's heart rate and respiration rate.

Warning



Precautions for use

⚠ Warning

1. Use only for indoor pets, and do not use in damp or dusty places.
2. This product is not waterproof. Do not contaminate the product with water or other foreign matter.
3. Dolittle is for dogs and cats only. Do not use it on other animals or people.
4. This product is worn on the chest of a pet. Don't put it on your pet's head, legs or neck.
5. Do not overtighten the band so that the breathing or movement of the pet is uncomfortable.
6. Put the device on the pet when the fur is completely dry. Do not touch the product with wet hands.
7. Use when your pet is in a stable state or sleeping.
8. Do not allow pets to chew or eat this product.
9. Consult a veterinarian after measurement.
10. Measured data cannot be used as a substitute for medical care.
11. This product is applied to pets. People with limited physical, sensory or mental abilities or those with less experience and knowledge (including children) can use with a guardian for safety reasons.
12. When there is a problem with the product, stop using it and do not disassemble, repair, or modify it arbitrarily, but ask the manufacturer for an inspection.

⚠ Caution

1. Do not use the device while charging. Turn the device OFF when not in use for a long time.
(When storing, keep it out of direct sunlight.)
2. Do not use it when the room temperature is too low or too high.
3. Charge the device according to the instruction manual.
4. Position the device correctly at the measurement spot.
(See the user guide video)

Regulatory Compliance Information

FCC Caution

Any changes or modifications to the equipment not expressly approved by the party responsible for compliance could void user's authority to operate the equipment.

6. If the diaphragm is aged, damaged, or contaminated, please visit our website.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation

Industry Canada Statement:

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence.

L'exploitation est autorisée aux deux conditions suivantes :

(1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

Consumer Damage Compensation Criteria

1. This product is of excellent quality through a strict inspection.
2. The warranty period is 1 year from the date of purchase.
3. If the following situations occur within the warranty period, a refund will be made.

- Malfunctions caused by incorrect use, such as an operation that is not under the instruction manual.

- Failure or damage caused by falling or external impact after purchasing the product.

- In case of natural disasters such as fire, earthquake, flood, lightning, etc.

- In case of damage occurs due to use and storage in unsuitable environments such as high temperature and humidity, compounds, microorganisms, etc.

- Failure caused by ingress of sand, soil, etc. into the product

- In case of failure or damage due to other consumer negligence.

RF du IC d'exposition aux radiations: Cet équipement est conforme à l'exposition de la IC rayonnements RF limites établies pour un environnement non contrôlé. L'antenne pour ce transmetteur ne doit pas être même endroit avec d'autres émetteur sauf conformément à la IC procédures de produits Multi-émetteur.

Cet équipement doit être installé et utilisé avec une distance minimale de 20 cm entre le radiateur et votre corps.

Warranty

Thank you for using Zentry product.

Please fill out the warranty form and keep it.

Warranty Policy

1. Shipping costs for returns or exchanges without quality problems borne by the user. Shipping costs for re-delivery borne by the seller. Relevant invoices and promotion goods (if any) must be returned together.

2. If there is a problem with the product, you can receive an exchange service from the seller within 10 days from the date of purchase.

3. Guarantee service is available under consumer damage compensation standards.

If there is no warranty

1. Arbitrarily repair, misuse, collision, negligence, abuse, accident, alteration, torn or altered label, anti-counterfeiting markings.

2. Change, replacement of the product serial number and other similar cases.

3. If the product serial number is missing, or the product mode and the serial number of the warranty voucher and the actual product are different.

4. When the warranty period has expired.

5. If the appearance of the product is damaged due to abnormal use.

6. In the case of offerings of external packages and attached free goods, services or additional configurations, or free goods.

7. Damage to function due to improper storage or long-term use, or exposure to abnormal humidity or temperature.

8. In case of damage, abuse or abnormal use, arbitrary modification, or force majeure caused by not following the instruction manual.

Broadcasting communication equipment, etc.
electromagnetic compatibility test
220500101-EMC
220500102-EMC
Broadcasting communication equipment, etc.
(wireless) test
22050101-EL1
-R-Z7D1

Office hour
am 10:0~12:00
pm 14:0~17:00

If any problems occur,
please visit our website.

Office hour
am 10:0~12:00
pm 14:0~17:00

Hotline : 070-84-04-655
Website : www.zentry.k
Address : 1399-4, flc-3-10, Dong-gu, Seo-gu, Seoul, Republic of Korea

zentry

ver.2020.20