



TABLE OF CONTENTS

Introduction	4
Intended Use	6
Important Safety Information	6
Package Contents	9
Requirements	9
Using QardioCore For The First Time	10
How To Turn On/Off QardioCore	12
Detailed Instructions For The Correct Fitting Of QardioCore	12

Taking An ECG/EKG Recording	15
Checklist For A Correct And Accurate ECG/EKG Recording	15
Visualizing Your ECG/EKG Data	16
Important Facts About ECG/EKG And Self-Measurement	17
Activity Tracking	18
Recharging QardioCore Battery	19
Resetting The Pairing	20

Accuracy Testing And Maintenance	20
Contraindications	22
User Responsibility	22
General Use, Safety And Precautions, Cleaning	23
Customer Service Contact	26
Limited Warranty	26
Error Messages And Troubleshooting	28
Qardiocore Technical Specifications	30

Disposal	34
Certifications	34
FCC Statement	36
RF Statement	40



INTRODUCTION

Qardio offers a better way of tracking heart health that fits effortlessly into your life. Our devices are powerful and smart, have a beautiful design with a delightful user experience so you can use them anytime, anywhere.

QardioCore is a clinical-quality wearable electrocardiogram recorder. It records electrocardiogram (ECG/EKG) and physical activity.

Health conscious individuals or patients with known or suspected heart conditions can use QardioCore to record every day their ECG, physical activity, sport performance and medical symptoms during daytime or whenever they feel like, and share their data with their doctors. Medical professionals can use the QardioCore to quickly assess rate and rhythm, screen for arrhythmias, and remotely monitor and manage patients who use QardioCore. Athletes

can use QardioCore to precisely track their sport cardiovascular performance.

QardioCore is the smart way to track your electrocardiogram with contextual information about your physical activity. This device was developed in collaboration with physicians and clinical tests were conducted to prove its measurement accuracy.

With its ease of use and accuracy, QardioCore is ideal for monitoring your ECG trace throughout your day. With this system, you can record your ECG over extensive periods of time. This allows your doctor to assess ECG traces accompanying intermittent or infrequent symptoms not usually captured by conventional ECG and holter systems.

Please read through these instructions carefully so you understand all functions and safety information.

We want you to be happy with your QardioCore. If you have any questions, problems or suggestions, please contact Qardio's Customer Service at support. getqardio.com, or visit our website www.getqardio.com for more information.

QardioCore has features that are only available to those users who are under the care of a physician. These features are available to prescription users only and may not be available in all regions or all languages.

INTENDED USE

The QardioCore ECG monitor is intended to record single-channel electrocardiogram (ECG) traces of an adult individual. The QardioCore ECG monitor also records ECG rhythms and detects the presence of atrial fibrillation and normal sinus rhythm (when prescribed and used under the care of a physician). It may be useful as a diagnostic tool for patients with chest pain, palpitations, neurologic symptoms, shortness of breath, or to monitor response to treatment of cardiac disorders.

The QardioCore ECG monitor is intended for use by healthcare professionals, individuals with known or suspected heart conditions and health conscious individuals.

INDICATIONS

The clinical indications for mini Holter Event electrocardiography include: assessment of symptoms that may be related to arrhythmias assessment of risk in patients with or without symptoms of arrhythmia assessment of efficacy of antiarrhythmic therapy

QardioCore is intended to record, store and transfer [single-channel] electrocardiogram (ECG) rhythms. QardioCore also displays ECG rhythms and is indicated for the assessment of symptoms that may be related to arrhythmias and the assessment of risk in patients with or without symptoms of arrhythmia (when prescribed or used under the care of a physician). QardioCore is intended for use by healthcare professionals, individuals with known or suspected heart conditions and health conscious individuals. The device has not been tested for and it is not intended for pediatric use.

QardioCore has features that are only available to those users who are under the care of a physician. These features are available to prescription users only and may not be available in all regions or all languages.

IMPORTANT SAFETY INFORMATION

- Please read the User Manual carefully before using the QardioCore ECG monitor.
- If you have a pacemaker or other internal electronic device, consult your physician before using QardioCore.
- This device is not designed or intended for complete diagnosis of cardiac conditions. This device should never be used as a basis for starting or modifying treatment without independent confirmation by professional medical examination.
- Always consult your physician before you begin or modify any exercise program.
- This device does not detect or measure all heart rate, heart rhythm and heart waveform changes, especially those related to ischemic heart conditions.
- Do not attempt self-diagnosis or self-treatment

based on the recording results and analysis.

- Interpretations made by this device are potential findings, not a complete diagnosis of cardiac conditions. All interpretations should be reviewed by a medical professional.
- Users should notify their doctor of a possible change in health: a labeling of the ECG as normal should not be relied on as a guarantee of absence of arrhythmias or other health conditions. Interpretations made by this device are potential findings, not a complete diagnosis of cardiac conditions. All interpretations should be reviewed by a medical professional for clinical decisionmaking.
- The heart rate analysis is only valid if there is a valid rhythm (QRS complex visible).
- Do not use this device during an MRI scan.

- QardioCore relies on sensors that track your movement and other metrics. The data and information provided by these devices is intended to be a close estimation of your activity and metrics tracked, but may not be completely accurate, including step, distance, and calorie data.
- ECG reports viewed or printed at any magnification other than 100% may appear distorted and could lead to misdiagnosis.
- Qardio does not recommend using on humans less than 20kg.
- Do not store in extremely hot, cold, or wet conditions.
- Do not expose to strong electromagnetic fields.
- Do not attempt to use QardioCore while it is charging.

- Do not use if the electrodes are dirty: clean them first.
- Do use this device to record heart rate and heart rhythm only.
- Review carefully the Cautions and the General Use, Safety and Precautions indications enclosed in this User Manual.
- Self-diagnosis of measurement results and selftreatment are potentially dangerous. You should always consult your doctor.
- If you suffer from an irregular heartbeat, measurements taken with this device should be evaluated with your doctor.
- This device may only be used for the purposes described in this User Manual. The manufacturer cannot be held liable for damage or injury caused by incorrect use. Always follow the operating

- procedures described in this User Manual.
- Always follow the operating and storage conditions described in this User Manual to use and store your QardioCore.
- DO NOT use during magnetic resonance imaging (MRI), external defibrillation, or cautery procedures.
- Radio Frequency (RF) interference between this device and any existing RF transmitting or receiving equipment at the installation site, including electrosurgical equipment, in close proximity to the cardiograph should be evaluated before the equipment is operated as they may seriously degrade performance.
- This device is susceptible to interference from RF energy sources (lowered RF immunity) which exceed the IEC 60601-1-2 limits, such as power line bursts, other medical devices, and certain cellular products, information technology equipment and

radio/television transmission.

- Artifacts on the ECG caused by electromagnetic interference should be evaluated by a physician or physician authorized personnel to determine if it will negatively impact patient diagnosis or treatment.
- Like all electronic devices, this cardiograph is susceptible to electrostatic discharge (ESD).
 QardioCore should be removed before conducting AED or Cardioverter defibrillator. Electrostatic discharge typically occurs when electrostatic energy is transferred to the patient, the electrodes, or the cardiograph. ESD may result in ECG artifact that may appear as narrow spikes on the cardiograph display or on the printed report. When ESD occurs, the cardiograph's ECG interpretation may be inconsistent with the physician's interpretation.

- QardioCore has features that are only available to those users who are under the care of a physician.
 These features are available to prescription users only.
- Do not drop or bump with excessive force.
- Do not use in the presence of flammable anesthetics, drugs or pressurized oxygen.
- Use exclusively the charging cable and the straps provided.

PACKAGE CONTENTS

You will find in your QardioCore package:

- QardioCore device
- Three chest straps, to facilitate your everyday use
- Charging Cable
- Quick Guide
- User Manual

REQUIREMENTS

QardioCore requires a device with Bluetooth 4.0 (or later) and iOS 9.0 (or later) in order to use your QardioCore device, you:

- Must download the free Qardio App from the Apple App Store, or go to www.getqardio.com
- The first time you use QardioCore you require an internet connection and you must sign-up for a free Qardio account, via the Qardio App.

QardioCore works with iPhone, iPod, iPad and Apple Watch.

USING QARDIOCORE FOR THE FIRST TIME


- Download the free Qardio App: On your mobile phone or tablet go to www.getqardio.com and when prompted, download the app. Alternatively, go on the iTunes App.
- Open the Qardio App on your phone or tablet.
 If requested, you should enable Bluetooth on your device. You can enable Bluetooth under the Settings menu on your smartphone or tablet.
- Create a new user login, or login with your existing user name and password. Follow the on-screen instructions to register and set up your personal account.
- 4. Connect QardioCore to one of the chest straps provided. Always connect the right side first.
- 5. Adjust the length of the strap to adjust it for your chest size

- 6. Fit the QardioCore to your chest. You should wear the heart rate monitor directly on your skin, just below your sternum. Wear your QardioCore with the Qardio logo to the left. All the electrodes should be touching your skin and QardioCore should be snug enough to stay in place during your movement.
- When you connect the strap to the left side of QardioCore you should see a green light blinking to indicate that QardioCore has switched on.
- 8. While wearing QardioCore with the Qardio App open, touch and hold your phone or tablet on your chest to perform the pairing of your QardioCore with your phone or tablet. When prompted, accept the pairing request.

- 9. After a few seconds, QardioCore will automatically start recording your electrocardiogram. Your ECG can be affected by the position of QardioCore on your chest and your physiological condition. It is very important that QardioCore is correctly placed. Please read the "Detailed instructions for the correct fitting of QardioCore" and the "Checklist For a Correct and Accurate ECG/EKG Recording" sections on the User Manual with particular care.
- 10. The ECG recording can be stopped at any time by detaching the left clip from QardioCore and removing it from the chest.



Adjust length of strap QardioCore should be snug enough to stay in place during your movement.



Insert the buckle into place Make sure QardioCore is fully touching your body and the strap fits comfortably.



For first time use only: With the Qardio App open hold your iOS device near your QardioCore to pair it.

HOW TO TURN ON/OFF QARDIOCORE

QardioCore turns itself on when you attach the chest strap to its left side and it turns off when you unclip the left side of the strap. When QardioCore is not being worn, if the strap is left attached to the device, QardioCore will remain in a low power mode to preserve battery life.

To check if your device is switching on, look for the short blink of a green light twice on the top of the device when connecting the chest strap to the left side of QardioCore. Always store QardioCore with the chest strap detached from the device.

DETAILED INSTRUCTIONS FOR THE CORRECT FITTING OF QARDIOCORE

Wearing your QardioCore with the right fit will keep you comfortable: the device should be snug but comfortable, and it should not move while running. You should wear the heart rate monitor directly on your skin, just below your sternum. An overly tight strap can cause skin irritation, while a strap that's too loose can cause rubbing and affect the performance of QardioCore.

- Connect QardioCore to one of the chest straps provided. Always connect the right side first.
- 2. Adjust the length of the strap to adjust it for your chest size.
- Wrap the QardioCore around your chest, and connect the strap to the left side of QardioCore. You should wear the QardioCore directly on your skin, just below your sternum. Wear your

QardioCore with the Qardio logo to the left. All the electrodes should be touching your skin and QardioCore should be snug enough to stay in place during your movement.

After you put on the QardioCore, it switches on automatically and starts recording the electrocardiogram after a few seconds. When you connect the strap to the left side of QardioCore, you should see a green light blinking twice to indicate that QardioCore has switched on.



NOTE: Make sure the care tag does not fold over.

If the ECG data is erratic or does not appear, you can try these tips:

- Clean and dry your chest before putting on the device.
- Avoid wearing sunscreen under the device.
- Wear the device below your sternum bone.
 The device should be snug but comfortable.

NOTE: You can try wearing the device lower on your chest.

NOTE: Rinse the device with fresh water after each workout.

Tips to avoid skin irritation, or if you have allergies or skin sensitivities.

Keeping your QardioCore and straps—as well as your skin—clean and dry will maximize comfort and prevent long-term damage to the device. This is especially important after workouts or exposure to liquids such as sweat, soap, sunscreen, and lotions that can cause skin irritations.

These irritants can make skin reactions more likely if found between QardioCore and your skin.

If you have known allergies or sensitivities to substances like metals or plastics, check the materials in QardioCore and its straps on support.getqardio. com.

If you experience redness, swelling, itchiness, or any other irritation, immediately remove QardioCore and consult your physician before you put QardioCore back on.

Ensure to allow your skin to breathe for a few hours a day by not wearing QardioCore.

Individuals with sensitive skin or in poor physical condition may experience skin irritation when wearing QardioCore. Such individuals should keep their QardioCore and straps particularly clean.

List of Materials

- QardioCore Cases: Acrylonitrile butadiene styrene (ABS) and Thermoplastic vulcanizates (TPV)
- 2. Connectors (Electrodes): Silver plated brass
- 3. Strap: 20% Nylon Wooly, 44% Spandex, and 36% Nylon
- 4. Buckles on the strap: Acrylonitrile butadiene styrene (ABS)
- 5. Rating Label on the strap: 100% Polyester
- 6. Logo labels on the strap: 100% Nylon.

TAKING AN ECG/EKG RECORDING

When you wear a fully charged QardioCore, it will automatically start recording your electrocardiogram after approximately one minute.

It is very important that QardioCore is correctly placed. Please read the "Detailed instructions for the correct fitting of QardioCore" and the "Checklist For a Correct and Accurate ECG/EKG Recording" sections on the User Manual with particular care.

The ECG recording can be stopped at any time by detaching the left clip from QardioCore and removing it from the chest.

CHECKLIST FOR A CORRECT AND ACCURATE ECG/EKG RECORDING

- Adjust the length of the strap to adjust it for your chest size.
- Wear your QardioCore with the Qardio logo to the left.
- Fit the QardioCore to your chest. When you
 connect the strap to the left side of QardioCore
 you should see a green light blinking twice to
 indicate that QardioCore has switched on.
- You should wear the QardioCore directly on your skin, just below your sternum.
- Ensure all the electrodes are touching your skin
- Ensure QardioCore is fitting snug enough to stay in place during your movement.
- You can stop the recording at any time by pressing the detaching QardioCore from your chest.

VISUALIZING YOUR ECG/EKG DATA

Press the History button on the QardioCore page to see your ECG trace and heart rate data.

The ECG trace is visualized by default at 25mm/s and 10mm/mV. You can pinch and zoom when viewing the recording.

The ECG visualization can also be set to 25mm/s, or 50mm/s, and 5mm/mV, 10mm/mV or 20mm/mV. These options are only available to health professionals.

NOTE: ECG recording viewed or printed at any magnification other than 100% may appear distorted and could lead to misdiagnosis.

NOTE: You will not be able to view your recordings or utilize any of OardioCore's analytical tools until you are under the care of a physician. To gain access to these features, you must create an account on the

Qardio App. The one minute of your first recording will then be automatically sent for a free analysis by a cardiologist. Once you have received the ECG analysis, you will have access to view that recording and subsequent recordings. This service is not intended to replace medical advice, please seek professional medical assistance if you are suffering from any medical problem.

IMPORTANT FACTS ABOUT ECG/EKG AND SELF-MEASUREMENT

What is an electrocardiogram?

An electrocardiogram – often abbreviated, as ECG or EKG – is a test that measures the electrical activity of the heartbeat. With each beat, heart an electrical impulse (or "wave") travels through the heart. This wave causes the muscle to squeeze and pump blood from the heart.

Why is it done?

An ECG gives two major kinds of information. First, by measuring time intervals on the ECG, a doctor can determine how long the electrical wave takes to pass through the heart. Finding out how long the wave takes to travel from one part of the heart to the next shows if the electrical activity is normal or slow, fast or irregular. Second, by measuring the amount of electrical activity passing through the heart muscle, a cardiologist may be able to find out if parts of the heart are too large or are overworked. During and

ECG, several sensors, called electrodes capture the electrical activity of the heart.

Wearing QardioCore during exercise

During exercise the body requires more oxygen. As the level of physical activity increases, the heart has to work harder to deliver more oxygen-rich blood to the exercising muscles, so the heart beats faster. By monitoring the electrical signals of the heart as it beats faster, it is often possible to see coronary problems that cannot be seen when the body is at rest. As the intensity of the exercise increases, your heart rate will increase. You should stop your physical exercise if you experience dizziness, fatigue, chest pain, or other symptoms. Normally you should maintain your heart rate below your target level (85% of a predicted maximum heart rate, based on your age, or lower depending on your any known medical conditions). You should consult your doctor about your target heart rate level.

Target heart rate

While using QardioCore, the Qardio App displays the current Heart rate and the target Heart rate during recording, in the overview area.

The target Heart rate can be adjusted at any time in the settings

You can either set the target Heart rate directly, or set it as a percentage of the predicted maximum Heart rate. By default, the target Heart rate value is set as 85% of the predicted maximum HR.

The predicted maximum Heart rate is calculated with the following formula: Predicted Maximum Heart Rate = (220 – user's age in years)

Heat rate detection method

- ECG signal is filtered with a set of filters to amplify QRS parts of the signal and detect QRS complex. Each detected QRS complex is tested for valid amplitude, duration and refractrory period from previous QRS complex.
- RR interval is calculated between two consecutive valid QRS complexes as time difference between its maximal amplitude peaks. Each detected RR intervals is tested for valid duration and heart rate is reported as reciprocal of each interval in minutes (bpm) if tests are passed.

ACTIVITY TRACKING

Turning On Activity Tracking

The activity tracking feature records your daily step count, distance traveled, intensity minutes, and calories burned for each recorded day. Your calories

burned includes your base metabolism plus activity calories.

Activity goal

QardioCore sets automatically a goal of 10,000 steps per day. You can change your daily steps goal in the Qardio App settings. The Qardio App shows your progress toward your daily steps goal as you move during the day.

Intensity Minutes

To improve your health, organizations such as the World Health Organization, the U.S. Centers for Disease Control and Prevention, the American Heart Association recommend at least 150 minutes per week of moderate intensity activity, such as brisk walking.

The Qardio App monitors your activity intensity and tracks your time spent moving in moderate to vigorous intensity activities. Move in at least 10

consecutive minutes of moderate to vigorous intensity activities to work toward achieving your weekly intensity minutes' goal.

RECHARGING QARDIOCORE BATTERY

First charge your QardioCore by connecting the charging cable to the QardioCore and connecting the cable to a USB power source. QardioCore will not work until the battery has enough power.

When you charge QardioCore, a small LED on the device will blink indicating the charge status:

- Blinking green light every two second indicates QardioCore is charging
- Solid green light indicate QardioCore is fully charged

Qardio recommends recharging your QardioCore battery daily.

RESETTING THE PAIRING

In order to reset the pairing, remove the buckle from the QardioCore and use a paper clip to press the button on the pinhole on the QardioCore. You should see a green light shining through.

If necessary, go into the Settings of your phone or tablet, select the QardioCore and select "Forget this device".



ACCURACY TESTING AND MAINTENANCE

QardioCore comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the "Technical Specifications" section.

If you can't fix the problem using the trouble shooting instructions, please contact Qardio customer service at support.getqardio.com.

We recommend the QardioCore be tested for accuracy every 2 years or after mechanical impact (e.g. being dropped). Please contact Qardio customer service at support.getqardio.com to arrange the test.

CONTRAINDICATIONS

There are no known contraindications for the recording of ECG tests with QardioCore.

USER RESPONSIBILITY

This product is designed to perform in conformity with the description thereof contained in this manual and accompanying documentation, when operated and maintained in compliance with the instructions provided. A defective product should not be used. Parts that are broken, plainly worn, distorted or contaminated should be replaced immediately. Should any repair or replacement become necessary, we recommend to contact Qardio Customer Support. The user of the product shall have the sole responsibility for any malfunction, which results from improper use, faulty maintenance, damage or alteration.

GENERAL USE, SAFETY AND PRECAUTIONS, CLEANING

- If you have a pacemaker or other internal electronic device, consult your physician before using the QardioCore.
- Always consult your physician before beginning or modifying any exercise program.
- QardioCore is not intended to diagnose, cure, or prevent any disease or medical condition.
- QardioCore is not intended to be a replacement for a 12- lead ECG.
- Do not wear QardioCore while recharging.
 QardioCore recharging connector has been purposely designed to prevent user wearing QardioCore while its battery is being recharged.
 Do not perform any modifications to QardioCore or its charging cable to allow wearing QardioCore while recharging.

- Do not use QardioCore for any purpose other than as specified in this User Manual.
- Do not leave QardioCore unattended around children or persons who cannot express their consent to use.
- Do not bend QardioCore.
- Do not apply strong shocks and vibrations to QardioCore, as this may result in damage to the device.
- Do not put QardioCore electrodes in contact with any other metal objects.
- Do not drop QardioCore.
- Do not expose QardioCore to temperatures outside the storage or operating range.
- Do not expose QardioCore to direct sunlight for

extensive periods of time.

- Do not disassemble, modify, remanufacture, puncture or damage QardioCore.
- Do not immerse or expose QardioCore to liquids, fire, explosion, or other hazard.
- Do not use a power cable that is not approved or supplied by Qardio.
- QardioCore has features that are only available to those users who are under the care of a physician.
 These features are available to prescription users only and may not be available in all regions or all languages.
- QardioCore should not be used in case of emergency.
- Certain services may give you ability to request professional clinical interpretation and analysis

of your ECG recordings. Your location may restrict your ability to use this service as certain telemedicine restrictions may apply to your area. Qardio does not know your location and it is your responsibility to ensure this service is legal according to your local telemedicine laws.

QardioCore dust and water resistance

Your QardioCore is dust resistant and splash and water resistant, but not waterproof (QardioCore has a water resistance rating of IP65 under IEC standard 60529). For example, you may wear and use your QardioCore during exercise (exposure to sweat is OK), in the rain.

Submerging your QardioCore isn't recommended. Water resistance isn't a permanent condition, and your QardioCore can't be rechecked or resealed for

water resistance. The following may affect the water resistance of your QardioCore and should be avoided:

- Dropping or bending your QardioCore or subjecting it to other impacts.
- Submerging your QardioCore in water for long periods of time.
- Swimming or bathing with your QardioCore.
- Exposing your QardioCore to pressurized water or high velocity water, for example, water skiing, wake boarding, surfing, jet skiing, and so on.
- Exposing your QardioCore to temperatures outside of its operating or storage range
- Wearing your QardioCore in a sauna or steam room.

Cleaning

- Remember that you wear your QardioCore
 in contact with your skin, just like a clothing
 item. Because of this, you should always keep
 your QardioCore clean. Even if the surface of
 QardioCore appears clean, sweat, soap, sunscreen,
 and lotions on the surface or in crevasses can
 cause skin irritation.
- Keep your QardioCore clean and dry. a soft toothbrush or a nonabrasive, lint-free cloth to scrub it with a weak solution of fresh water and a mild neutral detergent.
- The bacteria and odor resistant QardioCore straps protect against odor generated by the formation of bacteria from sweat, which ensures good comfort and hygiene. In order to ensure maximum bacteria and odor resistance, keep your QardioCore and the straps clean.

- Cleaning products, and abrasive materials shouldn't be used on your QardioCore.
 Compressed air and external heat sources like hair dryers can cause damage.
- Direct sunlight, high temperatures, and humid conditions can cause damage.

To clean QardioCore, follow these steps:

- Take off your QardioCore and keep the charging cable unplugged.
- Clean your QardioCore with a soft toothbrush or a nonabrasive, lint-free cloth to scrub it with a weak solution of fresh water and a mild neutral detergent.
- 3. Dry your QardioCore with a recommended cloth.
- 4. To clean the QardioCore strap, wash it by hand or

washing machine with a mild soap and cold water. Let the strap air dry.

CUSTOMER SERVICE CONTACT

Qardio customer service contact is available at support.getqardio.com.

LIMITED WARRANTY

This device is covered by a three-year limited warranty from the date of purchase. The strap has a functional warranty (elasticity) for two years or 365 sessions of use (whichever comes sooner), while battery and other wearing parts are not covered by the limited warranty.

The limited warranty is valid only on presentation of the purchase receipt confirming date of purchase. Opening or altering the device invalidates the limited warranty.

The guarantee does not cover damage caused by improper handling, discharged batteries, accidents or non-compliance with the operating instructions and normal wearing of parts.

If a defect arises during the warranty period, Qardio, at its option and to the extent permitted by law will (1) repair the product at no charge, using new parts or parts that are equivalent to new in performance and

reliability, (2) exchange the product with a functionally equivalent product that is new or equivalent to new in performance and reliability, or (3) refund the original purchase price. This warranty excludes damage resulting from abuse, accident, modifications or other causes that are not defects in materials and workmanship.

Other than the consumer law rights to which you are entitled, all warranties, conditions and other terms not set out in this warranty document are excluded from the limited warranty. Some countries do not allow limitations on how long such warranties, conditions and/or implied terms may last, so the limitations described above may not apply to you.

In no event shall Qardio be liable for (a) any losses that were not caused by our breach of this limited warranty; (b) any incidental, special, indirect or consequential damages, whether resulting from

use, misuse, or inability to use this product or from defects in the product, (c) losses relating to any business of yours, loss of profits, loss of data or loss of opportunity.

The provisions of this limited warranty shall not apply to any other liability, except those that cannot be limited or excluded as a matter of law. Depending on where you live, some of the above limitations or exclusions may not apply to you. To obtain warranty service, contact Qardio at www.getgardio.com.

Online Auction Purchases: Products purchased through online auctions are not eligible for rebates or other special offers from Qardio warranty coverage. Online auction confirmations are not accepted for warranty verification. To obtain warranty service, an original or copy of the sales receipt from the original retailer is required. Qardio will not replace missing components from any package purchased through an online auction.

Australian Purchases: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. The benefits under our Limited Warranty are in addition to other rights and remedies under applicable law in relation to the products.

ERROR MESSAGES AND TROUBLESHOOTING

Lead quality problems

The following table describes different lead quality problems and possible actions to try to resolve this:

Condition	Possible Causes	Actions
Bad signal	Poor electrode contact with skin	Make sure QardioCore is fitting tightly and that all electrodes are touching the skin
Wandering baseline (an upward and downward fluctuation of the waveforms)	 Electrodes that are dirty, loose, or positioned on a bony area. Oily skin or body lotions Rising and falling of chest during rapid or apprehensive breathing. 	 Clean skin with isopropyl alcohol or acetone Reposition QardioCore higher (or lower) on the chest User should relax If wandering baseline persists, turn the baseline filter on.
Condition	Possible Causes	Actions

Muscle tremor interference (random irregular voltage superimposed on the waveforms). May resemble or coincide with AC interference

- User is uncomfortable, tense, nervous.
- User is cold and shivering.
- QardioCore strap is too tight.

- User should get comfortable.
- Ensure that all electrodes are touching the skin.
- If interference still persists, the problem is probably electrical in nature. See the following suggestions for reducing AC interference.

AC interference (even-peaked, regular voltage superimposed on the waveforms). May resemble or coincide with muscle tremor interference

- Electrodes that are dirty, loose, or positioned on a bony area.
- Patient touching any metal objects.
- Electrical devices in the immediate area, lighting, concealed wiring in walls or floors.
- Improperly grounded electrical outlet nearby.
- Incorrect AC filter frequency setting or AC filter is turned off.

- Ensure that all electrodes are touching the skin.
- Verify that the user is not touching any metal.
- Verify that the proper AC filter is selected.
- Try moving to another room.
- If possible, unplug electrical devices in the immediate area.
- If interference still persists, the noise may be caused by other equipment in the room or by poorly grounded power lines.

QARDIOCORE TECHNICAL SPECIFICATIONS

Chest Size	26~43 in (66~109 cm)
Memory	Practically unlimited, when connected to smartphone or tablet. Up to 48 hours data storage when not connected to smartphone or tablet, depending on usage conditions.
ECG Data Compression	No data compression. All data are stored and transmitted in raw format.
ECG Recording type	Continuous
ECG Channel	Single channel
Input Dynamic range	50mV Peak-to-Peak
DC Dynamic span	±400mV
Gain accuracy	5%

ECG Amplitude resolution	0.8 μV
ECG Signal bandwidth	0.05 to 40 Hz
ECG External Sampling rate	external: 200 samples per second
ECG Sampling resolution	16 bit
ECG Common Mode Rejection	60dB
ECG Input Impedance	>100M Ω
ECG Calibration	Automatic
Movement sensor	3-axis accelerometer
Power source	The device is powered by a built-in, 3.7V lithium-ion battery that you can charge using the charging cable supplied in the product package

Battery life	Continual Working Time up to 24 hours of battery life, when fully charged.
Battery charge time	About 3 hours to 100%. Charge time depends on environmental factors; actual results will vary.
Water resistance	Water resistance rating of IP65 under IEC standard 60529
Weight (including strap)	130 gram including the battery
Dimensions	185 x 87 x 9 mm
Operating conditions	-20C to 45C temperature for discharging and 0C to 45C temperature for Charging, 25% to 90% (non-condensing) relative maximum humidity, atmospheric pressure 86~106kpa, maximum altitude: 2000m.
Storage conditions	-20C to 45C temperature, 45% to 85% (non condensing) relative maximum humidity, atmospheric pressure 86~106kpa, maximum altitude: 2000m.
Shelf life	Estimated 2 years.

Specifications and features are subject to change without prior notice or any obligation for the manufacturer, and may not be available in all regions or all languages. Certain features may require purchase of separate services.

DISPOSAL

Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal. The symbol applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste.

At the end of devices useful life, the user must deliver it to the able collecting centers for electric and electronic garbage, or give back to the retailer when purchasing a new device. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the environment and health. In case of abusive disposal of device by the user, will be applied administrative endorsements in

compliance with current standard. The device and its parts is made with regard to disposal, as appropriate, in accordance with national or regional regulations.

CERTIFICATIONS

This device complies with the following normative documents:

- EN ISO 13485:2003 /AC: 2009: Medical devices -Quality management stems – Requirements for regulatory purposes (ISO 13485:2003) Reference to standards contd.
- 2. IEC/EN 60601-1-11:2010 General requirements for basic safety and essential performance –Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- 3. FCC part B 15B:2013 Electromagnetic Compatibility
- 4. FCC Rule Part: 15.247 Cat: DSS (Bluetooth) FCC

- Rule Part: 15.247 Cat: DTS (BT4.0)
- EN ISO 10993-1:2009/AC:2010 Biological evaluation of medical devices-Part 1: Evaluation and testing within a risk management process (ISO 10993-1:2009)
- EN ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity (ISO 10993-5:2009)
- EN ISO 10993-10:2009 Biological evaluation of medical devices -- Part 10: Tests for irritation and skin sensitization
- EN 55011 Group 1 Class B:2009+A1:2010:
 Industrial, scientific and medical equipment
 Radio-frequency disturbance characteristics Limits and methods of measurement
- IEC60601-1-2: Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances

- 10. IEC 60601-2-47 Ambulatory electrocardiographic monitors
- 11. EN 300 328 V1.9.1: 2015 Electromagnetic compatibility and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the R&TTE Directive
- 12. EN 301 489-1 V1.9.2 (2011-09) Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements
- 13. EN301489-17 V2.2.1:2012 Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment; Part 17: Specific conditions for Broadband Data Transmission Systems
- 14. IEC 62133:2012 RLV: Secondary cells and batteries

containing alkaline or other non-acid electrolytes -Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications

- 15. UN38.3, Fifth Edition: Recommendations on transport of dangerous goods, manual of test and criteria, Section 38.3 – Lithium metal and lithium ion batteries
- EN 62366: Medical devices. Application of usability engineering to medical devices
- 17. IEC 60601-1-6: Medical electrical equipment -Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
- 18. IEC 62304:2006 Medical device software Software life cycle process

FCC STATEMENT

Federal Communications Commission (FCC)

Statement 15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be

determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna. -Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

A. This device complies with Part 15 of the FCC Rules/ Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference and
- 2) This device must accept any interference received, including interference that may cause undesired operation of the device.

B. This device and its antenna(s) must not be co-located or operating in conjunction with any other antenna or transmitter

C. Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user authority to operate this equipment

IMPORTANT NOTE (for portable device configuration):

Federal Communication Commission (FCC) Radation Exposure Statement. This EUT is in compliance with SAR for general population/uncontrolled exposure limits in ANSI/IEEE C95.1-1999 and has been tested in accordance with the measurement methods and procedures specified in OET Bulletin 65 Supplement C.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de license. L'explotation est autorisée aux deux conditions suivantes:

1) il ne doit pas produire de brouillage et

2) L'utilizateur du dispositif doit être prêt à recepter tout brouillage radioélectrique reçu, même si se brouillage est susceptible de compromettre le fonctionnement du dispositif.

ICES-003.

This Class B digital apparatus complies with Canadian ICES-003.

Cet appareil numérique de la classe B est conforme à la norme NMB-003 du Canada.

FCC RF Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions

for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

IC Radiation Exposure Statement / IC Déclaration sur la radioexposition.

This EUT is in compliance with SAR for general population-uncontrolled exposure limits in IC RSS-102 and has been tested in accordance with the measurement methods and procedures specified in IEEE 1528. This equipment should be installed and operated with minimum distance of 1.5cm between the ratiator and your body.

Cet appareil est conforme avec SAR pour la population générale/limites d'exposition abusive IC RSS-102 et a été testé en conformité avec les méthodes et procédures spécifiées dans la norme IEEE 1528 mesure. Cet équipement doit être installé et utilisé à une distance minimale de 1,5cm entre le radiateur et votre corps. La séparation de test SAR de

la distance de 10mm pour hotspot.

Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

En vertu de la réglementation de l'Industrie du Canada, cet émetteur de radio ne peuvent fonctionner en utilisant une antenne d'un type et maximum (ou moins) gain approuvé pour l'émetteur par Industrie du Canada. Pour réduire le risque de brouillage aux autres utilisateurs, le type d'antenne et son gain doivent être choisis de sorte que la puissance isotrope rayonnée équivalente (p.i.r.e.) ne dépasse pas ce qui est nécessaire pour la réussite de comunication.

RF STATEMENT

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following section.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. The use of accessories and cables other than those specified may result in increased emissions or decreased immunity.

The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

The device is suitable for use in all establishments,

including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Any other accessories, transducers and cables may result in increased emissions or decreased immunity and EMC performance.

The device should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, it should be observed in order to verify normal operation in the configuration in which it will

be used.

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following. Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

Information in this document is subject to change without notice. All changes will be in compliance with regulations governing manufacture of medical equipment. Qardio reserves the right to change or improve its products and to make changes in the content of its User Manuals without obligation to notify any person or organization of any such changes or improvements.

Visit the Qardio web site (www.getqardio.com) for current updates and supplemental information concerning the use and operation of this and other Qardio products.

Guidance and manufacturer's declaration-electromagnetic emissions

The QardioCore Wireless ECG Monitor system is intended for use in the electromagnetic environment specified below. The customer or the user of the QardioCore Wireless ECG Monitor system should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance	
CE emissions CISPR11	Group 1	The QardioCore Wireless ECG Monitor uses RF energy only for its internal function.	
RE emissions CISPR11	Class B	Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
Harmonic emissions IEC 61000-3-2	Not applicable	The QardioCore Wireless ECG Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Voltage fluctuations/			
Flicker emissions	Not applicable		
IEC 61000-3-3			

Designed by and manufactured for Qardio, Inc. California, USA. FOR US AND INTERNATIONAL

www.getqardio.com



Type BF Applied Part (cuff)



YA HORNG ELECTRONIC CO., LTD.
Tainan, Taiwan
Factory: ATTEN ELECTRONIC (DONGGUAN) CO., LTD.

188 Industrial District, Ping Shan Administrative District, Tang Shia Town, Dongguan, 190, CN, 518055



FCC ID: 2ABF2-888CORE IC: 11885A-888CORE



2016



Read this manual before use.

US Importer

Qardio, Inc. 340 S Lemon Ave #1104F, Walnut, California 91789, USA.



WEEE

C€0434

EC REP

Kahl Handelsvertretung Add.: Isarstr.33 40699 Erkrath, Germany

Designed by and manufactured for Qardio, Inc. California, USA. FOR CANADA

www.getqardio.com



Type BF Applied Part (cuff)



Factory: ATTEN ELECTRONIC (DONGGUAN) CO., LTD. 188 Industrial District, Ping Shan Administrative District, Tang Shia Town, Dongguan, 190, CN, 518055



FCC ID: 2ABF2-888CORE IC: 11885A-888CORE



2016



Read this manual before use.

US Importer

Qardio, Inc. 340 S Lemon Ave #1104F, Walnut, California 91789, USA.



WEEE

C E 0434

EC REP

Kahl Handelsvertretung Add.: Isarstr.33 40699 Erkrath, Germany

Declaration – electromagnetic emissions and immunity for equipment and systems that are not life-supporting and are specified for use only in a shielded location

The QardioCore Wireless ECG Monitor system declaration-electromagnetic immunity

The QardioCore Wireless ECG Monitor system is intended for use in the electromagnetic environment specified below. The customer or the user of the QardioCore Wireless ECG Monitor system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol ((**)*)
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 385 MHz 5 V/m 450 MHz 10 V/m 810 to 930 MHz 20 V/m	

Declaration – electromagnetic immunity

The QardioCore Wireless ECG Monitor system is intended for use in the electromagnetic environment specified below. The customer or the user of the QardioCore Wireless ECG Monitor system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	8 kV contact 15 kV air	8 kV contact 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	-5 % UT(95 % dip in UT) for 0.5 cycle -40 % UT(60 % dip in UT) for 5 cycles -70 % UT(30 % dip in UT) for 25 cycles -5 % UT(95 % dip in UT) for 5 sec	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Qardio, QardioArm, QardioBase, QardioCore, QardioMD and their respective logos are either registered trademarks or trademarks of Qardio. Inc. in the United States and/or other countries.

iPhone, iPad, iPod are trademarks of Apple Inc., registered in the U.S. and other countries.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Qardio, Inc. is under license. Other trademarks and trade names are those of their respective owners.

Qardio products should be only used in conjunction with professional medical advice for diagnosis or treatment, and not as a substitute or replacement of it. Patents assigned and pending.

Copyright © 2016 Qardio, Inc. All rights reserved. Version 0.9.0.1 - Eng