PowerDot PD-02 Muscle Stimulator

Operator's Instructions For Use

# INTENDED USE

PowerDot PD-02 uses electrical muscle stimulation (EMS) (also known as neuromuscular electrical stimulation or NMES) and transcutaneous electrical nerve stimulation (TENS) technologies to stimulate your muscles and nerves for therapeutic purposes.

For your convenience, you can operate PowerDot PD-02 wirelessly using PowerDot Doctor mobile application running from your iPad. Additionally, you can also prescribe/schedule certain stimulation programs for your patients to be executed by them in their home environment using PowerDot Patient mobile application from their iOS or Android mobile devices.

## Indications for Use

PowerDot PD-02 stimulator is a prescription device which is intended to be used following the directions of a healthcare provider. The device can be either used by the therapist in healthcare facility setting (when operated from PowerDot Doctor mobile application) or by patient/lay operator in a home environment (when operated from PowerDot Patient mobile application).

PowerDot PD-02 has the following indications for use:

### **NMES**

- Retarding or preventing disuse atrophy
- Maintaining or increasing range of motion
- Re-educating muscles
- Relaxation of muscle spasms
- Increasing local blood circulation
- Prevention of venous thrombosis of the calf muscles immediately after surgery
- Stimulation of healthy muscles in order to improve or facilitate muscle performance

# **TENS**

- Symptomatic relief and management of chronic, intractable pain
- Post-surgical and post-trauma acute pain
- Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities

Caution: Federal law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

### **Safety Guidelines**

#### 1. Counter-indications

*Do NOT use PowerDot PD-02 for the patients with the following conditions:* 

- Cardiac stimulator (pacemaker), implanted defibrillator or other implanted electronic device. Such use could cause electric shock, burns, electrical interference or death.
- Body-worn electro mechanical medical devices, i.e. insulin pump.
- Serious arterial circulation problems in lower limbs.
- Abdominal or inguinal hernia
- · Symptomatic local pain relief unless etiology is established or pain syndrome has been diagnosed
- Cardiac arrhythmia (do not use on chest)

#### 2. Additional Precautions

- Use caution for patients with suspected or diagnosed heart problems
- Use caution for patients with suspected or diagnosed epilepsy
- Use caution in the presence of the following:
  - when there is a tendency to hemorrhage following acute trauma or fracture
  - following recent surgical procedure when muscle contraction may disrupt the healing process
  - over a menstruating or pregnant uterus
  - over the skin areas which lack normal sensation
- Depending on the judgement of the responsible physician, the device may only be applied under supervision and with the parameters defined by the responsible physician. Otherwise the exercise may be too strenuous for the patients with:
  - hypertension (> stage 2), ischemic heart disease and cerebrovascular diseases
  - cardiovascular diseases
  - pregnancy. Safety of powered muscle stimulators for use during pregnancy has not been established.
  - under 16 years of age
- Proceed with caution after recent surgery
- Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium, or alternate electrode placement

# 3. Adverse Effects

- Skin irritation and burns beneath the electrodes have been reported with the use of powered muscle stimulators.
- Potential adverse effects with TENS are skin irritation and electrode burns.

# 4. General Safety Measures

- Read all instructions for operation before treating a patient.
- The long-term effects of chronic electrical stimulation are unknown.
- TENS is not effective for pain of central origin (including headache).

- TENS waveforms have no curative value.
- TENS is a symptomatic treatment, and as such, suppresses the sensation of pain which would otherwise serve as a
  protective mechanism.
- Stop treatment immediately if patient experiences discomfort or pain.
- The choice of the therapy parameters to program, of the stimulation programs and electrode placement is restricted to the responsible physician or therapist. It is the physician's or therapist's decision whether or not to use the device on a specific patient.
- Do not apply stimulation if you patient has progressive cancer or near any cancerous tumour. The increased metabolism, caused by certain modes of stimulation, is likely to encourage cancer cells to spread.
- Proceed with caution if stimulation is applied to areas of the skin whose level of sensation is lower than normal.
- Do not apply stimulation to a person who cannot express themselves.
- Stimulators should be kept out of the reach of children

### 5. Device Operating Safety Measures

- Do not charge PowerDot PD-02 device while the device remains placed on/attached to any part of the patient's body.
- The patient must be familiar with the ways to terminate their stimulation session, if needed. Patients unable to operate the emergency stop function (either by pressing button on the device or by unplugging lead cables from the device), e.g. paralytic patients, must never be left unattended during therapy.
- During the muscular contraction phase it is recommended to hold the extremities of the stimulated limbs to avoid any shortening of the muscle during contraction, which could cause cramps.
- Never begin an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which produces a feeling weakness and a tendency towards fainting. If this does occur, all that is required is to stop the stimulation and for the person to lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Do not wind cables around the neck. Tangled cables can cause strangulation.
- Do not apply stimulation during sleep.
- Do not use PowerDot PD-02 while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the patient at undue risk of injury.

# 6. Co-existence and Environmental Safety Measures

- Do not use PowerDot PD-02 in the areas where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents and disinfectants as well as in the oxygen-rich environments.
- Do not use device in water or in a humid atmosphere (sauna, bath, shower etc.) that would cause electrical failure.
- Do not use PowerDot PD-02 at altitudes of over 9,800 feet.

- Do not use PowerDot PD-02 if patient is connected to a high-frequency surgical instrument as this could cause skin
  irritations or burns under the electrodes.
- Do not use PowerDot PD-02 in the close proximity of medical devices such as MRI, CT, diathermy, X-Ray as well as RFID equipment (e.g. electromagnetic security systems) as those could alter the current generated by PowerDot PD-02, cause tissue damage, and can result in severe injury or death. If you have any doubts when using PowerDot PD-02 device in close proximity to another medical device, please contact the device manufacturer.
- Do not use the PowerDot PD-02 within 3 feet of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the stimulator in close proximity to another medical device, seek advice from the manufacturer.
- Do not apply stimulation with PowerDot PD-02 near the area of an implant, such as cochlear implants, pacemakers, skeletal anchorage or electric implants. This could cause an electrical shock, burns, electrical interference or death.
- Do not apply stimulation with PowerDot PD-02 close to metal. Remove jewelry, piercings, belt buckles or any other metallic product or device in the area of stimulation.
- Disconnect the electrode pads from the device before using electrosurgical equipment, or a defibrillator, to avoid cutaneous burns from the electrodes and destroying the device.
- Apply caution when using PowerDot PD-02 near electronic surveillance equipment (e.g. cardiac monitors, ECG alarms), as there is a risk they may not work properly whilst the electrical stimulation device is being used.
- Do not use PowerDot PD-02 in areas in which unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with the device.

### 7. Electrode Pads Placements Safety

# NEVER attach electrodes:

- Transcerebrally, or on the eyes
- On the front and sides of the patient's neck or mouth because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure
- Stimulation should not be applied over the carotid sinus nerve particularly in patients with a known sensitivity to the carotid sinus reflex.
- Counter-laterally, i.e. one stimulator should not be applied on opposite sides of the body.
- Transthoratically, to the patient's front torso (i.e. chest or abs) and back torso (i.e. upper back, lower back) simultaneously. Introduction of electrical current into the heart may cause cardiac arrhythmia.
- Over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins)
- · Over, or in proximity to, cancerous lesions

### 8. General Electrode Pads Safety Measures

- Always use electrode pads supplied by the manufacturer. Other electrodes may have electrical properties that are
  unsuitable for or may cause damage.
- For best results, wash and clean the skin of any oil and dry it before attaching the electrode pads.

- For reasons of hygiene, each patient must have their own set of electrode pads. Do not use the same electrodes on different patients.
- Contaminated electrode pads or hydrogel can lead to infection.
- Use of electrode pads with degraded hydrogel can result in burns to the skin.
- Do not place the electrode pads in water.
- Do not apply solvents of any kind to the electrode pads.
- Always stop the stimulation before moving or removing any electrode pads during a session, to avoid electrical shock to the patient.
- Attach the electrode pads in such a way that their entire surface is in contact with the skin.
- Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive
  medium. The irritation can usually be reduced by using an alternate conductive medium, or alternate electrode pad
  placement. If redness of the skin is observed under the electrode pad, do not start another stimulation session on the same
  area if the redness is still visible.
- Electrode pads usage and storage instructions are provided on the mylar bag with electrodes.

# PowerDot PD-02 Device & Accessories

- 1) Power/Stimulation indicator
- 2) Multifunction button
- 3) Micro USB connector (dual function: connects either leads or charging cable)
- 4) Charging indicator
- 5) Female magnetic pad connector
- 6) Lead cable with two female magnetic connectors

#### Power/Stimulation Indicator (1):

Can be either WHITE (when the PowerDot is on or on standby mode) or ORANGE (when it's being used for stimulation).

### Multifunction Button (2) Modes:

Multifunction Button carries out these PowerDot functions:

Multi-function Button Action	PowerDot Initial State	PowerDot Resulting State
Hold for ~1 second	PowerDot is OFF (no lights)	PowerDot is now ON (white light is on)
Hold for ~1 second	PowerDot is ON (white light is on)	PowerDot is now OFF (no lights)
Quick Click	PowerDot is in Stimulation (orange light is on)	Stimulation stops (orange light is ON )
Hold for 5 seconds	Device is ON	Full Factory Reset: activation lock is released, , Power LED turns OFF, after few seconds turns ON again, blinks several times and the PowerDot becomes OFF
Hold for 3 seconds	Device is OFF	Soft Reset: Activation lock is released, Power LED blinks several times

# Micro-USB Connector (3):

PowerDot uses the same high voltage micro-USB connector (3) for lead cable connection and for charging.

If the device is being used for stimulation, it will immediately stop once the lead cables are disconnected from the micro-USB connector.

# Charging Indicator (4):

While the device is still charging, you will see the ORANGE Charging Indicator (4) next to the micro-USB charging connector. Refer to the **Recharging PowerDot** section below for more information on PowerDot charging.

# Snap Connector (5):

Used to attach and hold the PowerDot device in place on the rectangular pad.

### Lead Cable (6):

For an easier, more comfortable PowerDot stimulation, two lead cables of different lengths (10 cm and 30 cm) (6) are provided with each PowerDot. These give you the option to choose the right cable to reach the muscle group you're targeting, depending on your physical measurements.

# Electrode Pads (7,8):

PowerDot supports 3 types of electrode pads -5.5 cm (2.2") diameter round pads, 4.5 cm (1.8") diameter round pads and 9x5 cm (3.5"x2") rectangular pads.

PowerDot pads use unique skin biocompatible hydrogel with superb conductive qualities and adhesiveness.

The lifetime of PowerDot pads depends a lot on your individual skin properties, level of hairlessness and the quality of maintenance. On average, each pad lasts for around 20 stimulation sessions. After that, adhesiveness and conductive properties of the pads may start deteriorating.



Store the pads on the safety film, in a dry environment (either in the original zipper plastic bag or inside the PowerDot carry case). Make sure you attach the pads to clean and dry skin!

#### **Directions For Use**

### 1. Installing/Launching PowerDot Doctor App

- 1) Make sure you iPad runs iOS 9.0 or later.
- Launch Apple App Store application, search for "PowerDot Doctor" Mobile Application and install it.
- 3) Launch the installed PowerDot Doctor application, create your profile and then follow instructions to locate and activate your PowerDot (or PowerDots) for the first use.

#### 2. Turning PowerDot ON/OFF

To turn on your PowerDot PD-02 unit, hold the Multifunction Button for approximately 1 second, until you see the power light turn on. When PowerDot is not in stimulation, you can turn it off by holding the Multifunction Button again for approximately 1 second.

When a stimulation is ongoing, click on the Multifunction Button once to stop the stimulation and then hold the Multifunction Button to turn off the PowerDot.

### 3. Activating PowerDot

Before PowerDot PD-02 can be used for muscle stimulation, it must be activated from within your PowerDot Doctor App.

When PowerDot is activated, it gets paired with your mobile phone securely. Neither you nor anyone else will be able to connect and use your PowerDot from any other mobile device.

Follow the App's onscreen advice to scan for and activate your PowerDots. Make sure your PowerDot (or PowerDots) is turned on before you start scanning.

Use the Devices screen (available from the PowerDot Doctor App menu) to activate another or additional PowerDot devices or to deactivate previously activated ones. You can activate and use up to 32 PowerDots from your PowerDot Doctor app.


### 4. Preparing and Placing PowerDot

1) Plug the lead cable you intend to use into the USB Type C connector on your PowerDot PD-02 device.



2) Detach both round pads, rectangular pad from their safety film and stick them to the targeted part of your body, according to the pad placement visual guidelines provided by PowerDot Doctor App or in this manual. Based on your experience and understanding of the principles of electrical muscle stimulation, you

may try to introduce some variations to the offered pad placement set ups for better efficiency. Make sure you create and store photos of the new setups once they prove efficient. PowerDot Doctor app provides you with the functionality to store and display custom pad placement photos for your patients.



- Attach female magnetic round pad connectors on the lead cable to the male snaps on the round pads.
   Attach PowerDot PD-02 device female magnetic connector to the male snap on the rectangular pad;
- 4) Make sure that the entire surface of the round and rectangle pads are completely and securely attached to your body;
- 5) It's a good time to turn on your PowerDot (or PowerDots), if you haven't done it yet.



Always thoroughly inspect the lead cable and PowerDot unit itself for any signs of damage BEFORE every stimulation session. Do not use damaged accessories or devices For your safety, you are strongly advised to replace them before using again.

# 5. Starting a Stimulation

PowerDot Doctor App offers multiple ways to initiate a stimulation session:

1) Initiate ad hoc stimulation session, which is not attributed to any specific patient.

You will be asked first to select stimulation program, relevant muscle groups or body areas and adjust stimulation parameters and optionally select number of PowerDot PD-02 devices to be used before launching your session.

We generally recommend using this option for testing or demo purposes.

- 2) Initiate ad hoc stimulation session for selected patient.
  - o First, the relevant patient is looked up and selected on the Patients screen.
  - o Secondly, you will be asked first to select stimulation program, relevant muscle groups or body areas and optionally adjust stimulation parameters before launching your session.
  - o Thirdly, if relevant, you will be asked to select number of PowerDot PD-02 devices you want to engage in the current session (one or two).
  - o Once this stimulation session is completed, relevant session statistics will be recorded on the patient's Sessions History screen.

We generally recommend using this option for casual/one off sessions or when testing patient for the tolerance and/or reaction to electrical stimulation.

- 3) Initiate scheduled stimulation session for the selected patient.
  - o First, you look up and select relevant patient from the Patients screen.

- Secondly, you will be asked to look up and select stimulation program, up to 3 relevant muscle groups or body areas and optionally select number of PowerDot devices to be used.
- Thirdly, you will be asked to select scheduling parameters, such as total duration of the program (in days), periodization of the sessions, starting date, notifications settings (about missing or late stimulation sessions).
- o Fourthly, you can optionally mark stimulation sessions to be accessible for the Remote (home use) stimulation from Patient Mobile App.

Be careful when confirming stimulation sessions for remote execution by the patient!

First, you have to make sure that your patients or their caregivers are prepared for and capable of running unsupervised electrical stimulation sessions.



Specifically, they should be aware on to assemble, apply and control PowerDot PD-02, are able to follow stimulation pad placements and stimulation positions recommendations and are able to adjust their stimulation intensities in the controlled manner.

They should be also aware of the product counter-indications, basic troubleshooting and maintenance guidelines. Remote stimulation sessions are not for everyone!

 Finally, you will be asked to adjust stimulation parameters for the session (in the relevant ranges only).

Use scheduled sessions functionality whenever a stimulation course is required for the specific patient. Scheduled sessions' parameters, listed above, can be always changed/re-adjusted if required.


Before commencing the stimulation, PowerDot Doctor App will be asking to turn on PowerDot PD-02 device (or devices) which you are planning to use for the current stimulation session (please make sure the device has been previously activated via Devices screen). Doctor App will automatically detect turned on devices, connect to them, check for activation, display serial numbers and blink Power LED to confirm that right PowerDot PD-02 device has been selected and associated with the session.

In addition, PowerDot Doctor App will automatically check for remaining battery level and lead cable connection status for all the PowerDot PD-02 devices associated with the current session.

After the checks, Get Started button will appear, which means you can already launch the stimulation session.

### 6. Controlling A Stimulation

Using the PowerDot Doctor App, you can control your stimulation session from the Stimulation Dashboard, which has the following controls available:



*Header Indicator* (1) displays the remaining time for the total stimulation session, current stimulation phase (if applicable), name of the stimulation program and stimulation phase (if applicable), and name of the patient (if applicable).

Contraction/Rest Indicator (2) – for the stimulation programs (and individual phases) involving muscle contractions, helps to understand time between and during contraction and rest intervals.

Intensity Adjustment controls – Tap on Intensity Increase (or "+") (3) and Intensity Decrease (or "-") (4) buttons an accurate adjustment, or hold either button down for a faster adjustment.

Device-level Intensities Lock (6) and Channel-level Intensities Lock controls (5) synchronize intensity changes either a) between PowerDot PD-02 devices (this applies only to the two-device Duo stimulation mode) or (b) between two stimulation channels of individual PowerDot PD-02 device.

Both controls are enabled by default and can be used to fine-tune the intensity of your current stimulation session. However, it's best to use them only when you really need them (e.g. to address muscle imbalance).

Use *Stop, Pause, Resume, Skip Phase* controls to control the flow of the stimulation session. Skip Phase control is only applicable to stimulation programs which are comprised of multiple and skippable stimulation phases (e.g. such as Performance *Strength* program).

Switch Session Control (7) - allows you to switch between actively running stimulation sessions (on a different patients). When switching between sessions, your iPad will be disconnecting from PowerDot PD-02(s) used in currently active session and re-connecting to PowerDot(s) being used in the newly selected session.

Home Screen Control (8) – helps you to navigate back to PowerDot Doctor App Main Menu screen, so you could select another patient and initiate a new stimulation session.

<u>/!\</u>

When applying Switch Session and Home Screen Controls, you will be leaving current stimulation session and disconnecting from currently used PowerDot PD-02 devices. Although the devices will be continuing maintaining your stimulation session parameters, you won't be able to control them programmatically (e.g. adjust intensities or stop/resume session using App controls) until your select your original stimulation session and re-connect to your devices again. Re-connection is not immediate and it can take up to 10 seconds, therefore, in all emergency situations, we strongly advice to get yourself familiar and follow manual Terminating Stimulation protocols 2 and 3, as provided below.



If your patient experiences major discomfort or pain – PAUSE your stimulation session and/or DECREASE intensities. For your patient's safety, after a PAUSE, your stimulation session will resume at only 80% of the previous intensity values.

### 7. Terminating Stimulation

It's important to know the quick ways to stop the stimulation session when something unexpected happens (e.g. pads get detached or partially detached; patient starts feeling pain; the stimulation area gets wet, etc.).

There are 3 main ways to immediately terminate stimulation:

- 1) The recommended/most commonly used option: Tap Pause or Stop on the Stimulation Dashboard
- 2) Alternatively press the Multifunctional Button on the PowerDot device
- 3) Or: quickly unplug the lead cable from PowerDot device

# 8. Carrying & Storing PowerDot

PowerDot PD-02 carry case is specifically designed for carrying and storing up to 2 PowerDot devices, their lead cables and electrode pads.

To keep electrode pads clean and make them last longer, always re-attach them to the safety film in between use, then store them in your carrying case pocket. You can use both sides of a single safety film to attach one set of PowerDot electrode pads (one side for the rectangular pad and the other for the two round pads).

In multi-patient environment, it is always a good idea to mark patient pads with the patient name or initials on the top surface or on the pads film, e.g. using resistant marker pen.

# 9. Deactivating PowerDot

PowerDots can be deactivated from your Doctor App either from within the Patient or or by manually resetting the PowerDot device.

For manual deactivation: when PowerDot is OFF, press and hold Multi-function button for 3 seconds until you see the power indicator blinking. Now you can activate your PowerDot.

Deactivation from the PowerDot Doctor App is performed from the Devices screen as per on-screen instructions.

Upon successful deactivation, the power indicator will blink several times.

# **Intensity Adjustment Guidelines**

For TENS programs or less demanding NMES programs, you are generally required to progressively increase the stimulation intensities until you see muscle twitches or feel comfortable sensation.

For the majority of NMES stimulation programs, that use powerful muscle contractions, the efficacy of the treatment can be proportional to the maximum number of fibers being recruited, and, therefore, it is often required to adjust your stimulation intensity to the maximum levels the patient can comfortably endure.

Keep in mind that the maximum intensity levels may vary, not only from one stimulation session to another, but also within the course of a single stimulation session. It's possible that patient's muscles will adapt to stimulation at a certain intensity level reasonably quickly. Various conditions, such as differences in skin's dampness or sweat level, or the rate of the electrode pad deterioration, may affect the intensity of stimulation.



During stimulation, it is important to maintain feedback with the patient to make sure that his/her stimulation intensity settings are always in balance with the level of comfort on one side and efficiency of stimulation on the other.

For efficiency in NMES and Performance sessions, based on patient feedback and feelings, try increasing slowly intensities during the course of the session.

# **Charging PowerDot**

PowerDot PD-02 device can be recharged from any reliable USB connection (e.g. laptop, mobile phone charger, wall USB charger, etc.). It takes around 60 minutes for the device to go from zero level to full charge.

NEVER charge PowerDot when it is attached to the patient's body.

Always use an original charging cable provided with PowerDot.



NEVER charge PowerDot from unreliable or problematic sources.

When using  $3^{rd}$  party party USB AC chargers, we recommend unplugging the AC plug from the wall before contacting PowerDot device.

PowerDot PD-02 uses UL 1642 and UN 38.3 certified built-in Lithium Polymer battery, which requires recharging after approximately 5-6 hours of continuous usage. The battery will last for at least 500 charging cycles.

If you plan to store PowerDot PD-02, unused, for longer than six months, charge it to at least 50% every six months.

In the PowerDot Doctor App, the current battery charge level is displayed at the Stimulation and is also available on the Devices screen.

When you see the ORANGE charging light next to micro-USB connector, this means PowerDot is currently charging.

Once PowerDot is fully charged, the ORANGE light will turn off.

# **Stimulation Protocols**

Stimulation protocols, accessible to doctors from PowerDot Doctor App, based on their intended use, can be divided into 3 main groups:

# 1. NMES Rehabilitation Protocols

# Disuse Atrophy

Protocol Name					Indication For Use		
Disuse Atrophy*					Retarding disuse atrophy		
Description							
Targeting slow to	witch (Type I) mus	scle fibers to facili	tate increases in mu	scle v	olume after long pe	riods of immobiliz	ation or
diminished move	ement.						
Contraction, sec	Contraction, Hz	Rest, sec	Rest, Hz	Ra	mp Up/Down, sec	Pulse Width, uS	Duration, min
(allowed range)	(allowed range)	(allowed range)	(allowed range)		(allowed range)	(allowed range)	(allowed range)
5	40	5	4		2/0.5	208	25
(3-8)	(30-50)	(3-8)	(0-6)		(0.5-3)	(160-320)	(15-30)
Reinforcement 1	Post Disuse Atrop	hy*			Stimulation of hea	2	der to improve
					or facilitate muscl		
		muscles that have	e gained the volume	. Incre	eases strength and st	tability by targetin	g fast twitch
(type II) muscle						T	
4	80	8	4		1.5/0.75	320	20
(3-6)	(60-100)	(5-12)	(0-6)		(0.5-3)	(208-432)	(15-30)
Reinforcement l	Post Disuse Atrop	hy (Agonist/Anta	ngonist)		Stimulation of hea	2	der to improve
					or facilitate muscl	1	
			e gained the volume	. Incre	eases strength and st	tability by targetin	g fast twitch
(type II) muscle	fibers. Agonist/An	tagonist sequencir				,	
4	70	4	4		1.5/0.75	320	15
(3-6)	(60-100)	(3-6)	(0-6)		(0.5-3)	(208-432)	(10-20)
Prevention of D					Preventing disuse		
			contractions compar	able t	to the normal level of	of activity during a	day.
	e immobilized lim		one fracture.			T T	
5	35	12	4		2/0.5	320	45
(3-8)	(30-50)	(9-18)	(0-6)		(0.5-3)	(208-432)	(30-60)
	isuse Atrophy (Ag				Preventing disuse		
			ontractions compara				day.
			one fracture. Agon	ist/An			
6	35	6	4		2/0.5	320	20
(4-8)	(30-50)	(4-8)	(0-6)		(0.5-3)	(208-432)	(15-30)

# Muscle Re-education and Increasing Range of Motion

Protocol Name					Indication For Use		
Muscle Re-education, Level 1 (Slow Twitch)*					Muscle re-educati	on	
Description							
Targets re-educa	tion and muscle fu	nction improveme	ents of slow twitch (	Type :	I) muscle fibers.		
Contraction, sec	Contraction, Hz	Rest, sec	Rest, Hz	Ra	ump Up/Down, sec	Pulse Width, uS	Duration, min
(allowed range)	(allowed range)	(allowed range)	(allowed range)		(allowed range)	(allowed range)	(allowed range)
5	40	5	4		2/0.5	208	25
(3-8)	(30-50)	(3-8)	(0-6)		(0.5-3)	(160-320)	(15-30)
Muscle Re-educ	ation, Level 2 (F	ast Twitch)*			Muscle re-educati	on	
Targets re-educa	tion and muscle fu	nction improveme	ents of fast twitch (T	ype I	I) muscle fibers. Ap	ply/schedule after	at least 2 weeks
of usage of Leve	l 1 program.						
4	80	8	4		1.5/0.75	320	20
(3-6) (50-100) (5-12) (0-6) (0.5-3) (208-432) (15-30)						(15-30)	
Neuro Rehab* Muscle re-education							
Targets facilitation	on and re-learning	of motor skills for	stroke patients.				

5	40	15	0		4/2	208	20
(3-8)	(30-50)	(9-24)			(1.5-6)	(160-320)	(15-30)
Increase Range of Motion*  Maintaining or increasing range of motion							
Targets increase	in range of motior	by creating contr	action and moveme	nt in t	he antagonist muscl	e, thus causing ma	aximum stretch
of the problemat	ic muscle.						
5	40	10	4		2/0.5	208	20
(3-8)	(30-50)	(6-16)	(0-6)		(0.5-3)	(160-320)	(15-30)
Spasticity Treat	ment*				Maintaining or inc	creasing range of a	notion
Targets reduction	n of spasticity by i	nhibition of motor	neurons of the spas	stic mu	scle via the recipro	cal inhibition refle	ex. To be used
on the antagonis	on the antagonist of the spastic muscle with the intensity sufficient to create movement to the maximum range of motion.						
5	40	10	4		2/0.5	208	20
(3-8)	(30-50)	(6-16)	(0-6)		(0.5-3)	(160-320)	(15-30)

# Specific Rehabilitation

Protocol Name				Indication For Us	e	
<b>Hip Prosthesis</b> ,	Level 1			Retarding disuse a	atrophy	
Description						
			Maximus and Medius	s muscles following the o	rthopedic surgery	of hip. Level 1 i
targeting slow tw	itch (Type I) muse	cle fibers.				
Contraction, sec		Rest, sec	Rest, Hz	Ramp Up/Down, sec	Pulse Width, uS	Duration, mir
(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allowed range
6	40	6	0	1.5/0.75	208	30
(4-8)	(30-50)	(4-8)		(0.5-3)	(160-320)	(20-40)
Hip Prosthesis,				Retarding disuse a		
				s muscles following the o		of hip. Level 1
is targeting fast	twitch (Type II) m	uscle fibers. Use/s	schedule Level 2 at 1	east 2 weeks after using l	Level 1 program.	
4	80	10	0	1.5/0.75	208	15
(3-6)	(60-100)	(7-15)		(0.5-3)	(160-320)	(10-20)
Patellar Syndro				Retarding disuse a		
				atrophy in Medial Vastus	and larger Quadri	ceps muscle.
Level 1 is targeti	ng slow twitch (Ty	pe I) muscle fiber	rs.			
6	40	6	0	1.5/0.75	208	30
(4-8)	(30-50)	(4-8)		(0.5-3)	(160-320)	(20-40)
Patellar Syndro	me, Level 2			Retarding disuse a	atrophy	
Targets increase	of the stability of t	he knee and retard	dation of the disuse	atrophy in Medial Vastus	and larger Quadri	ceps muscle.
Level 2 is targeti	ng fast twitch (Typ	e II) muscle fiber	s. Use/schedule Lev	rel 2 at least 2 weeks after	r using Level 1 pro	gram.
4	80	10	0	1.5/0.75	208	15
(3-6)	(60-100)	(7-15)		(0.5-3)	(160-320)	(10-20)
Hemiplegia*				Maintaining or inc		
Targets re-learni	ng of motor skills	and reduction of s	pasticity. Features lo	ong ramp up time o the co	ontraction and a loa	ng rest phase.
Can be used to co	omplement tradition	onal physiotherapy	after a stroke.			
10	45	20	0	4/2	208	20
(6-15)	(30-50)	(12-30)		(1.5-6)	(160-320)	(15-30)
Rotator Cuff, L				Muscle re-educati		
				of the rotator cuff. Use to		ge of motion,
shoulder tendopa	thies and pain aro	und the shoulder.	Training can be con	nbined with active mover		
6	40	6	4	1.5/0.75	208	20
(6-15)	(30-50)	(4-8)	(0-6)	(0.5-3)	(160-320)	(15-30)
Rotator Cuff, L				Muscle re-educati		
				of the rotator cuff. Use to		
		und the shoulder.	Training can be con	nbined with active mover	nents. Apply/sched	dule after at leas
	Level 1 program.					
4	80	10	4	1.5/0.75	208	15
(3-6)	(60-100)	(7-15)	(0-6)	(0.5-3)	(160-320)	(10-20)
	tment (Agonist/A			Muscle re-educati		
				ACL. This is a contraction	n offset session: st	timulation starts
on Hamstrings at	nd then continues of	on the Quadriceps	, preventing the risk	of the anterior drawer.		
3+6	40	8	4	1.5 + 3/0 + 0.75	320	25

	(30-50)					(240-432)	
Shoulder Subluxation Prevention and retardation of atrophy							
Targets stimulati	on of Deltoid and	Supraspinatus mus	scles as prevention	and tre	atment for atrophy	of the sub-disloca	ted shoulder.
8	40	8	0		3/1.5	208	25
(6-10)	(30-50)	(6-10)			(1-5)	(160-320)	(20-30)
Back-Trunk Sta	bilization*				Muscle re-educati	on	
			mulating abdomina				
in the bank and trunk due to long term pain or neurological disorder. This program can be also combined with active movements.							
6	40	12	4		2/1	208	30
(3-8)	(30-50)	(6-16)	(0-6)		(0.5-3)	(160-320)	(20-40)

# $Circulation, \ Relaxation \ of \ Muscle \ Spasms, \ Vascular$

Protocol Name				Indication For Us		
Heavy Legs				Increasing local b	lood circulation	
• 0				Relaxation of muscle spasms		
Description					-	
compress the dee		and eject venous	blood upwards. The	g of the legs being very h stimulation will also hel		
	Contraction, Hz	Rest. sec	Rest. Hz	Ramp Up/Down, sec	Pulse Width. uS	Duration, min
(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allowed range
continuous	7/5/3	0	0	1/0.5	208	21
	(7/7/7 min)				(160-320)	(15-30)
Capillarization				Increasing local b	lood circulation	
	uscle fibers more i	resistant to fatigue	. Capillarization dec	reases the amount of lact	ic acid being prod	uced and creates
a larger area of e	xchange and distri	bution of oxygen a	and metabolites.			
continuous	8	0	0	1/0.5	208	25
	(5-10)				(160-320)	(20-30)
De-contraction Relaxation of muscle spasms						
				Post-surgical and	post-trauma acute	pain
Targets reduction	of muscular tensi	on by using very l	ow stimulation freq	uency.		
continuous	1	0	0	1/0.5	208	20
					(160-320)	(10-40)
Arterial Insuffic	eiency			Increasing local b	lood circulation	
improves the tol		iding the tetaniza e.		tion increases the capacit ble fatigue. Electrodes		
15	10	15	3	1/1	208	15
(10-20)	(7-12)	(10-20)	(0-6)	(0.5-1.5)	(160-320)	(10-20)
	enous Thrombosi			Prevention of ven		
				Uses short tetanic contra eles (Gastrocnemius) and		y long active
4	50	20	8	1.5/1.5	208	20
(3-5)	(40-60)	(15-25)	(6-10)	(1-2)	(160-320)	(15-30)
Cramp Preventi		( /	(* */	Increasing local background	lood circulation	( /
Targeting cramps Each sequence re		quence increases	local blood circulation	on, second sequence prov		muscle tone.
continuous	8/3 x 4 (8 min/2 min) (2-10)	0	0	1/0.5	208 (160-320)	40

# 2. TENS Pain Management Protocols

# By Type of Pain

Protocol Name		Indication For Use		
Acute Pain		Post-surgical and post-trauma acute pain		
Description				
High frequency TENS is a popular	pains. Adjust pulse width to patient pain sen	sitivity level.		
Frequency, Hz	Burst Parameters	Pulse Width, uS	Duration, min	
(allowed range)	(allowed range)	(allowed range)	(allowed range)	
100	N/A	208	30	
(80-150)		(32-320)	(15-60)	
Fracture Pain		Post-surgical and post-trauma acute pain		
Medium frequency/medium pulse	width TENS is a good choice for treat	ing sensitive fracture pains. Pain manageme	ent helps to	
prevent various complications such	as pulmonary complications or thron	nbosis.		
75	N/A	160	30	
(60-80)		(128-240)	(15-60)	
Chronic Pain		Symptomatic relief and management of chronic,		
		intractable pain		
		Temporary relief of pain associated with s		
		muscles due to strain from exercise or nor	mal household	
		and work activities		
		timulation. Low frequency stimulation allev	viates pain by	
		es also increase local blood circulation.	Τ	
5	N/A	240	30	
(1-8)		(160-320)	(15-60)	
Radiating Pain		Symptomatic relief and management of chronic,		
Projected Pain		intractable pain		
		d) pains such as mononeuropathy, central p		
		y bursts alleviate pain by stimulating muscl	es to release the	
	itches also increase local blood circula		20	
80	5 pulses per burst, 2 bursts/sec	208	20	
(80-150)	(5-10 ppb, 2-4 bursts/sec)	(160-320)	(10-40)	

# By Body Part

Protocol Name		Indication For Use			
Cervical Pain		Symptomatic relief and management of	chronic,		
Neck Pain		intractable pain	r		
		Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities			
Description					
For treating neck pain resulting from chronically contracted of Levator scapulae and/or superior Trapezius muscles. Low fre stimulation alleviates pain by stimulating muscles to release the endorphins. In addition, muscle twitches also increase local circulation.					
Frequency, Hz	Burst Parameters	Pulse Width, uS	Duration, min		
(allowed range)	(allowed range)	(allowed range)	(allowed range)		
5	N/A	240	30		
(1-8)		(160-320)	(15-60)		
Lower Back Pain		Symptomatic relief and management of chronic, intractable pain			
		Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities			

		ar paravertebral muscles. Low frequency stition, muscle twitches also increase local blo	
5	N/A	240	30
(1-8)	14/11	(160-320)	(15-60)
Lumbosciatica Sciatica		Symptomatic relief and management of clintractable pain	
		Temporary relief of pain associated with a muscles due to strain from exercise or not and work activities	mal household
nerve root exit sites and on the lov	wer part of the buttock/posterior face	paravertebral muscles. Electrodes are placed to of the thigh. Low frequency stimulation al hes also increase local blood circulation.	
5	N/A	240	30
(1-8)		(160-320)	(15-60)
Mononeuropathy Central Pain Deep Muscular Pain		Symptomatic relief and management of clintractable pain	
Rhizopathy		Temporary relief of pain associated with a muscles due to strain from exercise or not and work activities	mal household
sensory of touch as well as for deep	p muscular pains. High frequency pu	(projected) pains, for conditions with reduce lses packaged into low frequency bursts alle- hes also increase local blood circulation.	
(80-150)	(5-10 ppb, 2-4 bursts/sec)	(160-320)	(10-40)
Arthralgia	(3-10 ppo, 2-4 bursts/sec)	Symptomatic relief and management of cl	
	n or affecting a joint, caused by eithe tion of stimulation frequency at very	Temporary relief of pain associated with a muscles due to strain from exercise or not and work activities  r degenerative or destructive processes. To a short pulses.	mal household
50-150 Hz at 2 sec	N/A	48	20
(50-150 Hz at 2-5 sec)	1771	(32-160)	(10-40)
Epicondylitis		Symptomatic relief and management of clintractable pain	nronic,
		Temporary relief of pain associated with a muscles due to strain from exercise or not and work activities	
For treating epicondylitis pain. Epi adaptation, this protocol uses conti	icondylitis is a tendinopathy which minuous variation of stimulation freque	ay result from repetitive gripping of the objective with short width pulses.	ects. To avoid
50-150 Hz at 2 sec	N/A	48	20
(50-150 Hz at 2-5 sec)		(32-160)	(10-40)
		Symptomatic relief and management of c	ronic.
Lower Back Muscle Pain		intractable pain	,
		Temporary relief of pain associated with muscles due to strain from exercise or nor and work activities	sore and aching mal household
For treating lower back muscle pair		Temporary relief of pain associated with muscles due to strain from exercise or not	sore and aching

Muscle Pain		Symptomatic relief and management of intractable pain	chronic,			
		Post-surgical and post-trauma acute pai	Post-surgical and post-trauma acute pain			
		Temporary relief of pain associated wit muscles due to strain from exercise or r and work activities				
		h modulation, which creates undulating sensat	ion and is generally			
more comfortable than constant pulse 80 Hz at 2 sec (60-150 Hz at 2-5 sec)	N/A	64-320 (32-432)	30 (15-60)			
Knee Pain		Symptomatic relief and management of intractable pain				
		Post-surgical and post-trauma acute pai	n			
		Temporary relief of pain associated wit muscles due to strain from exercise or r and work activities	normal household			
		nedical conditions such as osteoarthritis, rhe h creates undulating sensation and is generall				
80 Hz at 2 sec (60-150 Hz at 2-5 sec)	N/A	64-320 (32-432)	30 (15-60)			
Shoulder Pain		Symptomatic relief and management of intractable pain	chronic,			
		Post-surgical and post-trauma acute pai	n			
		Temporary relief of pain associated wit muscles due to strain from exercise or r and work activities				
	protocol uses pulse width mod	orts injuries or underlying medical conditions lulation, which creates undulating sensation ar				
80 Hz at 2 sec (60-150 Hz at 2-5 sec)	N/A	64-320 (32-432)	30 (15-60)			
Joint Pain		Symptomatic relief and management of intractable pain	chronic,			
		muscles due to strain from exercise or r and work activities	Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities			
For treating various kinds of joint pain undulating sensation and is generally r		knees. This protocol uses pulse width modulant pulse width stimulation.	tion, which creates			
60 Hz at 2 sec (50-100 Hz at 2-5 sec)	N/A	48-160 (32-240)	30 (15-60)			
Lumbago		Symptomatic relief and management of intractable pain	chronic,			
For treating lumbago pain, the result of	f sharp and intense contraction	Temporary relief of pain associated wit muscles due to strain from exercise or r and work activities ons of lower back muscle. This protocol produc	normal household			
frequency stimulation twitches to redu		of the stimulated muscles.	co very low			
1	N/A	208 (160-320)	20 (10-40)			
Torticollis		Symptomatic relief and management of intractable pain				

		Temporary relief of pain associated with sore and aching				
		muscles due to strain from exercise or normal household				
	and work activities					
For treating torticollis pain: a sharp pain in the neck, often accompanies by considerable reduction of the mobility of the cervical region. This protocol produces very low frequency stimulation twitches to reduce muscle tension during rest of the stimulated						
muscles.						
1	N/A	208	20			
		(160-320)	(10-40)			

# By Type Of TENS

Protocol Name		Indication For Use			
Conventional TENS, High Frequ	lency	Post-surgical and post-trauma acute pain			
		Temporary relief of pain associated with sore and ach muscles due to strain from exercise or normal household a work activities			
Description		1			
		e pains. Adjust pulse width to patient pain se- uency TENS/Hans Stimulation and Mediun			
Frequency, Hz	Burst Parameters	Pulse Width, uS	Duration, min		
(allowed range)	(allowed range)	(allowed range)	(allowed range)		
100	N/A	208	30		
(80-150)		(32-320)	(15-60)		
Conventional TENS, Medium Fi	requency	Post-surgical and post-trauma acute pain			
		muscles due to strain from exercise or no work activities	Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities		
	sed as a more comfortable alternative				
75	N/A	160	30		
(60-80)		(128-240)	(15-60)		
Acupuncture TENS Endorphin TENS		Symptomatic relief and management of chronic, intractable pain			
Low Frequency TENS		Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities			
Burst TENS or less adaptive Low	I is generally used to manage chroni Range Frequency Modulated TENS In addition, muscle twitches also inc	c pains, can be varied with or accompanied by protocols. Low frequency pulses alleviate parrease local blood circulation.	y more comfortable in by stimulating		
5	N/A	240	30		
(1-8)		(160-320)	(15-60)		
Conventional TENS, High Frequ	ency/Short Pulse	Post-surgical and post-trauma acute pain			
		Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities			
	choice for the management of acute e also varied with Mixed Frequency	e pains. Short pulse version should be used fo TENS/Hans Stimulation protocol.	r the patients with		
100	N/A	48	20		
(80-150)		(32-160)	(10-40)		
Frequency Modulated TENS, Hi	gh Range	Post-surgical and post-trauma acute pain			
This protocol can for treatment of	vorious chronia poin conditions such	Temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities arthralgia, epicondilitis, lower back muscle pains etc.			
		th to patient pain sensitivity level in the treat			

50-150 Hz at 2 sec	N/A	48	20	
(50-150 Hz at 2-5 sec)	N/A	(32-160)	(10-40)	
Frequency Modulated TENS, Me	dium Range	Symptomatic relief and management of chronic, intractable		
1.1.1.		pain	,	
		Temporary relief of pain associated wi	th save and eahing	
		muscles due to strain from exercise or no		
		work activities	ina noasenora ana	
		ns such as arthralgia, epicondilitis, lower bac		
	avoid adaptation. Medium range mo	dulation provides more comfortable stimula	tion than high range	
option. 40-100 Hz at 3 sec	NT/A	240	30	
(20-100 Hz at 3 sec (20-100 Hz at 2-5 sec)	N/A	(160-320)	(15-60)	
Frequency Modulated TENS, Lo	w Range	Symptomatic relief and management of		
14		pain		
		Temporary relief of pain associated wi		
		muscles due to strain from exercise or no work activities	rmai nousenoid and	
This protocol is a more effective ve	ersion of Low Frequency TENS and I	Burst TENS protocols for the long-term treat	ments as it avoids	
		in by stimulating muscles to release the end-		
muscle twitches also increase local	blood circulation.	•		
1-8 Hz at 5 sec	N/A	240	30	
(1-15 Hz at 3-10 sec)		(160-320)	(15-60)	
Pulse Modulated TENS, High Ra	nge	Symptomatic relief and management of pain	chronic, intractable	
		Temporary relief of pain associated with sore and aching		
		muscles due to strain from exercise or normal household and		
This protocol can used for treating	of chronic pain conditions requiring	work activities g higher levels of comfort such as muscle p	oin charn knoo and	
	age like effect on certain muscle grou		am, sharp knee and	
80 Hz at 2 sec	N/A	64-320	30	
(60-150 Hz at 2-5 sec)		(32-432)	(15-60)	
Pulse Modulated TENS, Low Rai	nge	Symptomatic relief and management of	chronic, intractable	
		pain		
		Temporary relief of pain associated wi	th sore and aching	
		muscles due to strain from exercise or no		
		work activities		
		ing higher levels of comfort. Low Range		
60 Hz at 2 sec	N/A	, more suitable for treating of the joint pains 48-160	30	
(50-100 Hz at 2-5 sec)	14/14	(32-240)	(15-60)	
Burst TENS		Symptomatic relief and management of	chronic, intractable	
		pain		
		Towns and a first of a six a second and	41 11	
		Temporary relief of pain associated wi muscles due to strain from exercise or no		
		work activities	illiai nouschoid and	
Burst TENS is usually effective ag	ainst radiating pain in the arms or leg	gs (rhizopathy), for conditions with reduced	or changed sensory	
of touch, for deep muscular pain or	when the post-treatment of TENS is	s too short. High frequency pulses packaged	l into low frequency	
bursts alleviate pain by stimulating	muscles to release the endorphins. In	addition, muscle twitches also increase loca	al blood circulation.	
Ruret TENS can be offeetively a	altered with more aggressive I ow	Fraguency/Endorphin TENS and loss add	entable Low Dance	
Frequency Modulated TENS protoc		Frequency/Endorphin TENS and less ada	ipianie Low Kallge	
80	5 pulses per burst, 2 bursts/sec	208	20	
(80-150)	(5-10 ppb, 2-4 bursts/sec)	(160-320)	(10-40)	
Mixed Frequency TENS		Post-surgical and post-trauma acute pain		
Hans Stimulation				

		Temporary relief of pain associated with sore and aching			
		muscles due to strain from exercise or normal household			
		work activities			
Mixed Frequency TENS is more co	omfortable and convenient form of Hi	gh Frequency Conventional TENS. Also co	mmonly referred as		
Hans Stimulation. Can be altered w	ith more aggressive High Frequency	Conventional TENS protocol.			
80/2, 3 sec each	N/A	208 30			
		(32-320) (15-60)			
1 Hz TENS		Post-surgical and post-trauma acute pain			
Relaxation of muscle spasms					
Targets reduction of muscular tension by using very low stimulation frequency.					
1	N/A	208 20			
		(160-320)	(10-40)		

# 3. NMES Performance and Massage Programs

Protocol Name					cation For U			
Muscle Enduran	ice*				Stimulation of healthy muscles in order to improve			
						le performance	-	
Description								
Targets improven	nents of muscle en	durance and fatigu	e resistance through	n activation	of slow twite	ch (Type I) muscle	fibers.	
Contraction, sec	Contraction, Hz	Rest, sec	Rest, Hz		p/Down, sec	Pulse Width, uS		
(allowed range)	(allowed range)	(allowed range)	(allowed range)	(allow	ved range)	(allowed range)	(allowed range	
8	15	4	3		1/0.5	208	45	
(6-10)	(30-25)	(3-5)	(0-6)		).5-2)	(160-432)	(35-60)	
Strength Endura	ınce*					althy muscles in o	rder to improve	
						le performance		
			durance), ability to					
		of slow twitch (T	Type I) and some fas	t twitch ox	idative (Type	IIa) muscle fibers	at medium	
stimulation freque								
7	40	7	3		5/0.75	240	35	
(5-10)	(30-50)	(5-10)	(0-6)		).5-3)	(160-432)	(25-45)	
Resistance *						althy muscles in o	rder to improve	
						le performance		
			se and prolonged ef		e hypertrophy	. Works through a	ctivation of fast	
	Type IIa) muscle f	ibers at medium s	timulation frequenc	ies.				
7	55	7	3		2/1	320	25	
(5-10)	(45-65)	(5-10)	(0-6)		).5-3)	(160-432)	(20-30)	
Strength*					Stimulation of healthy muscles in order to improve			
						le performance		
Targets improven stimulation freque		ength. Works thro	ough activation of fa	st twitch ox	idative (Type	e IIa) muscle fibers	s at high	
4	80	20	3		2/1	352	20	
(3-8)	(70-100)	(15-40)	(0-6)		).5-3)	(160-432)	(15-30)	
<b>Explosive Streng</b>	th*			Stin	nulation of he	althy muscles in o	rder to improve	
_				or fa	or facilitate muscle performance			
Targets improven very high stimular		strength and powe	er. Works through ac	ctivation of	fast twitch gl	ycolic (Type IIb)	muscle fibers at	
3	120	30	3		2/1	400	15	
(2-8)	(100-150)	(20-80)	(0-6)		).5-3)	(160-432)	(10-20)	
<b>Active Recovery</b>						lood circulation	·	
·					Relaxation of muscle spasms			
Targets accelerati muscle relaxation		very after intensive	e exercise. Increases				ns, provides	
210 (for each	10, 8, 5, 3, 2, 1	N/A	N/A		1/1	208	18	
phase) (120-360)	(1-10)	IV/A	IVA		1/1	(160-432)	(12-36)	
Relaxation Mass	age					plood circulation		
				Rela	Relaxation of muscle spasms			

Targets decrease	in muscular tensio	n. Provides effect	of wellbeing and re	laxati	on.		
420 (for each	7,5,3	N/A	N/A		1/1	208	21
phase)	(1-10)					(160-432)	(15-36)
(300-720)							
Circulation Mas	sage				Increasing local b	lood circulation	
		on and uncomfort	able body tensions	via inc	creasing blood circu		et areas.
120 (for each	1,3,5,7,5,3,1,	N/A	N/A		1/1	208	16
phase)	1-7@5 sec					(160-432)	(12-32)
(90-240)	(1-10)						
De-contraction N					Increasing local b	lood circulation	
De contraction i	1433450				Relaxation of mus		
Enables de-contra	ction of the muscl	les through the cor	nfortable low freque	ency v	vibrations that incre		ulation in the
	l improve oxygena		•	•			
60, 60, 60, (30,	3,2,1,(1-5@5,	N/A	N/A		1/1	208	17
30, 30)x5,	1-3@5, 1)x5,						
(60,30,60,30)x2,	(1-5@5, 1, 1-						
30	3@5, 1)x2,1						
Dual Wave Mass	sage				Increasing local b	lood circulation	
	O				C		
					Symptomatic relie	ef and managemen	nt of chronic,
					intractable pain		
Helps to get rid of heaviness sensation and uncomfortable body tensions via increasing blood circulation in the target areas.							
Combines massag	ge with TENS pain				by isolating stimul		
devices in Duo m	ode.						
60, 60, (30x4,	5, 100, (1-8x4,	N/A	N/A		1/1	208	22
30x4)x5	100x4)x5						

<sup>\*</sup> Marked programs are preceded with 2 minutes Warm Up phase at 5-10 Hz stimulation and 3 minutes of Recovery at 5-1 Hz stimulation using 208 uS pulses.

For each stimulation program, you can adjust stimulation parameters within allowed ranges. Contraction/rest duty cycles can be only adjusted proportionally to the proposed default values.

Adjusted stimulation programs will be saved and always available from the Sessions History screen (for ad hoc sessions) or from the Scheduled Sessions screen (if the session has been previously scheduled).

To identify stimulation parameters for recommended intended use, manufacturer has used stimulation parameters and programs from the number of substantially equivalent NMES and TENS devices, available on the market.

# **Basic Troubleshooting**

#### DEVICE DOESN'T TURN ON

PowerDot PD-02 is probably very low on the battery. Charge your PowerDot for a few hours.

# DURING PRE-STIMULATON SET UP, POWERDOT(S) CANNOT BE FOUND (OR ARE NOT CONNECTED).

Check out Devices screen and make sure that PowerDot device you're trying to use is in the list of active devices

Make sure your PowerDot is turned on (the white light is on in the device).

In some rare cases, you would be required to forget (remove) relevant PowerDot device from Bluetooth Settings (use device serial number to identify the right one).

#### STIMULATION DOES NOT PRODUCE THE USUAL SENSATION

Check that your electrode pads are firmly attached to your body and are correctly positioned (as advised on the Pre-Stimulation screen). Put the stimulation on pause, re-attach or reposition pads, then resume the stimulation.

#### THE STIMULATION CAUSES DISCOMFORT OR A BURNING FEELING

If you're using your standard intensity modes, than most probably your pads are worn out and and/or losing their bonding strength. Pause the stimulation and re-attach your pads firmly, then resume stimulation. If the same sensation continues, pause the stimulation again and replace your pads.

# ELECTRODE PADS DON'T STICK TO THE BODY OR STIMULATION IS SURPRISINGLY WEAK EVEN ON HIGHER INTENSITIES

Replace your electrode pads. Most probably they're worn out. If it doesn't work, check your lead cable for physical damage. If there is any damage, replace the lead cable.

### CAN'T ACTIVATE POWERDOT DEVICE OR IPAD STOPPED CONNECTING TO THE POWERDOT

- Remove relevant PowerDot(s) from Devices menu (use serial number to identify the right one(s)).
- Perform manual factory reset of your PowerDot (or PowerDots) by turning the device off and then holding the button on the device for around 5 seconds, until you see power light blinking several times.
- Go to Bluetooth Settings and remove/forget relevant PowerDot(s) using device's serial number.
- Terminate PowerDot Doctor App using Task Manager.
- Launch PowerDot Doctor App again and try to activate your device again.
- If after all steps above, you still experience connectivity problem, please send us your phone model, OS version and the list of actions you have performed to service@powerdot.com. We will respond within 24 hours.

#### **PowerDot Maintenance**

PowerDot PD-02 device, together with its accessories, should be kept in PowerDot carry case and carefully stored on a secure surface and protected conditions listed in the Warnings above.



Keep replacing your electrode pads after 20 uses as recommended. Deteriorated & worn out pads can cause major discomfort during stimulation, affect the effectiveness and even lead to minor injury.



Cleaning: only clean PowerDot device using a dry soft cloth.



Keep PowerDot device and electrode pads away from water. Store them in a dry place, in protective packaging or in the PowerDot carry case.

PowerDot PD-02 devices do not require calibration or verification of performance parameters. The characteristics are systemically verified and validated for each device manufactured. Those characteristics are stable and do not vary when used under normal conditions.

The manufacturer states that PowerDot cannot be repaired by personnel external to the company. Any work of this nature carried out by personnel not authorized by the manufacturer will be classified as tampering with the unit, releasing the manufacturer from any responsibility with regards to the warranty and hazards that the operator or user may be exposed to.

# **PowerDot Warranty**

PowerDot PD-02 is covered by a worldwide warranty of 2 years, which comes into effect on the date of purchase of the device (proof of purchase is required).

The warranty does not apply to the electrode pads and carry cases. Within the warranty period, manufacturer will replace your faulty PowerDot or accessories at no charge (except shipping & handling fees in some cases), provided that the product:

- Has been used for the intended purpose and in the manner described in this manual
- Has not been connected to an unsuitable power source
- Has not been subjected to misuse or neglect
- Has not been modified or repaired
- Has not been damaged further by shock

Legal rights are not affected by this warranty.

# **Technical Specifications**

All electrical specifications are given for an impedance of 1000  $\Omega$  per channel.

Battery: 2x Lithium Polymer (LiPo) rechargeable 3.7 V, 210 mAh

Charging Input: 5V through USB Type C connection (custom USB charging cable is provided as part of

the package), I/P rating: 5Vdc === 1-2.1A

Stimulation Channels: 2 independent, optically isolated

Stimulation Waveform: Bi-phasic rectangular with zero mean (under load)

Supported Stimulation Frequency Range: 1-150 Hz

**Supported Stimulation Pulse Width:** 32-416 µs (for main/positive phase)

Maximum output voltage/amperage: 130 V/130 mA (+-5%)

Bluetooth: Built-in Bluetooth Low Energy 4.0

Electro-compatibility (EMC): ETSI EN 301 489-1/EN 301 489-17/EN 50385/EN 55011/IEC 60601-1-2

#### C RF Data:

• Operating Frequency Range: 2402 MHz-2480 MHz (ISM range)

Modulation Type: GFSK with AFH

• Peak Transmit Power: -15.86dBm (0.026 mW)

• Channel Spacing/Number of Channels: 2 MHz, 40 channels (3 for advertising, 37 for data)

Antenna Type: PCB Antenna, 2 dBi (1.58 mW) gain

### **Mobile Application Compatibility:**

• Android 5.0 Lollipop (or later) powered smart phone with Bluetooth Smart Ready compatibility and High Definition (or better) touch screen

• Apple iPhone 4S/iPod 5<sup>th</sup> Gen/iPad 3<sup>rd</sup> Gen or newer smart phone/tablet powered by iOS 9.0 (or later)

**Device Dimensions:** 55x55x13.6 mm

Device Weight: 30 g

# **Environment Specifications:**

• **Operating/Storage/Transport:** Temperature from 0 C to +40 C

• **Humidity:** 10-90% RH

Atmospheric pressure: from 700 hPa to 1060 hPa

**Product Expected Lifetime:** 5 years

Housing: ABS & TPU

**Limitations:** product is not suitable for use in the environments with a high concentration of oxygen and/or flammable liquids and/or flammable gas; do not use with equipment for electro surgery or short-wave or microwave therapy; the device may be interfered by other equipment, even if that other equipment complies with CISPR EMISSION requirements.

PowerDot PD-01 has been tested to the compliance with the following Emission and Immunity standards:

### Emission:

STANDARD	ITEM	REMARKS
CISPR 11: 2011	Conducted	Class B
	Radiated	Class B
IEC 61000-3-2:2014	Harmonic current emissions	
IEC 61000-3-3:2013	Voltage fluctuations & flicker	

### IMMUNITY:

STANDARD	ITEM	IEC 60601-1-2 Test Levels for Home Healthcare Environment	PowerDot PD-01 Test Levels	REMARKS
IEC 61000-4-2:2008	ESD	± 8 kV contact; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±2 kV, ± 4 kV, ± 8 kV contact; ±2 kV, ± 4 kV, ±8 kV, ± 15 kV air	No performance degradation observed.
IEC 61000-4-3:2010	RS	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 5.785 GHz 80% AM at 1 kHz	No performance degradation observed.
IEC 61000-4-4:2012	EFT	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	No performance degradation observed.
IEC 61000-4-5:2014	Surge	$\pm 0.5 \text{ kV}, \pm 1 \text{ kV}$	$\pm 0.5 \text{ kV}, \pm 1 \text{ kV}$	No performance degradation observed.
IEC 61000-4-6:2013	CS	3V 0.15 MHz – 80 MHz 6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	10V 0.15 MHz – 80 MHz 80% AM at 1 kHz	No performance degradation observed.
IEC 61000-4-8:2009	PFMF	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz	No performance degradation observed.
IEC 61000-4-11:2004	Voltage dips & voltage variation s	Voltage Dips: 1) 0% U <sub>T</sub> ; 0,5 cycle at 0°. 45°, 90°, 135°, 225°, 270°, 315°  2) 0% U <sub>T</sub> ; 1 cycle; Single phase at 0°  3) 70% U <sub>T</sub> ; 25/30 cycles; Single phase at 0°	As on the previous column	Voltage Dips: 1) No performance degradation observed  2) No performance degradation observed

	<b>Voltage Interruptions:</b> 0% U <sub>T</sub> ; 250/300 cycle;	3) No performance degradation observed
		Voltage Interruptions:
		Performance degradation
		(device stopped charging) has
		been observed only during
		voltage interruption
		testing, but no degradation
		observed after the testing

# **Bluetooth Connectivity**

The unique feature of PowerDot system is that PowerDot stimulators are controlled through Bluetooth Low Energy wireless radio interface.

PowerDot PD-02 is specifically designed to be used together with PowerDot Doctor App, which is supported for the selected Apple iPad devices.

### Secure Pairing

Your PowerDot is paired with your iPad using secure 8 digit numeric code which, by design, prohibits any other mobile phones or wireless devices to connect to your PowerDot. Secure pairing takes place during PowerDot activation process (see *Activating PowerDot* above) and, once your PowerDot becomes active, numeric activation code is written into PowerDot's flash memory and gets verified by your PowerDot Doctor App after every PowerDot restart.

All Bluetooth commands sent from your mobile phone to PowerDot device are securely encrypted using Bluetooth AES-128 encryption protocol.

### Disconnections and Quality of Service

PowerDot Doctor App and PowerDot PD-02 Bluetooth communication interface are specifically designed to accommodate temporarily and permanent Bluetooth disconnections during a stimulation session.

PowerDot PD-02 device is capable of independent execution of a pre-loaded stimulation program with the latest intensity values as well as implements automatic Bluetooth re-connections.

In this regard, temporary radio frequency interference (e.g. caused by co-existence of multiple Bluetooth and/or Wi-Fi devices in your range) should not affect the overall efficiency and safety of your stimulation session.

Due to hardware-level emergency stop mechanisms (see *Directions For Use*), Bluetooth disconnections of more permanent nature should not affect the safety of stimulation, and can only cause temporary inconvenience by forcing you to postpone your planned stimulation session until a more favorable Bluetooth connectivity environment is established.

Like any wireless device, PowerDot PD-02 emits very low levels in the radio frequency (RF) interval, and, is therefore not likely to cause any interference with nearby electronic equipment (e.g. radios, computers, telephones, etc.).

PowerDot PD-02 is designed to withstand foreseeable disturbances originating from electrostatic discharges, mains supply magnetic fields, or radio frequency transmitters.

Despite this, it is not possible to guarantee that the stimulator will not be affected by strong RF (radio frequency) fields emitting from other sources (such as in the proximity of working microwave oven).



Try not to use PowerDot closer than 1.5 meters to the working microwave oven as radio interference from microwave is likely to cause disconnection between PowerDot and your mobile phone.

### Troubleshooting Wireless Connectivity

If you run into issues with Bluetooth wireless connectivity (e.g. your PowerDot PD-02 device becomes unresponsive to PowerDot Doctor App commands during stimulation session or you were not able to connect to your PowerDot and initiate stimulation), do not panic and consider terminating your stimulation session manually by shortly pressing Power button on your PowerDot device.

PowerDot Doctor App has built-in re-connection and disconnection detection mechanism and, in most cases, it will re-connect back to your PowerDot shortly and allow you to resume your stimulation using Resume button on the stimulation screen.

If you fail to re-connect and resume stimulation after several attempts, consider stopping your stimulation session using Stop button on the Stimulation Screen and postponing it for later.

# **Used Symbols**

SN Serial Number

(I) Stand by

/! Attention

Direct Current (DC)

# **IMPORTANT NOTE:** (For Portable Device Configuration)

FCC Radiation Exposure Statement:

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment for bodyworn configuration in direct contact to the phantom.

The device complies with Part 15 of FCC. Operation is subject to the following conditions:

a. This device may not cause harmful interference, and

b. This device must accept any interference received, including interference that may cause undesired operation.

FCC ID: 2AC2KPD001201400SMT

# **EU Only Symbols:**

( <del>(</del>

This symbol on your PowerDot unit is to indicate conformity with the requirements of the Medical Device Directive (94/42/EEC)



Manufacturer



EU Authorized representative



Internally powered device Class II with Type BF applied parts



Product subject to WEEE regulations concerning separate waste collection



Read the instructions for use carefully before using this device

IP22

IP Rating IP22

### **Manufacturer & After-Sale Service:**

Smartmissmo Technologies Pte Ltd

4 Shenton Way, #28-01 SGX Centre II

Singapore 068807

E-Mail: service@powerdot.com

Phone: +1-844-479-7368

Contact for any assistance in setting up, using, maintaining, or reporting unexpected operation or events.

#### EU AUTHORISED REPRESENTATIVE:

Medical Technology Promedt Consulting GbmH,

Altenhostrasse 80, 66386, St. Ingbert, Germany

# **Electromagnetic Compatibility (EMC)**

PowerDot PD-02 is designed to be used in home healthcare environments in accordance with the EMC safety standard IEC 60601-1-2 (4<sup>th</sup> Edition) and with limitations, defined by the warnings and precautions in this manual (e.g. operation near RFID emitters, working microwave ovens, etc.).

Examples of home healthcare environment include restaurants, cafes, shops, stores, markets, schools, churches, libraries, outdoors (streets, sidewalks, parks), domiciles (residences, homes, nursing homes), vehicles (cars, buses, trains, boats, planes, helicopters), train stations, bus stations, airports, hotels, hostels, pensions, museums, theatres.


PowerDot PD-02 is designed to support anticipated disturbance originating from electrostatic discharge, magnetic fields for the power supply or radiofrequency emitters.


However, the performance of PowerDot PD-02 device can still affected by radio frequency fields originating from other sources.

For more information about EMC emissions and immunity, contact the manufacturer.
















PowerDot PD-02 battery charging performance might degrade in the environments

with frequent voltage interruptions (when charging from the wall adapter). To support consistent and reliable charging, the usage of uninterrupted power supply (UPS) is highly recommended, if operating in such environments.

The following device function is considered essential to the safety of the user: ability to maintain consistent stimulation intensity (amplitude), pulse frequency and pulse waveform (both shape and width).

In case if the essential performance is lost or degraded due to electromagnetic disturbances, stimulation safety and effectiveness can be compromised. Whenever the patient realizes unexpected change in any of stimulation parameters, it's is advised to terminate the stimulation session immediately by using one of the methods provided in the Terminating Stimulation section.