



USER MANUAL HHM-A001

LIST OF REVISIONS

Version	Date	Status / Change description
01	2022-03-06	Initial version

1 ABOUT THIS MANUAL

Updates to this manual will be made available continuously, and news about and/or updated content will be posted on the Sonitor partner portal at <http://partner.sonitor.com>. Please contact Sonitor Technologies for access permission information.



Figure: Photo of the The SenseClean™ Hand Hygiene Module, HHM-A001

2 DISPENSER COMPATIBILITY

The SenseClean™ Hand Hygiene Module can be used with the following GOJO dispensers:

- 1930V3-04 – PURELL LTX-12 SMARTLINK Ready Wht/Wht
- 1930V3-04-CHR – PURELL LTX-12 SMARTLINK Ready Blk/Chr
- 1931V3-04 – PROVON LTX-12 SMARTLINK Ready Wht/Gry
- 1931V3-04-CHR – PROVON LTX-12 SMARTLINK Ready Blk/Chr

3 MOUNTING

Mounting of the SenseClean™ Hand Hygiene Module is to be performed by a representative of the dispenser manufacturer. The module shall be plugged into one of the compatible dispensers as described in this user manual. The plugging is a simple procedure that does not demand any training. The dispenser itself can remain mounted, e.g. on a wall, during the procedure.

The following step-by-step instruction should be followed.

3.1 Open the GOJO dispenser

Press the door latches inward to release the door and pull the door forward.

3.2 Orient the module and plug the cable

The module is to be mounted with the label facing the person doing the assembly. This will make the symbol with the hand and the droplet visible through the dispenser.

Plug the module into the mechanical guiding structures on the inside of the dispenser and plug the cable to the input of the dispenser electronics.

3.3 Close and inspect

Close the lid and check that the module is visible through the dispenser. Place a hand below the dispenser to check that the LED blinks.



Figure: Illustration of how the dispenser with a module will appear when correctly mounted.

4 FCC NOTE

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to the equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation.

5 ISSED NOTE

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'ISDE Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes :

- (1) l'appareil ne doit pas produire de brouillage, et
- (2) l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

6 RECYCLING

Please dispose the electronics in line with local environmental requirements.

For EU/EEA countries, all Sonitor products are subject to the European Waste Electrical and Electronic equipment (WEEE) Directive 2012/19/EU, the Battery Directive 2006/66/EC, and the Packaging Directive 94/62/EC in EU/EEA countries. The party responsible to comply with the disposal regulations of the Directive is the party specified as liable according to the legislation in the country where the end user is located.

7 LABELLING

It is a requirement that a sanitizer in which an HHM-A001 module is installed is labelled with the following information:

CONTAINS FCC ID: 2AD7T21122013101

CONTAINS IC: 20330-21122013101

Labelling must follow regulatory requirements for placement and integrity.

The content of this document is the property of Sonitor Technologies AS.
Copying and distribution is not permitted without written permission from Sonitor Technologies AS.