

Linkhealth Smart Thermometer user manual

Linkhealth Corporation



Linkhealth smart thermometer, a wearable health care equipment, is a collection of temperature, heart rate and sports measurement.

Double taps the equipment, begin new experience.



Download linkhealth APP

1. Download APP

Visit Apple store or Google Play , search "linkhealth" or land website: <http://app.linkhealth.me/> , download linkhealth APP.


2. Connect Smart Thermometer

Double taps the equipment to boot the thermometer, open APP, click "menu" , select "Connect Thermometer" to search for equipment, find "LINKHEALTH" in the available equipment list; click on icon or name for connection. (Once connect the equipment, will show the thermometer information.)

3. The component of product



4. Product Wear



Please pay attention to avoid wearing the thermometer in a wrong way. You need to make sure that the sensor should be in the armpit.

5. Operation

Thermometer boot: Double taps the thermometer screen.

Thermometer Power off: Not wear state and no operation,the thermometer will power off in 15 minutes.

Double taps will switch display:

Temperature, Time , Heart Rate ,Power.

(The screen will turn off in 3 seconds, double taps display.)

6. The Common Failures and Countermeasures

Failures	Reasons	Solution
Time error	Not connect App	Download and open Linkhealth APP, and connect the thermometer.
No response on double taps	The battery is too low	To charge the equipment.
Charging failure(charging icon is not shown)	Poor contact charging interface;	1.To connect the charging cable again; 2.To clean the charging terminal interface and connect again;

Warranty Card

Product Model:

Customer Name:

Purchase Date:

Adress:

Store Name:

Zip Code:

Personal Address:

Tel:

Seal of Supplier



FCC Warning

Caution: Any changes or modifications in construction of this device which are not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device meets the FCC requirements for RF exposure in public or uncontrolled environments.

Linkhealth Corporation

Address : 472 Sandhurst Dr Milpitas, CA 95035, USA

E-mail: service@linkhealth.me