

CC-418P TRANSMITTER

USERS MANUAL

Functional Description:

The CC-418P Radio Frequency (RF) transmitter is designed to operate with all Digital Care Systems, Inc. products. Its small size and long battery life make it an ideal product for Emergency Call (E-call) systems. The CC-418P can be worn as a pendant transmitter, wrist transmitter, belt clip transmitter or key fob transmitter. Activation of the transmitter is via the button on the front of the unit indicated with three brail bumps. The CC-418P is a water proof device and can be worn in the shower or bathtub. However, it will **not** work underwater. Each CC-418P transmits a unique code number and will only activate the device for which it is trained to.

Operation:

Pressing the button on the front of the unit, indicated by the three-brail bumps will cause the CC-418P to transmit. The button must be held for at least one half (1/2) second before transmission will begin. This helps prevent accidental activation. Transmission will be indicated by the RED light next to the activating button. The duration of the transmission is limited to a 5 second period, regardless of how long the button is held. The button must be released and pressed again to cause a second transmission. This feature saves battery life while still assuring adequate transmission time to be received by the receiving device.

Limited Warranty and Disclaimer:

Digital Care Systems Inc. warrants the CC-418P transmitter for a period of one (1) year from the date of purchase to be free from defects in materials and workmanship. Liability shall be limited to the repair or replacement of the equipment. Any consequential damages arising from the use of this product are the sole responsibility of the user. This warranty does not include batteries. Tampering with this equipment will void this warranty and the user's authority to operate it.

This warranty is in lieu of all other expressed warranties, obligations, or liabilities. There are no warranties, which extend beyond the description on the face hereof. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING, BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE WARRANTY PERIOD OF THIS WRITTEN LIMITED WARRANTY. Some states do not allow limitation on how long an implied warranty lasts, so this limitation may not apply to all purchasers. Digital Care Systems, Inc's LIABILITY SHALL BE LIMITED TO THE COST OF THE PRODUCT AND IN NO EVENT SHALL Digital Care Systems, Inc. BE LIABLE FOR ANY SPECIAL, INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR THE INABILITY TO USE THIS PRODUCT OR FOR THE BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESSED OR IMPLIED, WHATSOEVER. Some states do not allow this exclusion, so this limitation may not apply to all purchasers.

FCC/Canadian Information:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an output on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- A shielded cable must be used when connecting a peripheral to the serial ports.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference. And (2) this device must accept any interference received, including interference that may cause undesired operation.

This Class B digital apparatus complies with Canadian ICES-003.

Cet appareil numérique de la classe B est conforme à la norme NMB-003 du Canada.

This device complies with Industry Canada's license-exempt RSSs. Operation is subject to the following two conditions:

(1) This device may not cause interference; and

(2) This device must accept any interference, including interference that may cause undesired operation of the device.

Digital Care Systems, Inc., 2000C Hartel Ave. Levittown, PA 19057
(215) 946-7700