

User Manual for Yuejing

Introduction

"Yuejing" smart period tracker is used to measure women's basal body temperature. And the temperature data will be sent to the app "Yuejing" via the bluetooth. The app will predict women's period according to their basal body temperature.

Mobile version

Android: Android 4.3 and above

iOS: iPhone 4s and above

Using Steps

1. Connecting the device with the app

- (1) Download the app "Yuejing" through your mobile application market or through the website: www.yuejingapp.com.
- (2) Enter into the app "Yuejing", click the button "non-connected".
- (3) Press the switch on the Yuejing smart device, and press the button "connect right now" on the app. Make use your phone's bluetooth is opened. And then the app will search and connect the device in 30 seconds.

2. Measuring steps

- (1) Start to record temperature. After finished the connecting steps above, press the button "record the temperature". Put Yuejing smart device under your armpit and make sure the temperature sensor to stick close to your skin.
- (2) Finish the record process. It takes about 5 minutes to reach a final stable temperature, and the app will shake to indicate that the measuring process is finished. Then device light will turn off in 30 seconds.

3. Notes

- (1) Make sure to measure your body temperature when you wake up in the morning before drinking water or doing any other exercise.
- (2) Make sure to connect your smart device with the app "Yuejing".
- (3) Please keep the device clean and dry.

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

Note: The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate this equipment. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

RF exposure warning

This device is measured RF output power is less than the SAR exclusion threshold value for human head and body. The equipment must not be co-located or operating in conjunction with any other antenna or transmitter.