



USER MANUAL

Stealth Products

Alternative Drive Controls



i-Drive[®] 4.0

Stealth's User Manual and Maintenance Guide for i-Drive 4.0
Alternative Control with *Bluetooth*[®] Wireless Technology



Stealth Products strives for 100% customer satisfaction. Your complete satisfaction is important. Please contact us with feedback or suggested changes that will help improve the quality and usability of our products.

You may reach us at:



Stealth Products, LLC

104 John Kelly Drive, Burnet, TX 78611

Phone: (512) 715-9995


Toll Free: 1(800) 965-9229

Fax: (512) 715-9954

Toll Free: 1(800) 806-1225

info@stealthproducts.com

www.stealthproducts.com



MDSS GmbH

Schiffgraben 41

30175 Hannover, Germany



General

Read and understand all instructions prior to the use of the product. Failure to adhere to instructions and warnings in this document may result in property damage, injury, or death. Product misuse due to failure of the following instructions will void the warranty.

Immediately discontinue use if any function is compromised, parts are missing, loose, or shows signs of excessive wear. Consult with your supplier for repair, adjustment, or replacement.

If this document contains information you do not understand, or there are concerns about safety or operation, contact your supplier.

Google, the *Google Logo* and *Google Play Logo* are registered trademarks of Google Inc. *Apple Logo* and the *AppStore Logo* are registered trademarks of Apple Inc.

The *Bluetooth®* word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Stealth Products, LLC is under license. Other trademarks and trade names are those of their respective owners.

Important Information!

All persons responsible for fitting, adjustment, and daily use of the devices discussed in these instructions must be familiar with and understand all safety aspects of the devices mentioned. In order for our products to be used successfully, you must:

- Read and understand all instructions and warnings.
- Maintain our products according to our instructions on care and maintenance.
- Devices should be installed and adjusted by a trained technician.

Supplier Reference

Supplier: _____

Telephone: _____

Address: _____

Purchase Date: _____

Model: _____

Before you install or begin using this product, it is important that you read and understand the content of these installation and operating instructions. The installation instructions will guide you through the options and possibilities with the product.

Instructions are written with the expressed intent of use with standard configurations. They also contain important safety and maintenance information, as well as describe possible problems that can arise during use. For further assistance, or more advanced applications, please contact your supplier or Stealth Products at (512) 715-9995 or toll free at 1-800-965-9229.

Always keep the operating instructions in a safe place so they may be referenced as necessary.

All information, pictures, illustrations, and specifications are based on the product information that was available at the time of printing. Pictures and illustrations shown in these instructions are representative examples and are not intended to be exact depictions of the various parts of the product.

NOTICE

Instructions are designed with the intended purpose of use with the standard configurations of i-Drive. Further assistance or more advanced steps in setup, please contact Stealth Products, LLC.

Ordering Documentation

You can download additional copies of this user manual on the Stealth website:
www.stealthproducts.com

Or search:

i-Drive 4.0 User Manual in the search bar at the top of the page.

Our products are designed, manufactured, and produced to the highest of standards. If any defect in material or workmanship is found, Stealth Products will repair or replace the product at our discretion. Any implied warranty, including the implied warranties of merchantability and fitness for a particular purpose, shall not extend beyond the duration of this warranty. Stealth Products, LLC does not warrant damage due to, but not limited to:

- Misuse, abuse, or misapplication of products.
- Modification or product without written approval from Stealth Products, LLC.
- Any alteration or lack of serial number, where applicable, will automatically void this warranty.
- Stealth Products, LLC is liable for replacement parts only.
- Stealth Products, LLC is not liable for any incurred labor costs.

No person is authorized to alter, extend, or waive the warranties of Stealth Products, LLC.

Stealth Products warrants against failure due to defective materials or workmanship:

Covers: 180 days

Hardware: 5 years

Electronics: 3 years

In Case of Product Failure

In the event of product failure covered by our warranty, please follow the procedures outlined below:

1. Call Stealth at +1 (512) 715-9995 or toll free +1-800-965-9229.
2. Request the Returns Department or obtain an RA from the Returns Department and follow department or documentation instructions.

Federal Communications Commission (FCC) Statement (USA)

FCC ID: 2AJXV_2AJXVIDHBT500

FCC RF Exposure Statement

This device meets the FCC requirements for RF exposure in public or uncontrolled environments.

This device complies with FCC part 15 FCC part 15 FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

FCC Warning

Changes or modifications not approved by Stealth Products, LLC could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help

1.0 Customer Satisfaction.....	i
2.0 Important Information	ii
3.0 Introduction	iii
4.0 Warranty	iv
5.0 FCC Statement	v
6.0 Table Of Contents	vi
7.0 Warning Labels	ix
7.1 Warning Labels	ix
7.2 Limited Liability	ix
7.3 Testing.....	ix
8.0 Design and Function	1
8.1 Intended Use	1
8.2 Features	1
8.3 Factory Settings.....	2
9.0 Parts and Accessories	3
9.1 i-Drive Packages Available.....	3
9.2 Replacement/Additional Parts Available	3
10.0 R-Net™ Installation.....	4
10.1 R-Net™ Omni Setup.....	4
Set up	5
11.0 Q-Logic™ Installation	6
11.1 Q-Logic™ Setup	6
Setup.....	7
11.2 Q-Logic™ Drive Control System Further Configuration	8
Back Pad Toggle.....	8
Side Pad Toggle/Control	8
Egg Switch Function	9

12.0 Seating Functions	10
12.1 Programming Seating Functions	10
13.0 Advanced Programming Software	11
13.1 i-Drive Advanced Programming Software	11
13.2 Getting Started.....	11
13.3 Downloading the i-Drive Software to Your PC.....	12
13.4 Connecting to the i-Drive	13
13.5 Running The Setup Wizard	14
13.6 i-Drive Configuration Settings.....	15
13.7 Configuration: Sensor Channels	15
13.8 Configuration: Double Tap Settings.....	17
13.9 Running Real Time Diagnostics	18
13.1.0 Help Feature	19
13.1.1 Closing the i-Drive Programming Software	20
14.0 i-Drive Mobile App Installation	21
14.1 Downloading the i-Drive Bluetooth® Application.....	21
For iPhone Users	21
For Android Users.....	22
14.2 Opening the i-Drive Bluetooth® App	23
14.3 Config Mode	24
Channel Assignment.....	24
Channel Assignment Settings.....	25
Multi-Tap Setting.....	26
Diagnostics.....	27

15.0 Testing..... 28

15.1 Operational Test..... 28

15.2 Test Drive 28

15.3 Stop Test 28

16.0 First Time Use29

16.1 Dealer Assistance..... 29

16.2 User Testing29

16.3 Conditions Of Use29

17.0 Maintenance.....30

17.1 Care And Maintenance 30

17.2 Safety 30

17.3 Technical Data.....30

Warning Labels

7.1

Warnings are included for the safety of the user, client, operator and property. Please read and understand what the signal words **SAFETY, NOTICE, CAUTION, WARNING** and **DANGER** mean, how they could affect the user, those around the user, and property.

 DANGER	Identifies an imminent situation which (if not avoided) will result in severe injury, death, and property damage .
 WARNING	Identifies a potential situation which (if not avoided) will result in severe injury, death, and property damage .
 CAUTION	Identifies a potential situation which (if not avoided) will result in minor to moderate injury, and property damage .
NOTICE	Identifies important information not related to injury, but possible property damage .
SAFETY	Indicates steps or instructions for safe practices, reminders of safe procedures, or important safety equipment that may be necessary.

Limited Liability

7.2

Stealth Products, LLC accepts no liability for personal injury or damage to property that may arise from the failure of the user or other persons to follow the recommendations, warnings, and instructions in this manual.

Testing

7.3

Initial setup and driving should be done in an open area free of obstacles until the user is fully capable of driving safely.

The i-Drive should always be programmed and tested without any person sitting in the wheelchair until every alteration of the physical installation or adjustment is complete.

Intended Use


8.1

The i-Drive 4.0 with *Bluetooth*® wireless technology provides the ultimate head support and has a sophisticated setup within its technology that allows the end user to program the device in diverse ways. An integrated circuit board and microprocessor gives the new i-Drive 4.0 this *Bluetooth*® wireless technology capability.

Features

8.2

- *Bluetooth*® Programming App with full wireless programming capabilities.
- Smart Input Connections-sensors and mechanical switches can now be placed at any port.
- Easy interface exchange capability
- Waterproof interface box




NOTICE

Existing dongles will not work with the new 4.0 interface, however 4.0 dongles will work with the 3.0 interface.

Factory Settings

8.3



The i-Drive 4.0 has the capability of taking inputs from mechanical switches, proximity sensors, fiber optic sensors, joysticks, and Sip & Puff modules. With the new i-Drive interface, there is the capability of programming the interface with the associated Advanced Programming Software from your PC, and now also includes the option to download the i-Drive App and program the interface directly from your phone via *Bluetooth*[®] wireless technology.

NOTICE

The software will not read correctly if a non *Bluetooth*[®] dongle is plugged into a port . Any dongle plugged into the interface **MUST** be *Bluetooth*[®] .

NOTICE

All connectors have blue heat shrink on the cord for visual reference and easy recognition.

i-Drive 4.0 Packages Available**9.1**

Part Number	Description
IDHBT200-1	i-Drive Tri-Array Head Support, includes 3 Proximity Sensors, non-proportional, with egg switch for reset mode change and i-Drive Bluetooth® Interface
IDHBT210-1	i-Drive Tri-Array Head Support with Sip & Puff, includes 3 Proximity Sensors, with egg switch for reset mode and i-Drive Bluetooth® Interface
IDHBT250-1	i-Drive Pediatric Tri-Array Head Support, includes 3 Proximity Sensors, non-proportional with egg switch for reset mode change and i-Drive Bluetooth® Interface
IDHBT260-1	i-Drive Pediatric Tri-Array Head Support with Sip & Puff, includes 3 Proximity Sensors, non-proportional, with egg switch for reset mode change and i-Drive Bluetooth® Interface
IDHBT300-1	i-Drive Pro Series Head Support, includes 3 Proximity Sensors, occipital pad, swing away hardware for left and right side, 2 spot pads with egg switch for reset mode change and i-Drive Bluetooth® Interface
IDHBT310-1	i-Drive Pro Series Head Support with Sip & Puff, includes 3 Proximity Sensors, non-proportional, with egg switch for reset mode change and i-Drive Bluetooth® Interface
IDHBT350-1	i-Drive Pediatric Pro Series Head Support, includes 3 Proximity Sensors, non-proportional, with egg switch for reset mode change and i-Drive Bluetooth® Interface
IDHBT360-1	i-Drive Pediatric Pro Series Head Array with Sip & Puff, includes 3 Proximity Sensors, non-proportional, with egg switch for reset mode change and i-Drive Bluetooth® Interface

Replacement/Additional Parts Available**9.2**


Part Number	Description
IDHBT500-17MM*	i-Drive Bluetooth® Interface for 17mm Ball Mount
IDHBT500-1*	i-Drive Bluetooth® Interface for 1" Ball Mount
IDHBT500*	i-Drive Bluetooth® Interface

* These packages include 2 mono plug dongles

The i-Drive plugs into the 9 pin port on the power chair’s control interface. The location of the port depends on the manufacturer of the interface. Follow the instructions for your interface manufacturer in these next sections.

R-Net™ Omni Setup

10.1



Mode

Switch mode to drive, power seating, or OBP (onboard programming)

Navigation

Up & Down button to navigate the menus, Right button to open a menu item and Left button to return to previous menu.

Profile


Switch between preset drive profiles and activate device as control.

Power

Turns chair ON/OFF.

Plus & Minus

Toggle settings or changes value of the highlighted parameter.



NOTICE

An R-Net™ Programming Dongle is required.

Setup

1. Power off the chair with the power button on the joystick or display.



DANGER

Disengage the motors before programming the chair.

2. Connect the R-NET programming dongle into the system.
3. Ensure all the cords for the display are properly connected to the chair.
4. Connect i-Drive Alternative Control into *Port 1* or *Port 2* on the bottom of the display with 9 pin connector.
5. With the chair powered off, open up the app on a Bluetooth enabled device and log in using the email and password provided to you. Power the chair on and connect to the application.
6. Press the *Mode 'M'* button until you enter the *OBP* menu.
7. In the *OBP* menu, navigate with the *Up & Down* buttons  to the *Omni* menu option and select with the *Right* button. 
8. In the *Omni* menu, navigate with the *Up & Down* buttons  to *Global* on the menu and select with the *Right* button. 
9. In the *Global* menu, navigate with the *Up & Down* buttons  to *Sleep 12V* then toggle to *On* with the *Right* button. 
10. Press the *Left* button  until you return to the *Omni* menu.
11. In the *Omni* menu, navigate with the *Up & Down* buttons  to *Profiled* and select with the *Right* button.  Ensure your chair is programmed for the port your i-Drive is connected to.
12. In the *Omni* menu, navigate to the port you have the i-Drive connected to (*Port 1* or *Port 2*).
13. In the *Port* menu, navigate with the *Up & Down* buttons  to *SID*, toggle to *Prp*, and with the *Right* button,  select. *For 3 switch head array applications, set SID to 3-Prp.*
14. Press the *Left* button  to back out of the *OBP* menu. Scroll down to *Systems* and select with the *Right* button. 
15. Scroll down to *Joystick Calibration* and with the *Right* button,  select. After you select *Joystick Calibration*, toggle with the *Up & Down* buttons  to the port you have the i-Drive plugged in to. After you select, your chair will automatically go in to calibration mode. The chair display will ask you to move the chair joystick forward.
16. Go to the i-Drive application on the Phone/Tablet/PC Device. Select *Config Chair* tab. Select *Calibrate Chair*. Click *Okay*.
17. Wait until the calibration finalizes. Your chair display will show a green checkmark meaning calibrations was a success.
18. On the *Omni* display you will need to exit the *OBP* menu and power the chair off. Now you will be able to set up the chair as required.
19. Check your configuration on the *Diagnostics* tab. If all is working correctly, your i-Drive is ready for use. Re-engage your motors.

Q-Logic™ Setup

11.1

Navigation

Arrows navigate the main menu. *Up & Down* to navigate the menus, *Right* to open a menu item and *Left* to return to the previous menu.

Bookmark Buttons


Select *Options* in the main menu. Button actions are displayed on screen above corresponding buttons. In other menus, hold button to *Bookmark Settings*, and press button to quickly go to bookmark.

Plus & Minus


Toggle settings or changes values of the highlighted parameter.

Help Button

Displays information about options selected on screen.



Q-Logic Enhanced Display



NOTE: Q-Logic™ Enhanced Display and Q-Logic™ Handheld Programmer required.

Setup

1. Power off the chair with the power button on the joystick or display.



DANGER

Disengage the motors before programming the chair.

2. Ensure all the cords for the display are properly connected to the chair.
3. Connect the i-Drive Alternative Control into the port on the bottom of the display with the 9-pin connector.
4. With the chair powered off, open up the app on a Bluetooth enabled device and log in using the email and password provided to you. Power the chair on and connect to the application.
5. Plug in the Q-Logic Programmer to the back of the Q-Logic Display.
6. On the Programmer, select the *Program Adjustment* option.
7. In the *Program Adjustment* menu, navigate with the *Up & Down* buttons $\updownarrow$ to the *Specialty Control* option and with the *Right* button, $\rightarrow$ select.
8. In the *Specialty Control* menu, navigate with the *Up & Down* buttons $\updownarrow$ to *Active Device* and toggle with the *Plus & Minus* buttons $\pm$ to *Proportional Head 3-Direction*.
9. With the *Up & Down* buttons, $\updownarrow$ scroll down to *Proportional Head-3 Direction* folder. With the *Right* button, $\rightarrow$ select and enter folder. Scroll down to *Joystick Setup* folder and select.
10. In the *Joystick Setup* menu, with the *Up & Down* buttons, $\updownarrow$ scroll down and select *Calibrate*. Select *Okay*.
11. The Q-Logic programmer will read, *Roll around joystick in two full circles*.
12. Go to the i-Drive application on the Phone/Tablet/PC Device. Select *Config Chair* tab. Select *Calibrate Chair*. Click *Okay*.
13. Wait until the calibration finalizes, then select *Okay* with the Q-Logic programmer and verify the calibration was successful. (If calibration was not successful, power off the chair and repeat *Steps 5-12*.)
14. Return to the Main Menu and exit the Q-Logic programming.
15. Unplug the Q-Logic handheld programmer.
16. Power chair off.



WARNING

Powering off the chair completely before operation is mandatory.

17. The chair is now programmed to recognize the i-Drive. Power the chair back on.
18. Check your configuration on the diagnostics tab. If all is working correctly, your i-Drive is ready for use. Re-engage your motors.

Q-Logic™ Drive Control System Further Configuration

11.2

Back Pad Toggle

If you want to disable or enable the back pad's ability to toggle into forward and reverse:

1. Reconnect the Q-Logic Handheld Programmer. Follow *Steps 1-7* from *Setup*.
2. On the Q-Logic Programmer, select the *Program Adjustments* option.
3. In the *Program Adjustments* menu, navigate with the *Up & Down* buttons Δ to *Specialty Control*, and with the *Right* button, $\triangleright$ select.
4. In the *Specialty Control* menu, navigate with the *Up & Down* buttons Δ to *Proportional Head 3-Direction*, and with the *Right* button, $\triangleright$ select.
5. In the *3 Switch Head* menu, navigate with the *Up & Down* buttons Δ to *Back Switch Toggle* and select with the *Right* navigation button. $\triangleright$
6. In the *Back Switch Toggle* menu, navigate with the *Up & Down* buttons Δ to *m: Back Toggle* and with the *Plus & Minus* buttons $\pm$, toggle to *Enable*. This will implement the *forward/reverse* toggle with the back pad, or with the *Plus & Minus* buttons $\pm$, toggle to *Disable* to disable toggle.
7. After configuring the back pad, disconnect the Q-Logic™ Programmer and restart the chair's electronics (follow *Steps 10-13*), or continue to configure other operations.

Side Pad Toggle/Control

To disable or enable *Mode* switching with the side pad:

1. Reconnect the Q-Logic Handheld Programmer (follow *Steps 1-7* from *Setup*).
2. On the Q-Logic Programmer, select the *Program Adjustments* option.
3. In the *Program Adjustments* menu, navigate with the *Up & Down* buttons Δ to *Specialty Control* and with the *Right* button, $\triangleright$ select.
4. In the *Specialty Control* menu, navigate with the *Up & Down* buttons Δ to *Proportional Head 3-Direction* and with the *Right* navigation button, $\triangleright$ select.
5. In the *3 Switch Head* menu, navigate with the *Up & Down* buttons Δ to *Device Options/Timing* and select with the *Right* button. $\triangleright$
6. In the *Device Options/Timing* menu, with the *Up & Down* buttons, Δ navigate to *m:Device Double Command* and with the *Plus & Minus* buttons $\pm$ toggle to *Enable*. This will set up the mode switch with the side pad, or with the *Plus & Minus* buttons $\pm$ toggle to *Disable* to disable the mode switch.
7. After configuring the back pad, disconnect the Q-Logic™ Programmer and restart the chair's electronics (follow *Steps 10-13* from *Setup*), or continue to configure other operations.

Egg Switch Function

If you would like to change the function of the egg switch to *Mode* or *Toggle Forward* and *Reverse*:

1. Reconnect the Q-Logic Handheld Programmer (follow *Steps 1-7* from *Setup*).
2. On the Q-Logic Programmer, select *Program Adjustments* option.
3. In the *Program Adjustments* menu, with the *Up & Down* buttons, Δ navigate to *Specialty Control* and select with the *Right* button. $\triangleright$
4. In the *Specialty Control* menu, navigate with the *Up & Down* buttons Δ to *Proportional Head 3-Direction* and with the *Right* button, $\triangleright$ select.
5. In the *Proportional Head 3-Direction* menu, navigate to *Switch Options/Timing* with the *Up & Down* buttons Δ and select with the *Right* button. $\triangleright$
6. In the *Switch Options/Timing* menu, navigate with the *Up & Down* buttons Δ to *m:Mode Jack Switch Type* and toggle with the *Plus & Minus* buttons $\pm$ to *Mode* to enable mode switch with the egg switch, or *Toggle* for forward/reverse toggle.
7. After configuring the back pad, disconnect the Q-Logic™ Programmer and restart the chair's electronics (follow *Steps 10-13* from *Setup*) , or continue to configure other operations.

Programming Seating Functions

12.1

R-Net

To setup the chair's seating functions, follow the instructions below:

1. Reconnect the R-Net programming dongle.
2. On the Omni display, enter into the *OBP* menu.
3. In the *OBP* menu, navigate with the *UP & Down* buttons  and select *Omni* with the *Right* button. 
4. In the *Omni* menu, select *Port 1* or *Port 2* with the *Right* button. 
5. After selecting the port, select *Controls*.
6. Once in the *Controls* menu, select *Actuator Axis* and switch to either *L/R* or *R/L*.

NOTICE

The left/right pads on the headrest control the seating. One side will toggle between seat functions while the other side will activate the function and toggle directions.

Q-Logic

When the option of *Proportional Head-3 Direction* is chosen during the setup of the chair, seating functions will automatically default to *left/right* and the back pad will be deactivated. By double tapping the left pad, mode settings can be chosen.

No further programming will be needed.

i-Drive Advanced Programming Software

13.1

The i-Drive's Advanced Programming Software is not required to operate the drive control, but is available to ensure a customized fit and driving experience. Programming the i-Drive is simple. By connecting a tablet or smartphone and using the i-Drive Advanced Programming Software application you can:

- Dynamic Channel Assignment— any port can control any of the outputs
- Adjust Activation Proximity— Change the range of motion required for switch activation
- Provide an adjustable double tap feature that extends chair double tap limits
- Configuration— Programmed by Bluetooth or PC

Getting Started

13.2

The first step to begin use of the i-Drive 4.0 is downloading the software. The link to download the software is listed on pg. 11. If there are any difficulties obtaining the software, contact Stealth Products, LLC.

After the software has been installed onto your device, you can begin utilizing the Advanced Programming Software. Launch the "i-Drive Programmer" application on your device, and with a USB to Mini-USB cable, connect the i-Drive to your PC or Tablet. Allow time for the i-Drive to make a full connection with your computer.

NOTICE

Sections 12.3 through 12.1.1 will guide you on the setup and configuration processes on your PC or tablet for the i-Drive.

WARNING

Please manually disengage both drive motors before utilizing the *Setup Wizard*, *Diagnostics*, and *Config*. If motors are not disengaged, the chair could move during setup.

Downloading the i-Drive Software

13.3

To program the i-Drive from your PC, the Advanced Programming Software™ will need to be downloaded from the Stealth Products website:

<https://stealthproducts.com/idrive/download/>



 WARNING

To ensure security on your i-Drive, contact Stealth Products to obtain your user name and password for the software.

Connecting To The i-Drive

13.4

After adequate time is given to allow the i-Drive to connect to your computer and you have launched the application, click the *Connect* button.



NOTICE

If may be required to change the *COM* in the drop-down menu next to the *Connect* button.

NOTICE

When the i-Drive is connected to the software, the serial number of the interface, the firmware version of the software, and the system voltage should appear on the bottom of the screen.


Running The Setup Wizard

13.5

The *Setup Wizard* will guide you through the initial configuration and channel assignment of the i-Drive. To begin the *Setup Wizard*, click the *Setup Wizard* button on the left of the screen.

Setup Wizard Feature

The *Setup Wizard* will guide you through the initial port configuration and channel assignment of the i-Drive. A window will pop-up on the screen and walk you through a 7 step setup process. Follow the instructions on the screen or push *Cancel* to quit the setup.




i-Drive Configuration Settings

13.6

To begin the configuration, click the *Config Mode* button on the left of the screen.

Configuration Feature

Configuration allows you to adjust the way the i-Drive activates and allows you to fine tune the function and sensitivity of the sensor/joystick/chair.



- A. *Sensor Channels* Tab– Navigates to the Sensor Channels menu
- B. *Double Tap Settings* Tab– Navigates to the Double Tap Settings menu
- C. *Reboot Device* Button– Reboots the i-Drive
- D. *Factory Reset* Button– Resets all settings to factory default
- E. *Save Settings* Button– Saves your settings

Configuration: Sensor Channels

13.7


- 1-Diagram shows channel switches and channels currently in use
- 2-Drop down menus that allow you to assign function to each channel
- 3-Channel Deadband slider that allows sensitivity change of deadband
- 4-Changes the channel of which the deadband has been adjusted

NOTICE

Channels that are not being used will appear in red on the screen.

NOTICE


Be sure to save your settings with the *Save Settings* button after all the desired changes have been made.



Configuration: Double Tap Settings

13.8

- A. *Double Tap*–Enabling Double Tap allows for adjustment of Double Tap input and output timing
- B. *Double Tap Input Delay*– Amount of time for Double Tap command capture
- C. *Double Tap Output Speed*– Time period between Mode output pulses. This value should match the Double Tap Input Timing of the chair



NOTICE

It is recommended to use *Calibrate Double Tap* option in the *Diagnostics* to get a properly timed double tap setting for the user.

Running Real Time Diagnostics


13.9

To begin the *Diagnostics*, click the *Diagnostics* button on the left of the screen.

- A. Mode Switch activation
- B. Left, Right, Forward & Reverse activation
- C. Corresponding i-Drive port when switch is pressed
- D. Double Tap Calibration

Diagnostics Feature

With real time *Diagnostics*, you are able to observe pad behavior as the head array is activated by the user. This is useful to fine tune the i-Drive and confirm it is operating properly.



⚠ CAUTION

Double Tap needs to be enabled in the *Configuration* menu *Double Tap Settings* before *Double Tap Calibration* can be used.

Help Feature

13.1.0

Clicking the *Help* button will activate help. When *Help* is active, the help button will be green and a question mark (?) will be displayed next to your cursor. With *Help* active, clicking buttons and menu items will pop up with an explanation of what the clicked object is and how it can be utilized. To disable *Help*, click on the *Help* button again.



NOTICE

The *Help* button will appear black when inactive and green when active.

Closing the i-Drive Programming Software

13.1.1

When you are finished configuring the i-Drive, you can click on the *Quit* button on the bottom right of the screen to close the application.




Download the i-Drive *Bluetooth*® Application

14.1




For iPhone® Users:

- Search, “Trident Research, LLC” in the Apple® App Store®
- Install the App and then “Open” once App has been successfully installed on your phone.



For Android Users:

- Search “Trident Research, LLC” in the Play Store
- Select “Install” and then “Open” when the App has been successfully installed on your phone.



Opening the i-Drive App

14.2

- A battery, battery charger, dongle cable, and a 9-pin pig tail plugged into the USB port on the interface is required in order to connect the i-Drive to your phone.
1. To Log In, enter the User Name and Password assigned to you.
 2. When the Main Menu appears, *Connect* to the i-Drive device (*Figure 1*).
 3. Your Bluetooth® device will search for the i-Drive device.
 4. When you have successfully logged in and the App has located the i-Drive device, new options will appear in the *Main Menu: Config Mode, Diagnostics, and About* (*Figure 2*).

NOTICE

Ensure your phone or wireless device has the *Bluetooth®* option turned on.




Figure 1




Figure 2

Config Mode

14.3

1. In the *Main Menu*, Select *Config Mode* (Figure 3). Like the Advanced Programming Software on your PC or tablet, channel assignments, double tap settings, and other sensor settings are adjusted here.
2. Under the *Name* section, you have the option of renaming your i-Drive. Changing the name of your i-Drive will require a power cycle.




Figure 3

Channel Assignment

1. Select the *Channels* tab under the *Config Mode* option in the App (Figure 4).
2. The *Channels* tab view shows the standard channel assignments on the i-Drive device, and any channels that are disconnected from the device.
3. To assign a function to each channel, select the channel you are wishing to modify.




Figure 4

Channel Assignment Settings

- 1. Within the channel option chosen, you will see a list of modifications and assignments that can be made to the channel.
- 2. To select the assignment for the channel, click on the assignment. To go back to the channel options, click back on the upper left hand corner.
- 3. Under the *Sensor Settings* header, the sensitivity of the sensor connected to the device can be adjusted. Slide your finger to the left or to the right to increase or decrease the sensitivity.
- 4. After channel assignments have been made, in the *Channels* tab, all of the selections will appear to the right of the channel they have been assigned to (Figure 6). If the channels are to the correct assignments, select *Save* in the top right hand corner.
- 5. Select back in the top left hand corner to go back to the *Main Menu*.




Figure 5




Figure 6

Multi-Tap Setting

- 1. To activate the Multi-Tap function, select *Multi-Tap* in the *Config Mode* menu.
- 2. Under the *Multi-Tap* setting, *Input Delay* and *Output Speed* can be adjusted. Move the cursor to the right or left to make adjustments. (*Figure 7*)
- 3. When adjustments have been made, select the *Save* option in the upper right hand corner.




Figure 7

Diagnostics

- 1. After channel assignments and any other modifications necessary have been made, go back to the *Main Menu*.
- 2. Select *Diagnostics* in the *Main Menu*.
- 3. Test the assignments on your i-Drive device by activating each channel. When the direction assigned to the channel has been activated, the direction should light up in green. (Figure 8)
- 4. After the diagnostics have been tested, select back in the top left hand corner to go back to the *Main Menu*.
- 5. After all updates to the interface have been made, exit the App.




Figure 8

After the installation of the i-Drive, execute the following tests before the wheelchair is delivered or put into service:

- Operational Test
- Test Drive
- Stop Test

Operational Test

15.1

WARNING

Execute this test only on a level surface, in an open area free of obstacles.

1. Activate the wheelchair operating system.
2. Check for any error messages.
3. Apply pressure to the head pad/sensor/joystick, etc. designated for moving the chair in a forward motion until you hear the parking breaks switch off.



CAUTION

The wheelchair may start moving.

4. Immediately release the pressure from the head pad/sensor/joystick, etc. You should hear the parking break react within a few seconds.
5. Repeat *Steps 3 & 4* three times, while carefully applying pressure to the device.
6. Check whether the double tap features function properly.

Test Drive

15.2

- Perform a test drive with the wheelchair.
- Check whether the wheelchair and all of its operations are fully functioning.
- Check that no cabling or parts may get damaged or hindered in any possible position of the wheelchair.

Stop Test

15.3

- Drive full speed ahead and shut off the wheelchair.
- The wheelchair may not suddenly stop, but must slow down to a gradual stop.

Dealer Assistance

16.1

During first time use by the client, it is advised that the dealer or service technician assists and explains the different drive configurations to the customer (the user and/or the attendant). If needed, the dealer can make final adjustments.

User Testing

16.2

It is important that the customer is fully aware of the installation of the i-Drive, how to use it, and what settings in the software or App can be adjusted in order to gain as much mobility as possible. As a dealer, proceed as follows:

- Explain and show the customer how you have executed the installation, and explain the functions that have been assigned to the device.

Have the user test all configurations of the device:

- Is the device being used by the client within easy reach and in an optimal position for the user?
- Can the user safely operate the wheelchair with minimal effort?

Conditions Of Use

16.3

The i-Drive 4.0 with *Bluetooth*[®] wireless technology is intended for use as installed by the dealer, in accordance to the installation instructions in this manual.

- The foreseen conditions of use are communicated by the dealer or service technician to the user and/or attendant during the first time use.
- If the conditions of use change significantly, please contact your dealer or a qualified service technician to avoid excessive wear and tear or unintended damage.

CAUTION

These products are designed to be fitted, applied, and installed exclusively by a health care professional trained for these purposes. The fitting, application and installation by a non-qualified individual could result in serious injury.

Care And Maintenance

17.1

Ensure the hardware you are using stays in working order by keeping it cared for and maintained.

- Keep electronics dry and out of the water.

WARNING

Getting electronics wet will cause irreversible damage to the circuitry.

- Periodically check the hardware for loose screws or worn parts. Replace or repair the parts as needed.
- To clean the aluminum, use a mild, nonabrasive household cleaner.
- To wash covers, remove the cover from the pad, machine wash cold on a delicate cycle and drip dry afterwards.
- On the hardware, lightly tighten set screws until they are snug, then tighten an extra quarter turn.

WARNING

Do not overtighten set screws. Over tightening set screws will prevent hardware from functioning properly and could cause irreversible damage to the hardware.

Safety

17.2

Avoid getting electronics wet and working on electronics with wet or damp hands. This could cause minor to severe electrocution, resulting in personal injury and/or product damage.

Technical Data

17.3

The i-Drive Bluetooth Interface was tested and found to be in compliance with the required criteria and met the standards of IEC60601-1-2:2014 (Edition 4.0); RESNA WC-2:2009.

This image shows a full page of blank handwriting practice paper. It features approximately 20 evenly spaced, horizontal blue lines across the entire page. The lines are thin and consistent in color, providing a guide for letter height and placement. There are no margins, text, or other markings on the paper.

This image shows a full page of blank handwriting practice paper. It features a series of evenly spaced, horizontal blue lines across the entire width of the page. The lines are thin and light blue, set against a plain white background. There are no margins, text, or other markings on the page.

[illegible]



Stealth Products, LLC. • info@stealthproducts.com • www.stealthproducts.com
+1(800) 965-9229 | +1(512) 715-9995 | 104 John Kelly Drive, Burnet TX 78611