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


<h3>ECG/Blood glucose/Blood pressure Smart Bracelet</h3>  <p>User Manual</p>	<h4>Device's structure</h4>  <p>1</p>	<h4>Disassembly and wear</h4> <p><b>1. Disassembly</b> Align the metal pin on left end of the watch strap with the hole at the connection point of the bracelet, flip the metal ball on the right, and align the hole on the right to complete the installation and disassembly.</p> <p><b>2. Wear</b> Better worn the watch at 1-2CM behind the ulnar neck of the wrist. When wearing, make sure that the bottom of the watch is close to the skin to avoid light leakage of the sensor and affect the accuracy of the measurement.</p>  <p><b>Charge</b> This is magnetically charge, so only need to align it with the charging port and automatically attach it.</p>  <p><b>Power on / off</b></p> <p>2</p>	<p><b>Switch on:</b> Long press the screen to power on.</p> <p><b>Switch off:</b> On settings and find out switch off.</p>  <p>3</p>	<h4>Install APP and Bluetooth connection</h4> <p>1. Turn on the device, turn on the Bluetooth of the mobile phone, scan the QR code below or search "H Band" in the App store / Google play store to download and install.</p>  <p>4</p>	<h4>Operation guide of the device function interface</h4>  <p>5</p>	<h4>Device function introduction</h4> <p><b>1. Time, date &amp; week</b> For the first time, please connect to the APP to sync the correct time, date, and week.</p> <p><b>2. Dial switch</b> Long press the screen to switch the dial style. Note: To customize the watch face, please go to APP Settings ( Mine →EP01→Dial Settings).</p>  <p>6</p>	<h4>Wear and maintenance</h4> <p>The products are suitable for wear all day, so when you wear and care for the smart device, please follow the following brief instructions:</p> <ol style="list-style-type: none"> <li>Keep the smart device dry:       <ol style="list-style-type: none"> <li>Clean your wrist and the smart band/smart watch you are wearing regular, especially after exercise, sweating, or skin contact with soap or detergent and other substances, the substance may get inside of the product so should be cleaned regularly.</li> </ol> </li> <li>Keep wearing the product correctly:       <ol style="list-style-type: none"> <li>Please do not wear the device too tightly, should be able to move back and forth on the wrist.</li> <li>If use a heart rate device to record your heart rate throughout a day, the device should be flat one finger wide among the wrist bone. In order to better check the data during active exercises, wear the device among the wrist bone about 2-3 Finger-width position.</li> <li>Prolonged friction and restraint may irritate your skin, after a long time wearing, take off the device and give skin a rest period.</li> </ol> </li> <li>Keep smart devices clean:       <ol style="list-style-type: none"> <li>If device is damp, such as sweating or splashing, please clean and dry it, make sure the skin layer is in a suitable dry humidity, and then wear again.</li> </ol> </li> <li>Keep a correct sitting posture, make the hands flat and keep it parallel with heart.</li> </ol> <p>The data collected by this device is for reference only, not for medical purpose</p> <p>7</p>
<h4>8. ECG test</h4> <p>Switch to the interface of ECG. Click the button to the ECG test state, and put your finger against the right electrode pad. (Please make sure the device is worn on your wrist properly and the skin is moist if the result is out). The result by this manual testing on bracelet will not be saved on App. If you want to save the data, please start the ECG testing on App. The data will be shown on a graph.</p>  <p>9. Blood glucose monitoring</p> <p>Switch to this interface, Tap the triangle button below, The icon starts to flash, then shows the result.</p> <p>10. Blood oxygen monitoring</p> <p>Switch to this interface, Blood oxygen automatic measurement, The icon starts to flash, then shows the result.</p> <p>11. Body temperature monitoring</p> <p>Switch to this interface, Tap the triangle button below, The icon starts to flash, then shows the result.</p> <p>12. Sleep monitoring</p> <p>Switch to this interface, Check the duration of deep sleep, light sleep, duration of total sleep and sleep quality.</p> <p>13. Weather</p> <p>When the device connected to the APP, the weather will automatically be sync.</p> <p>8</p>	<h4>14. Notification</h4> <p>Synchronize incoming calls, text messages, QQ, WeChat and other application software messages, never missing any important messages anytime, anywhere.</p> <h4>15. Alarm</h4> <p>Switch to this function interface, to find clock, when the alarm not set, it is shown on the right side: APP find "My", Tap "EP01" to find "Alarm Settings", tap on the upper right corner "+", after setting completed will show on the device.</p> <h4>16. Stopwatch</h4> <p>Tap to enter the stopwatch interface, Tap the middle button to start or pause, To recalculate, press the reset button, and press the exit button to return.</p> <h4>17. Timer</h4> <p>Find out timer in the function list, adjust the time to start timing.</p> <h4>18. Music</h4> <p>When the device connected to the APP, The music player software has been installed on the phone just tap play to realize the synchronization of the device. Tap the middle play/pause button to control the phone to play music.</p> <h4>19. Find phone</h4> <p>After Bluetooth connected, tap on the device find phone function, the phone will vibrate and ringing.</p> <h4>20. Female</h4> <p>Save personal information as women. This function can be set according to personal circumstances.</p> <h4>21. Camera</h4> <p>After the Bluetooth connection, enter the watch camera interface, you can control the phone to take remote pictures.</p> <h4>22. Settings</h4> <p>Find out set in the function list, As follows: Brightness/Auto-lock/Switch/Date&amp;time/H Band/Dial/ Notification/Language/About/Erace Date/Power off.</p> <p>9</p>	<h4>APP functions and instructions</h4> <p><b>1. Dashboard</b> ①Steps and time;②More (Heart rate / Blood pressure / setting);③Steps, distance and calories data interface;④Sleep data interface;⑤Heart rate monitoring data interface;⑥Blood pressure monitoring data interface;⑦ECG monitoring data interface;⑧Blood oxygen monitoring data interface;⑨HRV data interface;⑩Body temperature data interface;⑪Women's menstrual cycle data interface;⑫Blood glucose data interface.</p>  <p>10</p>	<h4>2. Movement</h4> <p>①Select the GPS Movement button in the App. ②Select the map software installed on the phone. ③Select Start sport button. ④Go to the open area and select Pause button if the GPS single is weak.</p>  <p>11</p>	<p>When GPS exercise completed select movement of the device ⑤, tap start button, and select exercise mode to confirm. At the same time, device will automatically switch to the data interface in this sport mode the APP can be used to pause exercise or on device left side to pause exercise, while the exercise data will be automatically saved in the APP.</p>  <p>3. Mine</p> <p><b>Personal information:</b> Tap the avatar and nickname (3) to enter. The content are as follows: nickname, gender, date of birth, weight, height, BMI, skin color.</p> <p><b>My device:</b> the content are as follows 2: Alarm notification: receive/close the application message reminder. Alarm settings: check or create an alarm. Sedentary setting: Set a sedentary reminder time. Heart rate alert: Set the upper limit of resting heart rate, and automatically alarm when it reaches or exceeds it. Turn wrist detection: turn your wrist to brighten the screen. Pair device: binding the device. Find device: find the device. BP Private Mode: set the private mode to obtain more</p>  <p>12</p>	<p>accurate data for individually, normally use general mode.</p> <p><b>Screen-on Time:</b> Set screen time.</p> <p><b>Weather setting:</b> Set temperature unit/switch on weather push function.</p> <p><b>Switch Setting:</b> HR Automatic monitoring/BP Automatic monitoring/Blood Oxygen Nighttime monitoring/Disconnect Alert/Scientific sleep/Automatic temperature monitoring</p> <p><b>Take photo:</b> After connecting Bluetooth and clicking app to take photos, the watch can control the mobile phone to take photos.</p> <p><b>Reset device password:</b> reset the device password.</p> <p><b>Dial settings:</b> support choose to use the inherent dial (3) of the device, or choose the "photo dial", (customize to upload pictures as you like), also can choose exquisite watch face in "More Watch Faces".</p> <p><b>Firmware update:</b> current firmware information and update firmware.</p> <p><b>Clear Data:</b> Clear all data on the device and shut down.</p> <p><b>Reset:</b></p> <p><b>Disconnect:</b> Disconnect device.</p>  <p>13</p>	<p><b>Help:</b> check answers of "Cannot connect to the device", "Easy to disconnect the device", and "Message cannot be pushed".</p> <p><b>Permission management:</b> Obtain part of mobile phone permissions during APP use.</p> <p><b>Steps goal:</b> Set target steps.</p> <p><b>Sleep:</b> Set target sleep time.</p> <p><b>Unit setting:</b> Set imperial or metric units.</p> <p><b>APP style:</b> Set the APP background style.</p> <p><b>About us:</b> check current APP version information.</p> <p><b>WeChat:</b> bind WeChat sports.</p> <p><b>Log out:</b> Log out of the current account.</p> <p><b>Parameter</b></p> <p>Chip: GR5515 Display: 1.47-inch HD color screen Material: Tempered glass/ABS+PC Silicone Weight: 32g Charge time: 2hour</p> <p>Battery: 200mAh Waterproof: IP67 Standby: 15days Bluetooth version: 5.1 Support OS: iOS 8.0 or above Android 4.4 or above</p> <p><b>Note</b></p> <ol style="list-style-type: none"> <li>Before measuring, please make sure that you have saved correct personal information on the APP, such as height, age, gender, weight, etc., and synchronize data.</li> <li>Sit for 5-10 minutes before measurement to keep your heart beat calm and breathing stable.</li> <li>Please keep relaxed during measurement, panting, talking walking and exercise will affect the accuracy of the measurement results.</li> <li>Keep a correct sitting posture, make the hands flat and keep it parallel with heart.</li> </ol> <p>14</p>	<h4>Wear and maintenance</h4> <p>The products are suitable for wear all day, so when you wear and care for the smart device, please follow the following brief instructions:</p> <ol style="list-style-type: none"> <li>Keep the smart device dry:       <ol style="list-style-type: none"> <li>Clean your wrist and the smart band/smart watch you are wearing regular, especially after exercise, sweating, or skin contact with soap or detergent and other substances, the substance may get inside of the product so should be cleaned regularly.</li> </ol> </li> <li>Keep wearing the product correctly:       <ol style="list-style-type: none"> <li>Please do not wear the device too tightly, should be able to move back and forth on the wrist.</li> <li>If use a heart rate device to record your heart rate throughout a day, the device should be flat one finger wide among the wrist bone. In order to better check the data during active exercises, wear the device among the wrist bone about 2-3 Finger-width position.</li> <li>Prolonged friction and restraint may irritate your skin, after a long time wearing, take off the device and give skin a rest period.</li> </ol> </li> <li>Keep smart devices clean:       <ol style="list-style-type: none"> <li>If device is damp, such as sweating or splashing, please clean and dry it, make sure the skin layer is in a suitable dry humidity, and then wear again.</li> </ol> </li> <li>Keep a correct sitting posture, make the hands flat and keep it parallel with heart.</li> </ol> <p>The data collected by this device is for reference only, not for medical purpose</p> <p>15</p>

**FCC Caution:**

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.