CPRBAND

Model: CREDO-CB

User Manual

Index

Overview	3
Summation	4
Important Information	6
Warnings and Cautions	7
Items Included	8
CPR BAND (Model: CREDO-CB)	9
Appearance	9
Operation	10
Battery Charging	11
Battery Indicator	12
Wearing CPR BAND (Model: CREDO-CB)	13
Preparation for use	14
CPR Mode	15
Healthcare in CPR BAND (Model: CREDO-CB)	19
Application	20
Install and Using	20
Device Connection	21
Setting	22
Real Time CPR	23
CPR Report	24
Healthcare	27
Cyber security requirement	28
Maintenance	29
Troubleshooting tips / service life	29
Storage after use	30
Storage when not in use	32
Specifications	33
Frequently used function	35
Essential Performance (EP)	36
Symbol Glossary	37
Regulatory Information	39
Additional Information	43

Overview

- Manufacturer: CREDO
- Place of manufacture:
 - #302, 86, Baeul-ro, Wonju-si, Gangwon-do, Republic of Korea
- EC Representative: KTR Europe GmbH
- Address: Mergenthalerallee 77, Eschborn, Hessen, 65760, Germany
- Product name: Cardiopulmonary Resuscitation Assistance Product
- Device Classification
 - CE: Class I (Rule 13)
 - FDA: Class II (Product code: LIX)
- FCC ID: 2ANY8-CREDO-CB
- Model name: CREDO-CB
- The term of warranty: A separate mark
- Manufacturing number: A separate mark
- Date of manufacture: A separate mark
- User Manual Date of Issue: 2019-01-01
- Intended use: The CPR BAND (Model: CREDO-CB) is used as a guide in administering cardiopulmonary resuscitation (CPR) to a suspected SCA victim at least 8 years old.

This product is '**Medical device**'.

Summation

The CPR BAND (Model: CREDO-CB) is a small, lightweight device that uses rechargeable batteries and is worn on the wrist. This device is designed for first aid personnel trained in CPR and use of the CPR BAND (Model: CREDO-CB).

The CPR BAND (Model: CREDO-CB) provides real-time information about chest compression according to the latest CPR guideline when CPR is performed on patient who is estimated a sudden cardiac arrest (SCA). This device displays CPR feedback indicator for the chest compression depth. It guides chest compression rhythm through sound and indicates the total chest compression time and number of chest compression.

The CPR BAND (Model: CREDO-CB) is a device that provides Bluetooth function. It can transmit CPR record and confirm information on the smart device. In the smart device that provides the wireless communication function by receiving information provided by the CPR BAND (Model: CREDO-CB), it can display the pressure axis angle, pressure axis average value, compression depth average value, compression rhythm average value, and chest compression graph which are not displayed on the CPR BAND (Model: CREDO-CB).

The device provides other functions on the user's wrist while waiting without CPR. The time, step count, movement distance and calorie consumption function are displayed.

Operating Principle

The CPR BAND (Model: CREDO-CB) measures the chest compression depth during CPR via the built-in accelerometer and displays it on the LED. The measured data is temporarily stored in the CPR BAND (Model: CREDO-CB) and transmitted to the application via the Bluetooth interface. In the application, detailed data information such as compression depth, compression rhythm, and compression axis angle can be checked.

Important Information

Read this User Manual and become familiar with the operation of the product prior to use. Use the product only as described in this User Manual.

Warnings and Cautions

'Warning' states a condition, hazard, or unsafe behavior that can result in serious injury or death.

'Caution' states a condition, hazard, or unsafe behavior that can result in minor personal injury.

Notes

A note states important information about the product or its operation.

Warning

Do not use device for SCA patients under the age of 8 years.

Note

CPR cannot assure survival, no matter how well it is performed. In some patients, the underlying problem causing the cardiac arrest is not survivable despite any available care.

Warnings and Cautions

Warnings

- Do not use the device in conjunction with any mechanical or automated compression device.
- When using this device, make sure that there are no objects between the performer's palm and the patient's chest.
- Do not delay CPR. If there are any problems using the device, continue CPR without it.
- If the device appears to be damaged e.g. develops cracks or sharp edges, do not use.
- No modification of device is allowed.
- Do not modify device without authorization of the manufacturer.
- If device is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the device.


Cautions

- If difficulty is encountered in entering CPR mode, do not delay initiation of CPR.
- The band is made of natural latex, but if you are sensitive to latex, it may cause an allergic reaction.
- The '!' status indicates a technical error. If this occurs, stop using the CPR BAND (Model: CREDO-CB) and continue CPR. Contact CREDO technical support after the event.
- The module is small and can be swallowed by children or pets. It should be kept out of the reach of children or pets.

Contact a local CREDO representative for further assistance, or visit www.credo-ltd.com for more information and Frequently Asked Questions (FAQ).


Items Included

Charge Cable may differ in appearance and may be changed.

Contact www.credo-ltd.com for more information including latest product downloads, spare parts and accessories.

Appearance

- | | |
|----|---|
| 1 | Combining Part: Combined with the '8' part of the band |
| 2 | OLED: Displays information such as time, battery amount, CPR count, CPR time, and healthcare function |
| 3 | Buzzer: Compression rhythm guide at CPR |
| 4 | LED: Compression depth feedback at CPR |
| 5 | Button: Change mode such as Power On, CPR Mode entry |
| 6 | Charging Port: Connect the cable to the battery when charging |
| 7 | Band Head: Connected to '11' to wear on user's wrist |
| 8 | Band Clip: combined with the '1' part of the module. press this part when separating it. |
| 9 | Band Body: Combined with the bottom of the module |
| 10 | Band Body Rubber: Combined with the '6' charging port of the module. When combining modules and bands, these parts must be joined first. |
| 11 | Band Tail: Connected to '7' to wear on user's wrist |

CPR BAND (Model: CREDO-CB)

Operation

Power on	Press the button briefly to turn on the power
Sleep mode	If you do not press the button for 4 seconds, it becomes sleep mode.
Additional function	If you press the button briefly, additional functions including healthcare mode will be sequentially indicated.
CPR mode	Press button for 3 seconds to start CPR mode. CPR BAND (Model: CREDO-CB) guides the compression rhythm through the Buzzer and provides depth feedback through the LED. To exit CPR mode, press and hold the button for 2 seconds. OLED displays the CPR time and the total number of compression for 10 seconds and then goes into sleep mode.
Rebooting	When charging, press button to reboot.

CPR BAND (Model: CREDO-CB) LED

- The LED of the device operates for feedback of pressure depth in CPR mode.
- Module has 5 LEDs on each side, total 10 LEDs.
- The colors of LEDs are Orange, Blue, Blue, Blue, Orange in order.
- If the compression depth is lower than the correct range, only the Orange LED at the bottom blinks.
- If the compression depth is within the correct range, blue LEDs blink.
- If the compression depth exceeds the correct range, all Blue LEDs and Orange LED at the top blink.

Battery Charging

- Charge the battery by connecting the cable to the charging port of the module. Use the enclosed cable. Use a power source certified by IEC 60950 or IEC 60601-1.
- When charging, the OLED indicates the charging status.
- To reboot and reset, press the button during charging. The set correct depth, rhythm, and stored CPR data will be initialized

connect the cable

Battery charging

Warning


Use a power source certified by IEC 60950 or IEC 60601-1. The rating should be output 5V-1A (max 2A).

Do not be posited device to the area where it is difficult to operate the USB Cable

CPR BAND (Model: CREDO-CB)

Battery Indicator

The CPR BAND (Model: CREDO-CB) constantly monitors battery power. Check the battery status regularly by turning on the device. The battery amount is displayed in the top right corner of the OLED as follows


- Duration based on remaining battery power

Battery amount	CPR Mode Duration	Sleep Mode Duration
100%	5 hours	40 hours
80%	4 hours	32.5 hours
60%	3 hours	25 hours
40%	2 hours	17.5 hours
20%	1 hours	10 hours

- Low battery

If the low battery icon appears, the product might turn off during operation. Please use after charging enough

Warning

Do not interrupt CPR to charge the battery. Continue CPR without feedback from the CPR BAND (Model: CREDO-CB).

CPR BAND (Model: CREDO-CB)

Wearing CPR BAND (Model: CREDO-CB)

When combining the module and the band, the charging port of the module should be connected to the Band Body Rubber. To separate the module and the band, press the clip.

how to combine module and band		
		
1. Combine Band and Module charging port	2. Combine Band and Module combining part	3. Wear the band on wrist
how to separate module and band		
		
1. Press the Band Clip	2. Lift up the module combining part	3. Remove the module from the band


Cautions

- The band must be worn properly on the user's wrist. Incorrect wearing can result in inaccurate function.
- Failure to observe the combining method may damage the band.
- Please note that it may irritate the skin when worn for a long time.

Preparation for use

- The product must be powered on.
- Make sure the patient is on a solid surface.
- Remove the patient's upper garment and reveal the chest.

Caution

- Provide chest compressions according to local CPR protocol.
- Keep your wrist at 90° for correct CPR.

Warning

When performing CPR on a patient lying on a mattress, a backboard must be used to limit the amount of compressed depth which is absorbed by the mattress. Depending on characteristics of the mattress, backboard and patient, the depth compensation does not guarantee that the patient chest is pressed correctly.

CPR BAND (Model: CREDO-CB)

CPR Mode

- Press button for 3 seconds in Power on state to enter CPR mode.
- CPR BAND (Model: CREDO-CB) guides the compression rhythm through the Buzzer and provides depth feedback through the LED.
- When CPR is finished, press button for 2 seconds to exit CPR mode.
- Immediately after quitting CPR mode, OLED displays 'Total number of compression' and 'CPR Time' for 10 seconds and then turns off.

Depth feedback

Rhythm guide

CPR time

Total number of compression

CPR Mode

1) Compression Depth Feedback

The LED lights according to the correct compression depth set in the [Setting] tab. The maximum value is 10 cm and the minimum value is 1 cm of the correct compression depth you can set. The factory-set values of correct range is 4.5 cm ~ 6.0 cm.


* DC: Depth of Compression

* M: Maximum depth of setting

* m: Minimum depth of setting

$$* A1: m + (M - m) \times \frac{1}{3}$$

$$* A2: m + (M - m) \times \frac{2}{3}$$

①	$DC < m$	Lower 1 orange LED lights up
②	$m \leq DC < A1$	1 blue LED lights up
③	$A1 \leq DC < A2$	2 blue LEDs light up
④	$A2 \leq DC < M$	3 blue LEDs light up
⑤	$M \leq DC$	3 blue LEDs and upper Orange LED light up

CPR Mode

2) Compression Rhythm Guide

The CPR BAND (Model: CREDO-CB) guides the chest compression rhythm through the built-in Buzzer. Compression rhythm (BPM) can be set in the [Setting] tab of the application according to international standards. Factory-set is 110BPM.

Caution

- Change the setting according to local CPR protocol and provide chest compressions.
- There is no the time required for the CPR BAND (Model: CREDO-CB) to warm from the minimum storage temperature between uses until the device is ready for its intended use when the ambient temperature is 20 °C
- There is no the time required for the CPR BAND (Model: CREDO-CB) to cool from the maximum storage temperature between uses until the device is ready for its intended use when the ambient temperature is 20 °C.

How to process CPR data

Data process	How to process
Data collection	<p>- Data storage items:</p> <ul style="list-style-type: none">① Total CPR time (CPR mode start~CPR mode stop)② Compression time (2cm or more compression start ~ last compression)③ Each compression max depth / count to the next max depth④ Compression angle <p>- Storage capacity: 128KB / Approximately 120 minutes CPR available</p>
Data release	The data is automatically backed up when connecting with the application, and the data of the module is deleted.
Data reset	Data reset function does not exist separately, but data is automatically backed up when the product is connected to the application, and the data of the module is deleted (reset).

CPR BAND (Model: CREDO-CB)

Healthcare in CPR BAND (Model: CREDO-CB)

The usual CPR BAND (Model: CREDO-CB) provides time and healthcare functions. If the button is pressed briefly while the product is on, the following healthcare functions will appear in sequence on the OLED. If there is no touch for 4 seconds, it goes to sleep mode.

Time: Indicates the current time.

Walk Steps: Indicates the number of steps of the day.

Distance Moved: Indicates the distance moved of the day.

Calories Burned: Indicates the calorie consumption of the day.


Notes

- In order to use more accurate healthcare functions, please input your information in application. See 'Setting' on page 22 for more information.
- In all status, press button for 3 seconds to enter CPR mode.

Frequently of administration	Daily
Duration of application	Within 24 hours

Install and Using

1. Android products search for 'credo cpr band' in Play store.
IOS products search for 'credo cpr band' in App store.
2. Click the application icon to run it.


- ① **Setting:** User Information, Device Name, Correct Compression Depth, Correct Compression Rate, App Information, Company Information, Open Source License
- ② **Bluetooth:** CPR BAND Connecting management
- ③ **CPR DATA / PEDOMETER:** Click the desired tab from CPR DATA and PEDOMETER.
- ④ **Real Time CPR:** You can see real-time CPR information.
- ⑤ **CPR History:** CPR records. The latest record is displayed at the top and click the desired record to check the report.

Application

Device Connection


1. Turn on Bluetooth on smart device.
2. Click the Bluetooth icon in the top right corner of the application and click [Find CPR Band] to connect the device to the application.
3. If you keep Auto Connecting ON, my CPR BAND will connect automatically.

If the device is in sleep mode, it may not work. Press the button on the module to wake up the sleep mode.

Application

Setting

User Information

User information setting. Set user gender, age, height and weight for healthcare functions.

Device Name

The name of the connected device should be changed for identification with other devices. (Up to 8 characters)
Factory set: CPRBAND

Correct Compression Depth

Setting the CPR compression depth according to local regulations. Select Min/Max value of compression depth.

Correct Compression Rate

Setting the CPR compression rhythm according to local regulations. Select BPM of compression rhythm.

App Version

Version information of the application. The Android product is updated in the Play Store, and the iOS product is updated in the App Store.

App Name

Application name

Company Information


Company name, address, contact information, email information

Open Source License

Open source license information

Real Time CPR

If CPR Mode of CPR BAND is operated while the application is connected, the CPR Data is displayed on the Real Time CPR Screen of application at real time.

① A graph of compression depth is displayed.

② Total time / The number of correct CPR / Total number of CPR are displayed.

③ DC<M: displayed in blue.
M≤DC: displayed in red.


④ If the wrist angle is vertical, the arrow points to the center.

If the wrist angle is wrong, the arrow points to the left or right.

⑤ If the compression rate is correct, the bar is centered. The bar is on the left when the compression is slow, and on the right when it is fast.

Application

CPR Report

Case date: Date CPR was performed

Gender: Gender of patient

Episode Summary

Episode start time

: CPR start time

Episode end time

: CPR end time

Total length of episode

: Total time to CPR

Compression Data

Total number of compressions: Total number of chest compressions

Total number of correct compressions: Total number of correct compressions

Percentage of correct compression [%]: correct compression as a percentage

Average compression depth [mm]: Average of chest compression depth

Average compression rate [/min]: Average of chest compression rhythm

Wrist Angle accuracy [%]: Accuracy of the wrist angle as a percentage

Graph: Indicates peak value of compression depth


Comment / Note: Indicates a comment left by the user

Application

CPR Report

On the CPR report, click [].

Reports can be modified, exported, and deleted.

Modify

: You can enter patient gender and comment.

Export

: Reports can be converted to PDF files for export.


Delete

: Delete the report.

Application

CPR Report

Click [Export] to convert the report to a PDF file for e-mailing or printing.


Converts Report to PDF
file and sends it by e-mail.

Converts Report to PDF
file and prints it.

Healthcare

If you click [PEDOMETER] on the main screen, you can check the healthcare information. Data is stored for up to 30 days.

Cyber security requirement

Set the security function of the smart device itself so that non-users cannot access the application

When updating the application, update is performed by securing the unique identification ID and password of our update manager

- 1) When updating from the app store of iOS, update after the pre-check and verification process according to Apple's regulations
- 2) When updating from play store of android, update after pre-check and verification process according to Google's regulations
- 3) The security rules of the application are based on the internal security rules of Apple and Android.
- 4) When updating, a unique update number is given to prevent product connection through unauthorized updates.

Troubleshooting Tips

These troubleshooting tips provide immediate corrective actions for possible problems. If a problem persists, contact your local CREDO representative for assistance.

OBSERVATION	CORRECTIVE ACTION
The LED display is visible, but does not provide the correct compression feedback.	Confirm the [correct compression depth] in the application settings.
The product can't be connected to the application.	Confirm the smart device's Bluetooth feature is turned on. CPR BAND (Model: CREDO-CB) may be in sleep mode. Press the button once and try connecting again.
Module does not work.	Connect the charger and press the button to reboot.
The app doesn't work.	Do not connect the app and CPR BAND, please contact CREDO.

Expected Service Life

During product development, the CPR BAND (Model: CREDO-CB) is subjected to rigorous life testing. However, both rapid technological changes and the availability of replacement parts limit the useful life of all modern medical devices. The routine cleaning and disinfection program recommended in this manual will assist in reliable operation.

Item	Expected Service Life
CPR BAND (Model: CREDO-CB)	3 years
Battery	3 years; Recycling 500 times
Charge Cable	1 year

Storage after use

► Routine maintenance

If the button is pressed while the charger is connected, CPR BAND (Model: CREDO-CB) is rebooted and all setting values are initialized. The frequency of maintenance depends on usage, if it is used daily, it is once after 4300 hours.

► After each use

After use on a patient, the CPR BAND (Model: CREDO-CB) may be contaminated and should be handled appropriately.

1. Place the contaminated CPR BAND (Model: CREDO-CB) in a plastic bag until it can be cleaned.
2. If it is visibly soiled, wipe the CPR BAND (Model: CREDO-CB) with a soft cloth or paper towel to remove as much contamination as possible.
3. Separate the module and the band. The upper part of the module must be separated first.
4. Clean the CPR BAND (Model: CREDO-CB) as described under Cleaning and Disinfection. Proper cleaning is required to achieve disinfection.
5. Check the exterior of the CPR BAND (Model: CREDO-CB) for signs of damage. Contact CREDO for replacement parts or service if needed.

Maintenance

► Cleaning and Disinfection

Contaminated CPR BAND (Model: CREDO-CB) should be properly cleaned and disinfected.

Cleaning

1. Prepare a small brush (e.g. toothbrush) moistened with warm water (40 - 50 °C).
2. Scrub the device with a wet brush for at least 2 minutes. Do not clean the charging port directly.
3. If the port is clogged, use the brush to remove any obstruction.
4. Wipe the exterior with a soft cloth dampened in lukewarm water (22 - 40°C).

Disinfection

1. The exterior can be disinfected in a 0.55% solution of ortho-phthalaldehyde, or a 70% solution of isopropyl alcohol (isopropanol). Spray the exterior surfaces with the solution and ensure it remains wet for a minimum of 12 minutes. Repeat spraying as necessary to account for evaporation. Do not spray the solution directly on the charging port.
2. Wipe the exterior at least 3 times with a clean damp cloth to remove all traces of the disinfectant.
3. Allow to dry completely.

Warning

Do not immerse the CPR BAND (Model: CREDO-CB) in running water or allow moisture to penetrate it. Do not sterilize the CPR BAND

Storage when not in use

When not in use, keep the product free from damage.

When you use again, use it after charging.

Caution

If the Customer Service Indicator appears,

Contact www.credo-ltd.com for further support.

Specifications

CPR Target

Category	Specification
Compression Depth Target	>50mm Depth accuracy: $\pm 10\%$
Compression Rate Target	100 ~ 120/min ± 3 /min

CPR BAND: Module (Model: CREDO-CB-MO)

Category	Specification
Dimensions	47.1mm x 21.1mm x 12mm
Weight	11.0 g
Sleep mode hold time	40 hours at full charge in room temperature
CPR mode operating time	5 hours at full charge in room temperature
OLED display	Number of Pixels: 64 x 32
Electromagnetic compatibility	Meets IEC 60601-1-2
INPUT	5 V / 1 A
Bluetooth	BLE4.1

CPR BAND: Band

Category	Specification
Dimensions	260.9 mm x 17.7 mm x 10.9 mm
Weight	18.0 g

Environmental conditions

	Operation	Transport / Storage
Ambient temperature	5°C ~ 40°C	-25°C ~ 70°C
Ambient humidity	15% ~ 90% (Non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa)	0% ~ 90% (Non-condensing)
Atmospheric pressure	700hPa ~ 1060hPa	700hPa ~ 1060hPa (Water vapour pressure: <50hPa)

Specifications

CPR BAND (Model: CREDO-CB): consumable

consumable	Specification
Rechargeable battery	Lithium-ion Polymer Battery, DC 3.7 V, 120mAh service life: 500 times from full charge to full discharge

Software information

Component	Specification
Module	CREDO-CB-FW V1.0, 2018-01-01
Application	CREDO-CB-AP V1.0, 2018-01-01

Environmental Considerations

Product	Information
CPR BAND (Model: CREDO-CB)	The CPR BAND (Model: CREDO-CB) contains electronic components. Dispose of it at an appropriate recycling facility in accordance with local regulations.

Safety Classification

Category	Specification
protection against electric shock	CPR BAND (Model: CREDO-CB) is internally ME equipment. This device does not have applied parts.
Protection Against Solid Object and Liquid Ingress	IP 22 (IEC 60529 IP Code) Only when module and band combined
Method(s) of sterilization	No sterilization device
Suitability for use in an Oxygen rich environment	No intended for use in an Oxygen rich environment
Mode of Operation	Continuous Operation





Frequently used function

Function	Description
Charging Battery	Charge the battery by connecting the cable to the charging port of the module.
Combining Components	Connect Module and Band
Installing the Mobile Application	<ul style="list-style-type: none"> • Install the mobile application, which is called as the 'credo cpr band' or 'credo band', through the Play Store (Android Products) or APP Store (IOS Products) • Connect the application with CPR BAND (Model: CREDO-CB) through the Bluetooth
Power on	Power on the CPR BAND (Model: CREDO-CB) through pressing the button.
Setting	<ul style="list-style-type: none"> • Module Setting CPR Mode Setting : Press Button for 3 seconds (To escape the CPR Mode, Press button for 2 seconds on the CPR Mode) Healthcare Mode Setting : After power on, it is normally healthcare mode • Mobile Application Setting : Following information should be indicated on the mobile application <ul style="list-style-type: none"> - Auto connection - Auto Saving data - Correct Compression Depth - Correct Compression Rate - User Information - App Version


Essential Performance (EP)

EP ID	Essential Performance	Performance Limit
EP 1	CPR Mode Setting	1.CPR BAND (Module) 1) Move to the CPR mode To move the CPR Mode, Press button for three seconds. 2) Escape from the CPR mode To leave the CPR Mode, Press button for two seconds.
EP 2	CPR Mode	1.CPR BAND (Module) 1)Depth of Compression (DC) The LED lights up according to the DC. 2)Compression Rhythm Use the Built-in Buzzer

Symbol Glossary

Symbol	Description
	Manufacturer
	Date of manufacture
	Authorized representative in the European Community
	WEEE sign. Should not be disposed in general garbage.
IP22	Degree of protection provided by enclosures (IP Code)
	manufacture serial number
	Attached document number
	Warning / Caution
	Note
	Storage and transportation temperature : -25°C to 70°C
	Storage and transportation humidity : 0% to 90% (Non-condensing)
	Storage and transportation atmospheric pressure : 700 to 1060hPa (Water vapour pressure: <50hPa)
	Refer to user manual
	Not for use on children under 8 years old

Symbol Glossary

Symbol	Description
	Bluetooth symbol
	Contains or presence of natural rubber latex
	Do not use if package is damaged
	Stand-by
	Loudspeaker
	Keep dry
	Device includes RF transmitter.
	Korea Certification Mark
	Lithium ion Battery rechargeable
	CE mark according to (EU) 2017/745 for medical devices
	Korea GMP certification mark
	ISO 13485:2016 certification mark
Rx Only	Caution: Federal law restricts this device to sale by or on the order of a physician

Warnings

- When the CPR BAND (Model: CREDO-CB) is used together with a defibrillator, make sure to follow the defibrillator manufacturer's instructions. Stop the compression, remain clear of all patient contact during defibrillation or when otherwise required, in accordance with a proper defibrillation protocol.
- The CPR BAND (Model: CREDO-CB) is not intended for use in a moving environment, such as an ambulance. If used during patient transport, the device may provide inaccurate feedback. If CPR is indicated in a moving environment, do not rely on the depth feedback during such conditions. It is not necessary to remove the device.
- Do not practice CPR by using the CPR BAND (Model: CREDO-CB) on a people. It may be with a training manikin or simply on a compliant surface for practice.
- Properly performed CPR may result in fracturing of the patient's ribs.* If rib integrity has been compromised, continue to provide CPR in accordance with your local protocol.
- Properly performed CPR may result in chest injuries* e.g. external chest wall bruising or abrasion.
- Do not rely on CPR BAND (Model: CREDO-CB) feedback during aircraft ascent and descent, as its accuracy is reduced in such conditions.
- Not intended to be used within or outside the MR environment

Side Effect: Chest injuries from rib fractures

Cautions

- This device is designed to be used only with CREDO-approved accessories and may perform improperly if non-approved accessories are used. Do not attempt to modify the device in any way.
- If the lithium battery needs to be replaced, please contact CREDO.
- The CPR BAND (Model: CREDO-CB) is not a serviceable device. If experiencing technical difficulties contact a local CREDO representative for support.
- A device can be reused only if it is reconditioned under the responsibility of the manufacturer to comply with the general safety and performance requirements

* Black CJ, Busuttill A, Robertson C. Chest wall injuries following Cardiopulmonary resuscitation. Resuscitation. 2004;63:339-343

Regulatory Information

Note

- Changes or modifications not expressly approved by CREDO could void the user's authority to operate this equipment.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Recommendation

Responders should receive training, including regular refresher training, in use of the CPR BAND (Model: CREDO-CB). When training with the device on a CPR manikin, disable or ignore feedback from the manikin.

WEEE

This appliance is marked according to European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment.

Disposal must be carried out in accordance with local environmental regulations for waste disposal.

For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or the CREDO representative where you purchased the product.

The product is in compliance with the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Regulatory Information

Electromagnetic Conformity

- Electromagnetic environment Guidance

The CPR BAND (Model: CREDO-CB) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The CPR BAND (Model: CREDO-CB) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

There are no special requirements with respect to electrostatic discharge.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CPR BAND (Model: CREDO-CB) is used exceeds the applicable RF compliance level above, the CPR BAND (Model: CREDO-CB) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CPR BAND (Model: CREDO-CB).

Product complies SAR Exclusion threshold with 5mm distance.

- Cable List

Name	Model	Length	Shielded/ Unshielded	Remark
Charge Cable	CREDO-CB-CC	20cm	Shielded	USB - Micro 5pin cable

Regulatory Information

- Electromagnetic Emissions / Immunity Conformity

EMI: Group 1 Class B

EMS: ESD (Contact $\pm 8\text{kV}$, Air $\pm 15\text{kV}$), Burst (Power Line: $\pm 2\text{kV}$)

- Recommended separation distances between portable and mobile RF communication equipment and the CPR BAND (Model: CREDO-CB)

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150kHz~80MHz $d = 1.2\sqrt{P}$	80MHz~800MHz $d = 1.2\sqrt{P}$	800MHz~2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Notes

1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
2. The ISM (industrial, scientific, and medical) bands between 150kHz and 80MHz are 6.765MHz to 6.795MHz; 13.533MHz to 13.567MHz, 26.957MHz to 27.283MHz; and 40.660MHz to 40.700MHz.
3. An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.
4. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
5. Transmitters/antennas of this power-level are most likely mounted on an emergency vehicle chassis. The distances cited here are for open field. For an external antenna, the separation distance is most likely shorter.

Regulatory Information

FCC Warning Statement

FCC Part 15.19

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Part 15.21

Any changes or modifications (including the antennas) to this device that are not expressly approved by the manufacturer may void the user's authority to operate the equipment.

FCC Radiation Exposure Statement

This equipment complies with FCC Radiation exposure limits set forth for an uncontrolled environment. This device and its antenna must not be co-located or operating in conjunction with any other antenna or transmitter. This equipment should be installed and operated with a minimum distance of 10cm between the radiator and your body.

Additional Information

Warranty

The CREDO CPR BAND (Model: CREDO-CB) has a one-year limited Warranty.

Refer to the CREDO Global Warranty for terms and conditions.

About this edition

The information in this applies to the model CREDO-CB. This information is subject to change.

The CPR BAND (Model: CREDO-CB) is protected by Korea and international registered patents. The design of CPR BAND (Model: CREDO-CB) and its feedback symbols are protected by several jurisdictions under international design registration.

CPR BAND (Model: CREDO-CB) is a trademark or registered trademark of **CREDO**

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