





Quick Start Guide


eufy Breast Pump S1




Disassembling



1 Hold the pump hub and pull up to remove it.



2 Remove the diaphragm from the groove.




3 Unwrap the outer ring of the flange, pull the heating pin out of the oval notch, and then pull the entire flange out.



4 Flick the duckbill valve with your finger to remove it.

Cleaning and Sterilizing

All components except the pump hub are washable.




Do not use a microwave to sterilize.

- **Before first use:** Sanitize all washable parts. Sterilize the washable parts by boiling them in water for no more than 5 minutes. While cleaning, swirl the parts continuously to prevent them from being damaged by sticking to the container. You can also sterilize with steam.
- **For daily use:** Rinse all washable parts with warm, soapy water, then air dry.


Assembling




1 Hold the bottom of the container up, insert the duckbill valve into the small hole inside the container with your thumb and index finger and flick it to ensure it does not fall out.



2 Push the flange onto the container, making sure the heating pin fits in the oval notch, and the notch on the flange aligns with the milk pouring spout. Wrap and press the edges to ensure the seal is in place. Then check again that the outer ring is properly assembled to prevent milk leakage.




3 Insert the diaphragm into the groove of the container, with the concave side of the diaphragm facing the flange side.




4 Install the pump hub on the top of the container, pressing firmly to seal the gap between the pump hub and the diaphragm.

Wearing



Slightly lean forward to make sure your nipple is centered in the flange. The breast pump should be firmly against your breast.



Use a full-coverage nursing bra with a flap to create a tight seal between your breasts and the flange.

Two criteria will help you determine if the breast pump is positioned correctly:

- Whether you feel a noticeable suction that can pull out milk;
- Whether it is comfortable to wear.

If the wearing position doesn't meet either of these criteria, you need to adjust the position of the pump.

It takes a bit of time for your body to adapt to a new pump, so don't be discouraged if the results aren't good at the very beginning.


Setting Up the System

Download and install the eufy Baby app, then sign in or create an account.



Follow the in-app instructions to add this device and complete the setup.


Button Controls & Light Guide




Button Controls		
	Press and Hold for 2 Seconds	Turn Breast Pump On / Off (The breast pump will turn off automatically after pumping for 30 minutes or pausing for 10 minutes.)
	Press	Pause / Continue
	Press	Switch between Stimulation Mode and Expression Mode
	Press	Turn Heating On / Off
	Press	Increase / Decrease the Suction Level
	Press and Hold Simultaneously for 2 Seconds	L/R Pairing in the App

Light Guide		
●	White	Battery Level >20%
	Orange	Battery Level 20%- 40%
	Red	Battery Level <20%
	White	In Stimulation / Expression Mode
	White	In Heating Mode
●●●●●●●	1 to 7 Indicators Solid White	Suction Level: 1 to 7
●●●●●●	1 to 3 Indicators Solid Orange	Heat Level: 1 to 3
L R	Left / Right Indicator Solid White	Left / Right Side Paired in the App

Pouring Milk



1 Press and hold the power button to turn off, lean forward when taking out the pump to prevent leakage.




2 Remove the pump hub from the container.

3



Leave the container on the table ready for transfer.

4



Hold the container with the flange side facing up and turn it to pour the milk out. When the pouring spout faces down, turn the container further to empty every drop of milk.

⚠️ Avoid disassembling any parts other than the pump hub before transferring milk to prevent leakage.

Best to Know

1. To prevent leakage:

- When you complete assembling, check again that the outer ring of the flange is perfectly sealed.
- Make sure your nipple is at the center of the flange tunnel, and use a full-coverage nursing bra to create a tight seal between your breasts and the flange.
- Do not lean back while pumping. At the end of a session, lean forward to remove the pump.

2. Can't feel suction or have zero or very little output?

- Make sure the pump has power.
- Check that the assembly is correct and that the four pieces (duckbill valve, flange tunnel, diaphragm, and pump hub) are perfectly sealed.
- Adjust the wearing position, make sure the nipple is at the center of the flange tunnel, the pump is sealed to your breast, and you can feel

- noticeable suction.
- Try different suctions and rhythms.

3. The lifespan of the washable parts is 1-3 months; please replace them if you feel the suction and performance are decreasing, or if you have stored them for too long when planning to pump for your next child.

Customer Service



12-month limited warranty

(US) +1 (800) 988 7973
 (UK) +44 (0) 1604 936200
 (DE) +49 (0) 69 9579 7960
 (中国) +86 400 0550 036
 (日本) 03 4455 7823
 (AU) +61 3 8331 4800
 (Egypt) +20 8000000826
 (UAE) +971 8000320817
 (Kuwait) +965 22069086
 (KSA) +966 8008500030



support@eufy.com

Anker Innovations Limited
 Room 1318-19, Hollywood Plaza, 610 Nathan Road, Mongkok, Kowloon,
 Hong Kong



@EufyOfficial



@EufyOfficial



@EufyOfficial

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING

DANGER – To reduce the risk of electrocution:

- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for a product that has fallen into water. Unplug immediately.
- Always disconnect from the charger before using it.

WARNING – To reduce the risk of burns, electrocution, fire, or injury to persons.

- Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
- Keep the cord away from heated surfaces.
- Never use while sleeping or drowsy.
- Never drop or insert any object into any opening or hose.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- This is a personal care product for a single user and should not be shared between users.
- If irritation or discomfort occurs, use your finger to break the seal between the flange and the breast tissue, discontinue use and see a healthcare provider.
- Do not modify the product.
- This is not a toy, keep it out of the reach of children.
- Do not wear and use the product for over 35 minutes. Consult your healthcare practitioner to understand how long you should pump in a pumping session.
- Do not wash or sanitize the pump hub. Wipe with clean, damp cloth only.
- Do not place the breast pump in refrigerator or freezer. Do not use a microwave oven.
- Do not boil the milk container and all the silicon parts for more than 5 minutes every time.
- Do not use the product in an environment that is too hot or too cold, never expose the product to strong sunlight or a humid environment.
- The suitable temperature for this product and accessories is 0°C–40°C.
- It is recommended to charge the device in an environment with a temperature that ranges from 5°C–25°C.
- The expected service life of the washable parts is 1–3 months. It is

recommended to replace the washable parts after long-term storage or before pumping for your next child.

SAVE THESE INSTRUCTIONS

Notice

Notice (United States Of America)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: (1) Reorient or relocate the receiving antenna. (2) Increase the separation between the equipment and receiver. (3) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. (4) Consult the dealer or an experienced radio/ TV technician for help.

FCC Radio Frequency Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure conditions without restriction.

The following importer is the responsible party:
Company Name: POWER MOBILE LIFE, LLC
Address: 10900 NE 8th St, Ste 501, Bellevue, WA 98004
Telephone: 1-800-988-7973

Notice (Canada)

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes:

- (1) l'appareil ne doit pas produire de brouillage, et
- (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

This Class B digital apparatus complies with Canadian ICES-003. Cet appareil numérique de la classe B est conforme à la norme NMB-003 du Canada.

IC RF Exposure Statement:

This device meets the government's requirements for exposure to radio waves.

The exposure standard for wireless device employs a unit of measurement known as the Specific Absorption Rate, or SAR. The SAR limit set by the IC is 1.6 W/kg. Tests for SAR are conducted using standard operating positions (0mm) accepted by the IC with the device transmitting at its highest certified power level in all tested frequency bands.

The highest reported SAR value for your device is: 0.112 W/kg
Cet appareil est conforme aux exigences gouvernementales en matière d'exposition aux ondes radio.

La norme d'exposition pour les appareils sans fil utilise une unité de mesure connue sous le nom de débit d'absorption spécifique (DAS). La limite de DAS fixée par le CI est de 1,6 W/kg. Les tests de DAS sont effectués en utilisant les positions de fonctionnement standard (0 mm) acceptées par l'IC, l'appareil émettant à son niveau de puissance certifié le plus élevé dans toutes les bandes de fréquences testées.

La valeur SAR la plus élevée rapportée pour votre appareil est de 0,112 W/kg

Notice (Europe and United Kingdom)

 This product complies with the radio interference requirements of the European Community.

Declaration of Conformity

Hereby, Anker Innovations Limited declares that this device is in compliance with the essential requirements and other relevant provisions

of Directive 2014/53/EU. For the declaration of conformity, visit the website: <https://support.eufy.com/s/article/Recommend?type=Download&secondType=doc>.

Hereby, Anker Innovations Limited declares that this device is in compliance with The Product Security and Telecommunications Infrastructure (Security Requirements for Relevant Connectable Products) Regulations 2023. For the Statement of Compliance, visit the Web site <https://www.eufy.com/uk/psti-related>.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safety.

Children should be supervised to ensure that they do not play with the appliance.

RF exposure information:

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure conditions without restriction.

Bluetooth Operating Frequency Range: 2402~2480MHz;

Bluetooth Max Output Power: 2.11 dBm (EIRP)

Anker Innovations Deutschland GmbH | Georg-Muche-Strasse 3, 80807 Munich, Germany



This product is designed and manufactured with high quality materials and components, which can be recycled and reused.



This symbol means the product must not be discarded as household waste, and should be delivered to an appropriate collection facility for recycling. Proper disposal and recycling helps protect natural resources, human health and the environment. For more information on disposal and recycling of this product, contact your local municipality, disposal service, or the shop where you bought this product.

Warning:

- Replacing a battery with an incorrect type may bypass safeguards and create danger;
- Disposal of a battery in a fire or other source of significant heat, or by crushing, or cutting, may result in an explosion as well as leakage of flammable liquids and gases;
- Leaving a battery in an extremely hot environment may result in an explosion as well as leakage of flammable liquids and gases;
- A battery subjected to extremely low air pressure may result in an explosion as well as leakage of flammable liquids and gases.

Attention:

- Le remplacement d'une batterie par un type incorrect peut contourner les protections et créer un danger;
- L'élimination d'une pile dans un feu ou une autre source de chaleur importante, ou par écrasement ou découpage, peut entraîner une explosion ainsi qu'une fuite de liquides et de gaz inflammables;
- Le fait de laisser une batterie dans un environnement extrêmement chaud peut entraîner une explosion ainsi qu'une fuite de liquides et de gaz inflammables;
- Une batterie soumise à une pression d'air extrêmement faible peut provoquer une explosion ainsi qu'une fuite de liquides et de gaz inflammables.



CAUTION!
THERE IS A RISK OF EXPLOSION IF AN INCORRECT BATTERY TYPE IS USED.
DISPOSE OF USED BATTERIES ACCORDING TO THEIR INSTRUCTIONS.



ATTENTION!
IL EXISTE UN RISQUE D'EXPLOSION EN CAS D'UTILISATION D'UN TYPE DE PILE INCORRECT.
ÉLIMINEZ LES PILES USAGÉES CONFORMÉMENT AUX INSTRUCTIONS.

Anker Innovations Limited
Room 1318-19, Hollywood Plaza, 610 Nathan Road, Mongkok, Kowloon,
Hong Kong