
Anatel PAT700 - Safety Notice

Reference: D02736, D02737, D02740, D21313, D21315

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	59.7 x 22.9 x 25.4 cm (23.5 x 9 x 10 in.)
Enclosure	Quick-connect wiring model: IP46 Conduit wiring model: IP56
Weight	13.6 kg (30 lb)
Pollution degree	2
Installation category	II
Protection class	I
Power requirements	100-240 VAC, $\pm 10\%$, 65 W, 50/60 Hz
Altitude	4000 m (13,125 ft) maximum
Operating temperature	10 to 40 °C (50 to 104 °F)
Humidity	5–95%, non-condensing
Inlet temperature	1 to 95 °C (34 to 203 °F)
Sample cell temperature	1 to 65 °C
Sample flow rate	60 to 300 mL/min
Sample pressure	0.70 to 7.0 kg/cm ² (10 to 100 psi)
Temperature accuracy	± 0.5 °C
UV lamps	2x, with UV Detect technology
User interface	Color touch screen
TOC	Range: 0.5 to 2000 ppb as carbon
	Accuracy: ± 3 ppb or $\pm 5\%$, whichever is greater
	Repeatability: ± 0.3 ppb or $\pm 1\%$, whichever is greater
	Detection limit: 0.5 ppb
Conductivity	Range for conductivity mode: 0.05 to 150 $\mu\text{S}/\text{cm}$ at 25 °C
	Accuracy: $\pm 2\%$ over the full range
	Range for TOC mode: 0.05–5.0 $\mu\text{S}/\text{cm}$
Outputs	Analog: Three isolated
	Discrete: Four isolated, configurable
Inputs	Discrete: Two isolated
Certifications	CE Certified. Certified to US and Canadian standards by CSA.
Warranty	1 year

General Information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety Information

IMPORTANT The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of Hazard Information

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

Indicates a potentially hazardous situation that may result in minor or moderate injury.

IMPORTANT Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Chapter 3 Installation > Installation Guidelines

DANGER

Fire hazard. This product is not designed for use with flammable liquids.

WARNING

Electrical shock hazard. Externally connected equipment must meet applicable country safety requirements.

WARNING

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

IMPORTANT To prevent measurement error, make sure that the instrument is attached in an upright, level position.

CAUTION

Risk of instrument damage. This device is intended for indoor use only. To avoid device damage, do not install the instrument outdoors.

Chapter 3 Installation > Mechanical Installation

Attach the Instrument to a Wall

WARNING

Personal injury hazard. Make sure that the wall mounting is able to hold 4 times the weight of the equipment.

Chapter 3 Installation > Electrical Installation

Wiring Safety Information

 DANGER

Electrocution hazard. Always remove power to the instrument before making electrical connections.

Wiring for Power

 DANGER

Electrocution hazard. Protective Earth Ground (PE) connection is required.

 DANGER

Electrocution hazard. This device is intended for indoor use only. If used in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used for connecting the equipment to its main power source.

About Quick-Connect Wiring Connections

 DANGER

Electrocution and fire hazards. For instruments with a power cord, the power cord serves as the primary disconnecting device. Do not position the instrument so that it is difficult to access the disconnection device (power cord).

About Conduit Wiring Connections

 WARNING

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Electrical shock and fire hazards. For conduit installation, clearly identify the local power disconnect and make sure there is easy access to it.

Chapter 3 Installation > Plumbing

Plumb the Instrument

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Chapter 8 Maintenance

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.

Clean Leaks and Spills

WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

WARNING

Personal injury hazard. The sample line contains water under high water pressure that can burn skin if hot. Qualified personnel must remove the water pressure and wear personal protective equipment during this procedure.

Replace a Fuse

DANGER

Fire hazard. Use the same type and current rating to replace fuses.

Water Filter Replacement

WARNING

Burn and pressure hazard. Water is pressurized and may be hot (up to 95 °C [203 °F]). Protective eye wear is required while performing this task, and the sample line must first be depressurized before the water filter is disconnected.

Send the Instrument for Service

WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

Safety Precautions

WARNING

Risk of operator injury if:

- All doors covers and panels are not closed and secured in place prior to and during instrument operation.
- The integrity of safety interlocks and sensors is compromised.
- Instrument alarms and error messages are not acknowledged and acted upon.
- You contact moving parts.
- You mishandle broken parts.
- Doors, covers and panels are not opened, closed, removed and/or replaced with care.
- Improper tools are used for troubleshooting.

To avoid injury:

- Keep doors, covers and panels closed and secured in place while the instrument is in use.
- Take full advantage of the safety features of the instrument.
- Acknowledge and act upon instrument alarms and error messages.
- Keep away from moving parts.
- Report any broken parts to your Beckman Coulter Representative.
- Open/remove and close/replace doors, covers and panels with care.
- Use the proper tools when troubleshooting.

CAUTION

System integrity could be compromised and operational failures could occur if:

- This equipment is used in a manner other than specified. Operate the instrument as instructed in the product manuals.
- You introduce software that is not authorized by Beckman Coulter into your computer. Only operate your system's computer with software authorized by Beckman Coulter.
- You install software that is not an original copyrighted version. Only use software that is an original copyrighted version to prevent virus contamination.

 CAUTION

If you purchased this product from anyone other than Beckman Coulter or an authorized Beckman Coulter distributor, and, it is not presently under a Beckman Coulter service maintenance agreement, Beckman Coulter cannot guarantee that the product is fitted with the most current mandatory engineering revisions or that you will receive the most current information bulletins concerning the product. If you purchased this product from a third party and would like further information concerning this topic, contact your Beckman Coulter Service Representative.

 CAUTION

Risk of instrument damage. This device is intended for indoor use only. To avoid device damage, do not install the instrument outdoors.

 WARNING

Risk of personal injury. Safety protection can be impaired if used in a manner not specified by the manufacturer. To avoid personal injury, use the instrument according to the manufacturer's instructions only.

General Warnings and Cautions

 DANGER

Chemical or biological hazards. If this instrument is used to monitor a treatment process and/or chemical feed system for which there are regulatory limits and monitoring requirements related to public health, public safety, food or beverage manufacture or processing, it is the responsibility of the user of this instrument to know and abide by any applicable regulation and to have sufficient and appropriate mechanisms in place for compliance with applicable regulations in the event of malfunction of the instrument.

 DANGER

Fire hazard. Use the same type and current rating to replace fuses.

 DANGER

Fire hazard. This product is not designed for use with flammable liquids.

 DANGER

Electrocution hazard. This device is intended for indoor use only. If used in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used for connecting the equipment to its main power source.

 DANGER

Electrocution hazard. Always remove power to the instrument before making electrical connections.

 DANGER

Electrocution and fire hazards. For instruments with a power cord, the power cord serves as the primary disconnecting device. Do not position the instrument so that it is difficult to access the disconnection device (power cord).

 DANGER

Electrocution hazard. Protective Earth Ground (PE) connection is required.

 DANGER

Electrical shock and fire hazards. For conduit installation, clearly identify the local power disconnect and make sure there is easy access to it.

 DANGER

Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.

 WARNING

Electrical shock hazard. Externally connected equipment must meet applicable country safety requirements.

 WARNING

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

 WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

WARNING

Personal injury hazard. Make sure that the wall mounting is able to hold 4 times the weight of the equipment.

WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

WARNING

Personal injury hazard. The sample line contains water under high water pressure that can burn skin if hot. Qualified personnel must remove the water pressure and wear personal protective equipment during this procedure.

WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

WARNING

Burn and pressure hazard. Water is pressurized and may be hot (up to 95 °C [203 °F]). Protective eye wear is required while performing this task, and the sample line must first be depressurized before the water filter is disconnected.

EMC

This device complies with the emissions and immunity requirements as specified in the EN/IEC 61326 series of Product Family Standards for a “basic electromagnetic environment.” Such equipment is supplied directly at low voltage from public mains network. This equipment is not intended for residential use.

This device generates, uses, and can radiate unintentional radio-frequency (RF) energy. If this device is not installed and operated correctly, this RF energy can cause interference with other equipment. It is the responsibility of the end user to be sure that a compatible electromagnetic environment for the device can be maintained so that the device operates as intended.

This equipment is designed for use in a commercial environment. It is likely to perform incorrectly if used in a domestic environment. If it is suspected that performance is affected by electromagnetic interference, correct operation may be restored by increasing the distance between the equipment

and the source of the interference. In addition, other equipment can radiate RF energy to which this device is sensitive. If one suspects interference between this device and other equipment, Beckman Coulter recommends the following actions to correct the interference:

- Evaluate the electromagnetic environment before installation and operation of this device.
- Do not operate this device close to sources of strong electromagnetic radiation (for example: unshielded intentional RF sources), as these can interfere with proper operation. Examples of unshielded intentional radiators are hand-held radio transmitters, cordless phones, and cellular phones.
- Do not place this device near medical electrical equipment that can be susceptible to malfunctions caused by close-proximity to electromagnetic fields.
- This device has been designed and tested to CISPR 11, Class A emission limits. In a domestic environment, this device can cause radio interference, in which case, you need to take measures to mitigate the interference.

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.

2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

EU

Beckman Coulter, Inc. hereby ensures and declares that the product conforms to the current European Union directives for EMC, Low Voltage, RoHS, and Radio Equipment Device directives by meeting the requirements of the applicable harmonized standards. The Official E.U. Declaration of Conformity can be found at www.beckman.com under “/Support/Technical documents/CE Declaration-Hardware”.

Taiwan

Without permission granted by the NCC, any company, enterprise, or user is not allowed to change frequency, enhance transmitting power, or alter the original characteristic as well as performance of approved low power radio-frequency devices. Low power radio-frequency devices shall not affect aircraft safety nor interfere with legal communications. If interference is found, the user shall cease operation immediately until improvements eliminate the interference. The aforementioned legal communications refer to radio communications operated in accordance with the provisions of the Telecommunications Management Act. Low power radio-frequency devices must not be susceptible to interference from legal communications or ISM radio wave radiating devices.

Mexico

The operation of this equipment is subject to the following two conditions: (1) it is possible that this equipment or device may not cause harmful interference, and (2) this equipment or device must accept any interference, including interference that may cause undesired operation.

Brazil

This equipment is required to be connected to electricity mains socket that has grounded (three-prong), according to ABNT NBR 5410, seeking the safety of users against electric shock.

This equipment is not entitled to protection against harmful interference and may not cause interference in duly authorized systems.

For queries, visit: www.anatel.gov.br.

RFID Module

This instrument contains an internal radio frequency identification device (RFID) certified for the countries in which it will be marketed. Certification IDs are listed on the exterior label of the instrument.

Table 1 Registration Information

Parameter	Value (D02736, D02740)	Value (D21313, D21315)
USA (FCC)	FCC ID: 2AOSQM1RFID23	FCCID: 2AOSQM1MINI
Brazil (ANATEL)	TBD	TBD
Canada (ISED)	IC: 23864-M1RFID23	IC: 23864-M1RFID23
Dominican Republic	TBD	TBD
Egypt	TBD	TBD
India (WPC)	TBD	TBD
Israel (MOC)	TBD	TBD
Japan Radio Law (IFERW)	TBD	TBD
Malaysia (SIRIM)	TBD	TBD
Mexico (IFETEL)	TBD	TBD
Saudi Arabia (CITC)	TBD	TBD
Singapore (IMDA)	TBD	TBD
South Africa (ICASA)	TBD	TBD
South Korea (KCC)	TBD	TBD
Taiwan (NCC)	TBD	TBD
Frequency	13.56 MHz $\pm$ 7 kHz	13.56 MHz $\pm$ 7 kHz
RF output power	<180 mW	<180 mW

Precautionary Labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. NOTE For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal. The symbol of a crossed-out wheeled bin on the product is required in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive of the European Union. The presence of this marking on the product indicates: 1. that the device was put on the European Market after August 13, 2005 and 2. that the device is not to be disposed via the municipal waste collection system of any member state of the European Union. For products under the requirement of WEEE directive, please contact your dealer or local Beckman Coulter office for the proper decontamination information and take-back program which will facilitate the proper collection, treatment, recovery, recycling, and safe disposal of device. For the Japan market: This system is considered an industrial waste, subject to special controls for infectious waste. Before disposal of the system, refer to the Waste Disposal and Public Cleaning Law for compliance procedures.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	The "RCM" (Regulatory Compliance Mark) is depicted as a triangle with a partial circle and check. The mark is applied to products that comply with the EMC requirements of the Australian Communications Media Authority (ACMA) for use in Australia and New Zealand.
	This symbol indicates the need for protective eye wear.
	This symbol indicates that the marked item requires a protective earth connection. If the instrument is not supplied with a ground plug on a cord, make the protective earth connection to the protective conductor terminal.
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.

	This symbol indicates the presence of mercury. Components that contain mercury must be disposed in accordance with all applicable local and regional law. Storage Requirements for Used UV Lamp Assemblies Put used UV Lamp Assemblies into the replacement UV Lamp Assembly shipping container, or a similar container. It is preferable that the original packing materials be used where possible, or materials adequate to prevent breakage during storage and transportation. Boxes of used UV Lamp Assemblies should be labeled as such and stored in a location where the potential for accidental breakage is minimized. A UV lamp re-cycler may have specific procedures and UV lamp storage requirements. Consult with a UV lamp re-cycler to determine all applicable policies. This component contains Mercury. Dispose according to Local, State, or Federal Laws.
	This symbol indicates a puncture or pinch hazard. Keep hands and fingers away.
	This symbol indicates the presence of a UV light source that may cause eye and skin injury. Wear appropriate protective equipment and obey all safety protocols.
	A "CE" mark indicates that a product has been assessed before being placed on the market, and has been found to meet European Union safety, health, and/or environmental protection requirements.
	This label indicates that the electronic information product contains certain toxic or hazardous substances. The center number is the Environmentally Friendly Use Period (EFUP) date, and indicates the number of calendar years the product can be in operation. Upon the expiration of the EFUP, the product must be immediately recycled. The circling arrows indicate the product is recyclable. The date code on the label or product indicates the date of manufacture.
	A "UKCA" mark indicates that a product has been assessed before being placed in the UK market, and has been found to meet UK safety, health, and/or environmental protection requirements.
	This symbol indicates recognition by the Canadian Standards Association (Nationally Recognized Test Laboratory or NRTL) that the instrument has met the relevant product safety standards.
	The "RCM" (Regulatory Compliance Mark) is depicted as a triangle with a partial circle and check. The mark is applied to products that comply with the EMC requirements of the Australian Communications Media Authority (ACMA) for use in Australia and New Zealand.
	This product can expose you to chemicals known to the State of California to cause Cancer and Reproductive Harm. For more information go to https://www.P65Warnings.ca.gov.
	To indicate the presence of mercury. Components that contain mercury must be disposed in accordance with all applicable national, local, and regional law.

Product Overview

Chemical or biological hazards. If this instrument is used to monitor a treatment process and/or chemical feed system for which there are regulatory limits and monitoring requirements related to public health, public safety, food or beverage manufacture or processing, it is the responsibility of the user of this instrument to know and abide by any applicable regulation and to have sufficient and appropriate mechanisms in place for compliance with applicable regulations in the event of malfunction of the instrument.

This analyzer uses UV light to oxidize water samples for TOC analysis in pure and ultra-pure water process applications. The analyzer includes features that support compliance to 21 CFR Part 11 with password protection and an audit trail of user actions.

For power and data collection, the analyzer connects to the facility systems through quick-connect connections.

The analyzer has an OASIS bottle bay to measure individual grab samples or do a calibration/validation. Bottle information is read automatically with the radio frequency identification device (RFID) wireless system and RFID tags on the bottle (Not all models include the RFID wireless system).

Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 U.S.A.

Beckman Coulter Eurocenter S.A.
22, rue Juste-Olivier
Case Postale 1044
CH - 1260 Nyon 1, Switzerland
Tel: +41 (0) 22 365 36 11

Initial Issue AA, 7/2024

Issue AB, 4/2025

© 2025 Beckman Coulter, Inc.
All rights reserved.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

All other trademarks are the property of their respective owners.

May be covered by one or more pat. - see www.beckman.com/patents.

Find us on the web at: www.beckman.com.


Glossary of Symbols is available at beckman.com/techdocs (PN C24689).

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Versione AB, 4/2025