


## External Components



1	Charge connector	6	Recording LED
2	Conductors	7	Recording Button
3	Power button	8	Power LED
4	Cable connector	9	Transferring LED
5	Charging LED		

## Introduction

This device is able to calculate your heart rate, heart conditions, stress and nervous index. Provide instant health information for you.

## Contents

- device x1
- Charger cable x1
- Instruction(this page) x1

## Before using

Please charge the device with the USB cable for 3 hours. The charging LED will turn on when charging and turn off after the battery is full.

**Attention: you may not operate the device when charging.**

## Turn on Device

The device will turn on after pressing the power button for 2 seconds. The power LED will light up and the Transferring LED will blink slowly in standby mode.

## Before measuring

Put the device to your chest. The 4 conductors should contact your body solid.



Please also put on some cream on your chest and conductors before measuring.

## Start measuring

Press the recording button when the device is solid on your chest. The device will beep once. After several seconds, the device beeps constantly for a good signal recorded. After 30 seconds, the device beeps twice for end recording and turns into standby mode.

If you need terminate the measurement, press the recording button again and the device will turn into standby mode.

You can utilize the Cmate Gen2 APP for data analysis service.

## Download APP

Find the QR code below for APP download.

Android	iOS


## Bluetooth Pairing

Android user need to pairing the device first.

1. Turn on the device (standby mode).
2. Explore Bluetooth settings in your phone.
3. Scan the nearby Bluetooth device, find "CmateH" device and finish the pairing

## Data uploading

1. Explore the Cmate Gen2 APP
2. Point at [file uploading] button and select the file for analysis by recording time.
3. Wait until server responds the result. If it shows "poor signal" please take another measurement and upload again.
4. All result will stay in the [my report page]. The latest analysis result will show on the main page of the APP.



## Power Management

The device will automatically turn off after 3 minutes standby for the sake of power saving. You may push the power button for 2 seconds for manual turn off.

When the battery is low, the device will beep once by every 5 seconds. The Power LED also blinks every second. The device is still workable yet we advise you to charge the device immediately.

When the battery is extreme low, the device will shut down automatically. Please charge the device for 3 hours immediately.

### Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### FCC Caution:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### Non-modification Statement:

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## Error Messages

LED fast blinks and beeps twice every seconds	Device memory full, please upload or delete the file with APP. If the error still occur, please contact your sales channel for further help
The device is not able to charge	please contact your sales channel for further help

**Brand Name:**Cmate Model:CHI-100

一般設備(低功率電波輻射性電機管理辦法第 12、14 條)

—— 經型式認證合格之低功率射頻電機，非經許可，公司、商號或使用者均不得擅自變更頻率、加大功率或變更原設計之特性及功能。  
—— 低功率射頻電機之使用不得影響飛航安全及干擾合法通信；經發現有干擾現象時，應立即停用，並改善至無干擾時方得繼續使用。前項合法通信，指依電信法規定作業之無線電通信。低功率射頻電機需忍受合法通信或工業、科學及醫療用電波輻射性電機設備之干擾。

### Radiation Exposure Statement:

The device has been evaluated to meet general RF exposure requirement, The device can be used in portable exposure condition without restriction.

Revlis Biotech Co., Ltd.  
7F., No 108-1, Minquan Rd., Xindian Dist., New Taipei City, Taiwan,  
Tel : +886-2-22198878  
FAX : +886-2-22192587