BrainCo Inc. All rights reserved.

BrainCoTM

USER MANUAL

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① Front Indicator Light

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
② Micro USB Port

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③ Replaceable Hydrogel
Ear Electrode

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1

④ Replaceable Hydrogel
Forehead Electrode

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


Power Button 电源键

Power On

Remove the charging cable, press and hold the power button until all of the 5 white indicator lights turn on; Note: Headband cannot be used if the battery is depleted;

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Power Off

When in use, press and hold the power button until all of the 5 white indicator lights turn off;

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Battery Checking

While not charging, quickly press the power button and release, the side indicators glow according to the current battery level;

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WiFi Reset (for Retail version only)

When in use, press and hold the power button until all of the 5 white indicator lights turn OFF and turn ON AGAIN;

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
Notice: If the Headband does not work properly, remove the charging cable, press and hold the power button for 20 seconds, then release the power button and the Headband will go into Power OFF status;

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
Side Indicator Light
侧边指示灯

	power on 开机	turn on one by one from left to right 从左到右逐一亮起
	power off 关机	turn off one by one from right to left 从左到右逐一亮起
	pairing 配对	cycle back and forth repeatedly 往复循环
	working 正常工作	flash from left to right individually 从左到右逐个闪烁
	charge 充电	Each fully lit indicator represents 20% of the total power. A blinking light shows the current charging status. 每个稳定亮起的灯 表示20%，闪烁的 灯显示正在充电的进程
	reset 重置	Flash from middle to side 从中间至两边逐个亮起

PUTTING ON


Put the Focus 1 on your head and adjust the Headband to ensure the correct orientation for proper fit and correct electrode placement.

Notice: Make sure that both electrodes and your skin are clean and have direct contact with each other.

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CHARGING


To charge the Headband, plug the charging cable into the USB port on your computer or other power source. The LED indicator lights will progressively illuminate on your Focus 1 as it charges. Each steady light represents 20% of the maximum charge. Once all five lights glow solid white light, the Headband will be fully charged.

Notice: Please use the USB cable included in the package to avoid any charging problems.

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
ELECTRODE MAINTAINANCE

Please wipe your electrodes weekly with the sterile saline wipe included to keep the electrodes clean and extend their lifetime.

Notice:

1. Store at room temperature (25°C, 50%), and do not put the electrodes into water or fire.
2. Avoid contact acidic, basic, toxic, and erosive materials and chemicals.
3. Please follow the instruction in the picture below, and wipe the electrode gently; please contact the agency to purchase more wipes.

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
REPLACING ELECTRODES

Replace the hydrogel electrodes when necessary to ensure signal quality. The recommended time is 6 months.

Notice:

1. While replacing the electrodes, please pay attention to the groove on the bottom of the electrode that shows the proper placement direction. Please contact the agency to purchase more.
2. Do not put the used electrodes into water or fire.

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SPECIFICATIONS

Weight	95g
Head Circumference	55-75cm
Operating Temperature	0 -40°C(32-104°F)
Battery Capacity	400 mAh
Battery Life	Up to 4 hours
Battery Charge Time	1-1.5 hour
Wireless Connection	Wi-Fi IEEE 802.11 b/g/n 2.4GHz
Electrode Material	Hydrogel

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SAFETY INSTRUCTION

Before using your Focus 1 headband, please read the following safety instructions. This will ensure correct usage. After reading this manual, store it in a safe place for future reference.

IMPORTANT:

Please use the accessories included in your purchase. Due to voltage and current requirements (5V,200mA), please use the USB cable that was included in your Product Pouch.

WARNING:

1. Focus 1 headband contains hydrogel electrodes. Please see a doctor if swallowed.
2. Please store your Focus 1 headband in the Product Pouch when it is not in use. Keep children and pets away from the product.
3. Please follow the instructions below before using your Focus 1 headband:
 - Before using, please read the User Manual carefully.
 - Watch your child when she/he is using this product.
 - Do not expose this product to water, rain, or any liquid.
 - Use of a power supply or charger not recommended or sold by an authorized manufacturer may result in a risk of fire or injury.
 - Do not attempt to use a damaged or modified Focus 1 headband. A damaged or modified Focus 1 headband may cause a fire or unknown danger.
 - Do not dispose of your Focus 1 headband in a fire or heat. The product could explode in a fire or high temperatures over 100°C.
 - Servicing can only be performed by a qualified repair person using identical replacement parts.
4. This product contains small parts which may be a choking hazard. Not suitable for children under age 6.
5. Do not tamper with your Focus 1 headband without authorization.
6. Use this product only with an agency approved power supply which meets local regulatory requirements (e.g., UL, CSA, VDE, CCC).
7. Do not expose this product, which contains batteries, to excessive heat (e.g. storage in direct sunlight, fire or any similar environment).
8. If you are under 18 years old, please read this document with your guardians and pay attention to all warnings.

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DISCLAIMER

Since this manual lists safety instructions in detail, please read it carefully before using your Focus 1 headband.

Important warning:

1. Do not allow children under 6 to use this product.
2. Children above 6 are allowed to use this product with a parent's or teacher's supervision to prevent injuries.
3. The product and its accessories including the hydrogel electrodes are all inedible.

The responsibility for any risks and consequences arising from the violation of this Safety Manual are borne by the users and their guardians. The manufacturer, developer, and seller assume no responsibility for any of those risks and consequences.

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UNITED STATES:

FCC Regulator Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End user must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

The portable device is designed to meet the requirements for exposure to radio waves established by the Federal Communications Commission (USA). These requirements set a SAR limit of 1.6 W/kg averaged over one gram of tissue. The highest SAR value reported under this standard during product certification for use when properly worn on the body with the minimum distance 0mm.

CANADA:
IC Regulator Statement

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

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This equipment complies with ISSED radiation exposure limits set forth for an uncontrolled environment. End user must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

The portable device is designed to meet the requirements for exposure to radio waves established by the ISSED. These requirements set a SAR limit of 1.6 W/kg averaged over one gram of tissue. The highest SAR value reported under this standard during product certification for use when properly worn on the body with the minimum distance 0mm.

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EUROPE:

BrainCo Inc. hereby declares that this product is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU and all other EU directive requirements.

This symbol means the product must not be discarded as household waste, and should be delivered to an appropriate collection facility for recycling. Proper disposal and recycling helps protect natural resources, human health and the environment. For more information on disposal and recycling of this product, contact your local municipality, disposal service.