

GLUCONEXT - English - Welcome guide

DIABNEXT™
**GLUCONEXT™
INFRARED**
Welcome guide
www.diabnext.com

English

GLUCONEXT™ INFRARED description

Model No. GI1

Battery cover

LED (Red/Green)

Blue button
(Validation/Transmission)

Infrared transmitter

Connection method with BGM	Infrared
Wireless connection	Bluetooth
Device weight	9 g
Device size	L 39.2 x W 27.3 x H 13.6 mm
Battery type	Lithium coin cell battery (CR2032)
Battery life	Over 1000 transfers

Device LED signal description

Slow green blinkings (3 seconds): waiting for Bluetooth connection

Quick green blinkings (1 second): connected to mobile device

5 Red blinkings (0.25 second): battery level low

If you have any questions about the use of the GLUCONEXT™ INFRARED, do not hesitate to visit our website www.diabnext.com or to contact us by email: support@diabnext.com. Our team is here to help you.

www.diabnext.com

Set up

1

Download DIABNEXT™ APP

D

DIABNEXT

Download on the App Store

GET IT ON Google Play

f

facebook

G

Google

Install the battery in the device. Turn on the Bluetooth of your smartphone. To pair your device with the mobile app, open your DIABNEXT™ app, then press the **CENTRAL BLUE BUTTON** on your GLUCONEXT™ INFRARED. Then, follow the instructions in your DIABNEXT™ app.

折疊線

DIABNEXT

Automatic synchronization
of the GLUCONEXT™ INFRARED
with the app

After pairing your device, the data will be automatically transmitted if you keep your smartphone's Bluetooth turned on and the DIABNEXT™ app running in the background*.
*Even if not open, an app continues to run in the background of your smartphone. If you quit the app on the other hand, it won't receive the data automatically.

Using the GLUCONEXT™ INFRARED

2

Consult your blood glucose meter INSTRUCTION MANUAL and turn on your blood glucose meter transmission mode.

2

Then press the GLUCONEXT™ INFRARED **BLUE BUTTON** to validate and transfer the data automatically to your logbook*.

*The infrared transmitter of the GLUCONEXT™ INFRARED is in the bottom of the device. Please confirm the infrared transmitter location of your blood glucose meter, and be sure to have it close to the infrared transmitter (within 15cm) of your GLUCONEXT™ INFRARED for proper syncing.

4

Verify the Bluetooth of your smartphone is turned on.

Bluetooth

3

The green light will be quick blinking to signal that the GLUCONEXT™ INFRARED is connected to DIABNEXT™ APP. Your digital logbook will be auto-filled by GLUCONEXT™ INFRARED after a few seconds*.

*It is not necessary to re-open the DIABNEXT™ application to automatically transfer your data. If you have not quit it, your application is still running in the background and can receive the data transmitted by the GLUCONEXT™ INFRARED.

5

Enjoy your digital logbook & communities.

6

WGD-GI1E1-A-1901

FCC Caution:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.