PH304/PS304


Instructions for Use **EN**

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GYN-Pump

EN

This instructions for use contains information that is subject to copyright. All rights reserved. This instructions for use should not be photocopied, duplicated on microfilm or otherwise copied or distributed, completely or in part, without the approval of W.O.M. WORLD OF MEDICINE GmbH.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in this instructions for use bear registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or unclarity in this instructions for use.

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Copyright © W.O.M. WORLD OF MEDICINE GmbH

Manufacturer/Hersteller/Fabricant/Fabrikant/Производитель:

W.O.M. WORLD OF MEDICINE GmbH
Salzufer 8
10587 Berlin, Germany
Phone: +49 30 39981-550
Fax: +49 30 39981-545
E-mail: info.berlin@wom.group

W.O.M. WORLD OF MEDICINE GmbH
Alte Poststraße 11
96337 Ludwigsstadt, Germany
Phone: +49 9263 877-0
Fax: +49 9263 877-152
E-mail: info.ludwigsstadt@wom.group

CE0197

CE marking according to Directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
EG-markering conform Richtlijn 93/42/EEG
Маркировка CE согласно директиве 93/42/EEC

Type: PH304/1201563/10000020621 00/2017-09/Kubetzek

EN

DE

FR

NL

RU

	Description	Beschreibung	Description	Beschrijving	Описание
	Caution	Vorsicht; Achtung	Prudence; Attention	Voorzichtig; Let op	Осторожно; внимание
	Follow instructions for use (white image on a blue background)	Gebrauchsanweisung befolgen (weißes Bild auf blauem Grund)	Respecter le mode d'emploi (image blanche sur fond bleu)	Gebruiksaanwijzing opvolgen (witte afbeelding op blauwe ondergrond)	Соблюдать инструкцию по применению (белое изображение на синем фоне)
	Consult instructions for use	Gebrauchsanweisung beachten	Respecter les instructions d'utilisation	Gebruiksaanwijzing opvolgen	Соблюдать инструкцию по применению
	Type BF applied part	Anwendungsteil des Typs BF	Appareil avec les parties appliquées de type BF	Toepassingsgedeelte type BF	Рабочая часть типа BF
	Equipotentiality	Potentialausgleich	Connexion équipotentielle	Potentiaalvereffening	Выравнивание потенциалов
IP 21	Degrees of protection provided by enclosures (IP-Code)	Gehäuseschutzklasse (IP-Code)	Degrés de protection procurés par les enveloppes (code IP)	Beschermingsklasse (IP-code) behuizing	Класс защиты корпуса (IP-код)
	Alternating current	Wechselstrom	Courant alternatif	Wisselstroom	Переменный ток
	Service	Service	Service	Service	Сервис
	Connection for Canister Scale	Anschluss Behälterwaage	Raccord pour unité de pesage à récipient	Aansluiting reservoirweegstelsysteem	Подключение весов для взвешивания резервуаров
	Input/Output	Eingang/Ausgang	Entrée/Sortie	Ingang/uitgang	Вход/выход
	Catalogue number	Artikelnummer	Numéro de commande	Bestelnummer	Номер для заказа
	Do not resterilize	Nicht resterilisieren	Ne pas restériliser	Niet opnieuw steriliseren	Не стерилизовать
	Do not reuse	Nicht zur Wiederverwendung	Ne pas réutiliser	Niet voor hergebruik	Не для повторного использования
	Sterilized using ethylene oxide	Sterilisiert mit Ethylenoxid	Stérilisé à l'oxyde d'éthylène	Gesteriliseerd met ethyleenoxide	Стерилизовано этиленоксидом
	Batch code	Chargencode	Codes de lot	Chargencode	Обозначение партии
	Serial number	Seriennummer	Numéro de série	Serienummer	Серийный номер
	Date of manufacture (YYYY-MM-DD)	Herstellungsdatum (JJJ-MM-TT)	Date de fabrication (AAAA-MM-JJ)	Fabricagedatum (YYYY-MM-DD)	Дата изготовления (ГГГГ-ММ-ДД)



EN

DE

FR

NL

RU

	Description	Beschreibung	Description	Beschrijving	Описание
	Use by date (YYYY-MM-DD)	Verwendbar bis (JJJ-MM-TT)	Utilisable jusqu'à (AAAA-MM-JJ)	Te gebruiken tot (YYYY-MM-DD)	Использовать до (ГГГГ-ММ-ДД)
R_x ONLY	Authorized for Sale or use by Physician only	Nur für autorisiertes Vertriebspersonal oder Arzt	Uniquement destiné à un personnel commercial ou un médecin autorisé	Uitsluitend voor bevoegd personeel of arts	Только для врача или авторизованного персонала
	Number of autocla- ving cycles	Anzahl der Autoklavie- rungen	Nombre d'autoclavages	Aantal behandelingen in de autoclaaf	Число циклов автоклавирования
	Quantity	Menge	Quantité	Hoeveelheid	количество
	Not made with phthal- ates	Enthält keine Phthalate	Contient aucun phtalates	Bevat geen ftalaten	Не содержит фталатов
	Contains DEHP	Enthält DEHP	Contient du DEHP	Bevat diëthylhexylftalaat (DEHP)	Данное изделие содержит диэтилгексилфталат (DEHP)
	Keep dry	Vor Nässe schützen	Protéger de l'humidité	Beschermen tegen vocht	Беречь от влаги
	Top-Bottom	Oben-Unten	Haut-Bas	Boven-Beneden	Верх-низ
	Fragile	Zerbrechlich	Fragile	Breekbaar	Хрупкое
	Waste management	Entsorgung	Élimination	Verwijdering	Утилизация
	Manufacturer	Hersteller	Fabricant	Fabrikant	Изготовитель
	Do not use if package is damaged	Inhalt beschädigter Packung nicht verwen- den	Ne pas utiliser le con- tenu si l'emballage est endommagé	Inhoud van beschadigde verpakking niet gebrui- ken	Не использовать содержимое повреждённой упаковки
	Keep away from sun- light	Vor Sonnenlicht schützen	Protéger contre les ray- ons du soleil	Uit het zonlicht houden	Беречь от солнечного света
	Protect from heat and radioactive sources	Vor Hitze und radioaktiver Strahlung schützen	Protéger de la chaleur et des ray- onnements radioactifs	Beschermen tegen hitte en radioactieve straling	Оберегать от воздействия тепла и радиоактивного излучения
	Humidity limitation	Luftfeuchtigkeit, Begrenzung	Limitation de l'humidité de l'air	Vochtbe grenzing	Ограничение влажности
	Atmospheric pressure limitation	Luftdruck, Begrenzung	Limitation de la pression atmosphérique	Begrenzing van de atmosferische druk	Ограничение атмосферного давления
	Temperature limit	Temperaturbegrenzung	Limite de la température	Temperatuurbegrenzing	Ограничение температуры
	Non sterile	Unsteril	Non stérile	Niet steriel	Нестерильно

	Description	Beschreibung	Description	Beschrijving	Описание
	Max. use 30 days	max. 30 Tage verwenden	Utiliser en 30 jours maximum	max. 30 dagen gebruiken	Применять на протяжении не более 30 дней
	Non-ionizing electromagnetic radiation	Nicht ionisierende elektromagnetische Strahlung	Rayonnement électromagnétique non-ionisant	Niet-ioniserende elektromagnetische straling	Неионизирующее электромагнитное излучение
	Pushing prohibited	Schieben verboten	Interdiction de pousser	Duwen verboden	Не двигать
	On/Off (push button)	Ein/Aus (druckbetätigt)	Marche/arrêt (Interrup-teur actionné)	Aan/Uit (drukknop)	Вкл/выкл (активируется нажатием)

EN

DE

FR

NL

RU

Table of contents

1	Important User Notes	3
2	Safety Information	4
3	General Information	5
3.1	Device Description	5
3.2	Intended Use and Contraindications	5
3.2.1	Intended Use	5
3.2.2	Contraindications	5
3.3	General Warnings and Precautions	5
3.3.1	General Warnings	5
3.3.2	Precautions	9
3.4	Indication-Specific Warnings	12
3.4.1	General Hysteroscopy Warnings	12
3.4.2	Laparoscopy Warnings	15
4	Initial Startup	16
4.1	Scope of Delivery	16
4.2	Setting up and Connecting the Device	16
4.3	Set up and Connect the Pump with Fluid Monitoring Unit	17
5	Description of Device	20
5.1	Front of Device	20
5.2	Rear of the Device	20
5.3	Touchscreen	21
6	Operating the Device	23
6.1	Switching the Device On	23
6.1.1	Device Displays After Unsuccessful Completion Of Device Check	23
6.2	Using the Tube Sets	23
6.3	Usable Tube Sets Overview	24
6.4	Irrigation tube set connection	25
6.5	Removing a Tube Set	27
6.6	Attaching and Connecting the Irrigation Fluid Bags (Use without Fluid Monitoring Unit)	27
6.7	Suction Tube Set Connection	28
6.8	Changing a Fluid Bag	29
6.9	Switching Device Off	30
7	Using the Suction Function	31
7.1	Connecting a Suction Tube	31
7.2	Starting/Stopping Suction	32
7.3	Replacing the Secretion Container	32
7.4	Safety Functions	32
8	Using the Laparoscopy Indication	33
8.1	Device-Inherent Dangers: Laparoscopy Indication	33
8.2	Screen Displays of the Laparoscopy Indication	33
8.3	Inserting a Tube Set	34
8.4	Starting and Stopping Irrigation	34
8.5	Opening and Configuring the Laparoscopy Indication Menu	34
8.5.1	Fluid Bag Configuration	35
8.5.2	Resetting to Factory Settings	36
9	Using the Hysteroscopy Indication	37
9.1	Device-Inherent Dangers: Hysteroscopy Indication	37
9.2	Screen Displays of the Hysteroscopy Indication	39
9.3	Inserting a Tube Set	40
9.4	Surgical Mode	41
9.4.1	Max. flow Range Display (max. flow)	41
9.5	Automatic Instrument Recognition	41
9.6	Preselecting the Nominal Pressure	42
9.7	Safety Threshold Nominal Pressure Over 100 mmHg	42
9.8	Preselecting the max. flow range (max. flow)	42
9.9	Saving the Nominal Values	43
9.10	Resetting the Deficit	43
9.11	Starting and Stopping Irrigation	43
9.12	Changing Bag during Surgery	44
9.13	Changing Container during Surgery	44
9.14	Opening and Configuring the Hysteroscopy Indication Menu	45
9.14.1	Setting the Deficit Threshold	45
9.14.2	Fluid Bag Configuration	46
9.14.3	Scales Test	46
9.14.4	Setting the Instrument Recognition	46
9.14.5	Resetting to Factory Settings	47
9.15	Safety Functions	47
9.15.1	General Safety Functions	47
9.15.2	Exceeding the Deficit Rate	47

9.15.3	Reaching and Exceeding the Deficit Threshold	47
9.15.4	Exceeding the Nominal Pressure	48
10	Menu	49
10.1	Overview Menu - General Settings	49
10.2	General Settings	50
10.2.1	Selecting the Language	50
10.2.2	Setting the Volume of the Acoustic Signals and Audio Warnings	50
10.2.3	Setting the (Display) Brightness.....	50
10.2.4	Device Info.....	51
10.2.5	Demo.....	51
10.2.6	Factory settings	51
10.2.7	Service.....	51
11	Safety Functions.....	52
11.1	Device Self-Test.....	52
11.2	Sensor Error.....	52
11.3	Motor Error.....	52
11.4	Electronics Error	52
11.5	Calibration Error	52
12	Function Test	53
12.1	Preparing the Function Test of the Device	53
12.2	Performing Device Function Test (Only for Laparoscopy Indication).....	53
12.3	Performing Device Function Test (Only for Hysteroscopy Indication).....	54
12.4	Performing Function Test of Fluid Monitoring Unit (Only for Hysteroscopy Indication).....	54
13	Care and Maintenance.....	55
13.1	Cleaning the Device	55
13.2	Maintenance Intervals	55
13.3	Maintenance by Authorized Service Technician.....	55
13.4	Changing the Fuse	56
13.5	Cleaning the accessories	57
13.6	Care of the Reusable Tube Set	57
13.6.1	Cleaning the Reusable Tube Set	57
13.6.2	Disinfecting the Reusable Tube Set.....	59
13.6.3	Sterilizing the Reusable Tube Set.....	59
14	Annual Inspection	62
14.1	Electrical Safety Test.....	62
14.2	Basic Function Test.....	62
14.3	Pressure Measuring Test	63
14.4	Performing Scale Function Test (for Hysteroscopy Indication)	64
15	Electromagnetic Compatibility.....	65
15.1	Electrical Connections.....	65
15.2	Guidelines and Manufacturer's Statement – Electromagnetic Emissions	66
15.3	Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity	67
15.4	Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity.....	68
15.5	Recommended Safety Distances Between Portable and Mobile RF Telecommunications Devices and GYN-Pump	69
16	Error and Warning Messages.....	70
17	Technical Data Pump.....	72
18	Technical Data of Fluid Monitoring Unit	74
19	Accessory List	75
20	Glossary.....	76
21	Appendix.....	77
21.1	Test Log.....	77
21.2	Return Form	78
	Index	79

1 Important User Notes

Read the instructions for use carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing staff or service personnel, or
- to damage or malfunction of device and/or accessories.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of its products.

The words WARNING, CAUTION, and NOTE carry special meanings. Sections marked with these words must be given special attention.

Subject to technical changes

Please note

WARNING!
The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third party.

CAUTION!
These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.

NOTE!
These paragraphs contain information to clarify the instructions or provide additional useful information.

Federal Law (only for U.S. market)

Exclusion of liability

Authorized service technician

Care and maintenance

Contamination

Waste management

2 Safety Information

CAUTION: Federal Law restricts this device to sale by or on the order of a physician.

The manufacturer is not liable for direct or consequential damages, and the warranty becomes null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained;
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedule is not adhered to.

The handing over of technical documents does not constitute authorization to make repairs or alterations to the device or accessories.

WARNING: Modifying the GYN-Pump is not permitted.

Only an authorized service technician may perform repairs, adjustments, or alterations on the device or accessories and use the service menu. Any violation will void the manufacturer's warranty. Authorized service technicians are only trained and certified by the manufacturer.

The service and maintenance of the device and its accessories has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use.

Maintenance of the device may not be performed during the operation.

Before shipping, decontaminate device and accessories in order to protect the service personnel. Follow the instructions listed in this manual. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to reject contaminated products for repair.

This symbol indicates that the waste of electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately instead. Please contact the manufacturer or an accordingly authorized disposal or waste management company for further information.

3 General Information

3.1 Device Description

The GYN-Pump is a system providing a suction and irrigation function for two different indications:

- Laparoscopy
- Hysteroscopy

The device is non-invasive and designed for use in non-sterile areas. It pumps medically sterile irrigation fluids through a sterile tube. These fluids are used to distend and irrigate corresponding body cavities to provide space and improve visibility for the attending physician. The device can be used with electrolyte-free media (e.g., glycine 1.5% and sorbitol 3.0%) and with isotonic, electrolyte containing media (e.g., saline 0.9% and Lactated Ringer's). The device is equipped with a vacuum suction function used to suction off secreted fluids. The device functions only with the tube sets described in the accessory list (see chapter 19). The pump is equipped with a balancing system that monitors the volume differential between the irrigation fluid flowing in and out of the uterus during hysteroscopic procedures.

The device can be operated at max. flow rates of up to 2.0 l/min when used during laparoscopic procedures. During hysteroscopic procedures, a maximum flow rate of 800 ml/min can be selected. Pressure settings of up to 150 mmHg can be selected during hysteroscopic procedures.

The device operates with a completely non-contact pressure measurement of the irrigation medium. The contact-free pressure measurement is taken by integrating the pressure chamber into the tube system. The pressure membrane transfers the tube pressure to the electronics of the device via a pressure sensor. The pressure control circuit continuously compares the actual pressure with the nominal pressure to maintain the nominal pressure. If the nominal flow rate is set too low, the nominal pressure cannot be reached. Pay attention to possible leaks.

3.2 Intended Use and Contraindications

3.2.1 Intended Use

The GYN-Pump is designed for body cavity distension, fluid irrigation and fluid suctioning. The GYN-Pump monitors the volume differential between the irrigation fluid flowing into and out of the patient during hysteroscopic procedures. The GYN-Pump is designed for the following indications:

- Laparoscopy: Irrigation of body cavities and fluid aspiration during diagnostic and operative procedures.
- Hysteroscopy: Uterine distention, fluid aspiration and monitoring of fluid deficit during diagnostic and operative procedures.

3.2.2 Contraindications

- Use of the GYN-Pump is contraindicated for the distention of body cavities when endoscopic surgery is contraindicated as well.
- Use of the GYN-Pump is contraindicated if, in the opinion of the responsible physician, its use endangers the health of the patient.

3.3 General Warnings and Precautions

3.3.1 General Warnings

CAUTION!

The device may not be used in conjunction with a defibrillator since it is not equipped with corresponding safety elements. The manufacturer accepts no liability in this case for ensuing damage.

WARNING!**Replacement device and accessories**

In case the device or any of the accessories fail during surgery, a replacement device and replacement accessories should be kept within close proximity to be able to finish the operation with the replacement components.

WARNING!**Technique and procedures**

Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

WARNING!**Check all factory settings**

Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

WARNING!

Different default settings of the warning message for identical or similar devices in the operating room may cause a risk due to conflicting acoustic signals.

WARNING!**Original accessories**

For your own safety and that of your patient, use only original accessories.

WARNING!**Not explosion-proof**

The device is not explosion-proof. Do not operate the device in the vicinity of explosive anesthetic gases and not in the vicinity of oxygen-enriched environments.

WARNING!**Falls and crashes**

Place the device on a stable and level surface. Cables must be laid safely. Tubes between the device and the patient must not create any obstruction.

WARNING!**Risk of electrical shock**

To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

WARNING!**Risk of electrical shock**

To avoid the risk of electrical shock, only use this device when connected to a properly grounded power supply network.

WARNING!**Professional qualification**

The Instructions for Use do not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualifications working under the direction and supervision of a physician.

**WARNING!****Functional test**

The functional test must be performed prior to each surgery.

**WARNING!**

Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

**WARNING!****Device defect**

If a device defect is suspected or confirmed, do not use it. Make sure the device can no longer be used until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!****Sterile mediums and accessories**

Always work exclusively with sterile substances and mediums, sterile fluids, and sterile accessories if so indicated.

**WARNING!****Replacing fuse**

Replace the fuse only with a fuse of the same type and rating.

**WARNING!****Maintenance and calibration**

Do not open the device. The device may not be opened or calibrated by the user. Only authorized service technicians are permitted to repair, calibrate, or modify the device or its equipment.

WARNING!

Condensation / Water penetration

Protect device from moisture. Do not use if moisture has penetrated the device.

WARNING!

Cleaning the device

Do not sterilize the device.

WARNING!

Distention may be lost when resetting the nominal flow or nominal pressure settings.

WARNING!

The device is only intended for use with flexible fluid bags. Do not use glass containers as they might break. Fluid cannot flow quickly enough due to the vacuum being generated inside of the bottle. Risk of implosion.

WARNING!

Set up the device in such a way as to allow for easy monitoring of the display values, device functions, and access to the control elements.

WARNING!

Keep full fluid bags ready for use

Always keep a full fluid bag on hand to replace an empty one. This avoids having to interrupt surgery due to a lack of irrigation fluid.

WARNING!

Contamination

Do not use device and/or accessories if signs of contamination are detected. Make sure the device or/and accessories can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

WARNING!

Different default settings of the warning message for identical or similar devices in the operating room may cause a risk due to conflicting acoustic signals.

WARNING!

Sterilize reusable instruments and tubing before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

WARNING!**Reprocessing of sterile disposable products**

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

**WARNING!**

Sterilize reusable instruments and tubing before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

**WARNING!**

Do not touch the patient and the On/Off button at the same time.

**WARNING!****Disconnecting device from mains power supply**

However, the ON/OFF switch does not disconnect the device from the mains power supply. This requires pulling the plug located in the rear of the device.

**WARNING!**

Use only parts and/or devices from ME systems (see chapter 20 Glossary) in patient environments in compliance with the standard IEC60601-1 in the respective currently valid version.

**WARNING!**

Connect only items, that have been specified as part of the ME system or specified as being compatible with ME system.

**WARNING!****Additional equipment**

Additional equipment connected to medical electrical devices must be demonstrably compliant with their respective IEC or ISO standards (IEC 60601-1, IEC 60950 or IEC 62368 for data processing equipment). Furthermore, all configurations must comply with the normative requirements for medical systems (see section 16 of the last valid edition of IEC 60601-1). Anyone who connects additional devices to medical electrical equipment is a system configurator and as such is responsible for the system's compliance with the normative requirements for systems. Please contact the technical service if you have additional questions.

3.3.2 Precautions

CAUTION!

Patient group

There are no restrictions as to the specification of the patient group when using the device as intended and its use does not endanger the patient's health.

CAUTION!

Continuous operation

After 24 hours of continuous operation, switch the device off and then on again so that the device self test can be carried out.

CAUTION!

The device may not be used in conjunction with a defibrillator since it is not equipped with corresponding safety elements. The manufacturer accepts no liability in this case for ensuing damage.

CAUTION!

When a morcellator/shaver system is used, the combination of low nominal pressures and a too high vacuum pressures may result in a significant loss of intra-cavity distension pressure, which has the potential to affect the visibility of the surgical field. Conversely, when employing high distension pressures, the de-activation of the morcellator system can lead to pressure spikes exceeding 150 mmHg.

CAUTION!

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

CAUTION!

Any power cables employed by the user that are not provided by the manufacturer must meet the safety requirements of the national standards in the respective current valid version.

CAUTION!

Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!

Ventilation of the device

Avoid device overheating. Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).

CAUTION!

Endoscope

The device may only be connected with endoscopes designed for and featuring the technical specification permitting such a combined use. Any utilized endo-

scopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600. Combining/connecting with other devices generates a medical electrical system (MES). The system configurator is responsible for compliance with the standard IEC 60601-1 / EN 60601-1 in its latest version.

CAUTION!
Electrical Interference

(See chapter 15). Care was taken during the development and testing of this device that electrical interference of or from other devices or instruments was practically eliminated. However, if you still detect or suspect such interference, please follow these suggestions:

- Move this, the other or both devices to a different location
 - Increase distance between used devices
 - Consult an electro-medical expert
-

CAUTION!
Use of other accessories, other transducers and cables

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM (see chapter 20) as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

CAUTION!

The medical electrical (ME) device is suitable for integration in ME equipment systems (see chapter 15). Operation of the ME device in vicinity of non-ME devices may result in voiding the intended use of the ME device.

CAUTION!
Peripheral devices

Additional peripheral equipment connected to interfaces of the medical monitor has to meet the requirements of the following specifications in the respective current valid version: IEC 60601-2-18 / EN 60601-2-18 for endoscopic devices and IEC 60601-1 / EN 60601-1 for electrical medical devices. All configurations have to comply with IEC 60601-1 / EN 60601-1 specifications. Whoever connects additional equipment to signal output or signal input is considered the system configurator and as such is responsible for complying with requirements of the standard IEC 60601-1 / EN 60601-1.

CAUTION!

The device is transportable. The rollers of the cart/scale system are used for positioning at the usage site. To transport the device, remove all fluid bags from the hooks and make sure there are no canisters or only completely emptied canisters on the cart/scale. Inflow and outflow tubes must be completely removed. Make sure power supply line does not touch the ground and there are no other objects located on the GYN-Pump oder Bilanzierungseinheit. Always use the handle to move the system safely.

3.4 Indication-Specific Warnings

3.4.1 General Hysteroscopy Warnings

WARNING!

When performing monopolar electrosurgery, the distension medium must be electrically non-conductive. Examples include glycine, sorbitol and mannitol. Isotonic saline irrigation fluids may only be used when performing bipolar electrosurgical resective procedures.

WARNING!

Fluid Overload

There is a risk of irrigation fluid reaching the circulatory system of the patient's soft tissue. This can be affected by distention pressure, flow rate, perforation of the distended body cavity and duration of the endoscopic surgery. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!

Hypothermia (monitoring body temperature)

The continuous flow of the distension liquid may lower the body temperature of the patient. Lower body temperatures can cause coronary and cardiovascular problems. Always monitor the patient's body temperature during the entire surgical procedure. Make especially sure that the following, hypothermia promoting, operation conditions are avoided as best as possible:

- longer operating times
- use of cold irrigation fluid.

WARNING!

Air Embolisms

An air embolism can be the result of air contained in the tube set or connected instrument reaching the patient. Make sure that there is always fluid in the bag to prevent air from being aspirated into the patient.

WARNING!

When using the scale, follow the operating instructions in this manual.

WARNING!

The deficit and inflow values are lost in case of a power loss or "brownout".

WARNING!

Intrauterine distention is usually possible with pressure values between 35 to 70 mmHg. A pressure above 75 to 80 mmHg is required only in rare cases or if the patient has an excessively high blood pressure.

WARNING!

Fluid intake and output surveillance

Strict fluid intake and output surveillance should be maintained. If a low viscosity fluid distention medium is used, intrauterine instillation exceeding 2 liters

should be followed with great care due to the possibility of fluid overload.

WARNING!

Hyponatremia

Some distension fluids may lead to fluid overload and, consequently, hyponatremia with its attending sequelae. This can be affected by the distending pressure, flow rate, and duration of hysteroscopic procedure. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!

Pulmonary edema

A surgical procedure has the risk of pulmonary edema, which is formed due a "fluid overload" with isotonic fluid. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!

Cerebral edema

A surgical procedure has the risk of cerebral edema resulting from fluid overload and electrolyte disturbances with hypoosmolar (non-ionic) fluids such as glycine 1.5% and sorbitol 3.0%. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!

Idiosyncratic reactions

In rare cases, idiosyncratic reactions such as

- intravascular coagulopathy
- allergic reaction including anaphylaxis

may occur during a surgical procedure if a liquid distention medium is used. Specifically, idiosyncratic anaphylactic reactions have been reported when using Hyskon as an irrigation fluid during a surgical procedure. These should be managed like any other allergic reaction.

WARNING!

Rupture of the fallopian tube secondary to tubal obstruction

Distention of the uterus may lead to a tear of the fallopian tube should there be an obstruction or permanent occlusion. The rupture could lead to irrigation fluid flowing into the patient's peritoneal cavity, resulting in a fluid overload. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!

Instrument replacement

Stop the device using the START/STOP switch if replacing the instrument during surgery.

WARNING!
Deficit displays and warnings serve as a tool for the treating physician and do not replace the monitoring of the patient's condition.

WARNING!
Fluid volume/sodium concentration
The fluid left in the patient and the concentration of sodium in the blood serum must both be monitored. The deficit amount is the entire amount of fluid lost by or to the system. Take note of the measurement tolerance of the system. Estimating the fluid volume remaining in the patient is the physician's responsibility.

WARNING!
The pressure should be kept as low as possible to allow for a sufficient intrauterine distension and to reduce the forces that could allow fluid, ambient air, and/or gas into the circulatory system.

WARNING!
Filling the tubing with irrigation fluid and resetting the deficit display to zero are to be done at the physician's discretion.

3.4.2 Laparoscopy Warnings

CAUTION!
Connecting the tube
The tube outlet may only be connected to instruments which are intended for intra-abdominal fluid suction and irrigation.

4 Initial Startup

The operation of the device GYN-Pump is reserved for persons with the relevant professional qualifications and trained to use the device.

4.1 Scope of Delivery

- GYN-Pump PH304
- Fluid monitoring unit PS304 (scale)
- User Manual
- Power cable

Always check all parts and optional accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized service company.

If it becomes necessary to return the device, always use the original packaging. The manufacturer does not take responsibility for damage that has occurred during transportation if the damage was caused by inadequate transport packaging.

Please make sure that all required information has been supplied:

- Name of owner
- Address of owner
- Device type and model
- Serial number of the equipment (see identification plate)
- Detailed description of defect

Use the return shipping form included in chapter 21.2 Return Form.

4.2 Setting up and Connecting the Device

Place the device on a flat surface free of vibration located in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in chapter 17.

Delivery inspection

Returning the device

Setting up

WARNING!

Use only parts and/or devices from ME systems (see chapter 20 Glossary) in patient environments in compliance with the standard IEC60601-1 in the respective currently valid version.

CAUTION!

The medical electrical (ME) device is suitable for integration in ME equipment systems (see chapter 15). Operation of the ME device in vicinity of non-ME devices may result in voiding the intended use of the ME device.

CAUTION!

Equipment should be positioned such that power cord can be easily disconnected.

CAUTION!

Ventilation of the device

Avoid device overheating. Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).

CAUTION!
Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!
Position the device in such a way that it is easy to operate and switch off.

CAUTION!
Place the device outside the sterile field.

CAUTION!
Position of the user
To avoid a malfunction, the user must be positioned correctly

- within a display viewing angle of $\pm 50^\circ$ to operate the device
- up to 2 m/6.5 ft from the device front for monitoring the actual values

Mains connection

CAUTION!
Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

Make sure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The mains connection cable may be plugged only into a properly installed, grounded safety wall socket (shock-proof socket) (see DIN VDE 0100-710). Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

Grounding contact

The power connection must be equipped with a grounding contact. Use the original power cable (if included in scope of delivery) to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.

Only for U.S. operators

Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 and IEC 16320-C13. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

Potential equalization

The equipotential bonding is used as a protective measure against the failure of the protective conductor according to requirements of IEC 60601-1 in the respectively valid version. The installation must be according to the relevant local safety regulations.

4.3 Set up and Connect the Pump with Fluid Monitoring Unit

CAUTION!
The scale to be connected must be attached to the device before the device is switched on or the device does not detect the scale.

CAUTION!

The scale can be disinfected by wiping off. Do not sterilize the scale.

EN

The fluid monitoring unit is a cart/scale system consisting of a weighing unit for the irrigation fluid bags and containers as well as a roller foot.


Fig. 4-1 Pump with fluid monitoring unit


- ① Holder with bag hooks
- ② Irrigation fluid bags
- ③ Bag deflector
- ④ Pump
- ⑤ Pump tray
- ⑥ Bag scale/container scale
- ⑦ Container (Bemis, Medela, Serres, Abbott)
- ⑧ Container holder
- ⑨ Roller wheel base
- ⑩ Locking foot brake

1. Place pump on the pump tray of ⑤ the cart/scale system.
2. Run the power cable through the notch of the pump tray provided for this purpose ⑥. Connect the power cable ⑤ with the pump and a grounded safety wall socket (Fig. 4-2).
3. Run the fluid monitoring unit cable ⑦ through the notch of the pump tray provided for this purpose ⑥ (Fig. 4-2). Connect the fluid monitoring unit cable ⑦ with the connection in the rear of the pump (Fig. 4-2). Make a loop with the pump cable to prevent bending and kinking and position it above the pump support.

EN

Fig. 4-2 Connecting pump and fluid monitoring unit

- ① Pump
- ② Pump tray
- ③ Holder with bag hooks
- ④ Bag deflector
- ⑤ Power cable/pump connection
- ⑥ Pump tray notch
- ⑦ fluid monitoring unit cable/pump connection

CAUTION!
Do not lean against or on the fluid monitoring unit. The bag scale can be loaded with a maximum of 23 kg, the container scale a maximum of 16 kg. Make sure only bags with irrigation fluid are attached to the bag scale. Only place containers with connected tubes into the holder for containers.

5 Description of Device

5.1 Front of Device

Familiarize yourself with the control and display elements of the device.


Fig. 5-1 Front of device

- ① Touchscreen
- ② Tube retainer
- ③ Roller wheel
- ④ Vacuum pump connection
- ⑤ ON/OFF switch

5.2 Rear of the Device


Fig. 5-2 Rear of the device

- ⑥ Mains power connection with fuse holder
- ⑦ Potential equalization plug
- ⑧ USB service interface
- ⑨ Scale connection
- ⑩ Vacuum output

WARNING!

Additional equipment

Additional equipment connected to medical electrical devices must be demonstrably compliant with their respective IEC or ISO standards (IEC 60601-1, IEC 60950 or IEC 62368 for data processing equipment). Furthermore, all configurations must comply with the normative requirements for medical systems (see section 16 of the last valid edition of IEC 60601-1). Anyone who connects additional devices to medical electrical equipment is a system configurator and as such is responsible for the system's compliance with the normative requirements for systems. Please contact the technical service if you have additional questions.

5.3 Touchscreen

The touchscreen connects display and control elements and offers intuitive ease of use. Function fields or symbols are used to control the unit. An acoustic signal is emitted when pressing a function field or symbol. This signal confirms your input. The following different screen displays are depicted when operating the device:

- Indication selection
- Indication screen
- Menu
- Form of treatment (only hysteroscopy)

The following is an overview of the screen displays of the indications.

Fig. 5-3 Screen displays overview


Fig. 5-4 depicts the 6 areas with display and control elements.

Indication screen

Fig. 5-4 Indication screen overview


- ① Mode line
- ② Actual values
- ③ Nominal values
- ④ Status line
- ⑤ Footer
- ⑥ Suction level
- ⑦ Shift key

The indication screen is divided into the following 6 areas:

1. **Modus line:** Displays the indication type.
2. **Actual values:** Indicates the measured actual pressure on the left side. In Hysteroscopy mode, the right side depicts the deficit volume if a scale is connected.
3. **Nominal values:** Set the nominal pressure on the left side.
4. **Status line:** Depicts messages, warnings, and error messages.
5. **Footer:** This features only control elements: Function fields and symbols.
6. **Suction level:** Used to select one of two available suction pump levels. The respective suction level is depicted with little green boxes.
7. **Shift key:** The display can be switched between the inflow volume and the maximum flow range (max. flow) using the Shift key.

The indication screens are described in details in the following chapters:

- 8.2 Screen Displays of the Laparoscopy Indication
- 9.2 Screen Displays of the Hysteroscopy Indication

6 Operating the Device

6.1 Switching the Device On

WARNING!

Do not touch the patient and the On/Off button at the same time.

1. Plug the device into the power outlet.
2. Make sure a tube set is **not** inserted into the tube retainer. Press the ON/OFF key.
3. The device performs a device check after being switched on. Initially, the **company logo and a progress bar** are depicted on the touchscreen. After the successful device check, the screen depicts the indications in the form of function fields. The following indication types are available:
 - **LAPAROSCOPY**
 - **HYSTEROSCOPY**
4. If a tube set is in the tube retainer before the device is being switched on, the screen depicts the following message: **Remove tube set**.
If the device check is unsuccessful, the corresponding error messages are displayed on screen.

6.1.1 Device Displays After Unsuccessful Completion Of Device Check

If the device check was **not** successful, the touchscreen depicts the respective device error. The following device errors may be depicted:

- **Motor error**
- **Sensor error**
- **Electronic error**
- **Calibration error**

Please read chapter 16 Error and Warning Messages for additional information.

6.2 Using the Tube Sets

The transponder technology detects the tube type, the validity and reliability of a tube set automatically and a corresponding message is output in the status line of the touchscreen. This eliminates virtually all "operating errors" since non fitting, invalid, and not allowed tube sets are reliably detected. The corresponding transponder is located underneath the tube retainer in each allowed tube set (see chapter 19 Accessory List).

If an approved tube set is located in the device, the transponder technology automatically invalidates this tube set when the device is started. Depending on the tube set, this has the following consequences:

- **Disposable tube set:** After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If the irrigation cycle is stopped, the device can be restarted within 30 minutes. If the device is switched off or in case of a power failure, the tube set is invalidated. The irrigation cycle can no longer be started. If this is the case, you must insert a new, valid, and approved tube set.
- **Reusable tube set (NOT FOR SALE IN USA):** The reusable tube set can be used about 20 times. The status line on the touchscreen depicts the **Remaining cycles: X**. After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If irrigation is stopped during a use cycle, it can be restarted within 30 minutes without the set being invalidated. If the device is switched off or in case of a power failure, the current use cycle is invalidated.

Before the last cycle the status line on the touch screen depicts **Last cycle!**. Irrigation can no longer be started following the last use cycle. If this is the case, you must insert a new, valid, and approved tube set.

RFID technology (transponder technology)

Invalidating a tube set

The device stops and cannot be restarted if the transponder loses its signal during current use. If the signal is restored with 20 seconds, it is possible to continue using the tube set.

Transponder signal loss

6.3 Usable Tube Sets Overview

The following table lists the function of each type of tube set depending on the selected indication. Explanation: [x] = usable, [-] = not usable

Art. No.	Tube category	Tube type	LAPAROSCOPY	HYSTEROSCOPY
T0505-01	Irrigation tube set	Disposable (one-time use)	x	x
T0506-01	Irrigation tube set	Reusable (NOT FOR SALE IN USA)	x	x
T0504-01	Vacuum tube set, incl. filter	Can be used 30 days	x	x
T0502-01	Suction tube set (1 connection)	Disposable (one-time use)	x	x
T0503-01	Suction tube set (2 connections)	Disposable (one-time use)		x

Table 6-1

CAUTION!

The pressure chamber diaphragm is fragile and can be replaced if damaged. The reusable tube sets must be prepared before each use. Instructions are listed in chapter 13.6 Care of the Reusable Tube Set.

WARNING!

The vacuum tube sets for this device contain diethylhexylphthalate (DEHP), which is classified as toxic to reproduction according to the EU Directive 1272/2008/EEC on Classification, Labeling and Packaging of Dangerous Substances. DEHP may impair fertility and may cause harm to the unborn child. Therefore, this product must not be used for unauthorized applications. When applied within the intended use, the potential risk to pregnant or breastfeeding women as well as to children resulting from the DEHP contained in this product is not critical. In regard to the short exposure time and the physical characteristics, the eventuality of critical quantities of DEHP being dissolved from the tube sets is negligible.

WARNING!

Sterilize reusable instruments and tubing before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

WARNING!

Reprocessing of sterile disposable products
Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.


NOTE!
Tubes must be stored at room temperature. The shelf life for all tubes is 5 years.

NOTE!
Observe applicable hygiene regulations when disposing of the tube set.

Fig. 6-1 Irrigation tube set elements

- ① Protective caps
- ② Tap spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Ring
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑩ Luer lock connector

The irrigation tube set is available as a disposable (one-time use) or reusable (autoclavable) tube set (chapter 19 Accessory List).

The irrigation tube set consists of 3 tube components, (irrigation tube ⑤), roller tube ⑦ and instrument tube ⑨), a Y-connector ④ and 2 tap spikes ②). The tap spikes ② are used to connect the tube components with the fluid bags.

The Luer lock connector ⑩ connects the instrument tube with the instrument.


Fig. 6-2 Inserting the irrigation tube set

- ⑤ Irrigation tube
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑪ Instrument
- ⑫ Fluid bag
- ⑬ Tube retainer

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to Fig. 6-1 and Fig. 6-2.

1. To be carried out by non-sterile technician:

- ▶ If you are using a reusable tube set, open the sterile autoclavable container of the tube set and have a sterile technician remove the inner tube set.
- ▶ If you are using a disposable tube set, open the outer packaging of the tube set and have a sterile technician remove and open the inner tube set packaging.

2. To be carried out by sterile technician:

- ▶ Keep the Luer lock connector (10) in the sterile area and hand the tube end with the tap spikes (2) to the non-sterile technician.
- ▶ Connect the Luer lock connector (10) with the instrument (e.g. inflow cannula). Open inflow valve at instrument.


3. To be carried out by non-sterile technician:

- ▶ Turn the device on and select the desired indication. The indication screen shows the following: **Insert tubeset**.
- ▶ 1. Insert the unpressurized pressure chamber (8) carefully into the lower notch of the tube retainer (12) up to the stop.
- ▶ When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (8) only if chamber is not pressurized.
- ▶ 2. Place the roller tube (7) around the roller wheel (13). Pull the tube until the tube ring fits in the upper part of the tube retainer.
- ▶ The validity of the tube set is now checked. If the tube set is invalid, the status line depicts **Tubeset not valid** followed by **Remove tubeset**. A short acoustic signal is emitted 3x. The entire procedure must be repeated with a valid tube set if an invalid tube set has been inserted.
- ▶ Connect the tube end with the tap spikes (2) to the full fluid bags (11).
- ▶ If the **Fluid bag warning** function is enabled (see 8.5.1 Fluid Bag Configuration), press and hold the **[Fluid bag]** symbol for 1.5 s to reset the device to a full fluid bag. Make sure the bag volume is set correctly. The **Fluid bag warning** is disabled in the factory setting.
- ▶ The device is now ready for use.

Removing the tube set

Connect instrument

Insert tube set

Remove tube set

6.5 Removing a Tube Set

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to Fig. 6-1 and Fig. 6-2.

1. To be carried out by sterile technician:

- ▶ Disconnect the Luer lock connector (10) from the instrument (e.g. inflow cannula).
- ▶ Hand the Luer lock connector (10) to the non-sterile technician.

2. To be carried out by non-sterile technician:

- ▶ Disconnect the end with the tap spikes (2) from the fluid bags (11).
- ▶ Remove the tube from the tube retainer (12).

6.6 Attaching and Connecting the Irrigation Fluid Bags (Use without Fluid Monitoring Unit)


WARNING!

The physician must determine an irrigation fluid suitable for the application and medical procedure!

Fig. 6-3 Attaching the fluid bags

- ① Fluid bag
- ② Tap spikes
- ③ Tube clamps
- ④ Irrigation tube

To be carried out by non-sterile technician:

- ▶ The irrigation tube (4) can receive irrigation fluid from 2 fluid containers (1). Connect both tube clamps (3) to the branches of the irrigation tube (Fig. 6-3).
- ▶ Use the tap spikes (2) to connect the tube ends of the tube set with the fluid bags. Always grasp the tap spike (2) with the provided handle when connecting.
- ▶ Comply with sterile conditions when inserting the tap spike (2) into the fluid bags (1).
- ▶ Attach fluid bags (1) to the hooks of the fluid monitoring unit.
- ▶ Open **one** of the two tube clamps (3) of the irrigation tube.

6.7 Suction Tube Set Connection

EN


Fig. 6-4 Suction tube set elements

- ① Protective cap
- ② Tap spike
- ③ Tube clamp
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Roller tube
- ⑦ Pressure chamber with membrane and transponder
- ⑧ Instrument tube
- ⑨ Suction and irrigation instrument
- ⑩ Suction tube
- ⑪ Connection to the suction container


The suction tube set is available as a disposable tube set for one-time use (chapter 19 Accessory List).

The suction tube consists of irrigation tube (5), roller tube (6), instrument tube (8), suction and irrigation instrument (9), and suction tube (10), as well as a Y-connector (4) and 2 tap spikes (2). The tap spikes are used to connect the tube sections with the fluid bags.

EN

Fig. 6-5 Inserting the suction tube set

- ① Suction and irrigation instrument
- ② Apron
- ③ Suction instrument tube
- ④ Container
- ⑤ Vacuum tube with filter

WARNING!
Do not touch the patient and the On/Off button at the same time.

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to the Fig. 6-4 and Fig. 6-5.

Fluid bag warning

6.8 Changing a Fluid Bag

If the **Warning** function in the **Fluid bag** menu is set to **ON** (see 8.5.1 Fluid Bag Configuration), an acoustic signal (beep) will be triggered as soon as the fluid bag is almost empty (volume <15 %) and again when the fluid bag is completely empty. A short warning signal is emitted 3x and the status line depicts the following: **Change fluid bag**. If the **Warning** function is not enabled, the fill level of the irrigation fluid bags must be monitored by the medical personnel.

The **Fluid bag warning** is disabled in the factory setting.

Enabling the fluid bag warning is only practical if:

- the height between roller wheel and fluid bag remains constant,
- the same bag size is used each time, and
- only **one** of the two tube clamps of the irrigation tube is opened at any time.

If the fluid bag warning is indicated visually or acoustically, the empty bag must be replaced with a full one.

Change fluid bag

CAUTION!

Make sure the bag volume of the used fluid bags is set correctly in the menu (see 8.5.1 Fluid Bag Configuration).

CAUTION!

Make sure the bag height of the used fluid bags is set correctly in the menu (see 8.5.1 Fluid Bag Configuration).


To be carried out by non-sterile technician:

- ▶ Open tube clamp of the full fluid bag and
- ▶ close tube clamp of the empty fluid bag.
- ▶ Press **[Fluid bag]** on the touchscreen for 1.5 s to reset the device to a full fluid bag. The status line will display the message: **Fluid bag replaced.**
- ▶ Replace the almost empty fluid bag with a full one.
- ▶ Keep the tube clamp of the new, full fluid bag closed until this bag is also in need of replacement.

This procedure must be carried out each time a fluid bag is replaced.

6.9 Switching Device Off

- Press the OFF switch. The device shuts off.

7 Using the Suction Function

The device is equipped with a vacuum pump. The vacuum pump is used to suction off secretions of a patient by using a suction tube and a container.

The vacuum pump generates 2 different vacuum levels:


Suction level	Vacuum
LOW	30 kPa
HIGH	60 kPa

Table 7-1

The suction level is displayed in the lower right in the indication view. The following depicts the three switching states of the suction.

Fig. 7-1 Suction level display

- ① Vacuum off
- ② Vacuum LOW
- ③ Vacuum HIGH

7.1 Connecting a Suction Tube

The suction system consists of the following:

- Vacuum tube with filter
- Container
- Suction instrument tube
- Suction instrument

(Fig. 6-4 Suction tube set elements)

WARNING!

Filter

The vacuum tube with integrated filter is designed for max. 30 days. Replace the vacuum tube if it is obviously contaminated. The filter prevents that body fluids enter the interior of the device. Please note that the filter may reduce the suction capacity.

CAUTION!

Only if the tube system is pre-evacuated (building a vacuum), the full suction capacity is available. Pre-evacuation takes about 30 to 60 seconds depending on the volume of the container.

**CAUTION!**

Only use containers with overflow protection.

1. Connect the device with the vacuum tube with filter ①.
2. Connect the vacuum tube ① with the container ②.
3. Connect the container ② with the suction instrument tube ③.
4. Connect the suction instrument tube ③ with the suction instrument ④.
5. Close the inflow valve of the suction instrument.

7.2 Starting/Stopping Suction

On the touchscreen, press the **[VAC]** button to activate suction level LOW or HIGH. The vacuum pump generates a negative pressure. Once a negative pressure matching the respective suction level has been reached, the vacuum pump stops. As soon as the negative pressure falls below the respective suction level setting, the vacuum pump restarts.

Start suction

7.3 Replacing the Secretion Container**CAUTION!**

Full secretion containers must be replaced immediately without stopping surgery. If the overflow protection of the secretion containers is triggered, suction is stopped to prevent the ingress of fluids.

The full secretion container can be replaced with an empty one during surgery.

Container change during surgery

1. Stop the suction (see chapter 7.2 Starting/Stopping Suction).
2. Replace the full secretion container with an empty one.
3. Start the suction (see chapter 7.2 Starting/Stopping Suction).

7.4 Safety Functions

If the suction function or vacuum pump is malfunctioning, the status line depicts the following warning message: **Vacuum pump defective**. An acoustic signal (beep) is emitted 3x. Surgery can be continued. After surgery, safeguard the device until checked by an authorized service technician.

Defective vacuum pump

Intended use: Laparoscopy

Contraindications

Technical application scope: Laparoscopy

8 Using the Laparoscopy Indication

Irrigation of body cavities and fluid aspiration during diagnostic and operative procedures.

The device may not be used in laparoscopy indication for surgical procedures where a specific pressure must be preset and reached as required during hysteroscopy procedures, for example.

The following characteristics apply to the laparoscopy indication:

- The device GYN-Pump works with a permanently preset pressure value when used in the laparoscopy indication.

8.1 Device-Inherent Dangers: Laparoscopy Indication

CAUTION!


Connecting the tube

The tube outlet may only be connected to instruments which are intended for intra-abdominal fluid suction and irrigation.

8.2 Screen Displays of the Laparoscopy Indication

Fig. 8-1 Laparoscopy indication screen

- ① Indication
- ② [Low] function field
- ③ [Medium] function field
- ④ [High] function field
- ⑤ [Menu] symbol
- ⑥ Suction function display
- ⑦ [Fluid bag] symbol
- ⑧ [START]/[STOP] function field
- ⑨ [Save] symbol
- ⑩ Status line for messages
- ⑪ [Back] symbol

The [High] flow level (factory setting) is set when activating the laparoscopy indication. The user can switch between [Low] flow rate and the higher flow rates [Medium] and [High] at any time.

Low	Medium	High
1.0 l/min	1.5 l/min	2.0 l/min

8.3 Inserting a Tube Set

The following tube sets can be used with the laparoscopy indication (ordering information can be found in chapter 19 Accessory List):

Laparoscopy tube types

Article no.	Tube category	Tube type
T0505-01	Irrigation tube set	Disposable (one-time use)
T0506-01	Irrigation tube set	Reusable (NOT FOR SALE IN USA)
T0504-01	Vacuum tube set, incl. filter	Can be used 30 days
T0502-01	Suction tube set (1 connection)	Disposable (one-time use)

Table 8-1

CAUTION!

When inserting a tube set already filled with irrigation fluid into the tube retainer, it is possible to damage the membranes. Only insert unfilled tube sets to avoid such damage.

For additional information, see 6.2 Using the Tube Sets and 6.4 Irrigation tube set connection.

8.4 Starting and Stopping Irrigation

During the laparoscopy indication, the pump works with a permanently preset pressure. The nominal flow depends on the nominal pressure and may reach a max. value of **2.0 l/min.**


1. Insert the tube set (see chapter 8.3 Inserting a Tube Set).
2. Open tube clamps at irrigation tube.
3. Open inflow valve at instrument.
4. Press the **[START]** function field to start the irrigation process. The roller wheel starts to turn.
5. Wait until the tube set is completely filled with fluid.
6. Press the **[STOP]** function field. The roller wheel stops. The irrigation process is canceled or ended.
7. Relieve any remaining pressure in the irrigation tube by opening the instrument.

Starting the irrigation process

8.5 Opening and Configuring the Laparoscopy Indication Menu

Device parameters are displayed and changed in the [LAPAROSCOPY] menu. The following chart provides an overview over the possible settings.

Menu overview

1. Press the **[Menu]** symbol while in the LAPARO indication to open the settings window. Press **[LAPAROSCOPY]**.
2. Press one of the two function fields: **[Fluid bag]** or **[Factory settings]**.

The following is a detailed description of how to change the indication-specific device parameters.

8.5.1 Fluid Bag Configuration

Configure warnings

1. In the menu, press the **[Fluid bag]** function field.
2. Press the **[ON]** function field to enable or the **[OFF]** function field to disable the warning function. The active function field lights up.
3. Press the **[◀]** symbol to return to the previous menu level. The setting is automatically saved.

Factory setting: Off

Set height

1. In the menu, press the **[Fluid bag]** function field.
2. Press the **▲** or **▼** symbol to set the height between device and fluid bag to a value between 0 and 1.5 m.
Pressing the **▲** or **▼** symbol longer than 1.5 seconds activates scrolling of the values in increments of 0.5 m.
3. Press the **[◀]** symbol to return to the previous menu level. The setting is automatically saved.

Factory setting: 1.00 m

1. In the menu, press the **[Fluid bag]** function field.
2. Press the ▲ or ▼ symbol to set the volume of the fluid bags to a value between 0.5 and 5 l.
3. Press the  symbol to return to the previous menu level. The setting is automatically saved.
Factory setting: 3 l

Setting the volume

(see Fig. 6-3 Attaching the fluid bags)

8.5.2 Resetting to Factory Settings

1. Press the **[Factory settings]** function field to reset the LAPAROSCOPY device parameters to factory settings.
You can now do the following:
2. Press the **[OK]** symbol to save the setting, or
3. Press the symbol  to return to the previous menu level without saving.

Reset to factory settings

Intended use: Hysteroscopy

Contraindications

Technical application scope: Hysteroscopy

9 Using the Hysteroscopy Indication

Uterine distention, fluid aspiration and monitoring of fluid deficit during diagnostic and operative procedures.

The device may not be used to introduce fluids into the cavum uteri when a hysteroscopy is contraindicated.

The following characteristics apply to the hysteroscopy indication:

- The nominal pressure can be set within a range of 15 to 150 mmHg.
- The maximum flow range (max. flow) can be set in the range of 50 to 800 ml/min.
- The device features an optional instrument recognition function.
- The pump reacts to pressures greater than the set nominal pressure by emitting warnings. The device features safety thresholds.
- The pump reacts by emitting warnings when the deficit limit is reached or exceeded.
- The pump reacts to perforation or high fluid loss (exceeding the deficit threshold) by emitting warnings.

9.1 Device-Inherent Dangers: Hysteroscopy Indication

WARNING!
When using the fluid monitoring unit, follow the operating instructions in this user manual.

WARNING!
The deficit and inflow values are lost in case of a power loss or "brownout".

WARNING!
Intrauterine distention is usually possible with pressure values between 35 to 70 mmHg. A pressure above 75 to 80 mmHg is required only in rare cases or if the patient has an excessively high blood pressure.

WARNING!
Fluid intake and output surveillance
Strict fluid intake and output surveillance should be maintained. If a low viscosity fluid distention medium is used, intrauterine instillation exceeding 2 liters should be followed with great care due to the possibility of fluid overload.

WARNING!
Hyponatremia
Some distension fluids may lead to fluid overload and, consequently, hyponatremia with its attending sequelae. This can be affected by the distending pressure, flow rate, and duration of hysteroscopic procedure. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!
Pulmonary edema
A surgical procedure has the risk of pulmonary edema, which is formed due a "fluid overload" with isotonic fluid. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!**Cerebral edema**

A surgical procedure has the risk of cerebral edema resulting from fluid overload and electrolyte disturbances with hypoosmolar (non-ionic) fluids such as glycine 1.5% and sorbitol 3.0%. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Idiosyncratic reactions**

In rare cases, idiosyncratic reactions such as

- intravascular coagulopathy
- allergic reaction including anaphylaxis

may occur during a surgical procedure if a liquid distention medium is used. Specifically, idiosyncratic anaphylactic reactions have been reported when using Hyskon as an irrigation fluid during a surgical procedure. These should be managed like any other allergic reaction.

**WARNING!****Rupture of the fallopian tube secondary to tubal obstruction**

Distention of the uterus may lead to a tear of the fallopian tube should there be an obstruction or permanent occlusion. The rupture could lead to irrigation fluid flowing into the patient's peritoneal cavity, resulting in a fluid overload. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Instrument replacement**

Stop the device using the START/STOP switch if replacing the instrument during surgery.

**WARNING!**

Deficit displays and warnings serve as a tool for the treating physician and do not replace the monitoring of the patient's condition.

**WARNING!****Fluid volume/sodium concentration**

The fluid left in the patient and the concentration of sodium in the blood serum must both be monitored. The deficit amount is the entire amount of fluid lost by or to the system. Take note of the measurement tolerance of the system. Estimating the fluid volume remaining in the patient is the physician's responsibility.

**WARNING!**

The pressure should be kept as low as possible to allow for a sufficient intrauterine distension and to reduce the forces that could allow fluid, ambient air, and/or gas into the circulatory system.

WARNING!

Fluid Overload

There is a risk of irrigation fluid reaching the circulatory system of the patient's soft tissue. This can be affected by distention pressure, flow rate, perforation of the distended body cavity and duration of the endoscopic surgery. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!

Filling the tubing with irrigation fluid and resetting the deficit display to zero are to be done at the physician's discretion.

9.2 Screen Displays of the Hysteroscopy Indication


The hysteroscopy indication should always be used with a connected scale to measure the differential volume. In the event that a scale is connected, this screen appears:

In case a scale is not connected, the screen display changes directly to the indication screen (diagnosis).

Fig. 9-1 Indication screen

- ① Indication
- ② [Scales connected] symbol
- ③ [Instrument recognition] symbol
- ④ [Menu] symbol
- ⑤ Deficit
- ⑥ [Deficit Reset] symbol
- ⑦ Shift key
- ⑧ Inflow volume/maximum flow range (max. flow)
- ⑨ Suction function display
- ⑩ [Fluid bag] symbol
- ⑪ [START]/[STOP] function field
- ⑫ [PAUSE]/[CONTINUE] function field
- ⑬ [Save] symbol
- ⑭ Status line for messages
- ⑮ Nominal pressure setting button
- ⑯ Nominal pressure display (numeric)
- ⑰ Actual Pressure Display
- ⑱ [Back] symbol


Fig. 9-2 Hysteroscopy screen displays with no fluid monitoring unit connected

- ① Indication
- ② **[Scales not connected]** symbol
- ③ **[Instrument recognition]** symbol
- ④ **[Menu]** symbol
- ⑤ Inflow volume
- ⑥ **[Deficit Reset]** symbol
- ⑦ Max. flow range (max. flow)
- ⑧ Setting buttons, max. flow range (max. flow)
- ⑨ Max. flow range (max. flow) of the suction function
- ⑩ **[Fluid bag]** symbol
- ⑪ **[START]/[STOP]** function field
- ⑫ **[PAUSE]/[CONTINUE]** function field
- ⑬ **[Save]** symbol
- ⑭ Status line for messages
- ⑮ Nominal pressure setting button
- ⑯ Nominal pressure display (numeric)
- ⑰ Actual pressure display
- ⑱ **[Back]** symbol

If the fluid monitoring unit is not connected to the device, then the device will no longer display the deficit volume but instead the consumed amount of fluid (inflow volume (4)).

9.3 Inserting a Tube Set

The following tube sets can be used for the hysteroscopy indication (chapter 19 Accessory List):

Hysteroscopy tube types

Deficit amount with no fluid monitoring unit connected

Article no.	Tube category	Tube type
T0505-01	Irrigation tube set	Disposable (one-time use)
T0506-01	Irrigation tube set	Reusable (NOT FOR SALE IN USA)
T0504-01	Vacuum tube set, incl. filter	Can be used 30 days
T0502-01	Suction tube set (1 connection)	Disposable (one-time use)
T0503-01	Suction tube set (2 connections)	Disposable (one-time use)

Table 9-1

A standard tube set is not inserted into the tube retainer until the device has been switched on and the device check has been completed.

CAUTION!

When inserting a tube set already filled with irrigation fluid into the tube retainer, it is possible to damage the membranes. Only insert unfilled tube sets to avoid such damage.

For additional information, see 6.2 Using the Tube Sets and 6.4 Irrigation tube set connection.

The deficit volume is displayed if a scale is connected.

Max. flow range display (max. flow)

Improved pressure measuring due to the instrument detection function

9.4 Surgical Mode

In surgical mode, the volume display of the device depicts the amount of fluid remaining in the patient during surgery and is not collected. To determine this value, the fluid irrigation consumption and the volume of the collected irrigation fluid are measured. The difference between both values yield the fluid volume remaining in the patient or lost during surgery. This differential (deficit) volume is depicted on the screen: see deficit display ⑥ in Fig. 9-1 Indication screen.

CAUTION!
Try to collect all the fluid that runs out of the cavum uteri during the procedure in order to achieve the most exact balancing possible.

9.4.1 Max. flow Range Display (max. flow)

In surgical mode, the inflow volume is displayed below the deficit volume. Use the Shift key ⑦ (Fig. 9-1) to toggle between the inflow volume and max. flow rate range (max. flow) display.

9.5 Automatic Instrument Recognition

The device features the automatic instrument detection function when using the [HYSTEROSCOPY]. The automatic instrument detection function runs in the back after each start of the irrigation process and serves to compensate pressure losses due to the flow through the narrow working channel. The pressure drop at the instrument is here included in the measured value of the actual pressure.

The indication menu has an option to enable or disable the instrument recognition function (chapter 9.14.4).

9.6 Preselecting the Nominal Pressure

WARNING!

If the current pressure does not react to an increase of the flow value during surgery, a perforation of the cavum uteri might be the cause. This results in an increased risk for bacteria entering the body. Examine the uterus for injuries.


The nominal pressure can be increased or decreased while the device is being used or not used. Values may be in the range of 15 to 150 mmHg.

- Briefly press the ▲/▼ symbol to increase or decrease the nominal pressure in increments of 5 mmHg. Pressing the ▲/▼ symbol longer than 1.5 seconds activates scrolling of the values in increments of 10 mmHg.

Factory setting: 60 mmHg

Increase/decrease nominal pressure

9.7 Safety Threshold Nominal Pressure Over 100 mmHg

To avoid dangerous potential pressures, the device is equipped with a safety threshold when setting nominal pressures over 100 mmHg.

An acoustic warning sound is emitted and **Pressure over 100 mmHg?** is depicted in the status line. After about 2 seconds, the pressure can be set to a value exceeding the 100 mmHg threshold.

9.8 Preselecting the max. flow range (max. flow)

CAUTION!

If the max. flow rate is set too low, the nominal pressure cannot be reached.

The max. flow range (max. flow) can be set while the device is being used or while it is not used. Values may be in the range of 50-800 ml/min


- Briefly press the ▲/▼ symbol to increase or decrease the nominal flow range (max. flow) in increments of 50 ml/min. Pressing the ▲/▼ symbol longer than 1.5 seconds activates scrolling of the values in increments of 100 ml/min.
- Press the [Save] symbol for 1.5 s to save the entered value.

Max. flow factory setting: 800 ml/min

Increase/decrease nominal flow

Saving the nominal pressure value

Reset differential volume

9.9 Saving the Nominal Values

The set nominal pressure value can be saved. Next time the hysteroscopy indication is opened, the last saved value for the nominal pressure (for values 15 to 80 mmHg) is used automatically. If the nominal pressure value to be saved exceeds 80 mmHg, the nominal pressure value is reset automatically to 80 mmHg next time the hysteroscopy indication is opened.

- Press the **[Save]** symbol for 1.5 s to save the entered nominal pressure value.

9.10 Resetting the Deficit

You can reset the deficit to the value 0.

- Press the **[->0<-]** symbol (reset) for 1.5 s to reset the deficit to the value 0. The status line depicts the following: **Fluid deficit reset**.

9.11 Starting and Stopping Irrigation

- Insert the tube set as described in chapter 9.3.
- Open tube clamps at irrigation tube.
- Fully open the inflow valve of the instrument and close the outflow valve.
- Press the **[START]** function field to start the irrigation process. The actual pressure display shows the current measured value. The roller wheel starts to turn. The device now carries out an **instrument recognition**.

CAUTION!

The instrument recognition must be performed outside and at the level of the patient.

- The status line depicts **Instrument recognition**. Wait until the instrument recognition has finished.
- Once the instrument recognition is finished, the status line depicts **Instrument recognition completed** and 1x informational beep is emitted.
- Check the deficit value and reset it if necessary.
- You may now start the surgical procedure.
- Press the **[STOP]** function field or **[II]** (PAUSE) to stop the irrigation process.
- The instrument recognition is restarted after the pump is stopped with the **[STOP]** function field and started again.
- Press the **[II]** (CONTINUE) function field again to continue the OR without a renewed instrument recognition.

WARNING!**Instrument replacement**

Stop the device using the START/STOP switch if replacing the instrument during surgery.

NOTE!

The plausibility of the instrument recognition is always checked during operation. If the plausibility check yields an error, the instrument characteristic is rejected. The pump then uses a safe standard curve.

9.12 Changing Bag during Surgery

If a connected fluid monitoring unit is used, the device system recognizes a bag change automatically. Brief fluctuations in the deficit calculation (< 10 s) may occur when a bag is replaced. The change is indicated with the message **Bag Change**.

Bag change during surgery

WARNING!

Touching the bags and the bag hooks as well as vibrations of the balancing system should be avoided during surgery to prevent triggering a false detection of the bag change and not negatively affect the accuracy of the deficit calculation.

WARNING!

The empty bags should remain on the hooks so that the measurement accuracy is not diminished.

WARNING!

Bags should be changed quickly to avoid affecting the accuracy of the deficit calculation.

9.13 Changing Container during Surgery

If a connected fluid monitoring unit is used, the device system recognizes a container change automatically. Brief fluctuations in the deficit calculation (< 10 s) may occur when a container is replaced. The change is indicated with the message **Container change**.

Container change during surgery

WARNING!

Touching the containers and their holders as well as vibrations of the balancing system should be avoided during surgery to prevent triggering a false detection of the container change and not negatively affect the accuracy of the deficit calculation.

WARNING!

Containers should be changed quickly to avoid affecting the accuracy of the deficit calculation.


CAUTION!

Full containers must be replaced immediately without stopping surgery. If the overflow protection of the containers is triggered, suction is stopped to prevent the ingress of fluids.

9.14 Opening and Configuring the Hysteroscopy Indication Menu

Device parameters are displayed and changed in the **[HYSTEROSCOPY]** menu. The following chart provides an overview over the possible settings.

1. Press the **[Menu]** symbol while in the hysteroscopy indication to access the settings. Press **[HYSTEROSCOPY]**.
2. Press one of the 4 function fields: **[Deficit]**, **[Fluid bag]**, **[Scale test]**, **[Instrument recognition]** or **[Factory settings]**.

The following is a detailed description of how to change the indication-specific device parameters.

9.14.1 Setting the Deficit Threshold

Set deficit threshold

The deficit threshold defines the threshold value for triggering the warning messages.

1. Press the **[Deficit]** function field. The deficit threshold can be set within a range of 0.1 to 2.5 l.
2. Press the **▲** or **▼** symbol to increase or decrease the deficit threshold. Pressing this symbol briefly increases or decreases the value in increments of 0.1 l, pressing it longer in increments of 0.2 l/min.

3. Press the  symbol to return to the previous menu level while saving the values.

Default setting: 1.0 l

9.14.2 Fluid Bag Configuration

1. In the menu, press the **[Fluid bag]** function field.
2. Press the **[ON]** function field to enable or the **[OFF]** function field to disable the bag change warning function. The active function field is depicted in green.
3. Press the  symbol to return to the previous menu level. The setting is automatically saved.

Default setting: Off

Configure warnings

1. Press the  or  symbol to increase or decrease the height between device and fluid bag. Pressing this symbol briefly increases or decreases the value in increments of 0.1 m, pressing it longer in increments of 0.5 m.
2. Press the  symbol to return to the previous menu level while saving the values.

Default setting: 1.00 m

Adjusting the Height (Use without fluid monitoring unit)

1. Press the  or  symbol to set the volume of the fluid bags to a value between 0.5 and 5 l.
2. Press the  symbol to return to the previous menu level. The setting is automatically saved.

Default setting: 1 l

Set volume

(see Fig. 6-3 Attaching the fluid bags)

NOTE!

The adjustment of the bag height is inactive when using a fluid monitoring unit since the height of the bags is predetermined by the fluid monitoring unit.

9.14.3 Scales Test

See chapter 12.4

9.14.4 Setting the Instrument Recognition

1. In the menu, press the **[Instrument recognition]** function field.
2. Press the **[ON]** function field to enable, or the **[OFF]** function field to disable the instrument recognition function. The active function field is depicted in green.
3. Press the  symbol to return to the previous menu level. The setting is automatically saved.

Factory setting: AN (responsivity)

9.14.5 Resetting to Factory Settings

Use the hysteroscopy menu to reset the hysteroscopy device parameters to the factory settings.

Hysteroscopy device parameters	Hysteroscopy factory settings
Nominal pressure	60 mmHg
Max. flow range (max. flow)	800 ml/min
Deficit limit	1.5 l
Fluid bag height	1.0 m
Fluid bag size	1 l
Fluid bag warning	OFF

Table 9-2

Reset to factory settings

1. Press the **[Factory settings]** function field to reset the HYS device parameters to the factory settings.
You can now do the following:
2. Press the **[OK]** symbol to save the setting, or
3. Press  to return to the previous menu level without saving.

9.15 Safety Functions

9.15.1 General Safety Functions

During device operation, the status line will display safety and warning messages. These messages refer to the handling and safety of the device during use. Please see chapter 11 Safety Functions for additional information concerning the general safety functions.

Deficit rate exceeded

9.15.2 Exceeding the Deficit Rate

If the deficit rate of 300 ml/min is exceeded, 5 warning beeps are emitted twice and the screen displays the following in red letters: **High Fluid Loss!** The warning message is repeated after 15 s.

Differential volume = threshold value

9.15.3 Reaching and Exceeding the Deficit Threshold

The deficit threshold defines the threshold value for triggering the warning messages (see chapter 9.14.1 Setting the Deficit Threshold).

If the differential volume reaches the set threshold (0.1-2.5 l), 3 warning signals (beeps) are emitted and the screen depicts the following: **Deficit reached!**

The deficit value on the screen is now displayed with a red progress bar.

Differential volume > threshold value

If the deficit volume exceeds the set threshold value, 5 short warning beeps are emitted and the screen depicts the following in red letters: **Deficit exceeded!**

Every additional 100 ml the value exceeds the threshold value is indicated with the same warning signal.

9.15.4 Exceeding the Nominal Pressure

Overpressures may occur during device operation. The device features corresponding safety mechanisms that support the physician performing the surgery in deciding on appropriate responses.

If the actual pressure exceeds the nominal pressure by more than 10 mmHg or 150 mmHg, the following safety mechanism are automatically triggered by the device:

Actual pressure > nominal pressure + 10 mmHg or 150 mmHg

- Immediately: The status line depicts **Overpressure!** Three beeps are emitted.
- A pressure relief function is activated that consist of turning the roller wheel backwards.
- If the actual pressure falls below the nominal pressure, the roller wheel starts turning forward again and irrigation is continued. The warning message in the status line is no longer displayed.

If the actual pressure exceeds a value of 200 mmHg, the following safety mechanism are automatically triggered by the device:

Actual pressure > 200 mmHg

- Immediately: The status line depicts **Critical overpressure!** in red letters. Five warning beeps are emitted.
- A pressure relief function is activated that consist of turning the roller wheel backwards.
- If no pressure relief is achieved after a further 5 seconds, the roller wheel stops.
- If the actual pressure falls below the nominal pressure, the roller wheel starts turning forward again and irrigation is continued. The warning message in the status line is no longer displayed.

CAUTION!

The message "Critical Overpressure!" overwrites all other possible messages. Overpressure warnings serve as information for the attending physician to respond appropriately to the existing excess pressure. The intrauterine pressure can be reduced by opening the outflow cannula, for example. If possible, the attending physician should perform this task.


10 Menu

Use the menu to do the following:

- Change the general device parameters
- Change the indication-specific parameters
- Call the service menu (only for authorized service personnel!)

The following chart provides an overview over the possible general settings.


10.1 Overview Menu - General Settings

10.2 General Settings

1. Switch device on.
2. Press the **[Menu]** symbol in the indication selection

Opening the general settings menu

3. or in an indication screen, press the **[Menu]** symbol and then the **[General]** function field.

4. Select one of the following function fields on the first page: **[Language]**, **[Volume]**, **[Brightness]**, **[Factory settings]**, **[Device info]**, **[Demo]** or **[Service]**.


The following is a detailed description of the general device parameters.

10.2.1 Selecting the Language

1. In the **General** menu, press the **[Language]** function field.
2. Press the function field of the desired language (e.g., **[English]**). The enabled function field is depicted in black.
3. Press the  symbol to return to the previous menu level. The setting is automatically saved.
Factory setting: English

Select Language

10.2.2 Setting the Volume of the Acoustic Signals and Audio Warnings

1. In the **General** menu, press the **[Volume]** function field.
2. Press the  or  symbol to set the desired volume to a value between 1 and 7.
3. Press the  symbol to return to the previous menu level. The setting is automatically saved.
Factory setting: 4

Setting the volume

10.2.3 Setting the (Display) Brightness

1. In the **General** menu, press the **[Brightness]** function field.
2. Press the  or  symbol to adjust the desired brightness to a value between 1 and 11.
3. Press the  symbol to return to the previous menu level. The setting is automatically saved.
Factory setting: 6

Setting the brightness

10.2.4 Device Info

In the **General** menu, press the **[Device Info]** function field. The serial and software version are displayed.

10.2.5 Demo

In the **General** menu, press the **[Demo]** function field.

The **Demo** mode serves to demonstrate the pump functionality. The functionality of the irrigation stops after 2 minutes and is blocked for another 30 seconds. If the pump is stopped before the 2 minutes have expired, you have to wait 30 seconds before restarting.

Tubes are not devalued during operation in demo mode.

10.2.6 Factory settings

In the **General** menu, press the **[Factory Settings]** function field.

You can now do the following:

1. Press the **[OK]** symbol to reset all general and all indication-specific device parameters, or
2. Press the  symbol to return to the previous menu level without saving.

10.2.7 Service

The **Service** menu is password protected and access is only intended for trained service technicians.

The following functions are available in the **Service** menu:

- Sensor calibration
- Motor calibration
- Touchscreen calibration

Please consult the service manual for additional information.

Opening the service menu

11 Safety Functions

The electronic components continuously monitor the proper function of the device. Device malfunctions are indicated with audible warning signals (beeps), error messages, and/or by blocking device functions. A table listing a summary of possible error and warning messages is provided in chapter 16 Error and Warning Messages.

11.1 Device Self-Test

After being switched on, the device performs a self-test of the sensors, the motor, and electronic components. The following describes the messages for defects of the individual modules or components.

11.2 Sensor Error

The offset values are verified as part of the device check process that is executed when the device is switched on. If an impermissible deviation or error is detected in the electronic pressure measurement system, the following message will appear on the display: **Sensor error** and 5 warning signals are emitted. Restart the device. If the error recurs, please contact the Service Department. The **Service** key will appear for access to the Service menu.

Sensor errors

A **Sensor error** may also occur due to pressure on the sensor while a tube set is not inserted.

Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

11.3 Motor Error

If the motor drive is defective, the following message will appear on the display: **Motor error** and 5 warning signals are emitted. Restart the device. If the error recurs, please contact the Service Department. The **Service** key will appear for access to the Service menu.

Motor error

Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

11.4 Electronics Error

If the electronic system is defective, the following message will appear on the display: **Electronics error** and 5 warning signals are emitted. Restart the device. If the error recurs, please contact the Service Department. The **Service** key will appear for access to the Service menu.

Electronics error

Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

11.5 Calibration Error

If a calibration error occurs, the following message will appear on the display: **Calibration error** and 5 warning signals are emitted. Restart the device. If the error recurs, please contact the Service Department. The **Service** key will appear for access to the Service menu.

Calibration error

Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

12 **Function Test**

WARNING!
Functional test
The functional test must be performed prior to each surgery.

WARNING!
Checking the warning signals
The warning signals must be checked prior to each operation. The system is to be set up so that all warning signals can be perceived.

WARNING!
Sterilize reusable instruments and tubing before surgery to prevent infections.
Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

WARNING!
Original accessories
For your own safety and that of your patient, use only original accessories.

WARNING!
Device Errors
Do not use this device if a defect is suspected or detected during the function check. It is prohibited to use the device in the case of obvious defects.

12.1 Preparing the Function Test of the Device

Prepare device as outlined by the following function test instructions.

Preparing the function test

- ▶ Use your hand to control the roller wheel and slightly turn it to check for ease of movement.
- ▶ Hang the fluid bags as described in chapter 6.6 Attaching and Connecting the Irrigation Fluid Bags (Use without Fluid Monitoring Unit).
- ▶ Connect the tube and an instrument as described in chapter 6.4 or Fig. 6-5.
- ▶ Check to make sure all tube connections are free of mechanical stresses and are routed without snagging.

12.2 Performing Device Function Test (Only for Laparoscopy Indication)

1. Set the standard flow rate.
2. Press the **[START]** function field.
3. Open inflow valve at instrument. The roller wheel is turning forward and delivers irrigation fluid.
4. Close the inflow valve of the instrument. The roller wheel stops.
5. Press the **[STOP]** function field. The roller wheel turns approx. one rotation backwards.

The function test of the device is successfully completed when the roller wheel behaves as described in item 3., 4., and 5.

12.3 Performing Device Function Test (Only for Hysteroscopy Indication)

1. Set the nominal pressure to 45 mmHg and the nominal flow to 500 ml/min.
2. Press the **[START]** function field.
3. Close the inflow valve of the instrument about halfway so that device can still deliver fluid through the instrument.

The function test of the device is successfully completed when the progress bar of the actual pressure depicts a value of approx. 45 mmHg.

12.4 Performing Function Test of Fluid Monitoring Unit (Only for Hysteroscopy Indication)

1. The switched on device must be connected to the scale (chapter 4.3). There must be no weight on the bag hooks or in the receptacle for the containers.
2. First, select the **[Hysteroscopy]** indication and then the **[Settings]** function field.
3. Press the **[Scales test]** function field.
4. The current weights of the two scales modules are displayed. Without loading the scales, a value of 0 g should be displayed.
5. Place a defined weight between 1 to 5 kg on the bag scale (e.g., 3 kg).
6. The attached weight (e.g., 3000 g) must be displayed in the **Fluid bag value** field. The permissible tolerance is +/- 100 g. If a greater difference is detected, an authorized service technician must re-calibrate the scale.
7. Place a defined weight between 1 to 5 kg into the receptacle for the containers (e.g., 3 kg).
8. The attached weight (e.g., 3000 g) must be displayed in the **Container value** field. The permissible tolerance is +/- 100 g. If a greater difference is detected, an authorized service technician must re-calibrate the scale.

Test connected scale

The function check of the balancing module is successfully completed when the value does not exceed the permissible tolerance of max. 100 g.

13 Care and Maintenance

Special care is necessary when servicing, maintaining, and storing the device and its accessories to maintain the functionality of the device and its accessories.

13.1 Cleaning the Device

- 1. Use the **ON/OFF** button to turn the device off.
- 2. Remove the power cable.
- 3. Wipe the surface of the device with a soft cloth moistened with the surface disinfectant based on alcohol or aldehyde (for example Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the device.

NOTE!
Do not sterilize the device.

13.2 Maintenance Intervals

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. The device needs to be serviced at least every two years. The tests are described in chapter 14.

Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

NOTE!
Service or maintenance work may not be carried out during surgery.

13.3 Maintenance by Authorized Service Technician

An authorized service technician has to inspect and service the device at appropriate intervals to ensure the safety and functionality of the unit. The minimum service interval is two years, depending on frequency and duration of use. If the service interval is not maintained, the manufacturer does not assume any liability for the functional safety of the device.

A sticker located on the rear panel of the device will remind you of the latest date for the next service or maintenance check.

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians.

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

Ask the service technician for a certificate after he or she has inspected the unit or performed any service tasks. This certificate lists the type and scope of the service as well as the date and name of the servicing company together with the signature of the service technician.

Manufacturer's specifications

Two-year maintenance interval

Authorized trained personnel

Unauthorized personnel

Liability

Technical documents

Certification

13.4 Changing the Fuse

CAUTION!

Before replacing the fuse, check the values of the fuse to be inserted acc. to chapter 17 Technical Data Pump.

The fuse may be defective and is in need of replacement if:

- displays and LEDs (if available on your equipment) do not light up,
- the device does not function.

Check to make sure

- the main power supply cable is properly connected to the power supply input and to a safety socket,
- the house power supply fuse is functioning.

WARNING!


Unplug the power cable from the device before checking the fuse.

The device does **not** have to be opened to replace the fuse.

1. Switch device off.
2. Disconnect device from power supply by pulling the mains plug from the mains socket.
3. Remove power connection cable from mains socket.
4. The fuse holder is located next to the mains socket.
5. Remove fuse holder as depicted in Fig. 13-1.
6. **A** Undo the latch of the fuse holder with a small screwdriver.
7. **B** Remove the fuse holder.
8. **C** Check fuse.
9. Insert a new fuse. Use only the specified type of fuse (see chapter 17 Technical Data Pump).
10. Insert the fuse holder until it can be heard snapping into place.
11. Use the power cable to reconnect the shockproof safety socket with the rear mains socket.

Fig. 13-1 Open the fuse holder

13.5 Cleaning the accessories

- 1. Switch device off.
- 2. Disconnect device from power supply.

The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Wipe the surface of the device with a soft cloth moistened with the disinfectant. Make sure moisture does not enter the device.

The manufacturer recommends **Meliseptol® rapid** as disinfectant.

WARNING!
Cleaning the fluid monitoring unit
Do not sterilize the fluid monitoring unit.

13.6 Care of the Reusable Tube Set

The device can be used with a reusable tube set that is listed as an accessory (see chapter 19 Accessory List) Please comply with the following information when using a reusable tube set.

13.6.1 Cleaning the Reusable Tube Set

WARNING!
The reusable tube set is made of silicone, polysulfone (PSU), and stainless steel. For cleaning and disinfection, use only pH-neutral or slightly alkaline cleaning agents (e.g. neodisher MediClean 2.0 %), disinfectants (e.g. Lysetol V 8 %), and drying and rinsing agents approved for the tube set materials. Using unsuitable agents (e.g. neodisher MediKlar rinsing agent) may damage the tube system and especially the PSU connectors.

Reprocessing restrictions

The manufacturer has tested and approved the reusable tube set for a specific number of reprocessing cycles. Always comply with the instructions on the label.

Never exceed the number of uses indicated by the manufacturer.

WARNING!

Please check the reusable tube for signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation.

The reusable tube set should be reprocessed as soon as possible following use (recommended: max 30 minutes). In case of heavy contamination (blood, tissue) it is recommended to remove these residues with a soft cloth while still at the site of the surgery. This cloth may be moistened with a slightly alkaline cleaner. The tube set is then transported, preferably in a dry basket, to the actual reprocessing location where the complete process is carried out.

Cleaning preparations**1. Preparatory cleaning at the reprocessing location:**

Disassemble the individual parts of the tube set.

Notes concerning the pressure chamber membrane (for irrigation tube sets only).

WARNING!

Remove the membrane before cleaning.

- Before cleaning: Carefully remove the membrane (2) from the pressure chamber (1). Make sure not to damage the membrane. To do so, pull upwards on the membrane strap (6).
- Reinstalling membranes before sterilization: Place the lip of the membrane (3) into the ring groove of the pressure chamber (1). The strap (6) must be positioned in the notch (5) provided for this purpose.
- Press the lip of the membrane into the groove (4). A properly inserted membrane is flush with the pressure chamber and exhibits no wrinkles.

2. Preparatory cleaning:

Use a soft cloth to carefully wipe off each component. Remove any remaining blood and tissue residues. The cloth to be used may be moistened with a slightly alkaline cleaner.

1. Automatic cleaning and disinfecting:


The cleaning and subsequent disinfection of the individual components of the tube set can be carried out with the appropriate programmable cleaning and disinfecting systems for medical devices. The manufacturer has validated the applicability of this method in terms of hygiene success and material compatibility using the "Vario" program of the disinfectant manufactured by the company Miele.

Suitability of the cleaning device used in the context of automated cleaning and disinfection, the selected program, and the used cleaning or disinfecting agent, as well as a guarantee with respect to the attainable hygiene status and material compatibility are the sole responsibility of the operator!

Individual automatic processing methods must be independently validated by the operator.

2. Manual cleaning:

- Place the prepared individual components of the tube set in demineralized water (room temperature, 20 °C to 30 °C) for a period of at least 3 minutes to 5 minutes and then rinse them off under running water. Any remaining residues such as adherent clotted blood should be removed with a soft brush.
- Use a soft cloth to thoroughly dry the components.
- The completely dried tube set components must now soak for a period of 25 minutes to 30 minutes in a suitable cleaning agent (e.g., neodisher Med-iClean 2.0 % or equivalent). Comply with the instructions and notes of the manufacturer of the used cleaning agent. Make sure the components are fully immersed in the cleaning agent.
- Place the components once more in demineralized water (room temperature, 20 °C to 30 °C) for a period of at least 3 minutes to 5 minutes and then

**Cleaning**

Disinfection

rinse them off under running water.

- Flush the silicone tube for 15 seconds with demineralized water. Then use a high-pressure cleaning gun to emit 10 pressure surges, each one with a length of one second (1 second) into the silicone tube.
- After cleaning, allow the components to dry. For this purpose, they are to be placed in an appropriate tray and allowed to drip dry or use an air gun to dry the components. It is also possible to use, where appropriate, a sterile, soft cloth. Or, the components can be dried in an air dryer (10 minutes at 100 °C).

13.6.2 Disinfecting the Reusable Tube Set

After cleaning, the individual tube set components must be disinfected. Only disinfect a thoroughly cleaned tube set.

1. Automatic disinfection:

The cleaning and subsequent disinfection of the individual components of the tube set can be carried out with the appropriate programmable cleaning and disinfecting systems for medical devices. The manufacturer has validated the applicability of this method in terms of hygiene success and material compatibility using the "Vario" program of the disinfectant manufactured by the company Miele.

Suitability of the cleaning device used in the context of automated cleaning and disinfection, the selected program, and the used cleaning or disinfecting agent, as well as a guarantee with respect to the attainable hygiene status and material compatibility are the sole responsibility of the operator!

Individual automatic processing methods must be independently validated by the operator.

2. Manual disinfection:

- The cleaned and completely dried tube set components must now soak for a period of 25 minutes to 30 minutes or as specified by the manufacturer in a suitable disinfectant (e.g. Lysetol 8 % or equivalent). Make sure the components are fully immersed in the cleaning agent and are not stacked on top of one another. Comply with the instructions and notes of the manufacturer of the used disinfecting agent. The tube set can be damaged if the concentration is too high or the components are soaked too long.
- After the soaking time, the individual tube set components should be removed from the disinfectant liquid with pliers with smooth jaws.
- Place the components once more in demineralized water (room temperature, 20 °C to 30 °C) for a period of at least 3 minutes to 5 minutes and then rinse them off under running water.
- Flush the silicone tube for 15 seconds with demineralized water. Then use a high-pressure cleaning gun to emit 10 pressure surges, each one with a length of one second (1 second) into the silicone tube.
- After cleaning, allow the components to dry. For this purpose, they are to be placed in an appropriate tray and allowed to drip dry or use an air gun to dry the components. It is also possible to use, where appropriate, a sterile, soft cloth. Or, the components can be dried in an air dryer (10 minutes at 100 °C).

13.6.3 Sterilizing the Reusable Tube Set

Reassemble the tube set from the cleaned and disinfected parts as described in the manual.

Preparing for sterilization

Inspection and function test

Check the individual components and the assembled tube set for signs of damage. Never sterilize a tube set that shows signs of damage, especially brittleness and perforation. Tube sets with evidence of or suspected of damage must be blocked from further use.

Maintenance

Maintenance of the tube set is not required. In case of loss or damage to the pressure chamber membrane, this component must be replaced (see instructions in the manual).

Only clean, dry, disinfected, and assembled tube sets should be sterilized.

Sterilization

- **Packaging**

Pack the dry and reassembled tube set in a suitable sterile packaging for sterilization (cloth, bag, or tray).

- **Sterilization**

The manufacturer has validated saturated steam sterilization as a suitable method of sterilization for the tube set based on the parameters listed below using the pre-vacuum process. Sterilize the tube set according to these specifications.

	Pre-vacuum
Sterilization method	Pre-vacuum sterilizer
Packaging	Packed
Cycle time	5 min (273 °F/134 °C, 3 bar)
Drying time	10 min

CAUTION!

The indicated drying time depends on several variables, including the following: Altitude, humidity, type of packaging, preconditioning, size of chamber, mass of load, and placement in chamber. Users must verify that the drying time set in their autoclave yields dry surgical equipment when using the method of saturated steam sterilization described here.

Please comply with the instructions of the operating manual included with your autoclave.

- **Storage**

Always store the tube set after sterilization in a suitable sterile packaging or in a suitable environment to avoid contamination.

- **Additional information**

When sterilizing several tube sets, instruments etc. at once using one sterilization cycle, make sure not to exceed the max. load of the sterilizer.

Only clean, dry, disinfected, and assembled tube sets should be sterilized.

Sterilization

- **Packaging**

Pack the dry and reassembled tube set in a suitable sterile packaging for sterilization (cloth, bag, or tray).

- **Sterilization**

The manufacturer has validated saturated steam sterilization as a suitable method of sterilization for the tube set based on the parameters listed below using the pre-vacuum process. Sterilize the tube set according to these specifications.

	Pre-vacuum
Sterilization method	Pre-vacuum sterilizer
Packaging	Packed
Cycle time	5 min (273 °F/134 °C, 3 bar)
Drying time	10 min

CAUTION!

The indicated drying time depends on several variables, including the following: Altitude, humidity, type of packaging, preconditioning, size of chamber, mass of load, and placement in chamber. Users must verify that the drying time set in their autoclave yields dry surgical equipment when using the method of saturated steam sterilization described here.

Please comply with the instructions of the operating manual included with your autoclave.

- **Storage**

Always store the tube set after sterilization in a suitable sterile packaging or in a suitable environment to avoid contamination.

- **Additional information**

When sterilizing several tube sets, instruments etc. at once using one sterilization cycle, make sure not to exceed the max. load of the sterilizer.

Sterilization with ETO

Sterilization with ethylene oxide (ETO) is possible but has not been approved by the manufacturer.

Sterilization with gamma rays

Do not sterilize with gamma rays.

14 Annual Inspection

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. These inspections have to be carried out on an annual basis. Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

The tests described in this chapter are designed specifically for trained personnel or a hospital technician. The operation of the device as well as its functionality and serviceability are easily checked. Each test conducted has to be documented with date and signature in chapter 21.1.

The following measuring tools and resources were used by the manufacturer to determine the listed measurements and tolerances:

- Original tube set
- Fluid container (3 l)
- Graduated measuring cup with 1l scaling
- Stopwatch
- Time measuring device

WARNING!

An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

Manufacturer's specification

Inspection tests

Measured values and tolerances

14.1 Electrical Safety Test

1. Perform a visual inspection. Make sure that
 - the fuse corresponds with the specifications indicated by the manufacturer,
 - labels and stickers on device are legible,
 - the mechanical condition of the device allows for its safe use,
 - the device is clean to ensure proper and safe functionality.
2. Carry out the measurements for the ground leakage current, short-circuit current/housing leakage current, and the protective conductor resistance as per IEC 62353 in the current version or according to the applicable national standard.

14.2 Basic Function Test

A basic function test checks the displays, keys, and performance of the device. For this test, you will need the following:

- One disposable tube set
- Fluid container (3 l)
- Stopwatch
- One 1-liter measuring cup (100 ml scaling)

The layout of the test and its setup are depicted in Fig. 14-1 Setup of pressure measuring test.

1. Switch device on. Wait until the device check has finished.
2. Select one of the possible indications.
3. Insert a disposable tube set into the tube retainer.
4. Hang fluid bag at a height of $h=100$ cm as depicted in Fig. 14-1 Setup of pressure measuring test and connect bag with irrigation tube.
5. Place the end with the Luer lock into the measuring cup.
6. Set the following nominal values:

Hysteroscopy Indication:

- Nominal pressure = 100 mmHg
- Nominal flow = 500 ml/min

Laparoscopy indication:

- Flow rate low = 1000 ml/min

- Press the **[START]** function field.
- Wait until the tube set is completely filled with irrigation fluid.

Hysteroscopy indication:

Allow the device to pump fluid for at least 1 minute (instrument recognition function).

- Clamp off the instrument tube end in the measuring cup without stopping the pump.
- Empty the measuring cup and place the end of the tube back into the measuring cup.
- Release the tube end. The irrigation process runs. Press the start button of the stopwatch.

12. Hysteroscopy indication:

After 2 minutes, press the **[STOP]** function field and the stop button of the stopwatch at the same time. The measuring cup must contain approx. 1 l ($\pm 10\%$) of fluid.

Laparoscopy indication


After 1 minute, press the **[STOP]** function field and the stop button of the stopwatch at the same time. The measuring cup must contain approx. 1 l ($\pm 10\%$) of fluid.

The basic function test has been successfully completed once these values have been reached.

14.3 Pressure Measuring Test

The layout of the test and its setup are depicted in Fig. 14-1 Setup of pressure measuring test.

Fig. 14-1 Setup of pressure measuring test

The pressure measurement test checks the proper functioning of the pressure measurement. This test requires a complete tube set and a container filled with **water**. The height of the water column (hydrostatic pressure) is used for measuring and is then converted to mm mercury column (mmHg). The height of the water column above the pressure chamber has to match the value of the actual pressure display after conversion.

Conversion formula: $p \text{ (cm H}_2\text{O)} \times 0.74 = p \text{ (mmHg)}$

- Switch device on. Wait until the device check has finished.
- Select one of the possible indications.
- Press the **[START]** function field. Wait until the tube set is completely filled and bubble-free.

Important: The tube set behind the roller wheel must be filled completely with fluid.

4. Press the **[STOP]** function field.
5. Hold the end of the completely filled instrument tube at the height of the sensor.
6. The actual pressure display should show 0 (+5 mmHg).
7. Hold the instrument tube at a level of approx. 67.5 cm above the roller wheel. Make sure the tube segment between this point and the cartridge is completely filled with fluid. Also make sure that no point of the tube is higher than the indicated value.
8. The height of the water column is used to calculate the generated hydrostatic (water) pressure as follows: $67.5 \text{ cm H}_2\text{O} \times 0.74 = 50 \text{ mmHg}$.
9. The actual pressure display now should depict approx. 50 mmHg ($\pm 10 \text{ mmHg}$).
10. When varying the height of the tube set above the roller wheel, the actual pressure display also varies accordingly.

The test of the pressure measurement is successfully completed when the values of the actual pressure display correspond to those of the converted water column height. Please enter each completed test into the test log (see chapter 21.1 Test Log).

14.4 Performing Scale Function Test (for Hysteroscopy Indication)

1. A scale without weight is connected to the switched on device (see chapter 4.3 Set up and Connect the Pump with Fluid Monitoring Unit).
2. Select the hysteroscopy indication.
3. Select the **[Surgical]** function field.
4. Select the **[Settings]** function field.
5. Press the **[HYSTEROSCOPY]** function field to access the settings of the hysteroscopy mode.
6. Press the **[Scales test]** function field.
7. The weight display depicts **0 g**.
8. Place a defined weight between 1 to 5 kg on the scale (e.g., 3 kg).
9. The weight display of the device must depict the weight of the load in grams (e.g., 3000 g). The permissible tolerance is $< 100 \text{ g}$. If a greater difference is detected, an authorized service technician must re-calibrate the scale.

Testing the connected scale

The function test of the scale is successfully completed when the permissible tolerance of $< 100 \text{ g}$ is maintained.

Precautionary measures

ESD (Electrostatic Discharge) precautionary measures

15 Electromagnetic Compatibility

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

This device is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.

CAUTION!

Accessories

To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device GYN-Pump must be used only with the accessories listed in chapter 19.

To ensure the basic safety and essential functionality in relation to electromagnetic interference over the life of the device, the device must be restarted after 24 hours so that a diagnostic self-test can be performed. Furthermore, the service interval described in chapter 13.2 must also be adhered to.

15.1 Electrical Connections

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

Hospital employees should be informed about and trained in ESD precautionary measures.

15.2 Guidelines and Manufacturer’s Statement – Electromagnetic Emissions

The device GYN-Pump is intended for use in the electromagnetic environment specified below. The user/operator of the GYN-Pump should make sure the device is operated within such an environment.

Emitted interference measurements	Compliance	Electromagnetic environment guidelines
HF emission according to CISPR 11	Group 1	The device GYN-Pump uses HF energy solely for its internal functions. Therefore, the HF emission is very low and it is unlikely that devices in close proximity will experience interference.
HF emission according to CISPR 11	Class B	The device GYN-Pump is suitable for use in all facilities including those in residential areas and those directly connected to a public utility network supplying buildings used for residential purposes as well.
Emission of harmonic oscillations according to IEC 61000-3-2	Class A	
Emission of voltage fluctuations / flickers according to IEC 61000-3-3	In compliance	

15.3 Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity

The GYN-Pump is intended for use in an electromagnetic environment as described below. The user/operator of the GYN-Pump must make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	IEC 60601-1-2 test level	Compliance Levels	Electromagnetic environment/guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30%.
Electrical fast transients/bursts according to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines Modulation 100 KHz	± 2 kV for power lines ± 1 kV for input and output lines Modulation 100 KHz	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Voltage surges according to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brown-outs, and fluctuations of the power supply according to IEC 61000-4-11	$< 5\% U_T^*$ ($> 95\%$ dip in the U_T) for $1/2$ cycle $< 40\% U_T$ ($> 60\%$ dip in the U_T) for 5 cycles $< 70\% U_T$ ($> 30\%$ dip in the U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in the U_T) for 5 cycles	$< 5\% U_T^*$ ($> 95\%$ dip in the U_T) for $1/2$ cycle $40\% U_T$ (60% dip in the U_T) for 5 cycles $70\% U_T$ (30% dip in the U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in the U_T) for 5 seconds $0\% U_T$; $1/2$ cycle at 0° , 45° , 90° , 135° , 180° , 225° , 270° , and 315° $0\% U_T$; 1 cycle and $70\% U_T$; 25/30 cycles Single phase: at 0° $0\% U_T$; 250/300 cycles	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment. If the user/operator of the system requires the continuation of functionality after power interruptions/disruptions, it is recommended to supply the device with power from an uninterruptible power supply.
Supply frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.

*Note: U_T is the mains alternating voltage before applying the test levels.

15.4 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity

Electromagnetic interference immunity tests	Test level	Compliance	Electromagnetic environment guidelines
Conducted HF interference quantities according to IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V _{eff} 150 kHz to 80 MHz	<p>Portable and mobile wireless devices should not be used in closer proximity to the GYN-Pump (including cables/lines) than the recommended safety distance calculated based on the transmitting frequency and the applicable formula. Recommended safety distance:</p> <p>$d = 1.2\sqrt{P}$ for 150 KHz to 80 MHz $d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ for 800 MHz to 2.7 GHz</p> <p>With P as the rated output of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d as recommended safety distance in meters [m].</p> <p>The field strength of stationary transmitters for all frequencies tested on site ^a should be lower than the concordance level. ^b</p> <p>Interference is possible in the proximity of devices featuring the following pictograph.</p> 
Radiated HF interference quantities according to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	

Note 1: The higher frequency range applies for 80 and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

^a The field strength of stationary transmitters such as base stations of wireless phones and cell phones, ham radio operators, AM and FM radio and TV stations can theoretically not always determined in advance. A study of the installation site should be considered to determine the electromagnetic environment concerning the stationary transmitter. If the measured field strength at the proposed GYN-Pump installation and operation site exceeds the concordance levels listed above, the GYN-Pump should be monitored to document proper functionality and operation as intended. If unusual performance characteristics are observed, additional measures may be required such as changing orientation or the location of the GYN-Pump.

^b The field strength should be less than 3 V/m for the frequency range of 150 kHz to 80 MHz.

15.5 Recommended Safety Distances Between Portable and Mobile RF Telecommunications Devices and GYN-Pump

The GYN-Pump is intended for use in an electromagnetic environment where HF interferences are controlled. The user/operator of the GYN-Pump can contribute to lowering electromagnetic emissions by complying with the minimum distance between portable and mobile HF telecommunications devices (transmitters) and the GYN-Pump, depending on the output power of the communication equipment listed below.

Rated output of the transmitter [W]	Safety distance based on the transmitting frequency [m]		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The safety distance d in meters [m] for transmitters with a max. rated output not listed in the table above can be calculated by applying the corresponding formula in the respective column. P is the max. rated output of the transmitter in watts [W] as specified by the transmitter manufacturer.

Note 1: The higher frequency range applies at 80 MHz and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

WARNING!

Portable HF communication equipment can affect the performance characteristics of the GYN-Pump. Such equipment must therefore comply with a minimum distance of 30 cm (regardless of all calculations) from the GYN-Pump, its components and cables.

16 Error and Warning Messages

Error and warning messages are indicated by the device acoustically and visually. The visual indication is usually a message in the status line of the touchscreen or a full screen display. Depending on message or alert type, a certain number of beeps is emitted.

NOTE!

Each warning signal is to be assigned one of the corresponding priorities. If multiple signals of equal priority are called at the same time, they are processed according to an internal priority ranking.

There are four different message types for operating notes and errors:

- High priority errors: indicated with red letters and exclamation mark and 5 warning beeps.
- Medium priority errors: indicated with black letters and exclamation mark and 3 warning beeps.
- Low priority errors, input warnings or negative information signal: indicated with black letters and 1 warning beep.
- Operating notes: indicated with black letters and 1 informational beep.

Message	Cause	Troubleshooting
High priority		
Electronic error! Call service	Electronic error	<ul style="list-style-type: none"> • Restart the device. • If the malfunction or error occurs again, please contact the service department.
Sensor error! Call service	Pressure sensor defective	<ul style="list-style-type: none"> • Restart the device. • If the malfunction or error occurs again, please contact the service department.
Motor error! Call service	Motor error	<ul style="list-style-type: none"> • Restart the device. • If the malfunction or error occurs again, please contact the service department.
Calibration error! Call service	The device is not calibrated properly.	<ul style="list-style-type: none"> • The device must be recalibrated. Please contact the service department.
Critical overpressure!	Hysteroscopy: Actual pressure above 200 mmHg	<ul style="list-style-type: none"> • The physician must reduce the pressure, for example, by opening the outflow valve.
High fluid loss	Deficit change > 150 ml in 30 s	<ul style="list-style-type: none"> • Physician must respond appropriately.
Deficit limit exceeded!	In increments of 100 ml above the preset value	<ul style="list-style-type: none"> • Physician must respond appropriately.
Medium priority		
Overpressure!	Hysteroscopy: Actual pressure > nominal pressure + 10 mmHg and >150 mmHg	<ul style="list-style-type: none"> • The physician must reduce the pressure, for example, by opening the outflow valve.
Deficit limit reached!	Deficit exceeds the preset value	<ul style="list-style-type: none"> • Physician must respond appropriately.
Scales defective!	Incorrect scale values or communication with scale interrupted	<ul style="list-style-type: none"> • Connect new scale, restart device.
Container change!	Container change during running irrigation	<ul style="list-style-type: none"> • Reset container or press [Pause] button.
Changing Bags!	Bag change during running irrigation	<ul style="list-style-type: none"> • Finish bag change speedily.
Scales overloaded!	Bag scale overloaded with more than 23 kg Container scale overloaded with more than 16 kg	<ul style="list-style-type: none"> • The weight must be reduced.
Instrument recognition failed!	A valid instrument was not detected	<ul style="list-style-type: none"> • The instrument recognition must be carried out again.
Defective vacuum pump	Defective vacuum pump	<ul style="list-style-type: none"> • Restart the device. • If the malfunction or error occurs again, please contact the service department.

Tube set incorrectly inserted	Incorrect position of tube set	• Insert tube set correctly
Low priority		
Change fluid bag	Following consumption of 85% to 95% of the fluid bag	• Replace fluid bag
Pressure over 100 mmHg?	Nominal pressure setting	• Physician must respond appropriately.
Tuberset not valid	Tube set consumed or not valid	• Insert new tube set.

17 Technical Data Pump

Type designation:	GYN-Pump (PH304)	
Manufacturer information:	W.O.M. WORLD OF MEDICINE GmbH Salzufer 8, 10587 Berlin	
Software version	See user manual/operator instructions to determine software version (chapter 10.2.4).	
Mains voltage range	100-240 V~	
Supply frequency range	50/60 Hz	
Fuse designation	2 x T 3.15 AH, 250 V, UL-recognized	
Internal voltage supply	No	
Power consumption:	Current	Power consumption
Normal operation 100 V/60 Hz	0.85 A	85 VA
Normal operation 240 V/50 Hz	0.50 A	120 VA
Protection class (I, II)	I	
Application part type (B, BF, CF)	BF ¹	
Defibrillator protection (yes/no)	No	
Protection class (IP code)	IP21	
Classification according to Medical Device Directive 93/42/EEC, (I, IIa, IIb, III)	IIb	
Conformity with the following standards: (in the currently valid version)	IEC 60601-1/EN 60601-1 IEC 60601-1-2/EN 60601-1-2	
Operating conditions [°C] [°F], [%], [kPa]	10 to 40 °C/50 to 104 °F 30 to 70 % rel. humidity 70 to 106 kPa air pressure 3000 m max. altitude above sea level for device use	
Possible use with explosive anesthetic gases	This device is not designed for use with flammable anesthetic gases (Class AP) or flammable anesthetic gases with oxygen (Class APG).	
Storage conditions [°C] [°F], [%], [kPa]	5 to +40 °C/41 to 104 °F 5 to 85% rel. humidity 70 to 106 kPa air pressure	
Transport conditions [°C] [°F], [%], [kPa]	-20 to 60 °C/-4 to 140 °F 5 to 90% rel. humidity at 30 °C/86 °F 70 to 106 kPa air pressure	
Max. sound level	< 80 dB(A) (with acoustic signals)	
Maximum negative suction pressure	-60 kPa	
Adjustable values/ranges		
	Max. flow	Hysteroscopy: 50-800 ml/min Laparoscopy: 1.0-2.0 l/min
	Suction power	High: -60 kPa Low: -30 kPa
	Nominal pressure range	Hysteroscopy: 15-150 mmHg Laparoscopy: not adjustable

Measurement range		
	Throughput	0-2.0 l/min
	Pressure	0-300 mmHg
	Deficit/inflow volume	9995 ml
Precision (repeatability)		
	Flow	±10%
	Pressure	±10 mmHg
	Deficit	± 10% (relative to pumped inflow volume)
Accuracy		
	Flow	±10%
	Pressure	± 5% (of final value)
	Deficit/inflow volume	±10% (relative to pumped inflow volume)
Max. dimensions	Width x Height x Depth	380 mm x 148 mm x 388 mm
Weight		6.3 kg
Interfaces:		
	IN/OUT signal for components	1x scale connection (D-sub connector, DE-9) 1x service connection (USB port, USB 2.0)
	Mains power socket	IEC 60320-1 C14
	RFID Transponder Technology	Working frequency: 13.5609 MHz, Transmitting power: -7.51 dBμA/m at 10 m
Essential Performance		<ul style="list-style-type: none"> • Hysteroscopy Pressure build-up in the body cavity, control and measurement, limit value: max. 150 mmHg (normal condition), First error: 200 mmHg for max. 5 seconds • Laparoscopy Generate irrigation flow, 1.0 l/min, 1.5 l/min, 2.0 l/min, tolerance ±10% First error: 0 l/min (no function) or 2.0 l/min (limit due to design constraints)

¹ A tube set is not an application part in terms of the standard. However, it meets all the technical requirements for an application part.

18 Technical Data of Fluid Monitoring Unit

Type designation	Fluid monitoring unit (PS304)	
Software version	Can be viewed via connected GYN-Pump (PH304) (chapter 10.2.4)	
Protection class (I, II)	I	
Enclosure protection (IP code)	IP21	
Tested as per following standards (in the respectively valid version)	EN 60601-1/IEC 60601-1	
Max. load of container scale	16 kg/35.3 lbs, 4 container	
Max. load of bag scale	23 kg/ 50.7 lbs, 4 bags with 5 l each	
Weight	22 kg	
Max. dimensions (Width x Height x Depth)	634 mm x 1159 mm x 634 mm	
Interfaces/ports		
	Signal input/output:	1x data connection (D-Sub-connector, DE-9) for connection to pump.

19 Accessory List

EN

Table 1: Tube sets

Article No.	Packing units	Description
T0505-01	10	Irrigation tube set, disposable (one time use)
T0506-01	1	Irrigation tube set, reusable (NOT FOR SALE IN USA)
T0504-01	10	Vacuum tube set, incl. filter, can be used 30 days
T0502-01	10	Suction tube set (1 connection), disposable (one time use)
T0503-01	10	Suction tube set (2 connections), one time use

Table 2: Additional accessories

Article No.	Packing units	Description
2206000	1	Fluid monitoring unit (PS304)
Z0521-01	10	Pressure chamber membrane
Z0155-01	1	Cannula for fluid bag, reusable (1 sales unit = 1 individual units)
Z0708-01	1	Tube connection Luer Lock male, 8 mm (1 sales unit = 1 individual unit)
Z4101-01	1	Power cable EU, 5 m
Z4176-01	1	Power cable UK, 5 m
Z4102-01	1	Power cable US, 5 m
Z4178-01	1	Power cable CN, 5 m
Z4560-01	1	Power cable BR, 5 m

20 Glossary

Term	Explanation
Transponder	A RFID transponder is a wireless communication and monitoring device that receives and automatically responds to signals. The term transponder is a mixture of "transmitter" and "responder." Transponders can be active or passive.
RFID	<p>RFID (Radio Frequency Identification) is an automated identification process that makes it possible to detect objects without touching them and to log or acquire and save data automatically.</p> <p>A RFID system usually consists of a transponder and a reading device. While the transponder marks objects, the reading device makes it possible to read out data, which are then further processed.</p>
Actual value	The actual value is the measured value of a parameter or unit.
Target or nominal value	The nominal value is the specified or target value to be reached or maintained.
Pressure sensor	The pressure sensor measures the pressure applied to the pressure membrane.
Pressure membrane	The pressure membrane is located in the tube retainer of a tube set and transfers the pressure to the pressure sensor.
EMC	The electromagnetic compatibility describes the ability of a device to function satisfactorily within an electromagnetic environment without adding unacceptable electromagnetic interferences/disturbances to the environment that may cause problems for other devices or equipment located nearby.
ME device	Medical electrical device used for therapy, monitoring or diagnosis of patients, equipped with no more than one connection to a supply network and which necessarily comes into physical or electrical contact with the patient or which transfers energy to or from the patient or which records or captures such energy transfer to or from the patient.
ME system	Medical electrical system consisting of a combination of devices, of which at least one is classified as a medical electrical device and specified by the manufacturer as such and which are joined together by a functional connection or by using a power strip.
Touchscreen	A touchscreen (sometimes also spelled touch screen) makes it possible to control a device directly by touching its screen or display.
HF surgery	High-frequency surgery uses high-frequency alternating current conducted through the human body to cut tissue or achieve coagulation.
Function test	The function test must be carried out before each surgery to ensure the respective device is fully functional.
Safety test	The safety test is part of the annual inspection.
Basic function test	The basic function test checks the basic function of the device and is part of the annual inspection.

21 Appendix

21.1 Test Log

[illegible]

21.2 Return Form

Please fill out this form when returning the device:

Name of owner:

Sales partner:

Address of person returning unit:

Street:

House number:

ZIP/Postal code:

City:

Country:

IMPORTANT!

Serial number (see identification plate):

Device type:

Description of defect:

Contact

Signature

Date

Index

A

Authorized service technician 4
Authorized trained personnel 55

B

Bag change during surgery 44

C

Calibration error 52
Care and maintenance 4
Certification 55
Change fluid bag 30
Connect instrument 26
Container change during surgery 32, 44
Contamination 4
Contraindications 33, 37

D

Defective vacuum pump 32
Deficit amount with no balancing unit connected 40
Delivery inspection 16

E

Electronics error 52
ESD (Electrostatic Discharge) precautionary measures 65
Exclusion of liability 4

F

Federal Law 4
Fluid bag warning 29

G

Grounding contact 17

H

Hysteroscopy tube types 40

I

Improved pressure measuring 41
Indication screen 22
Insert tube set 26
Intended use
 Laparoscopy 33
Invalidating a tube set 23

L

Laparoscopy tube types 34
Liability 55

M

Mains connection 17
Manufacturer's specification 62
Manufacturer's specifications 55
Measured values and tolerances 62
Menu overview 34, 45
Motor error 52

O

Only for U.S. operators 17

P

Potential equalization 17
Precautionary measures 65

R

Remove tube set 27
Removing the tube set 26
Reset to factory settings 47
Returning the device 16
RFID technology (transponder technology) 23

S

Sensor error 52
Setting up 16
Start suction 32
Subject to technical changes 3
T
Technical application scope
 Laparoscopy 33
Technical documents 55
Transponder signal loss 24
Two-year maintenance interval 55
U
Unauthorized personnel 55
W
Waste management 4

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Deze gebruiksaanwijzing bevat auteursrechtelijk beschermde informatie waar copyright op bestaat. Alle rechten voorbehouden. Het is verboden om deze gebruiksaanwijzing zonder uitdrukkelijke, schriftelijke toestemming van W.O.M. WORLD OF MEDICINE GmbH geheel of gedeeltelijk door middel van fotokopieën, microfilm of met andere middelen te vermenigvuldigen of te verspreiden.