PP120


Instructions for Use **en**

Gebruiksaanweisung **de**

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LAP-Pump

en

These instructions for use contain information that is subject to copyright. All rights reserved. These instructions for use may not be photocopied, duplicated on microfilm, or otherwise copied or distributed, completely or in part, without the approval of W.O.M. WORLD OF MEDICINE GmbH.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in these instructions for use are associated with registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or issues concerning these instructions for use.

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Copyright © W.O.M. WORLD OF MEDICINE GmbH

CE0197

CE marking according to Directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
EG-markering conform Richtlijn 93/42/EEG
Маркировка CE согласно директиве 93/42/ЕЭС

	Symbols	Bildzeichen	Symboles	Symbolen	Символ
	Authorized for Sale or use by Physician only	Nur für autorisiertes Vertriebspersonal oder Arzt	Autorisé seulement pour la vente par un personnel autorisé ou l'utilisation par un médecin	Uitsluitend voor bevoegd personeel of arts	Только для врача или авторизованного персонала
	Consult instructions for use	Gebrauchsanweisung beachten	Respecter les instructions d'utilisation	Gebruiksaanwijzing opvolgen	Соблюдать инструкцию по применению
	Follow instructions for use (white image on a blue background)	Gebrauchsanweisung befolgen (weißes Bild auf blauem Grund)	Respecter le mode d'emploi (image blanche sur fond bleu)	Gebruiksaanwijzing opvolgen (witte afbeelding op blauwe ondergrond)	Соблюдать инструкцию по применению (белое изображение на синем фоне)
	Caution	Vorsicht; Achtung	Prudence; Attention	Voorzichtig; Let op	Осторожно; внимание
	Type BF applied part	Gerät des Typs BF	Symbole pour un appareil type BF	Apparaat van het type BF	Прибор типа BF
IP 21	Degrees of protection provided by enclosures (IP-Code)	Gehäuseschutzklasse (IP-Code)	Degrés de protection procurés par les enveloppes (Code IP)	Beschermingsklasse (IP-code) behuizing	Класс защиты корпуса (IP-код)
	Equipotentiality	Potentialausgleich	Connexion équipotentielle	Potentiaalvereffening	Выравнивание потенциалов
	Alternating current	Wechselstrom	Courant alternatif	Wisselstroom	Переменный ток
	Service	Service	Service	Service	Сервис
REF	Catalogue number	Artikelnummer	Numéro d'article	Artikelnummer	Артикул
	Do not reuse	Nicht wiederverwenden	Ne pas réutiliser	Niet voor hergebruik	Не для повторного применения
	Do not resterilize	Nicht resterilisieren	Ne pas restériliser	Niet opnieuw steriliseren	Не стерилизовать
	Contains DEHP	Enthält DEHP	Contient du phtalate de diéthylhexyle (DEHP)	Bevat DEHP	Содержит DEHP
	Not made with phthalates	Enthält keine Phthalate	Ne contient pas de phtalate	Bevat geen ftalaten	Не содержит фталатов
STERILE EO	Sterilized using ethylene oxide	Sterilisiert mit Ethylenoxid	Stérilisé à l'oxyde d'éthylène	Gesteriliseerd met ethyleenoxide	Стерилизовано этиленоксидом
LOT	Batch code	Chargencode	Codes de lot	Chargecode	Код партии

en

de

fr

nl

ru


en

de

fr

nl

ru

	Symbols	Bildzeichen	Symboles	Symbolen	Символ
	Serial number	Seriennummer	Numéro de série	Serienummer	Серийный номер
	Date of manufacture (YYYY-MM-DD)	Herstellungsdatum (JJJ-MM-TT)	Date de fabrication (AAAA-MM-JJ)	Productiedatum (JJJ-MM-DD)	Дата изготовления (ГГГГ-ММ-ДД)
	Use by date (YYYY-MM-DD)	Verwendbar bis (JJJ-MM-TT)	Utilisable jusqu'à (AAAA-MM-JJ)	Uiterste gebruiksdatum (JJJ-MM-DD)	Использовать до (ГГГГ-ММ-ДД)
	Manufacturer	Hersteller	Fabricant	Fabrikant	Изготовитель
	Quantity	Menge	Quantité	Hoeveelheid	Количество
	Number of autoclav- ing cycles	Anzahl der Auto- klavierungen	Nombre d'autoclavages	Aantal auto- claveringen	Число авто- клавирований
	Do not use if package is damaged	Inhalt beschädigter Packung nicht verwen- den	Ne pas utiliser le con- tenu si l'emballage est endommagé	Inhoud van beschadigde verpakking niet gebrui- ken	Не использовать содержимое повреждённой упаковки
	Keep dry	Vor Nässe schützen	Protéger de l'humidité	Beschermen tegen vocht	Беречь от влаги
	Keep away from sun- light	Vor Sonnenlicht schützen	Protéger contre les ray- ons du soleil	Niet aan zonlicht blijven	Беречь от солнечного света
	Protect from heat and radioactive sources	Vor Hitze und radioakti- ver Strahlung schützen	Protéger des sources chaleur et radioactives	Beschermen tegen hitte en radioactieve straling	Оберегать от воздействия тепла и радиоактивного излучения
	Top-Bottom	Oben-Unten	Haut-Bas	Boven-Beneden	Верх-низ
	Fragile	Zerbrechlich	Fragile	Breekbaar	Хрупкое
	Waste management	Entsorgung	Élimination	Verwijdering	Утилизация
	Temperature limit	Temperaturbegrenzung	Limitation de la température	Temperatuurbegrenzing	Ограничение температуры
	Humidity limitation	Luftfeuchtigkeit, Begrenzung	Limitation de l'humidité	Luchtvochtigheid, begrenzing	Влажность воздуха, ограничение
	Atmospheric pressure limitation	Luftdruck, Begrenzung	Limitation de la pression atmosphérique	Luchtdruk, begrenzing	Давление воздуха, ограничение
	Fluid bag	Spülbeutel	Poche de lavage	Vloeistofzak	Промывочный пакет

	Symbols	Bildzeichen	Symboles	Symbolen	Символ
	Non sterile	Nicht sterile	Non stérile	Niet steriel	Нестерильно
	Electrostatic sensitive devices	Elektrostatisch empfindliche Komponente	DES appliquées à la manipulation	Elektrostatisch bedreigd onderdeel	Компонент, чувствительный к статическому электричеству
	Suction	Saugen	Aspirer	Zuigen	Всасывание
	Irrigation	Spülen	Irriguer	Spoelen	Промывание
	Tube status	Schlauchstatus	Statut tubulure	Slangstatus	Состояние шланга
	Max. use 30 days	max. 30 Tage verwenden	Utiliser en 30 jours maximum	max. 30 dagen gebruiken	макс. срок использования 30 дней
	Max. use 24 hours/10 applications	max. 24 Stunden verwenden/10 Applikationen	Utiliser en 24 heures maximum/10 applications	max. 24 uur gebruiken/10 applicaties	макс. срок использования 24 часа/10 применений

en

de

fr

nl

ru

Table of contents

1	Important User Notes	3
2	Safety Information.....	4
3	General Information	5
3.1	Device Description	5
3.2	Intended Use and Contraindications	5
3.2.1	Intended Use.....	5
3.2.2	Contraindications.....	5
3.3	General Warnings and Precautions	5
3.3.1	General Warnings.....	5
3.3.2	Precautions.....	11
4	Initial Device Startup	14
4.1	Scope of Delivery	14
4.2	Setting up and Connecting the Device.....	14
5	Operating the Device	17
5.1	Front of the Device.....	17
5.2	Rear of the Device	17
5.3	Usable Tube Sets Overview	18
5.4	Using the Tube Sets	19
5.5	Inserting a Standard Tube Set.....	20
5.6	Inserting the Suction/Irrigation Set.....	22
5.7	Inserting a Day Patient Tube Set (NOT FOR SALE IN USA AND CHINA).....	23
5.7.1	Replacing the Patient Tube after Surgery.....	26
5.8	Connecting the Irrigation Fluid Bags.....	26
5.9	Irrigation	27
5.10	Suction.....	28
5.11	Replacing the Suction Container	29
5.12	Removing a Vacuum Tube.....	29
5.13	Switching Device Off	29
6	Safety Functions.....	30
6.1	Pressure Limitation at 450 mmHg	30
6.2	Motor Cutout at 600 mmHg	30
6.3	Device Self-Test	30
7	Functional Test	31
7.1	Preparing the Function Test	31
7.2	Testing Irrigation	31
7.3	Suction Test.....	32
8	Using the Device during Surgery	33
9	Care and Maintenance.....	34
9.1	Cleaning the Device	34
9.2	Maintenance Intervals.....	34
9.3	Maintenance by Authorized Service Technician.....	34
9.4	Replacing the Fuse	35
9.5	Cleaning the Device.....	36
9.6	Care of the Reusable Tube Set (NOT FOR SALE IN USA)	36
9.6.1	Cleaning the Reusable Tube Set (NOT FOR SALE IN USA)	36
9.6.2	Disinfecting the Reusable Tube Set (NOT FOR SALE IN USA).....	38
9.6.3	Sterilizing the Reusable Tube Set (NOT FOR SALE IN USA).....	38
10	Annual Inspection	40
10.1	Electrical Safety Test	40
10.2	Irrigation Test	41
10.3	Suction Test.....	41
11	Technical Data	42
12	Accessory List	44
13	Electromagnetic Compatibility.....	45
13.1	Electrical Connections.....	45
13.2	Guidelines and Manufacturer's Statement – Electromagnetic Emissions	45
13.3	Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity	46
13.4	Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity	47
14	Error and Warning Messages.....	49
15	Glossary.....	51
16	Test Log.....	52
16.1	Test Log.....	52
16.2	Return Form	53
	Index.....	54

1 Important User Notes

Read the instructions for use carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing staff or service personnel, or
- to damage or malfunction of device and/or accessories.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of its products.

The words WARNING, CAUTION, and NOTE carry special meanings. Sections marked with these words must be given special attention.

Subject to technical changes

Please note

WARNING!
The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third party.

CAUTION!
These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.

NOTE!
These paragraphs contain information to clarify the instructions or provide additional useful information.

Federal Law (only for U.S. market)

Exclusion of liability

Authorized service technician

Care and maintenance

Contamination

Disposal

2 Safety Information

CAUTION: Federal Law restricts this device to sale by or on the order of a physician.

The manufacturer is not liable for direct or consequential damages, and the warranty becomes null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained;
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedule is not adhered to.

The handing over of technical documents does not constitute authorization to make repairs or alterations to the device or accessories.

WARNING: Modifying the LAP-Pump is not permitted.

Only authorized service technicians are permitted to repair, change settings, or modify the device or associated accessories and equipment. Any violation will void the manufacturer's warranty. Authorized service technicians are only trained and certified by the manufacturer.

The service and maintenance of the device and its accessories has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use.

Maintenance of the device may not be performed during the operation.

Before shipping, decontaminate device and accessories in order to protect the service personnel. Follow the instructions listed in this manual. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to reject contaminated products for repair.

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. For disposal of the device and its accessories, please consult the manufacturer or an authorized disposal company, in compliance with legal or national regulations.

3 General Information

3.1 Device Description

The device **LAP-Pump** is a pump with an irrigation and suction function for laparoscopy.

The device itself is non-invasive and designed for use in multi-unit racks (non-sterile area). It pumps medically sterile irrigation fluids through a sterile tube. These fluids are used to irrigate corresponding body cavities. The device can be used with electrolyte-free media (e.g., glycine 1.5 % and sorbitol 3.0 %) and with isotonic, electrolyte containing media (e.g., saline 0.9 % and Lactated Ringer's). The device is equipped with a vacuum suction function used to suction off secreted fluids. The secretion liquid is sucked through a tube into the container. The device functions only with the tube sets described in the accessory list (see chapter 12 "Accessory List").

The essential performance of the device in laparoscopy mode is to provide a flushing flow. Another essential performance of the device is to provide suction pressure to suction the fluid from the patient.

Essential performance

3.2 Intended Use and Contraindications

3.2.1 Intended Use

The device **LAP-Pump** is a suction and irrigation pump for minimally invasive surgical procedures. It pumps fluid into and aspirates fluid from the abdomen.

3.2.2 Contraindications

Do not use this device **LAP-Pump** for procedures where a specific pressure has to be preset and reached such as a hysteroscopy or an arthroscopy.

3.3 General Warnings and Precautions

3.3.1 General Warnings

WARNING!

The product may be used only by doctors or health care professional in accordance with the intended use. Personnel must have training that qualifies them to operate laparoscopy pumps. In addition, personnel must read the user manual and operating instructions and become familiar with handling the device before first use. Personnel must have mastered handling and operation of the product before first use during a surgery. This user manual does not include descriptions or instructions for surgical procedures/techniques. It is also not intended to introduce to or train physicians in the use of surgical techniques.

WARNING!

Technique and procedures

Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

WARNING!

Check all factory settings

Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

WARNING!**Original accessories**

For your own safety and that of your patient, use only original accessories.

WARNING!**Not explosion-proof**

The device is not explosion-proof. Do not operate the device in the vicinity of explosive anesthetic gases and not in the vicinity of oxygen-enriched environments.

WARNING!**Falls and crashes**

Place the device on a stable and level surface. Cables must be laid safely. Tubes between the device and the patient must not create any obstruction.

WARNING!**Risk of electrical shock**

To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

WARNING!**Risk of electrical shock**

To avoid the risk of electrical shock, only use this device when connected to a properly grounded power supply network.

WARNING!**Professional qualification**

The Instructions for Use do not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualifications working under the direction and supervision of a physician.

WARNING!**Functional test**

The functional test must be performed prior to each surgery.

WARNING!

Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

WARNING!

Unplug the power cable from the device before checking the fuse.

WARNING!**Device defect**

If a device defect is suspected or confirmed, do not use it. Make sure the device can no longer be used until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!****Sterile mediums and accessories**

Always work exclusively with sterile substances and mediums, sterile fluids, and sterile accessories if so indicated.

**WARNING!****Replacing fuse**

Replace the fuse only with a fuse of the same type and rating.

**WARNING!****Maintenance and calibration**

Do not open the device. The device may not be opened or calibrated by the user. Only authorized service technicians are permitted to repair, calibrate, or modify the device or its equipment.

**WARNING!****Condensation / Water penetration**

Protect device from moisture. Do not use if moisture has penetrated the device.

**WARNING!**

The irrigation fluid bags must not be stored on the device or hung up directly above the device.

**WARNING!****Cleaning the device**

Do not sterilize the device.

**WARNING!**

The device is only intended for use with flexible fluid bags. Do not use glass containers as they might break. Fluid cannot flow quickly enough due to the vacuum being generated inside of the bottle. Risk of implosion.

**WARNING!****Keep full fluid bags ready for use**

Always keep a full fluid bag on hand to replace an empty one. This avoids having to interrupt surgery due to a lack of irrigation fluid.

**WARNING!****Contamination**

Do not use device or vacuum tube if signs of contamination are detected. Fluid that has penetrated can be detected, for example, by dirt encrustations on the connections for the suction or at the exit of the vacuum pump. Make sure the device or/and accessories can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!**

Different default settings of the warning message for identical or similar devices in the operating room may cause a risk due to conflicting acoustic signals.

**WARNING!**

Sterilize reusable instruments and tubing before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

**WARNING!****Reprocessing of sterile disposable products**

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

**WARNING!**

The vacuum tube sets for this device contain diethylhexylphthalate (DEHP), which is classified as toxic to reproduction according to the EU Directive 1272/2008/EEC on Classification, Labeling and Packaging of Dangerous Substances. DEHP may impair fertility and may cause harm to the unborn child. Therefore, this product must not be used for unauthorized applications. When applied within the intended use, the potential risk to pregnant or breastfeeding women as well as to children resulting from the DEHP contained in this product is not critical.

**WARNING!**

The reusable tube set (NOT FOR SALE IN USA) is made of silicone, polysulfone (PSU), and stainless steel. For cleaning and disinfection, use only pH-neutral or slightly alkaline cleaning agents (e.g. neodisher MediClean 2.0 %), disinfectants (e.g. Lysetol V 8 %), and drying and rinsing agents approved for the tube set materials.

Using unsuitable agents (e.g. neodisher MediKlar rinsing agent) may damage the tube system and especially the PSU connectors.

**WARNING!**

Please check the reusable tube (NOT FOR SALE IN USA) for signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation.

WARNING!

The indicated drying time depends on several variables, including the following: Altitude, humidity, type of packaging, preconditioning, size of chamber, mass of load, and placement in chamber. Users must verify that the drying time set in their autoclave yields dry surgical equipment when using the method of saturated steam sterilization described here.

**WARNING!**

Remove the membrane before cleaning.

**WARNING!**

Make sure to connect tubes as indicated by the matching connectors.

**WARNING!**

The day tube remains connected to the device after each surgery. The patient tube must be disposed of immediately after surgery. A new sterile cap (included with the patient tube) must be screwed onto the hygiene connector at the day tube after each surgery. This sterile protector remains connected to the day tube until the next surgery.

**WARNING!**

The device is only intended for use with flexible fluid containers. Do not use glass containers as they might break. Fluid cannot flow quickly enough due to the vacuum being generated inside of the bottle. Risk of implosion.

**WARNING!**

Select an irrigation fluid suitable for the subsequent medical procedure.

**WARNING!**

Do not touch patient and the ON/OFF switch at the same time.

**WARNING!**

ME System (Medical Electrical System)

Use only parts and/or devices from ME systems (see Chapter 15 Glossary) in patient environments in compliance with the standard IEC60601-1 in the respective currently valid version.

**WARNING!**

ME System (Medical Electrical System)

Connect only items that have been specified as part of the ME system or specified as being compatible with ME system.

WARNING!

Make sure the tubes for the suction-irrigation instrument do not have any kinks when installed. Make sure the tubes are not crimped or clamped off.

WARNING!**Filter**

The vacuum tube with integrated filter is designed for max. 30 days. Replace the vacuum tube if it is obviously contaminated. The filter prevents that body fluids enter the interior of the device. Please note that the filter may reduce the suction capacity.

WARNING!**Do not touch**

Avoid touching the roller wheel. Risk of injury!

WARNING!**Serum sodium concentration**

It is also necessary to monitor the concentration of sodium in the blood of the patient to prevent electrolyte disturbances. Monitoring of the concentration of sodium in the blood must be performed by the physician and is not performed or supported by the system.

WARNING!**Air Embolisms**

An air embolism can be the result of air contained in the tube set or connected instrument reaching the patient. Make sure that there is always fluid in the bag to prevent air from being aspirated into the patient.

WARNING!

Position the device in such a way that it is easy to operate and switch off.

WARNING!

An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

3.3.2 Precautions

CAUTION!

Electrical Interference

(See chapter 13 "Electromagnetic Compatibility"). Electrical interference with other devices or instruments was practically eliminated when developing this device and none was detected during testing. However, if you still detect or suspect such interference, please follow these suggestions:

- Move this, the other or both devices to a different location
- Increase distance between used devices
- Consult an electro-medical expert

CAUTION!

Endoscope

The device may only be connected with endoscopes designed for and featuring the technical specification permitting such a combined use. Any utilized endoscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600. Combining/connecting with other devices generates a medical electrical system (MES). The system configurator is responsible for compliance with the standard IEC 60601-1 / EN 60601-1 in its latest version.

CAUTION!

Mains voltage

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

CAUTION!

Indoor climate

Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!

Power cables

Any power cables employed by the user that are not provided by the manufacturer must meet the safety requirements of the national standards in the respective current valid version.

CAUTION!

Ventilation of the device

Avoid device overheating. Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).

CAUTION!

The pressure chamber diaphragm is fragile and can be replaced if damaged. Replacement diaphragms are supplied with reusable tubes. The reusable tube sets must be prepared before each use. See Chapter 9.6 Care of the Reusable Tube Set (NOT FOR SALE IN USA) for further instructions.

CAUTION!

Position the device in such a way that it is easy to operate and switch off.

CAUTION!

The indicated drying time depends on several variables, including the following: Altitude, humidity, type of packaging, preconditioning, size of chamber, mass of load, and placement in chamber. Users must verify that the drying time set in their autoclave yields dry surgical equipment when using the method of saturated steam sterilization described here.

CAUTION!

The day tube must be replaced at the end of the day of surgery but not later than after 24 hrs.

CAUTION!

Full suction containers must be replaced immediately without stopping surgery.

CAUTION!

The full suction output/volume is available only, if the suction system is pre-evacuated (vacuum generation). Pre-evacuation takes about 20 to 60 seconds depending on the volume of the secretion container.

CAUTION!

The suction system has a leak if the vacuum pump does not stop. Check the suction system for leaks and repeat function test. Do not use device, if vacuum pump is not operating or not operating correctly.

CAUTION!

Peripheral devices

Additional peripheral equipment connected to interfaces of the medical monitor has to meet the requirements of the following specifications in the respective current valid version: IEC 60601-2-18 / EN 60601-2-18 for endoscopic devices and IEC 60601-1 / EN 60601-1 for electrical medical devices. All configurations have to comply with IEC 60601-1 / EN 60601-1 specifications. Whoever connects additional equipment to signal output or signal input is considered the system configurator and as such is responsible for complying with requirements of the standard IEC 60601-1 / EN 60601-1.

CAUTION!

ME Device in Rack

The ME device may not be used when stacked or stored directly adjacent to or with other devices. If it should be necessary to operate the device close to stacked with other devices, the ME device or ME system should be monitored to ensure it works properly as configured.

CAUTION!**Use of other accessories, other transducers and cables**

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM (see Chapter 15 "Glossary") as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

Delivery inspection

Returning the device

4 Initial Device Startup

WARNING!
The product may be used only by doctors or health care professional in accordance with the intended use. Personnel must have training that qualifies them to operate laparoscopy pumps. In addition, personnel must read the user manual and operating instructions and become familiar with handling the device before first use. Personnel must have mastered handling and operation of the product before first use during a surgery. This user manual does not include descriptions or instructions for surgical procedures/techniques. It is also not intended to introduce to or train physicians in the use of surgical techniques.

4.1 Scope of Delivery

- Device LAP-Pump
- Instructions for Use
- Power cord

Always check all parts and optional accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized service company.

If it becomes necessary to return the device, always use the original packaging. The manufacturer does not take responsibility for damage that has occurred during transportation if the damage was caused by inadequate transport packaging.

When returning the device please make sure that all the required information has been supplied:

- Name of owner
- Address of owner
- Device type and model
- Serial number of the equipment (see identification plate)
- Detailed description of defect

Use the return shipping form included in Chapter 16.2 Return Form.

4.2 Setting up and Connecting the Device

NOTE!
Locations
The device may be used only in the hospital or operating rooms.

Place the device on a flat surface free of vibration located in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in Chapter 11 "Technical Data".

WARNING!
ME System (Medical Electrical System)
Use only parts and/or devices from ME systems (see Chapter 15 Glossary) in patient environments in compliance with the standard IEC60601-1 in the respective currently valid version.

Setting up

CAUTION!

ME Device in Rack

The ME device may not be used when stacked or stored directly adjacent to or with other devices. If it should be necessary to operate the device close to stacked with other devices, the ME device or ME system should be monitored to ensure it works properly as configured.

CAUTION!

Equipment should be positioned such that power cord can be easily disconnected.

CAUTION!

Ventilation of the device

Avoid device overheating. Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).

CAUTION!

Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!

Position the device in such a way that it is easy to operate and switch off.

CAUTION!

Place the device outside the sterile field.

CAUTION!

Position of the user

To avoid a malfunction, the user must be positioned correctly

- within a display viewing angle of $\pm 50^\circ$ to operate the device
- up to 2 m/6.5 ft from the device front for monitoring the actual values

Mains connection

CAUTION!

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

Make sure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The mains connection cable may be plugged only into a properly installed, grounded safety wall socket (shock-proof socket) (see DIN VDE 0100-710). Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

The power connection must be equipped with a grounding contact. Use the original power cable (if included in scope of delivery) to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.

Grounding contact

Initial Device Startup

en

Only for U.S. operators

Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 and IEC 16320-C13. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

Potential equalization

Integrate the device into the potential equalization system as specified by local safety rules and regulations.

5 Operating the Device

5.1 Front of the Device

Familiarize yourself with the control and function elements at the front of the device.


Fig. 5-1 Device front

- ① Irrigation start/stop key
- ② Irrigation status LED
- ③ Suction start/stop key
- ④ Suction status LED
- ⑤ Roller wheel
- ⑥ Tube retainer
- ⑦ Pressure sensor
- ⑧ Connection for hygiene filter and vacuum tube
- ⑨ LED tube status green
- ⑩ LED tube status red
- ⑪ ON/OFF power switch
- ⑫ Power switch status LED

5.2 Rear of the Device

Familiarize yourself with the connection elements at the rear of the device.


Fig. 5-2 Device rear

- ⑬ Device labels
- ⑭ Vacuum outlet
- ⑮ Connection for potential equalization
- ⑯ Fuse holder
- ⑰ Mains socket
- ⑱ Service interface

5.3 Usable Tube Sets Overview

The following table lists the function of each type of tube set depending on the selected indication.

Tube type	Item No.	Tube category
Disposable (one-time use) (NOT FOR SALE IN USA AND CHINA)	T0500-01	Irrigation day patient tube set
Disposable (one-time use) (NOT FOR SALE IN USA OR CHINA)	T0501-01	Irrigation day tube
Disposable (one-time use)	T0502-01	Suction tube set (1 connection)
Can be used 30 days	T0504-01	Vacuum tube set, incl. filter
Disposable (one-time use)	T0505-01	Irrigation tube set
Reusable (NOT FOR SALE IN USA)	T0506-01	Irrigation tube set
Disposable (one-time use)	T0508-01	High flow suction and irrigation set with handle

Table 5-1

CAUTION!
The pressure chamber diaphragm is fragile and can be replaced if damaged. Replacement diaphragms are supplied with reusable tubes. The reusable tube sets must be prepared before each use. See Chapter 9.6 Care of the Reusable Tube Set **(NOT FOR SALE IN USA)** for further instructions.

WARNING!
The vacuum tube sets for this device contain diethylhexylphthalate (DEHP), which is classified as toxic to reproduction according to the EU Directive 1272/2008/EEC on Classification, Labeling and Packaging of Dangerous Substances. DEHP may impair fertility and may cause harm to the unborn child. Therefore, this product must not be used for unauthorized applications. When applied within the intended use, the potential risk to pregnant or breastfeeding women as well as to children resulting from the DEHP contained in this product is not critical.

WARNING!
Sterilize reusable instruments and tubing before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

WARNING!
Reprocessing of sterile disposable products
Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

NOTE!
Tubes must be stored at room temperature. The shelf life for all tubes is 5 years.

NOTE!

Observe applicable hygiene regulations when disposing of the tube set.

en

5.4 Using the Tube Sets

Transponder technology detects the tube type, the validity and reliability of a tube set automatically. This eliminates virtually all "operating errors" since non fitting, invalid, and not approved tube sets are reliably detected. In any approved tube set (see chapter 12 Accessory List), the corresponding transponder is located under the tube holder.

If an approved tube set is located in the device, the transponder technology automatically invalidates this tube set when the device is started. The tube status LED lights up green. Depending on the tube set, this has the following consequences:

- **Disposable tube set:** After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If the irrigation cycle is stopped, the device can be restarted within 30 minutes. If the device is switched off or in case of a power failure, the tube set is invalidated. The irrigation cycle can no longer be started. The tube status LED lights up red. If this is the case, you must insert a new, valid, and approved tube set.
- **Reusable tube set (NOT FOR SALE IN USA):** The reusable tube set can be used about 20 times. After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If irrigation is stopped during a use cycle, it can be restarted within 30 minutes without the set being invalidated. If the device is switched off or in case of a power failure, the current use cycle is invalidated. If an invalid reusable tube set has been inserted, a warning signal (beep) is emitted. A red LED will be displayed for 8 seconds.
- **Day patient tube set also day tube set (NOT FOR SALE IN USA AND CHINA):** The following description refers only to the day tube (for the patient tube, please see chapter 5.7).
The reusable day tube set can be used for a total of 10 times. After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If irrigation is stopped, it can be restarted within 30 minutes without the set being invalidated.
After the device is switched off or in the case of a power failure, the current use cycle is invalidated. Irrigation can no longer be started following the last use cycle. The tube status LED lights up red. If this is the case, you must insert a new, valid, and approved day tube.

The device stops and cannot be restarted if the transponder loses its signal during current use. If the signal is restored within 20 seconds, it is possible to continue using the tube set.

RFID technology (transponder technology)


Invalidating a tube set

Transponder signal loss

5.5 Inserting a Standard Tube Set

Fig. 5-3 Tube set elements

- ① Protective caps
- ② Tap spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Ring
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑩ Luer lock connector


The **standard tube set** is available as a disposable or reusable tube set (autoclavable) (see chapter 12 Accessory List).

The tube set consists of 3 tube components, (irrigation tube (5), roller tube (7), and instrument tube (9)), a Y-connector (4) and 2 tap spikes (2). The tap spikes (2) are used to connect the tube components with the fluid bags.

The Luer lock connector (10) connects the instrument tube with the instrument.

Fig. 5-4 Inserting the tube set

- ② Tap spikes
- ⑤ Irrigation tube
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑪ Fluid bag
- ⑫ Tube retainer
- ⑬ Roller wheel

WARNING!
Do not touch patient and the ON/OFF switch at the same time.

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers

refer to the Fig. 5-3 "Tube set elements" and Fig. 5-4 "Inserting the tube set".

1. To be carried out by non-sterile technician:

- ▶ If you are using a reusable tube set, open the sterile autoclavable container of the tube set and have a sterile technician remove the inner tube set.
- ▶ If you are using a disposable tube set, open the outer packaging of the tube set and have a sterile technician remove and open the inner tube set packaging.

2. To be carried out by sterile technician:

- ▶ Keep the Luer lock connector (10) in the sterile area and hand the tube end with the tap spikes (2) to the non-sterile technician.
- ▶ Connect the Luer lock connector (10) with the instrument (e.g. inflow cannula). Open inflow valve at instrument.

3. To be carried out by non-sterile technician:

- ▶ Before inserting, close both tube clamps (3) (Fig. 5-3 Tube set elements) on the bags.
 - ▶ Switch device on. The power status LED lights up. Make sure all tube sets are removed from the pump before switching the device on. A long beep is emitted to signal the successful completion of the device test and the tube can now be inserted.
- Inserting the roller tube is depicted in Fig. 5-5 "Inserting the roller tube set".
- ▶ When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (8) only if not pressurized.
 - ▶ Insert the unpressurized pressure chamber (8) carefully into the right notch of the tube retainer (12) up to the stop.
 - ▶ Place the roller tube (7) around the roller wheel.
 - ▶ Insert the tube adapter into the left section of the tube retainer by tightening the tube.
 - ▶ Make sure the tube adapter is firmly snapped into place! Irrigation and drain tube must point away straight towards the front from the tube retainer.

Connect an instrument and the corresponding bags (see also chapter 5.8 "Connecting the Irrigation Fluid Bags").

Removing the tube set

Connect instrument

Insert tube set



Fig. 5-5 Inserting the roller tube set

5.6 Inserting the Suction/Irrigation Set

Fig. 5-6 Elements of the suction/irrigation tube set

- ① Protective caps
- ② Tap spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Roller tube
- ⑦ Pressure chamber with membrane and transponder
- ⑧ Instrument tube
- ⑨ Suction valve (blue)
- ⑩ Irrigation valve (white)
- ⑪ Suction and irrigation pipe
- ⑫ Trumpet valve
- ⑬ Suction and irrigation instrument
- ⑭ Suction tube
- ⑮ Connection to the suction container


The suction/irrigation tube set is available as a disposable tube set for one-time use (see Chapter 12 Accessory List).

The suction/irrigation hose set consists of ⑤, roll tubing ⑥, instrument tube ⑧, suction and rinsing instrument ⑬ and suction tube ⑭, as well as a Y connector ④ and 2 piercing thorns ② tube. The tap spikes are used to connect the tube components with the fluid bags.

Fig. 5-7 Inserting the suction/irrigation set

- ② Tap spikes
- ③ Tube clamps
- ⑤ Irrigation tube
- ⑥ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑬ Suction and irrigation instrument
- ⑮ Fluid bag
- ⑯ Tube retainer
- ⑰ Suction container

WARNING!
Do not touch patient and the ON/OFF switch at the same time.

In order to separate sterile from non-sterile areas, assign the following tasks to

the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to the Fig. 5-6 "Elements of the suction/irrigation tube set" and Fig. 5-7 "Inserting the suction/irrigation set".

1. To be carried out by non-sterile technician:

- ▶ Open the outer packaging of the tube set and let you remove the inner tube packaging by sterile personnel and open.

2. To be carried out by sterile technician:

- ▶ Keep the instrument in the sterile area and hand the tube end with the tap spikes (2) to the non-sterile technician.
- ▶ Open the trumpet valve (12) of the suction and rinsing instrument.

3. To be carried out by non-sterile technician:

- ▶ Before inserting, close both tube clamps (3) (Fig. 5-6 Elements of the suction/irrigation tube set) on the bags.
- ▶ Switch device on. The power status LED on the ON switch lights up. Make sure all tube sets are removed from the pump before switching the device on. A long beep is emitted to signal the successful completion of the device test and the tube can now be inserted.
- ▶ Inserting the roller tube is depicted in Fig. 5-8 "Inserting the roller tube set".
- ▶ When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (8) only if not pressurized.
- ▶ Insert the unpressurized pressure chamber (8) carefully into the right notch of the tube retainer (16) up to the stop.
- ▶ Place the roller tube (7) around the roller wheel.
- ▶ Insert the tube adapter into the left section of the tube retainer by tightening the tube.
- ▶ Make sure the tube adapter is firmly snapped into place! Irrigation and drain tube must point away straight towards the front from the tube retainer.

Connect the suction tube to the suction tank (17) (see Fig. 5-7 "Inserting the suction/irrigation set").

Remove the tube set

Connect the instrument

Insert tube set


Fig. 5-8 Inserting the roller tube set

5.7 Inserting a Day Patient Tube Set (NOT FOR SALE IN USA AND CHINA)

The **day patient tube set** consists of two tube sections comprised of the **day tube** and **patient tube**. The day tube establishes the tube connection between the irrigation fluid bags and the patient tube. The day tube is inserted into the tube retainer of the device and not replaced for the duration of one day between surgeries and/or surgical procedures. Only the patient tube (tube connection between day tube and the instrument) must be replaced between surgical procedures. Efficiency is increased with a day patient tube set if several surgeries are to be performed on one day. Make sure to connect tubes as indicated by the matching connectors.

CAUTION!


Use only original day patient tube sets of the manufacturer. Other day patient tube set are not compatible with the pump!

en

Fig. 5-9 Day tube set elements

- ① Protective caps
- ② Tap spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Ring
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑩ Hygiene connector
- ⑪ Protective cap

The day tube

The day tube includes the connectors to the fluid bags and to the hygiene connector (Fig. 5-10 "Patient tube elements" see (12)) of the patient tube.

The day tube may be used the entire day of surgery. However, it is designed for a maximum of **10 uses**.

CAUTION!

The day tube must be replaced at the end of the day of surgery but not later than after 24 hrs.


A new sterile cap (included with the patient tube in the delivery) must be screwed onto the hygiene connector (Fig. 5-9 (10)) of the day set immediately after the surgical procedure to protect against contamination. This sterile protector remains connected to the day tube until the next surgery.

The patient tube

The patient tube has a connector on one side to connect with the day tube (hygiene connector (12)). Luer lock connector (14) to connect to an instrument such as a trocar, for example, is fitted to the other side. . Connect the hygiene connector (12) on the patient tube with the matching piece on the day tube (Fig. 5-9 (10)).

Fig. 5-10 Patient tube elements

- ⑪ Protective cap for hygiene connector
- ⑫ Hygiene connector
- ⑬ PVC tube, 285 cm
- ⑭ Luer lock connector with screw cap
- ⑮ Protective cap for Luer lock

The patient tube set has to be replaced for every patient to prevent cross-contamination. The day tube remains connected to the device after each surgery. The used patient tube must be disposed of after each surgery.

WARNING!

Make sure to connect tubes as indicated by the matching connectors.

en


Fig. 5-11 Inserting the day patient tube set

- (2) Tap spikes
- (3) Tube clamps
- (5) Irrigation tube
- (8) Pressure chamber with membrane and transponder
- (9) Instrument tube
- (10) Hygiene connector
- (16) Fluid bag
- (17) Tube retainer
- (18) Roller wheel

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to Fig. 5-9 Day tube set elements, Fig. 5-10 Patient tube elements and Fig. 5-11 Inserting the day patient tube set.

1. To be carried out by non-sterile technician:

- Open the packaging of the day tube and keep it with you.
- Open the packaging of the patient tube and have a sterile technician remove it for you.

2. To be carried out by sterile technician:


- Keep the Luer lock connector (14) in the sterile area and hand the tube end with the hygiene connector (12) to the non-sterile technician.
- Connect the Luer lock connector (14) with the instrument (e.g., inflow cannula). Open inflow valve at instrument.

Remove day patient tube set

Connect patient tube

Insert new day tube

Fig. 5-12 Inserting the roller tube set

3. To be carried out by non-sterile technician:

- ▶ Connect the hygiene connector (12) of the patient tube (Fig. 5-10 "Patient tube elements") with the hygiene connector (10) of the day tube (Fig. 5-9 "Day tube set elements").
- ▶ Switch device on. The power status LED on the ON switch lights up.
- ▶ 1. Insert the unpressurized pressure chamber (8) carefully into the lower notch of the tube retainer (17) (Fig. 5-11 "Inserting the day patient tube set") up to the stop.
- ▶ When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber only if not pressurized.
- ▶ 2. Place the roller tube (7) around the roller wheel (18). Pull the tube until the tube ring fits in the upper part of the tube retainer.
- ▶ The validity of the tube set is now checked. When an invalid tube set is inserted, the LED showing the tube status turns red. A short acoustic signal is emitted 3x. The entire procedure must be repeated with a valid tube set if an invalid tube set has been inserted.
- ▶ Connect the tube end with the tap spikes to the full fluid bags.
- ▶ Press the Irrigation Start/Stop key and wait until the day patient tube set is completely filled with irrigation fluid.
- ▶ 3. The device is now ready for use.

5.7.1 Replacing the Patient Tube after Surgery

WARNING!

The day tube remains connected to the device after each surgery. The patient tube must be disposed of immediately after surgery. A new sterile cap (included with the patient tube) must be screwed onto the hygiene connector at the day tube after each surgery. This sterile protector remains connected to the day tube until the next surgery.

1. To be carried out by non-sterile technician

- ▶ Open packaging of the patient tube.
- ▶ A sterile technician then removes the patient tube from the inside of the package.

2. To be carried out by sterile technician

- ▶ Keep the luer lock connector (14) and hand the other end of the patient tube to the non-sterile technician.
- ▶ Connect the Luer lock connector (14) on the patient tube with the instrument (e.g., inflow cannula). B. .

3. To be carried out by non-sterile technician

- ▶ Remove the sterile protector of the day tube and connect the new patient tube to the day set without delay.

5.8 Connecting the Irrigation Fluid Bags

WARNING!

The device is only intended for use with flexible fluid containers. Do not use glass containers as they might break. Fluid cannot flow quickly enough due to the vacuum being generated inside of the bottle. Risk of implosion.

WARNING!

The irrigation fluid bags must not be stored on the device or hung up directly above the device.


Fig. 5-13 Hanging the irrigation fluid bags

- ① Irrigation fluid bags
- ② Tap spikes
- ③ Tube clamps
- ④ Irrigation tube

The fluid bags must be suspended at a height of **h** between 0.7 to 1.0 m above the pump (Fig. 5-13 "Hanging the irrigation fluid bags") (pump inlet to the middle of the fluid filling in the bag).

1. The inflow tube can take irrigation fluid from two fluid bags ①. Close both clamps ③ at the branches of the inflow tube.
2. Grasp the tap spike ② at the provided handle when connecting or disconnecting.
3. Comply with sterile precautions when inserting the tap spike into the bag.
4. Open at least one tube clamp ③ at the irrigation tube.

5.9 Irrigation

The device is designed to work with an instrument featuring a trumpet valve or stopcock. Make sure tube set and instrument are completely filled with fluid. Remove all air from tube and instrument.

1. Insert the tube set (see chapter 5.5 "Inserting a Standard Tube Set").
2. Connect the irrigation fluid bag (see chapter 5.8 "Connecting the Irrigation Fluid Bags").
3. Connect the instrument tube to the instrument.
Close stopcock if using instrument with stopcock.
4. Press the irrigation start/stop key (LED lights up).
An admission pressure of 450 mmHg is generated in the tube. This admission pressure is necessary to provide the full irrigation volume/output.
5. Start the irrigation process:
Actuate the trumpet valve or open the stopcock.
6. Stop the irrigation process:
Release the trumpet valve or open the stopcock.
The admission pressure begins to be rebuilt so that the full irrigation volume/output is available for the next irrigation procedure.
7. Press the irrigation start/stop key again (LED goes out).
The tube is now relieved of pressure. The roller wheel spins backwards

5.10 Suction

Fig. 5-14 Suction system with suction-/ Irrigation instrument

- ① Vacuum tube set with hygiene fil-
ter
- ② Suction container
- ③ Suction tube
- ④ Irrigation tube
- ⑤ Instrument



Fig. 5-15 Suction system with standard instrument


- ① Vacuum tube set with hygiene fil-
ter
- ② Suction container
- ③ Suction tube
- ④ Standard instrument

The suction system consists of the following:

- Vacuum tube
- Hygiene filter
- Suction container
- Drain tube
- Standard instrument or suction/irrigation instrument

1./2.

1. Connect device with the suction container (using vacuum tube and hygiene filter). The vacuum tube set with hygiene filter can remain connected to the device for max. 30 days.
2. Connect suction container with instrument (using drain tube).
Close stopcock if using instrument with stopcock.
3. Press the suction start/stop key (LED lights up).
The device generates a negative pressure inside of the suction system. The vacuum pump stops when a negative pressure of max. -60 kpa (-0.60 bar) has been reached.
The vacuum pump automatically reactivates, if the min. negative pressure is reached.
4. Starting suction:
Actuate the trumpet valve or open the stopcock.
5. Stopping suction:
Release the trumpet valve or open the stopcock.
6. Press the suction start/stop key again (LED goes out).

CAUTION!

The full suction output/volume is available only, if the suction system is pre-evacuated (vacuum generation). Pre-evacuation takes about 20 to 60 seconds depending on the volume of the secretion container.

NOTE!

Only use suction containers with overflow protection.

5.11 Replacing the Suction Container

CAUTION!

Full suction containers must be replaced immediately without stopping surgery.

1. Stop the suction (see chapter 5.11 Replacing the Suction Container).
Press the suction start/stop key (LED goes out).
2. Actuate the trumpet valve or Open the stopcock to vent the system.
3. Remove tubes from filled suction container.
4. Reattach tubes to empty suction container.
Close stopcock if using instrument with stopcock.
5. Press the suction start/stop key (LED lights up).
6. Starting suction:
Actuate the trumpet valve or open the stopcock.

5.12 Removing a Vacuum Tube

1. Remove the drain tube (4) from the instrument (5) (see Fig. 5-14 Suction system with suction-/Irrigation instrument or Fig. 5-15 Suction system with standard instrument).
2. Remove the suction container (3) from the drain tube (4).
3. Remove the vacuum tube (1) from the suction container (3).
4. Remove the hygiene filter from the vacuum tube (1).
5. Remove the vacuum tube from the device.

5.13 Switching Device Off

Ending operation of the pump safely:

1. Press the ON/OFF switch (see Fig. 5-1 "Device front" (11))
2. Disconnect device from power supply. For this, pull the mains plug from the mains socket (see Fig. 5-2 "Device rear" (17)).

Pressure limitation when using a day patient set

6 Safety Functions

The electronic components continuously monitor the proper function of the device. Device malfunctions are indicated with audible warning signals (beeps), error messages, and/or by blocking device functions.

6.1 Pressure Limitation at 450 mmHg

The pump generates a max. pressure of 450 mmHg to reach the max. flow capacity during the irrigation mode. If external factors, for example, cause the pressure to exceed the pressure limits, the roller wheel reverses in an attempt to reduce the pressure.

The wheel does not rotate backward to reduce the pressure if a day patient set is used. If the actual pressure exceeds 450 mmHg, the motor stops and irrigation is interrupted to prevent contamination.

6.2 Motor Cutout at 600 mmHg

The roller wheel is stopped if the pressure exceeds 600 mmHg for more than 5 seconds. The suction pump is also stopped in this state. Further irrigation is not possible. The roller wheel and suction pump automatically restart when the pressure has fallen below 600 mmHg.

6.3 Device Self-Test

After being switched on, the device performs a self-test of the sensors, the motor, and electronic components. All LEDs light up. A long beep is emitted after the successful device test.

A table listing a summary of possible error and warning messages is provided in Chapter 14 Error and Warning Messages.

7 Functional Test

WARNING!

Functional test

The functional test must be performed prior to each surgery.

WARNING!

Sterilize reusable instruments and tubing before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

WARNING!

Original accessories

For your own safety and that of your patient, use only original accessories.

7.1 Preparing the Function Test

1. The device is switched off and without an inserted tube set.
2. Switch device on using the power switch.
3. The device carries out a self-check. A long beep is emitted after the successful device check.
4. Insert tube set (see Fig. 5-4 "Inserting the tube set"). Connect the suction system (see Fig. 5-7 "Inserting the suction/irrigation set").

7.2 Testing Irrigation


Fig. 7-1 Testing irrigation

1. Connect the instrument.
Close stopcock if using instrument with stopcock.
2. Press the irrigation start/stop key. The irrigation LED lights up.
3. Actuate the trumpet valve or open the stopcock. Make sure tube set and instrument are completely filled with fluid. Remove all air from tube and instrument.
4. Release the trumpet valve or Close the valve on the instrument.

7.3 Suction Test

1. Close stopcock if using instrument with stopcock.
2. Press the suction start/stop key (LED lights up).
The device generates a negative pressure inside of the suction system. The vacuum pump is deactivated once a low pressure of max. – 60 kpa is reached.

CAUTION!

The suction system has a leak if the vacuum pump does not stop. Check the suction system for leaks and repeat function test. Do not use device, if vacuum pump is not operating or not operating correctly.

3. Open the instrument. Suction starts. The vacuum pump is automatically re-activated if the min. negative pressure is exceeded.
4. Press the suction start/stop key again (LED goes out).

The function test is complete. The device is tested and ready for surgery.

WARNING!

Device defect

Do not use this device if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage of the power plug and power cable.

8 Using the Device during Surgery

WARNING!

Functional test

The functional test must be performed prior to each surgery.

WARNING!

Select an irrigation fluid suitable for the subsequent medical procedure.

1. Tube set and instrument are filled with irrigation fluid. The suction system is connected. The stopcocks are closed.
2. Press the suction start/stop key (LED lights up).
The suction system is pre-evacuated.
Only then can the suction procedure be started with the instrument.
1. Press the irrigation start/stop key (LED lights up).
2. Insert instrument into abdomen.
Start the irrigation process:
Actuate the trumpet valve or open the stopcock.
3. Starting the suction process:
Actuate the trumpet valve or open the stopcock.

Before surgery

During Surgery

After surgery

NOTE!

The manufacturer recommends emptying the tube set using the pump when employing saccharine irrigation fluids (HF use). Remove irrigation fluid bags and instrument first.

1. Press the irrigation start/stop key and the suction start/stop key (LED goes out).
2. Use the power switch to turn device off.
3. Remove the tube sets.

NOTE!

Comply with hygiene rules and regulations when disposing of the tube set, collected fluid, and the suction container.

9 Care and Maintenance

Special care is necessary when servicing, maintaining, and storing the device and its accessories to maintain the functionality of the device and its accessories.

9.1 Cleaning the Device

- 1. Use the **ON/OFF** button to turn the device off.
- 2. Remove the power cable.
- 3. Wipe the surface of the device with a soft cloth moistened with the surface disinfectant based on alcohol (for example Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the device.

NOTE!
Do not sterilize the device.

9.2 Maintenance Intervals

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. The device needs to be serviced at least every two years. The tests are described in chapter 10 "Annual Inspection".

Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

NOTE!
Service or maintenance work may not be carried out during surgery.

9.3 Maintenance by Authorized Service Technician

An authorized service technician has to inspect and service the device at appropriate intervals to ensure the safety and functionality of the unit. The minimum service interval is two years, depending on frequency and duration of use. If the service interval is not maintained, the manufacturer does not assume any liability for the functional safety of the device.

A sticker located on the rear panel of the device will remind you of the latest date for the next service or maintenance check.

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians.

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

Ask the service technician for a certificate after he or she has inspected the unit or performed any service tasks. This certificate lists the type and scope of the service as well as the date and name of the servicing company together with the signature of the service technician.

Manufacturer's specifications

Two-year maintenance interval

Authorized trained personnel

Unauthorized personnel

Liability

Technical documents

Certification

9.4 Replacing the Fuse

WARNING!

Replacing fuse

Replace the fuse only with a fuse of the same type and rating.

CAUTION!

Before replacing the fuse, check the values of the fuse to be inserted acc. to Chapter 11 Technical Data.

The fuse may be defective and is in need of replacement if:

- displays and LEDs (if available on your equipment) do not light up,
- the device does not function.

Check to make sure

- the main power supply cable is properly connected to the power supply input and to a safety socket,
- the house power supply fuse is functioning.

WARNING!


Unplug the power cable from the device before checking the fuse.

The device does **not** have to be opened to replace the fuse.

1. Switch device off.
2. Disconnect device from power supply by pulling the mains plug from the mains socket.
3. Remove power connection cable from mains socket.
4. The fuse holder is located next to the mains socket.
5. Remove fuse holder as depicted in Fig. 9-1 "Opening the fuse holder".
6. **A** Undo the latch of the fuse holder with a small screwdriver.
7. **B** Remove the fuse holder.
8. **C** Check fuse.
9. Insert a new fuse. Use only the specified type of fuse (see Chapter 11 Technical Data).
10. Insert the fuse holder until it can be heard snapping into place.
11. Use the power cable to reconnect the shockproof safety socket with the rear mains socket.

Fig. 9-1 Opening the fuse holder

9.5 Cleaning the Device

- 1. Switch device off.
- 2. Disconnect device from power supply. For this, pull the mains plug from the mains socket.

The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Wipe the surface of the device with a soft cloth moistened with the disinfectant. Make sure moisture does not enter the device.

The manufacturer recommends **Meliseptol® rapid** as disinfectant.

WARNING!
Cleaning the device
Do not sterilize the device.

9.6 Care of the Reusable Tube Set (NOT FOR SALE IN USA)

The device can be used with a reusable tube set that is listed as an accessory (see chapter 12 Accessory List) Please comply with the following information when using a reusable tube set.

9.6.1 Cleaning the Reusable Tube Set (NOT FOR SALE IN USA)

WARNING!
The reusable tube set (NOT FOR SALE IN USA) is made of silicone, polysulfone (PSU), and stainless steel. For cleaning and disinfection, use only pH-neutral or slightly alkaline cleaning agents (e.g. neodisher MediClean 2.0 %), disinfectants (e.g. Lysetol V 8 %), and drying and rinsing agents approved for the tube set materials.
Using unsuitable agents (e.g. neodisher MediKlar rinsing agent) may damage the tube system and especially the PSU connectors.

Reprocessing restrictions

The manufacturer has tested and approved the reusable tube set **(NOT FOR SALE IN USA)** for a specific number of reprocessing cycles. Always comply with the in-

structions on the label.

Never exceed the number of uses indicated by the manufacturer.

WARNING!

Please check the reusable tube (NOT FOR SALE IN USA) for signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation.

The reusable tube set (**NOT FOR SALE IN USA**) should be reprocessed as soon as possible following use (recommended: max 30 minutes). In case of heavy contamination (blood, tissue) it is recommended to remove these residues with a soft cloth while still at the site of the surgery. This cloth may be moistened with a slightly alkaline cleaner. The tube set is then transported, preferably in a dry basket, to the actual reprocessing location where the complete process is carried out.

Cleaning preparations

1. Preparatory cleaning at the reprocessing location:

Disassemble the individual parts of the tube set.

Notes concerning the pressure chamber membrane (for irrigation tube sets only).

WARNING!

Remove the membrane before cleaning.

- Before cleaning: Carefully remove the membrane (2) from the pressure chamber (1). Make sure not to damage the membrane. To do so, pull upwards on the membrane strap (6).
- Reinstalling membranes before sterilization: Place the lip of the membrane (3) into the ring groove of the pressure chamber (1). The strap (6) must be positioned in the notch (5) provided for this purpose.
- Press the lip of the membrane into the groove (4). A properly inserted membrane is flush with the pressure chamber and exhibits no wrinkles.

2. Preparatory cleaning:

Use a soft cloth to carefully wipe off each component. Remove any remaining blood and tissue residues. The cloth to be used may be moistened with a slightly alkaline cleaner.

1. Automatic cleaning and disinfecting:


The cleaning and subsequent disinfection of the individual components of the tube set can be carried out with the appropriate programmable cleaning and disinfecting systems for medical devices. The manufacturer has validated the applicability of this method in terms of hygiene success and material compatibility using the "Vario" program of the disinfectant manufactured by the company Miele.

Suitability of the cleaning device used in the context of automated cleaning and disinfection, the selected program, and the used cleaning or disinfecting agent, as well as a guarantee with respect to the attainable hygiene status and material compatibility are the sole responsibility of the operator!

Individual automatic processing methods must be independently validated by the operator.

2. Manual cleaning:

- Place the prepared individual components of the tube set in demineralized water (room temperature, 20 °C to 30 °C) for a period of at least 3 minutes to 5 minutes and then rinse them off under running water. Any remaining residues such as adherent clotted blood should be removed with a soft brush.
- Use a soft cloth to thoroughly dry the components.
- The completely dried tube set components must now soak for a period of 25 minutes to 30 minutes in a suitable cleaning agent (e.g., neodisher Med-

Cleaning

iClean 2.0 % or equivalent). Comply with the instructions and notes of the manufacturer of the used cleaning agent. Make sure the components are fully immersed in the cleaning agent.

- Place the components once more in demineralized water (room temperature, 20 °C to 30 °C) for a period of at least 3 minutes to 5 minutes and then rinse them off under running water.
- Flush the silicone tube for 15 seconds with demineralized water. Then use a high-pressure cleaning gun to emit 10 pressure surges, each one with a length of one second (1 second) into the silicone tube.
- After cleaning, allow the components to dry. For this purpose, they are to be placed in an appropriate tray and allowed to drip dry or use an air gun to dry the components. It is also possible to use, where appropriate, a sterile, soft cloth. Or, the components can be dried in an air dryer (10 minutes at 100 °C).

Disinfection

9.6.2 Disinfecting the Reusable Tube Set (NOT FOR SALE IN USA)

After cleaning, the individual tube set components must be disinfected. Only disinfect a thoroughly cleaned tube set.

1. Automatic disinfection:

The cleaning and subsequent disinfection of the individual components of the tube set can be carried out with the appropriate programmable cleaning and disinfecting systems for medical devices. The manufacturer has validated the applicability of this method in terms of hygiene success and material compatibility using the "Vario" program of the disinfectant manufactured by the company Miele.

Suitability of the cleaning device used in the context of automated cleaning and disinfection, the selected program, and the used cleaning or disinfecting agent, as well as a guarantee with respect to the attainable hygiene status and material compatibility are the sole responsibility of the operator!

Individual automatic processing methods must be independently validated by the operator.

2. Manual disinfection:

- The cleaned and completely dried tube set components must now soak for a period of 25 minutes to 30 minutes or as specified by the manufacturer in a suitable disinfectant (e.g. Lysetol 8 % or equivalent). Make sure the components are fully immersed in the cleaning agent and are not stacked on top of one another. Comply with the instructions and notes of the manufacturer of the used disinfecting agent. The tube set can be damaged if the concentration is too high or the components are soaked too long.
- After the soaking time, the individual tube set components should be removed from the disinfectant liquid with pliers with smooth jaws.
- Place the components once more in demineralized water (room temperature, 20 °C to 30 °C) for a period of at least 3 minutes to 5 minutes and then rinse them off under running water.
- Flush the silicone tube for 15 seconds with demineralized water. Then use a high-pressure cleaning gun to emit 10 pressure surges, each one with a length of one second (1 second) into the silicone tube.
- After cleaning, allow the components to dry. For this purpose, they are to be placed in an appropriate tray and allowed to drip dry or use an air gun to dry the components. It is also possible to use, where appropriate, a sterile, soft cloth. Or, the components can be dried in an air dryer (10 minutes at 100 °C).

9.6.3 Sterilizing the Reusable Tube Set (NOT FOR SALE IN USA)

Preparing for sterilization

Reassemble the tube set from the cleaned and disinfected parts as described in the manual.

Inspection and function test

Check the individual components and the assembled tube set for signs of damage. Never sterilize a tube set that shows signs of damage, especially brittleness and perforation. Tube sets with evidence of or suspected of damage must be blocked from further use.

Maintenance of the tube set is not required. In case of loss or damage to the pressure chamber membrane, this component must be replaced (see instructions in the manual).

Maintenance

Only clean, dry, disinfected, and assembled tube sets should be sterilized.

Sterilization

• Packaging

Pack the dry and reassembled tube set in a suitable sterile packaging for sterilization (cloth, bag, or tray).

• Sterilization

The manufacturer has validated saturated steam sterilization as a suitable method of sterilization for the tube set based on the parameters listed below using the pre-vacuum process. Sterilize the tube set according to these specifications.

	Pre-vacuum
Sterilization method	Pre-vacuum sterilizer
Packaging	Packed
Cycle time	5 min (273 °F/134 °C, 3 bar)
Drying time	10 min

CAUTION!

The indicated drying time depends on several variables, including the following: Altitude, humidity, type of packaging, preconditioning, size of chamber, mass of load, and placement in chamber. Users must verify that the drying time set in their autoclave yields dry surgical equipment when using the method of saturated steam sterilization described here.

Please comply with the instructions of the operating manual included with your autoclave.

• Storage

Always store the tube set after sterilization in a suitable sterile packaging or in a suitable environment to avoid contamination.

• Additional information

When sterilizing several tube sets, instruments etc. at once using one sterilization cycle, make sure not to exceed the max. load of the sterilizer.

Sterilization with ethylene oxide is possible but has not been approved by the manufacturer.

Sterilization with ETO

Do not sterilize with gamma rays.

Sterilization with gamma rays

Manufacturer’s specification

Inspection tests

Measured Values and Tolerances

10 Annual Inspection

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. These inspections have to be carried out on an annual basis. Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

The tests described in this chapter are designed specifically for trained personnel or a hospital technician. The operation of the device as well as its functionality and serviceability are easily checked. Each test conducted must be documented with date and signature in Chapter 16.1 "Test Log".

Each test conducted must be documented with date and signature on the test log (chapter 16 Test Log).

The following measuring tools and resources were used by the manufacturer to determine the measurements and tolerance levels:

Tube set	
Beaker/meas. cup	with scaling up to 2.0 l
2 irrigation fluid bags	
Suction tubes	Vacuum tube, drain tube
Suction container	3.0 l
instrument	Trumpet valve or stopcock
Stop-watch	


WARNING!
An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

10.1 Electrical Safety Test

1. Perform a visual inspection. Make sure that
 - the fuse corresponds with the specifications indicated by the manufacturer,
 - labels and stickers on device are legible,
 - the mechanical condition of the device allows for its safe use,
 - the device is clean to ensure proper and safe functionality.
2. Carry out the measurements for the ground leakage current, short-circuit current/housing leakage current, and the protective conductor resistance as per IEC 62353 in the current version or according to the applicable national standard.

10.2 Irrigation Test

Preparation

1. Switch device on using the power switch.
2. Insert the tube set.
3. Connect the irrigation tube with both irrigation fluid bags. The fluid bags have to be suspended 1.5 m above the device.
4. Start the irrigation process:
Press the irrigation start/stop key (LED lights up).
5. Open the flushing valve to fill the tube set completely with liquid.
6. Close the flushing valve.


Testing

1. Start test
Empty measuring cup.
2. Hold tube into measuring cup, open outlet, and stop at 30 seconds.
3. End test:
Press the irrigation start/stop key (LED goes out). The measuring cup must contain approx. 1.5 l ($\pm 10\%$) of water.

10.3 Suction Test

1. Connect the suction system including suction container and instrument.
Close stopcock if using instrument with stopcock.
2. Press the suction start/stop key (LED lights up). The device generates a negative pressure inside of the suction system.
The vacuum pump stops when a negative pressure of max. -60 kPa (-0.60 bar) has been reached.
The suction system has a leak, if the vacuum pump is not deactivated. Check the suction system for leaks and repeat test.
3. Starting suction:
Open the instrument.
The vacuum pump automatically reactivates, if the min. negative pressure is exceeded.
4. End suction:
Press the suction start/stop key again (LED goes out).

Each successfully conducted test must be documented with the test log.

11 Technical Data

Type or model designation:	PP120	
Manufacturer information:	W.O.M. WORLD OF MEDICINE GmbH Salzufer 8, 10587 Berlin	
Mains voltage range	AC 100-240 V ~	
Supply frequency range	50/60 Hz	
Fuse designation	2 x T 3.15 A, H, 250 V~, UL-recognized	
Power consumption:	Current (A)	Power consumption (VA)
Normal operation 100V AC/60 Hz	0.72	72
Normal operation 240V AC/50 Hz	0.38	93
Max. Power consumption	-	95
Protection class (I, II, III)	I	
Application part of the type (B, BF, CF)	BF	
Defibrillator protection (yes/no)	No	
Protection class (IP-Code)	IP21	
Classification (I, IIa, IIb, III) acc. to Appendix IX of the Europ. MDD .	IIa	
Conformity with the following standards: (in the currently valid version)	EN 60601-1 / IEC 60601-1 EN 60601-1-2 / IEC 60601-1-2	
Operating conditions	10 to 40 °C / 50 to 104 °F 30 to 75 % rel. humidity 70 to 106 kPa air pressure 3000 m max. altitude above sea level for device use	
Possible use with explosive anesthetic gases	This device is not designed for use with flammable anesthetic gases (Class AP) or flammable anesthetic gases with oxygen (Class APG).	
Storage and transportation conditions	-25 to +55 °C / -13 to +131 °F 10 to 85 % rel. humidity 70 to 106 kPa air pressure	
Max. Sound level	< 80 dB(a)	
Maximum negative suction pressure	-60 kPa	
Max. pump pressure	450 mmHg	
Maximum irrigation performance	2.0 l/min (with T0505-01/T0506-01/T0500-01 and T0502-01) 3.0 l/min (with T0508-01)	
Accuracy		
	flow	±10 % (for max. flow)
	Pressure	±50 mmHg (for maximum pump pressure)
	Deficit	n/a
Accuracy		
	flow	±10 % (for max. flow)
	Pressure	±50 mmHg (for maximum pump pressure)
	Deficit	n/a
Dimensions	Width x Height x Depth	380 mm x 148 mm x 388 mm
Weight	6.0 kg	
Interfaces/ports/connections		

IN/OUT signal for components	Service interface (D-Sub connector 9-pin, only for service, not ESD protected.)
Transponder technology (RFID)	Operating frequency: 13.56 MHz Modulation method: Amplitude Shift Keying (ASK) Carrier power: -7.51 dBμA/m at 10 m/32.8 ft
Mains power connection	IEC-60320-1 C14

12 Accessory List

Table 1: Tubings

Item No.	Packing units	Description
T0500-01	10 + 1	Irrigation day patient tube set, disposable/one-time use (NOT FOR SALE IN USA AND CHINA)
T0501-01	10	Irrigation day tube set, disposable/one-time use (NOT FOR SALE IN USA AND CHINA)
T0502-01	10	Suction tube set (1 connection), disposable
T0504-01	10	Vacuum tube set, incl. filter, can be used 30 days
T0505-01	10	Irrigation tube set, disposable
T0506-01	1	Irrigation tube set, reusable (NOT FOR SALE IN USA)
T0508-01	10	High flow suction and irrigation set with handle

Table 2: Additional accessories

Item No.	Packing units	Description
Z0521-01	10	Pressure chamber membrane
Z0155-01	1	Cannula for fluid bag, reusable (1 sales unit = 1 individual units)
Z0708-01	1	Tube connection luer lock male, 8 mm (1 sales unit = 1 individual units)
Z0101-01	1	Power supply cord EU, 2 m
Z0176-00	1	Power supply cord UK, 2 m
Z0102-01	1	Power cable US 2.5 m
Z0178-00	1	Power supply cord CN, 1.8 m
Z0560-01	1	Power supply cord BR, 1.8 m

13 Electromagnetic Compatibility

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

This device is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.

CAUTION!

Accessories

To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device LAP-Pump must be used only with the accessories listed in Chapter 12 "Accessory List".

To ensure the basic safety and essential functionality concerning electromagnetic interference over the life of the device, the device must be restarted no later than after 24 hours so that a diagnostic self-test can be performed. In addition, the specified maintenance intervals must be observed (see Chapter 9.2 "Maintenance Intervals").

This device complies with the electromagnetic compatibility (EMC) requirements for medical electrical devices as defined by IEC 60601-1-2. The limits used in testing provide a basic level of safety against typical electromagnetic interference likely to occur in professional health care facilities. Nevertheless, it can happen that individual performance features are no longer available or only to a limited extent due to the presence of EM interference.

13.1 Electrical Connections

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

Hospital employees should be informed about and trained in ESD precautionary measures.

ESD (Electrostatic Discharge) precautionary measures

13.2 Guidelines and Manufacturer's Statement – Electromagnetic Emissions

The LAP-Pump is intended for use in an environment as described below. The user/operator of the LAP-Pump should make sure the device is operated within such an environment.

Emitted interference measurements	Compliance	Electromagnetic environment guidelines
HF emission according to CISPR 11	Group 1	The LAP-Pump uses HF energy solely for its internal functions. Therefore, the camera's HF emission is very low and it is unlikely that devices in close proximity will experience interference.
HF emission according to CISPR 11	Class B	The LAP-Pump is suitable for use in all facilities including those in residential areas and those directly connected to a public utility network supplying buildings used for residential purposes as well.
Emission of harmonic oscillations according to IEC 61000-3-2	Class A	
Emission of voltage fluctuations / flickers according to IEC 61000-3-3	In compliance	

13.3 Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity

The LAP-Pump is intended for use in an electromagnetic environment as described below. The user/operator of the LAP-Pump must make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	IEC 60601-1-2 test level	Compliance Level	Electromagnetic environment/guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/ bursts according to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines Modulation 100 KHz	± 2 kV for power lines ± 1 kV for input and output lines Modulation 100 KHz	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Voltage surges according to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brownouts, and fluctuations of the power supply according to IEC 61000-4-11	$< 5\% U_T^*$ ($> 95\%$ dip in the U_T) for $1/2$ cycle $< 40\% U_T$ ($> 60\%$ dip in the U_T) for 5 cycles $< 70\% U_T$ ($> 30\%$ dip in the U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in the U_T) for 5 cycles	$< 5\% U_T^*$ ($> 95\%$ dip in the U_T) for $1/2$ cycle $40\% U_T$ (60 % dip in the U_T) for 5 cycles $70\% U_T$ (30 % dip in the U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in the U_T) for 5 seconds $0\% U_T$; $1/2$ cycle at 0° , 45° , 90° , 135° , 180° , 225° , 270° , and 315° $0\% U_T$; 1 cycle and $70\% U_T$; 25/30 cycles Single phase: at 0° $0\% U_T$; 250/300 cycles	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment. If the user/operator of the system requires the continuation of functionality after power interruptions/disruptions, it is recommended to supply the device with power from an uninterruptible power supply.
Supply frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.

*Note: U_T is the mains alternating voltage before applying the test levels.

13.4 Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity

The device LAP-Pump is intended for operation in an electromagnetic environment as described below. The user of the device LAP-Pump should make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	IEC 60601-1-2 test level	Compliance levels	Electromagnetic environment/guidelines
Conducted HF interference quantities according to IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V _{eff} 80 MHz to 800 MHz	<p>Portable and mobile wireless devices should not be used in closer proximity to the device LAP-Pump (including cables/lines) than the recommended safety distance calculated based on the transmitting frequency and the applicable formula. Recommended safety distances:</p> $d = 1.2\sqrt{P} \text{ for 150 KHz to 80 MHz}$ $d = 1.2\sqrt{P} \text{ for 80 MHz to 800 MHz}$ $d = 2.3\sqrt{P} \text{ for 800 MHz to 2.7 GHz}$ <p>P is the max. rated output of the transmitter in watt (W) according to the information provided by the manufacturer of the transmitter.</p> <p>The recommended separation distance in meters (m) is d.</p> <p>The field strength of stationary transmitters for all frequencies tested on site^a should be lower than the concordance level.^b</p> <p>Interference is possible in the proximity of devices featuring the following pictograph.</p> <div data-bbox="971 1113 1230 1345" data-label="Image"> </div>
Radiated HF interference quantities according to IEC 61000-4-3	6 V _{eff} in ISM bands, 150 kHz bis 80 MHz 3 V/m 80 MHz to 2.7 GHz	6 V _{eff} in ISM bands, 150 kHz bis 80 MHz 3 V/m 80 MHz to 2.7 GHz	

Note 1: The higher frequency range applies at 80 MHz and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

^a The field strength of stationary transmitters, such as base stations of wireless phones and cell phones, ham radio operators, AM and FM radio and TV stations theoretically cannot always be determined in advance. A study of the installation site should be considered to determine the electromagnetic environment concerning the stationary transmitter. If the field strength measured at the usage site of the device LAP-Pump exceeds the compliance levels listed above, the device LAP-Pump should be monitored for proper function. If unusual performance characteristics are observed, additional measures may be required such as changing orientation or the location of the device LAP-Pump.

^b The field strength should be less than 3 V/m for the frequency range of 150 kHz to 80 MHz.

IMMUNITY TESTS

Transients HF IEC 61000-4-3

Test frequency (MHz)	Band (MHz)	Service	Modulation	Max. power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ±5 kHz deviation	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

WARNING!
Portable HF communication devices can have an effect on the performance characteristics of the device LAP-Pump. Therefore, such devices must be kept at a minimum distance of 30 cm (independent of any calculation) from the device LAP-Pump, its accessories, and the cables.

14 Error and Warning Messages

Error and Warning Messages	Cause	Troubleshooting
Device error		
All LEDs flash rapidly Short warning beeps, 1x or 2x, are continuously repeated	Electronics defective	Turn the device off and back on after approx. 10 seconds. If the error message is displayed again, the device may not be used any longer. Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.
All LEDs flash rapidly Short 4 warning beeps are continuously repeated	Sensor errors: Impermissible deviation or error in the electronics components measuring the pressure. Pressure on sensor while no tube is inserted.	
All LEDs flash rapidly Short 5 warning beeps are continuously repeated	Motor error: Defective motor control. The motor error warning can also occur during the operation of the device due to a jammed roller wheel.	
All LEDs flash rapidly Short 7x warning beeps are continuously repeated	RFID error: Errors in the transponder electronics	
All LEDs flash rapidly Short 10x warning beeps are continuously repeated	Calibration error: The device is not calibrated properly.	
Suction does not start Suction status LED flashes slowly Short 3x warning beeps repeated 3 times	Vacuum pump defective	Turn the device off and back on after approx. 10 seconds. If the error message is displayed again, the device can be used without the integrated vacuum pump. In this case, use external suction. After the surgical procedure, have the device checked by an authorized service technician.
Tube warning		
Irrigation cannot be started The red tube status LEDs lights for 8 seconds Warning beeps, 1x short beeps	An expired tube was inserted.	Remove tube. Insert a new tube. Start the irrigation.
Irrigation cannot be started Red and green tube status LEDs flash slowly Irrigation and suction status LEDs light up Warning beeps, 3x short beeps	The pump was switched on while a tube was inserted.	Remove tube. Reinsert tube after a successful device test. A tube may not be inserted while starting.
The red tube status LED lights up after inserting the tube Warning beeps, 3x short beeps	The tube has no remaining uses left or the tube is not approved for this pump.	Reinsert the tube. If the tube is not detected again, remove tube and insert a new one into the pump.
The red LED flashes slowly after pressing the start key while a tube is inserted Warning beeps, 3x short beeps Irrigation cannot be started	The tube has no remaining uses left or the tube is not approved for this pump. Start key was pressed.	Reinsert the tube. If the tube is not detected again, remove tube and insert a new one into the pump.
Roller wheel stops Red and green tube status LEDs light up Warning beeps, 3x short beeps	The tube was detected by the transponder but the roller tube is not correctly tensioned over the roller wheel. Then the pump was started.	Remove tube and reinsert. Check the correct fit of the tube on the left side of the roller wheel. The roller wheel must be tensioned across the roller wheel.

Error and Warning Messages

en

Error and Warning Messages	Cause	Troubleshooting
Green tube status LED is flashing Warning beep, 1x short	Last cycle for reusable tube After the surgical procedure is finished, the inserted tube cannot be used any longer.	Make sure a new tube is available for the subsequent surgical procedure.
After pressing the start key: Red and green tube status LEDs flash slowly Warning beep, 1x short	The pump was started without an inserted tube or the tube was not detected by the pump, e.g. defective transponder.	Insert a new tube into pump before starting.
Device information		
Green tube status LED lights up	Tube OK - irrigation can be started	
All LEDs light up briefly and then go out Warning beep, 1x long	Self-test concluded successfully. The pump is ready	
Suction does not automatically shut off when instrument is closed	Leakage in tube system or defective vacuum suction	Check the tube system and suction container for leaks. If there is no leakage, suction is defective and must be switched off. It is possible to continue using the device without integrated vacuum pump. In this case, use external suction. After the surgical procedure, have the device checked by an authorized service technician.

15 Glossary

en

Term	Explanation
Day Patient Tube Set (NOT FOR SALE IN USA OR CHINA)	The day patient set is a tube set for arthroscopy and hysteroscopy and consists of a day and patient tube. While the patient tube must be replaced for each surgery, the day tube can remain connected to the device for up to 24 hours.
EMC	The electromagnetic compatibility describes the ability of a device to function satisfactorily within an electromagnetic environment without adding unacceptable electromagnetic interferences/disturbances to the environment that may cause problems for other devices or equipment located nearby.
ME Device	Medical electrical device used for therapy, monitoring or diagnosis of patients, equipped with no more than one connection to a supply network and which necessarily comes into physical or electrical contact with the patient or which transfers energy to or from the patient or which records or captures such energy transfer to or from the patient.
ME system	Medical electrical system consisting of a combination of devices, of which at least one is classified as a medical electrical device and specified by the manufacturer as such and which are joined together by a functional connection or by using a power strip.
HF surgery	High-frequency surgery uses high-frequency alternating current conducted through the human body to cut tissue or achieve coagulation.
Functional Test	The function test must be carried out before each surgery to ensure the respective device is fully functional.
Safety test	The safety test is part of the annual inspection.
Basic Function Test	The basic function test checks the basic function of the device and is part of the annual inspection.

16 Test Log

16.1 Test Log

[illegible]

16.2 Return Form

en

Please fill out this form when returning the device:

Name of owner:

Sales partner:

Address of person returning unit:

Street:	House number:
<div></div>	<div></div>

ZIP/Postal code:	City:
<div></div>	<div></div>

Country:

IMPORTANT!

Serial number (see identification plate):

Device type:

Description of defect:

Contact	Signature	Date
---------	-----------	------

Index**A**

Authorized service technician 4
Authorized trained personnel 34

C

Care and maintenance 4
Certification 34
Connect instrument 21
Connect patient tube 25
Connect the instrument 23
Contamination 4

D

Delivery inspection 14
Disposal 4

E

ESD (Electrostatic Discharge) precautionary measures 44
Essential performance 5
Exclusion of liability 4

F

Federal Law 4

G

Grounding contact 15

I

Insert new day tube 26
Insert tube set 21, 23
Inspection tests 40
Invalidating a tube set 19

L

Liability 34

M

Mains connection 15
Manufacturer's specification 40
Manufacturer's specifications 34
Measured Values and Tolerances 40

O

Only for U.S. operators 16

P

Potential equalization 16
Precautionary measures 45

R

Remove day patient tube set 25, 25
Remove the tube set 23
Removing the tube set 21
Returning the device 14
RFID technology (transponder technology) 19

S

Setting up 14
Sterilization with ETO 39
Subject to technical changes 3

T

Technical documents 34
The day tube 24
The patient tube 24
Transponder signal loss 19
Two-year maintenance interval 34

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Deze gebruiksaanwijzing bevat auteursrechtelijk beschermde informatie waar copyright op bestaat. Alle rechten voorbehouden. Het is verboden om deze gebruiksaanwijzing zonder uitdrukkelijke, schriftelijke toestemming van W.O.M. WORLD OF MEDICINE GmbH geheel of gedeeltelijk door middel van fotokopieën, microfilm of met andere middelen te vermenigvuldigen of te verspreiden.