PP110


Instructions for Use en

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Technical Writing Change Request

en

These instructions for use contain information that is subject to copyright. All rights reserved. These instructions for use may not be photocopied, duplicated on microfilm, or otherwise copied or distributed, completely or in part, without the approval of W.O.M. WORLD OF MEDICINE GmbH.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in these instructions for use are associated with registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or issues concerning these instructions for use.

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Copyright © W.O.M. WORLD OF MEDICINE GmbH

CE0197

CE marking according to Directive 93/42/EEC

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
Table of Contents

1	Important User Notes	3
2	Safety Information.....	4
3	General Information	5
3.1	Device Description	5
3.2	Intended Use and Contraindications	5
3.2.1	Intended Use.....	5
3.2.2	Contraindications	5
3.3	General Warnings and Precautions	5
3.3.1	General Warnings	5
3.3.2	Precautions.....	9
4	Initial Device Startup	11
4.1	Scope of Delivery	11
4.2	Setting up and Connecting the Device.....	12
4.2.1	Mounting on Stand.....	14
4.2.2	Positioning in Device Rack	15
4.3	Switching Device Off	15
5	Operating the Device	16
5.1	Device Front	16
5.2	Device rear.....	16
5.3	Overview of Usable Tube Sets	17
5.4	Using the Tube Sets	18
5.5	Irrigation Tube Set Connection.....	19
5.6	Inserting the Suction/Irrigation Set.....	21
5.7	Inserting a Day Patient Tube Set (NOT FOR SALE IN USA AND CHINA).....	22
5.7.1	Replacing the Patient Tube after Surgery.....	25
5.8	Connecting the Fluid Bags.....	26
5.9	Irrigation	26
5.10	Suction.....	27
5.11	Replacing the Suction Container	28
5.12	Removing a Vacuum Tube.....	28
6	Safety Functions.....	29
6.1	Pressure Limitation at 450 mmHg	29
6.2	Motor Cutout at 600 mmHg	29
6.3	Device Self-Test	29
7	Functional Test	30
7.1	Preparing the Function Test	30
7.2	Irrigation Test	30
7.3	Suction Test.....	31
8	Using the Device during Surgery	32
9	Care and Maintenance.....	33
9.1	Cleaning the Device	33
9.2	Maintenance Intervals.....	33
9.3	Maintenance by Authorized Service Technician.....	33
9.4	Replacing the Fuse	34
9.5	Care of the Reusable Tube Set (NOT FOR SALE IN USA)	35
9.5.1	General Notes.....	35
9.5.2	Reprocessing	37
10	Annual Inspection	43
10.1	Electrical Safety Test	43
10.2	Irrigation Test	44
10.3	Suction Test.....	44
11	Technical Data	45
12	Accessory list.....	47
13	Electromagnetic Compatibility.....	48
13.1	Electrical Connections.....	48
13.2	Guidelines and Manufacturer's Declaration – Electromagnetic Emissions	49
13.3	Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity	50
14	Error and Warning Messages.....	52
15	Test Log.....	54
15.1	Test Log.....	54
15.2	Return Form	55
	Index.....	56

1 Important User Notes

Read the instructions for use carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing staff or service personnel, or
- to damage or malfunction of device and/or accessories.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of its products.

The words WARNING, CAUTION, and NOTE carry special meanings. Sections marked with these words must be given special attention.

Subject to technical changes

Please note

WARNING!
The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third party.

CAUTION!
These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.

NOTE!
These paragraphs contain information to clarify the instructions or provide additional useful information.

en

Federal Law (only for U.S. market)

Exclusion of liability

Authorized service technician

Care and maintenance

Contamination

Waste management

2 Safety Information

CAUTION: Federal Law restricts this device to sale by or on the order of a physician.

The manufacturer is not liable for direct or consequential damages, and the warranty becomes null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained;
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedule is not adhered to.

The handing over of technical documents does not constitute authorization to make repairs or alterations to the device or accessories.

WARNING: Modifying the PP110 is not permitted.

Only authorized service technicians are permitted to repair, change settings, or modify the device or associated accessories and equipment. Any violation will void the manufacturer's warranty. Authorized service technicians are only trained and certified by the manufacturer.

The service and maintenance of the device and its accessories has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use.

NOTE: Service or maintenance work may not be carried out during surgery.

Before shipping, decontaminate device and accessories in order to protect the service personnel. Follow the instructions listed in this manual. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to reject contaminated products for repair.

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. For disposal of the device and its accessories, please consult the manufacturer or an authorized disposal company, in compliance with legal or national regulations.

3 General Information

3.1 Device Description

The device **PP110** is a pump providing an irrigation and suction function for laparoscopy.

The device itself is non-invasive and designed for use in multi-unit racks (non-sterile area). It pumps medical sterile irrigation fluid through a sterile tube. This flushes the corresponding fluids into the body cavity. The device can be used with electrolyte-free media (such as glycine 1.5 % or sorbitol 3.0 %) and with isotonic electrolyte-containing media (such as saline 0.9 % or Ringer's Lactate solution). The device is equipped with a vacuum suction function used to suction off secreted fluids. The secretion liquid is sucked through a tube into the container. The device functions only with the tube sets described in the accessory list (see chapter 12 "Accessory list").

3.2 Intended Use and Contraindications

3.2.1 Intended Use

WARNING!
To be used only with appropriate training
The device may only be used by surgeons or operating room personnel with the necessary training in the corresponding indication.

The device PP110 is a suction and irrigation pump intended for use during diagnostic and/or therapeutic laparoscopic procedures to irrigate fluid into and remove fluid from the abdominal cavity.

Intended use

3.2.2 Contraindications

Do not use this device **PP110** for procedures where a specific pressure has to be preset and reached such as a hysteroscopy or an arthroscopy.

Contraindications

3.3 General Warnings and Precautions

3.3.1 General Warnings

WARNING!
Technique and procedures
Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

WARNING!
Original accessories
For your own safety and that of your patient, use only original accessories.

WARNING!
Not explosion-proof
The device is not explosion-proof. Do not operate the device in the vicinity of explosive anesthetic gases and not in the vicinity of oxygen-enriched environments.

WARNING!
Falls and crashes
Place the device on a stable and level surface. Cables must be laid safely. Tubes between the device and the patient must not create any obstruction.

WARNING!
Device defect
Do not use this device if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage of the power plug and power cable.

WARNING!
Risk of electrical shock
To avoid the risk of electrical shock, only use this device when connected to a properly grounded power supply network.

WARNING!
Risk of electrical shock
To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

WARNING!
Function test
The function test must be performed prior to each device use.

WARNING!
Sterile media and accessories
Always work exclusively with sterile substances and mediums, sterile fluids, and sterile accessories if so indicated.

WARNING!
Replacing fuse
Replace the fuse only with a fuse of the same type and rating see Chapter 12 "Accessory list".

WARNING!
Unplug the power cable from the device before checking the fuse.

WARNING!
Maintenance and calibration
Do not open the device. The device may not be opened or calibrated by the user. Only authorized service technicians are permitted to repair, calibrate, or modify the device or its equipment.

WARNING!**Condensation / Water penetration**

Protect device from moisture. Do not use if moisture has penetrated the device.

**WARNING!**

The device is only intended for use with flexible fluid bags. Do not use glass containers as they might break. Fluid cannot flow quickly enough due to the vacuum being generated inside of the bottle. Risk of implosion.

**WARNING!****Keep full fluid bags ready for use**

Always keep a full fluid bag on hand to replace an empty one. This avoids having to interrupt surgery due to a lack of irrigation fluid.

**WARNING!****Contamination**

Do not use device or vacuum tube if signs of contamination are detected. Fluid that has penetrated can be detected, for example, by dirt encrustations on the connections for the suction or at the exit of the vacuum pump. Make sure the device or/and accessories can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!**

Different default settings of the warning message for identical or similar devices in the operating room may cause a risk due to conflicting acoustic signals.

**WARNING!****Preventing infections**

Sterilize reusable instruments before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

**WARNING!****Reprocessing of sterile disposable products**

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

**WARNING!****ME System (Medical Electrical System)**

Use only parts and/or devices from ME systems (see chapter 13 Electromagnetic Compatibility) in patient environments in compliance with the standard IEC60601-1 in the respective currently valid version.

**WARNING!**

Make sure the tubes for the suction-irrigation instrument do not have any kinks when installed. Make sure the tubes are not crimped or clamped off.

WARNING!

Filter

The vacuum tube with integrated filter is designed for max. 30 days. Replace the vacuum tube if it is obviously contaminated. The filter prevents that body fluids enter the interior of the device. Please note that the filter may reduce the suction capacity.

WARNING!

The vacuum tube sets for this device contain diethylhexylphthalate (DEHP), which is classified as toxic to reproduction according to the EU Directive 1272/2008/EEC on Classification, Labeling and Packaging of Dangerous Substances. DEHP may impair fertility and may cause harm to the unborn child. Therefore, this product must not be used for unauthorized applications. When applied within the intended use, the potential risk to pregnant or breastfeeding women as well as to children resulting from the DEHP contained in this product is not critical.

WARNING!

Cleaning of reusable tube set

The reusable tube set (NOT FOR SALE IN USA) is made of silicone, polysulfone (PSU), and stainless steel. For cleaning and disinfection, use only pH-neutral or slightly alkaline cleaning agents (e.g. neodisher MediClean 2.0 %), disinfectants (e.g. Lysetol V 8 %), and drying and rinsing agents approved for the tube set materials.

Using unsuitable agents (e.g. neodisher MediKlar rinsing agent) may damage the tube system and especially the PSU connectors.

WARNING!

Remove the membrane before cleaning.

WARNING!

Please check the reusable tube (NOT FOR SALE IN THE USA) for any signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation.

WARNING!

Make sure to connect tubes as indicated by the matching connectors.

WARNING!

The day tube remains connected to the device after each surgery. The patient tube must be disposed of immediately after surgery. A new sterile cap (included with the patient tube) must be screwed onto the hygiene connector at the day tube after each surgery. This sterile protector remains connected to the day tube until the next surgery.

WARNING!

Do not touch

Avoid touching the roller wheel. Risk of injury!

WARNING!

Do not touch patient and the ON/OFF switch at the same time.

**WARNING!**

The fluid bags must not be stored on the device or hung up directly above the device.

**WARNING!**

Replacement device and accessories

In case the device or any of the accessories fail during surgery, a replacement device and replacement accessories should be kept within close proximity to be able to finish the operation with the replacement components.

3.3.2 Precautions

CAUTION!

Cleaning the device

Do not sterilize the device.

**CAUTION!**

Endoscope

The device may only be connected with endoscopes designed for and featuring the technical specification permitting such a combined use. Any utilized endoscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600. Combining/connecting with other devices generates a medical electrical system (MES). The system configurator is responsible for compliance with the standard IEC 60601-1 / EN 60601-1 in its latest version.

**CAUTION!**

Electrical Interference

(See chapter 13 "Electromagnetic Compatibility"). Care was taken during the development and testing of this device that electrical interference from this device or other devices or instruments was practically eliminated. However, if you still detect or suspect such interference, please follow these suggestions:

- Move this, the other or both devices to a different location
- Increase distance between used devices
- Consult an electro-medical expert

**CAUTION!**

Use of other accessories, other transducers and cables

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM (see chapter 13 "Electromagnetic Compatibility") as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

General Information

CAUTION!
Ventilation of the device

- Avoid device overheating.
- Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).

CAUTION!
Incorrect voltage
Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

CAUTION!
Indoor climate
Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!
Power cables
Any power cables employed by the user that are not provided by the manufacturer must meet the safety requirements of the national standards in the respective current valid version.

4 Initial Device Startup

WARNING!

The product may be used only by doctors or health care professional in accordance with the intended use. Personnel must have training that qualifies them to operate laparoscopy pumps. In addition, personnel must read the instruction for use and operating instructions and become familiar with handling the device before first use. Personnel must have mastered handling and operation of the product before first use during a surgery. These instruction for use does not include descriptions or instructions for surgical procedures/techniques. It is also not intended to introduce to or train physicians in the use of surgical techniques.

4.1 Scope of Delivery

- Device PP110
- Instructions for Use
- Power cord

WARNING!

Check for defects

Check the product and the packaging for defects prior to use. Do not use if package is damaged.

CAUTION!

Premature unpacking

Do not premature unpack the devices, sufficient time must have passed to adapt to the room climate.

Always check all parts and accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized service company.

Delivery inspection

If it becomes necessary to return the device, always use the original packaging. The manufacturer does not take responsibility for damage that has occurred during transportation if the damage was caused by inadequate transport packaging.

Returning the device

Please fill out the return form enclosed at the end of the instructions. Enclose the manual with the device.

Please make sure that all required information has been supplied:

- Name of owner
- Address of owner
- Device type and model
- Serial number (see identification plate)
- Description of defect

4.2 Setting up and Connecting the Device

WARNING!
Always hold on to the power plug when disconnecting the device from the power supply. Never pull on the cable itself.

NOTE!
Locations
The device may be used only in the hospital or operating rooms.

Setting up

Place the device on a flat surface free of vibration located in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in chapter 11 "Technical Data" .

CAUTION!
ME systems (Medical Electrical System)
The medical electrical device is suitable for integration in ME equipment systems (see chapter 13 "Electromagnetic Compatibility"). Operation of the ME device in vicinity of non-ME devices may result in voiding the intended use of the ME device.

CAUTION!
ME Device in Rack
The ME device may not be used when stacked or stored directly adjacent to or with other devices. If it should be necessary to operate the device close to stacked with other devices, the ME device or ME system should be monitored to ensure it works properly as configured.

CAUTION!
Device Setup

- Equipment should be positioned such that power cord can be easily disconnected.
- Position the device in such a way that it is easy to use and switch off and on.
- Place the device outside the sterile area.

CAUTION!
Ventilation of the device

- Avoid device overheating.
- Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).

CAUTION!
Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!

Position of the user

- To avoid a malfunction, the user must be positioned correctly
- within a display viewing angle of $\pm 50^\circ$ to operate the device
 - up to 2 m/6.5 ft from the device front for monitoring the actual values

CAUTION!

Disconnection from the power supply is only guaranteed if the mains plug is pulled from the mains wall socket.

Mains connection

CAUTION!

Mains connection

- Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.
- Make sure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The mains connection cable may be plugged only into a properly installed, grounded safety wall socket (shockproof socket) (see DIN VDE 0100-710).
- Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

The power connection must be equipped with a grounding contact. Use the original power cable (if included in scope of delivery) to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.

Grounding contact

Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 and IEC 60320-C13. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.


Only for U.S. operators

Integrate the device into the potential equalization system as specified by local safety rules and regulations.

Potential equalization

4.2.1 Mounting on Stand

Fig. 4-1 Mounting on Stand

**WARNING!**

Please read the instruction manual for the stand before mounting the pump on the stand.

**WARNING!**

The position of the device on the stand can affect the stand's stability. If the device is mounted too high on the stand, the stand may tilt or topple. The user/operator of the device is responsible for making sure the device cannot be tipped over or otherwise become unstable.


1. Position the pump on the rod of the stand at the desired height, usually between 80 and 120 cm (example: Base of stand $\varnothing = 65$ cm, max. height of $h < 100$ cm recommended).
2. Tighten the toggle screw on the rear of the pump.

**NOTE!**

Make sure the toggle screw is tightened sufficiently to fasten the pump securely in place. Also check that the toggle screw remains tightened.

4.2.2 Positioning in Device Rack

en

4.3 Switching Device Off

Ending operation of the pump safely:

1. Press the ON/OFF switch (see Fig. 5-1 "Device front" [11](#))
2. Disconnect device from power supply. For this, pull the mains plug from the mains socket (see Fig. 5-2 "Device rear" [17](#)).


5 Operating the Device

5.1 Device Front

Familiarize yourself with the control and function elements at the device front.

Fig. 5-1 Device front

- ① Irrigation START/STOP key
- ② Irrigation status LED
- ③ Suction START/STOP key
- ④ Suction status LED
- ⑤ Roller wheel
- ⑥ Tube retainer
- ⑦ Irrigation pressure sensor
- ⑧ LED tube status red
- ⑨ LED tube status green
- ⑩ Connection for hygiene filter and vacuum tube
- ⑪ ON/OFF power switch
- ⑫ Power switch status LED
- ⑬ Power LED


5.2 Device Rear

Familiarize yourself with the connection elements at the device rear.

Fig. 5-2 Device rear

- ⑭ Service interface
- ⑮ Non-heating device plug
- ⑯ Fuse holder(s)
- ⑰ Metal edge
- ⑱ Pole mount
- ⑲ Screw
- ⑳ Vacuum output

5.3 Overview of Usable Tube Sets

The following table lists the function of each type of tube set depending on the selected indication.

Tube type	Art. No.	Tube category
Disposable (one-time use) (NOT FOR SALE IN USA AND CHINA)	T0500-01	Irrigation day patient tube set
Disposable (one-time use) (NOT FOR SALE IN USA AND CHINA)	T0501-01	Irrigation day tube
Disposable (one-time use)	T0502-01	Suction tube set (1 connection)
Can be used 30 days	T0504-01	Vacuum tube set, incl. filter
Disposable (one-time use)	T0505-01	Irrigation tube set
Reusable (NOT FOR SALE IN USA)	T0506-01	Irrigation tube set

Table 5-1

WARNING!
This product contains phthalates!
Animal experiments have shown that phthalates are potentially endangering reproduction. Based on the latest findings in science, in the case of prolonged exposure, risks for male babies and children cannot be excluded. On and in pregnant women, nursing mothers, babies, and children, medical products containing phthalates should only be used for a short time.
(DEHP) Diethylphthalate

WARNING!
Preventing infections
Sterilize reusable instruments before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

WARNING!
Reprocessing of sterile disposable products
Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

CAUTION!
The pressure chamber diaphragm is fragile and can be replaced if damaged. Replacement diaphragms are supplied with reusable tubes. The reusable tube sets must be prepared before each use. See chapter 9.5 Care of the Reusable Tube Set (NOT FOR SALE IN USA) for further instructions.

NOTE!
Tubes must be stored at room temperature. The shelf life for all tubes is 5 years.

NOTE!
Observe applicable hygiene regulations when disposing of the tube set.

RFID technology (transponder technology)

5.4 Using the Tube Sets

Transponder technology detects the tube type, the validity and reliability of a tube set automatically. This eliminates virtually all "operating errors" since non fitting, invalid, and not approved tube sets are reliably detected. In any approved tube set (see chapter 12 "Accessory list"), the corresponding transponder is located under the tube holder.

Invalidating a tube set

If an approved tube set is located in the device, the transponder technology automatically invalidates this tube set when the device is started. The tube status LED lights up green. Depending on the tube set, this has the following consequences:

- **Single-use tube set:** After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If the irrigation cycle is stopped, the device can be restarted within 30 minutes. If the device is switched off or in case of a power failure, the tube set is invalidated. The irrigation cycle can no longer be started. The tube status LED lights up red. If this is the case, you must insert a new, valid, and approved tube set.
- **Reusable tube set (NOT FOR SALE IN USA):** The reusable tube set can be used about 20 times. After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If irrigation is stopped during a use cycle, it can be restarted within 30 minutes without the set being invalidated. If the device is switched off or in case of a power failure, the current use cycle is invalidated. If an invalid reusable tube set has been inserted, a warning signal (beep) is emitted. A red LED will be displayed for 8 seconds.
- **Day patient tube set also day tube set (NOT FOR SALE IN USA AND CHINA):** The following description refers only to the day tube (for the patient tube, please see chapter 5.7 Inserting a Day Patient Tube Set (NOT FOR SALE IN USA AND CHINA)).

The reusable day tube set can be used for a total of 10 times. After inserting the tube set and starting the irrigation cycle, the tube set is invalidated after 10 minutes. If irrigation is stopped, it can be restarted within 30 minutes without the set being invalidated.

After the device is switched off or in the case of a power failure, the current use cycle is invalidated. Irrigation can no longer be started following the last use cycle. The tube status LED lights up red. If this is the case, you must insert a new, valid, and approved day tube.

Transponder signal loss

The device stops and cannot be restarted if the transponder loses its signal during current use. If the signal is restored within 20 seconds, it is possible to continue using the tube set.

5.5 Irrigation Tube Set Connection


Fig. 5-3 Irrigation tube set elements

- ① Protective caps
- ② Spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Ring
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑩ Luer lock connector

The irrigation tube set is available as a disposable or as a reusable tube set (autoclavable) (NOT FOR SALE IN THE USA) (see chapter 12 "Accessory list").

The tube set consists of 3 tube components, (irrigation tube (5), roller tube (7), and instrument tube (9)), a Y-connector (4) and 2 spikes (2). The spikes (2) are used to connect the tube components with the fluid bags.

The Luer lock connector (10) connects the instrument tube with the instrument.



Fig. 5-4 Inserting the tube set

- ② Spikes
- ⑤ Irrigation tube
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑪ Fluid bag
- ⑫ Tube retainer

Removing the tube set

Connect the instrument

Insert tube set

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to the Fig. 5-3 "Irrigation tube set elements" and Fig. 5-4 "Inserting the tube set".

1. To be carried out by non-sterile technician:

- ▶ If you are using a reusable tube set, open the sterile autoclavable container of the tube set and have a sterile technician remove the inner tube set.
- ▶ If you are using a disposable tube set, open the outer packaging of the tube set and have a sterile technician remove and open the inner tube set packaging.

2. To be carried out by sterile technician:

- ▶ Keep the Luer lock connector (10) in the sterile area and hand the tube end with the spikes (2) to the non-sterile technician.
- ▶ Connect the Luer lock connector (10) with the instrument (e.g. inflow cannula). Open inflow valve at instrument.

3. To be carried out by non-sterile technician:

- ▶ Before inserting, close both tube clamps (3) (Fig. 5-4 Inserting the tube set) on the bags.
 - ▶ Switch device on. The power status LED lights up. Make sure all tube sets are removed from the pump before switching the device on. A long beep is emitted to signal the successful completion of the device test and the tube can now be inserted.
- Inserting the roller tube is depicted in Fig. 5-5 "Inserting the roller tube set".
- ▶ When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (8) only if not pressurized.
 - ▶ Insert the unpressurized pressure chamber (8) carefully into the right notch of the tube retainer (12) up to the stop.
 - ▶ Place the roller tube (7) around the roller wheel.
 - ▶ Insert the tube adapter into the left section of the tube retainer by tightening the tube.
 - ▶ Make sure the tube adapter is firmly snapped into place! Irrigation and drain tube must point away straight towards the front from the tube retainer.

Connect an instrument and the corresponding bags (see also chapter 5.8 "Connecting the Fluid Bags").

Fig. 5-5 Inserting the roller tube set

5.6 Inserting the Suction/Irrigation Set

en


Fig. 5-6 Elements of the suction/irrigation tube set

- ① Protective caps
- ② Spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Roller tube
- ⑦ Pressure chamber with membrane and transponder
- ⑧ Instrument tube
- ⑨ Suction valve (blue)
- ⑩ Irrigation valve (white)
- ⑪ Suction and irrigation pipe
- ⑫ Trumpet valve
- ⑬ Suction and irrigation instrument
- ⑭ Suction tube
- ⑮ Connection to the suction container

The **suction/irrigation tube set** is available as a single use tube set (see chapter 12 Accessory list).

The tube set consists of two spikes provided with protective caps that are connected to the tube clamp (3) and the Y-connector (4) by a connecting piece. This connects the supply part with the irrigation tube (5) connected to the roller tube (6) by a ring. The instrument tube (8) is located on the pressure chamber with membrane and transponder (7).


Fig. 5-7 Inserting the suction/irrigation set

- ② Spikes
- ③ Tube clamps
- ⑤ Irrigation tube
- ⑥ Roller tube
- ⑦ Pressure chamber with membrane and transponder
- ⑧ Instrument tube
- ⑬ Suction and irrigation instrument
- ⑮ Fluid bag
- ⑯ Tube retainer
- ⑰ Suction container

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to the Fig. 5-6 "Elements of the suction/irrigation tube set" and Fig. 5-7 "Inserting the suction/irrigation set".

Remove the tube set

Connect the instrument

Insert tube set


Fig. 5-8 Inserting the roller tube set

1. To be carried out by non-sterile technician:

- Open the outer packaging of the tube set and let you remove the inner tube packaging by sterile personnel and open.

2. To be carried out by sterile technician:

- Keep the instrument in the sterile area and hand the tube end with the spikes (2) to the non-sterile technician.
- Open the trumpet valve (12) of the suction and irrigation instrument.

3. To be carried out by non-sterile technician:

- Before inserting, close both tube clamps (3) (Fig. 5-3 "Irrigation tube set elements") on the bags.
- Switch the device on. The power status LED on the ON switch lights up. Make sure all tube sets are removed from the pump before switching the device on. A long beep is emitted to signal the successful completion of the device test and the tube can now be inserted.
- Inserting the roller tube is depicted in Fig. 5-8 "Inserting the roller tube set".
- When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (8) only if chamber is not pressurized.
- Insert the unpressurized pressure chamber (8) carefully into the right notch of the tube retainer (12) up to the stop.
- Place the roller tube (7) around the roller wheel.
- Insert the tube adapter into the left section of the tube retainer by tightening the tube.
- Make sure the tube adapter is firmly snapped into place! Irrigation and suction tube must point away straight towards the front from the tube retainer.
- Connect the suction tube with the suction container (17) (see Fig. 5-7 "Inserting the suction/irrigation set").

5.7 Inserting a Day Patient Tube Set (NOT FOR SALE IN USA AND CHINA)

The **day patient tube set** consists of two tube sections comprised of the **day tube** and **patient tube**. The day tube establishes the tube connection between the fluid bags and the patient tube. The day tube is inserted into the tube retainer of the device and not replaced for the duration of one day between surgeries and/or surgical procedures. Only the patient tube (tube connection between day tube and the instrument) must be replaced between surgical procedures. Efficiency is increased with a day patient tube set if several surgeries are to be performed on one day. Make sure to connect tubes as indicated by the matching connectors.

CAUTION!

Use only original day patient tube sets of the manufacturer. Other day patient tube set are not compatible with the pump!


Fig. 5-9 Day tube set elements

- ① Protective caps
- ② Spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Irrigation tube
- ⑥ Ring
- ⑦ Roller tube
- ⑧ Pressure chamber with membrane and transponder
- ⑨ Instrument tube
- ⑩ Hygiene connector
- ⑪ Protective cap

The day tube includes the connectors to the fluid bags and to the hygiene connector (Fig. 5-10 "Patient tube elements" see ⑫) of the patient tube.

The day tube may be used the entire day of surgery. However, it is designed for a maximum of **10 uses**.

CAUTION!
The day tube must be replaced at the end of the day of surgery but not later than after 24 hrs.

The new sterile cap delivered together with the patient tube must be screwed onto the hygiene connector (Fig. 5-9 "Day tube set elements" ⑩) of the day tube set immediately after the surgical procedure to provide protection from contamination. This sterile protector remains connected to the day tube until the next surgery.

The patient tube has a connector on one side to connect with the day tube (hygiene connector ⑫). A luer lock connector (⑮) to connect to an instrument such as a trocar, for example, is fitted to the other side. Connect the hygiene connector ⑫ of the patient tube with the matching piece on the day tube (Fig. 5-9 "Day tube set elements" ⑩).

The patient tube


Fig. 5-10 Patient tube elements

- ⑪ Protective cap for hygiene connector
- ⑫ Hygiene connector
- ⑬ PVC tube, 285 cm
- ⑭ Luer lock connector with screw cap
- ⑮ Protective cap for Luer lock

The patient tube set has to be replaced for every patient to prevent cross-contamination. The day tube remains connected to the device after each surgery. The used patient tube must be disposed of after each surgery.


WARNING!

Make sure to connect tubes as indicated by the matching connectors.

Fig. 5-11 Inserting the day patient tube set

- (2) Fluid bag
- (3) Spikes
- (5) Tube clamps
- (8) Irrigation tube
- (9) Pressure chamber with membrane and transponder
- (10) Instrument tube
- (16) Hygiene connector
- (17) Tube retainer
- (18) Roller wheel

In order to separate sterile from non-sterile areas, assign the following tasks to the "sterile" or "non-sterile" nurse or medical technician. The indicated numbers refer to Fig. 5-9 "Day tube set elements", Fig. 5-10 "Patient tube elements" and Fig. 5-11 "Inserting the day patient tube set".

Remove day patient tube set

1. To be carried out by non-sterile technician:

- Open the packaging of the day tube and keep it with you.
- Open the packaging of the patient tube and have a sterile technician remove it for you.

2. To be carried out by sterile technician:

- ▶ Keep the Luer lock connector (14) in the sterile area and hand the tube end with the hygiene connector (12) to the non-sterile technician.
- ▶ Connect the Luer lock connector (14) with the instrument (e.g., inflow cannula). Open inflow valve at instrument.


3. To be carried out by non-sterile technician:

- ▶ Connect the hygiene connector (12) of the patient tube (Fig. 5-10 "Patient tube elements") with the hygiene connector (10) of the day tube (Fig. 5-9 "Day tube set elements").
- ▶ Switch device on. The power status LED on the ON switch lights up.
- ▶ 1. Insert the unpressurized pressure chamber (8) carefully into the lower notch of the tube retainer (17) (Fig. 5-11 "Inserting the day patient tube set") up to the stop.
- ▶ When inserting the roller tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber only if not pressurized.
- ▶ 2. Place the roller tube (7) around the roller wheel (18). Pull the tube until the tube ring fits in the upper part of the tube retainer.
- ▶ The validity of the tube set is now checked. When an invalid tube set is inserted, the LED showing the tube status turns red. A short acoustic signal is emitted 3x. The entire procedure must be repeated with a valid tube set if an invalid tube set has been inserted.
- ▶ Connect the tube end with the spikes to the full fluid bags.
- ▶ Press the Irrigation START/STOP key and wait until the day patient tube set is completely filled with irrigation fluid.
- ▶ 3. The device is now ready for use.

Connect patient tube

Insert new day tube

Fig. 5-12 Inserting the roller tube set

5.7.1 Replacing the Patient Tube after Surgery

WARNING!

The day tube remains connected to the device after each surgery. The patient tube must be disposed of immediately after surgery. A new sterile cap (included with the patient tube) must be screwed onto the hygiene connector at the day tube after each surgery. This sterile protector remains connected to the day tube until the next surgery.

1. To be carried out by non-sterile technician

- ▶ Open packaging of the patient tube.
- ▶ A sterile technician then removes the patient tube from the inside of the package.

2. To be carried out by sterile technician

- ▶ Keep the luer lock connector (14) and hand the other end of the patient tube to the non-sterile technician.
- ▶ Connect the luer lock connector (14) of the patient tube with the instrument (e.g., inflow cannula).

3. To be carried out by non-sterile technician

- ▶ Remove the sterile protector of the day tube and connect the new patient tube to the day set without delay.


5.8 Connecting the Fluid Bags

WARNING!

The device is only intended for use with flexible fluid bags. Do not use glass containers as they might break. Fluid cannot flow quickly enough due to the vacuum being generated inside of the bottle. Risk of implosion.

Fig. 5-13 Hanging the fluid bags

- ① Fluid bags
- ② Spikes
- ③ Tube clamps
- ④ Irrigation tube

The fluid bags must be suspended at a height of **h** between 0.7 to 1.0 m above the pump (Fig. 5-13 "Hanging the fluid bags") (pump inlet to the middle of the fluid filling in the bag).

1. The inflow tube can take irrigation fluid from two fluid bags ①. Close both clamps ③ at the branches of the inflow tube.
2. Grasp the spike ② at the provided handle when connecting or disconnecting.
3. Comply with sterile precautions when inserting the spike into the bag.
4. Open at least one tube clamp ③ at the irrigation tube.

5.9 Irrigation

The device is designed to work with an instrument featuring a trumpet valve or stopcock. Make sure tube set and instrument are completely filled with fluid. Remove all air from tube and instrument.

1. Insert the tube set (see chapter 5.5 "Irrigation Tube Set Connection").
2. Connect the irrigation fluid bag (see chapter 5.8 "Connecting the Fluid Bags").
3. Connect the instrument tube to the instrument.
Close stopcock if using instrument with stopcock.
4. Press the irrigation START/STOP key (LED lights up).
An admission pressure of 450 mmHg is generated in the tube. This admission pressure is necessary to provide the full irrigation volume/output.
5. Start the irrigation process:
Operate the trumpet valve or open the stopcock.
6. Stopping the irrigation process:
Release the trumpet valve or close the stopcock.
The admission pressure begins to be rebuilt so that the full irrigation volume/

output is available for the next irrigation procedure.

7. Press the irrigation START/STOP key again (LED goes out).

The tube is now relieved of pressure. The roller wheel spins backwards

5.10 Suction


Fig. 5-14 Suction system with suction-/Irrigation instrument

- ① Vacuum tube set with hygiene filter
- ② Suction container
- ③ Suction tube
- ④ Irrigation tube
- ⑤ Instrument


Fig. 5-15 Suction system with standard instrument

- ① Vacuum tube set with hygiene filter
- ② Suction container
- ③ Suction tube
- ④ Standard instrument

The suction system consists of the following:

- Vacuum tube
- Hygiene filter
- Suction container
- Irrigation tube
- Standard instrument or suction/irrigation instrument

1. Connect device with the suction container (using vacuum tube and hygiene filter). The vacuum tube set with hygiene filter can remain connected to the device for max. 30 days (Fig. 5-16 "Connect device with the suction system with suction-/Irrigation instrument").
2. Connect suction container with instrument (using irrigation tube Fig. 5-15 "Connect device with the suction system with suction-/Irrigation instrument").
Close stopcock if using instrument with stopcock.
3. Press the suction START/STOP key (LED lights up).
The device generates a negative pressure inside of the suction system. The vacuum pump stops when a negative pressure of max. -60 kpa (-0.60 bar) has been reached.
The vacuum pump automatically reactivates, if the min. negative pressure is reached.
4. Starting suction:
Actuate the trumpet valve or open the stopcock.
5. Stopping suction:
Release the trumpet valve or open the stopcock.
6. Press the suction START/STOP key again (LED goes out).


Fig. 5-16 Connect device with the suction system with suction-/Irrigation instrument

CAUTION!

The full suction output/volume is available only, if the suction system is pre-evacuated (vacuum generation). Pre-evacuation takes about 20 to 60 seconds

depending on the volume of the secretion container.

NOTE!

Only use suction containers with overflow protection.

5.11 Replacing the Suction Container

CAUTION!

Full suction containers must be replaced immediately without stopping surgery.

1. Stop the suction (see chapter 5.11 "Replacing the Suction Container").
Press the suction START/STOP key (LED goes out).
2. Actuate the trumpet valve or open the stopcock to vent the system.
3. Remove tubes from filled suction container.
4. Reattach tubes to empty suction container.
Close stopcock if using instrument with stopcock.
5. Press the suction START/STOP key (LED lights up).
6. Starting suction:
Operate the trumpet valve or open the stopcock.

5.12 Removing a Vacuum Tube

1. Remove the irrigation tube (4) from the instrument (5) (see Fig. 5-14 "Suction system with suction-/Irrigation instrument" or Fig. 5-15 "Suction system with standard instrument").
2. Remove the suction container (3) from the irrigation tube (4).
3. Remove the vacuum tube (1) from the suction container (3).
4. Remove the hygiene filter from the vacuum tube (1).
5. Remove the vacuum tube from the device.

6 Safety Functions

The electronic components continuously monitor the proper function of the device. Device malfunctions are indicated with audible warning signals (beeps), error messages, and/or by blocking device functions.

6.1 Pressure Limitation at 450 mmHg

The pump generates a max. pressure of 450 mmHg to reach the max. flow capacity during the irrigation mode. If external factors, for example, cause the pressure to exceed the pressure limits, the roller wheel reverses in an attempt to reduce the pressure.

The wheel does not rotate backward to reduce the pressure if a day patient set is used. If the actual pressure exceeds 450 mmHg, the motor stops and irrigation is interrupted to prevent contamination.

Pressure limitation when using a day patient set

6.2 Motor Cutout at 600 mmHg

The roller wheel is stopped if the pressure exceeds 600 mmHg for more than 5 seconds. The suction pump is also stopped in this state. Further irrigation is not possible. The roller wheel and suction pump automatically restart when the pressure has fallen below 600 mmHg.

6.3 Device Self-Test

After being switched on, the device performs a self-test of the sensors, the motor, and electronic components. All LEDs light up. A long beep is emitted after the successful device test.

A table listing a summary of possible error and warning messages is provided in chapter 14 "Error and Warning Messages".

7 Functional Test

WARNING!
Function test
The function test must be performed prior to each device use.

WARNING!
Preventing infections
Sterilize reusable instruments before surgery to prevent infections. Check all the single-use/disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.


WARNING!
Original accessories
For your own safety and that of your patient, use only original accessories.

7.1 Preparing the Function Test

1. The device is switched off and without an inserted tube set.
2. Switch device on using the power switch.
3. The device carries out a self-check. A long beep is emitted after the successful device check.
4. Insert tube set (see Fig. 5-4 "Inserting the tube set"). Connect the suction system (see Fig. 5-7 "Inserting the suction/irrigation set").

7.2 Irrigation Test

Fig. 7-1 Irrigation test

1. Connect the instrument.
Close stopcock if using instrument with stopcock.
2. Press the irrigation START/STOP key. The irrigation LED lights up.
3. Operate the trumpet valve or open the stopcock. Make sure tube set and instrument are completely filled with fluid. Remove all air from tube and instrument.
4. Release the trumpet valve or close the instrument stopcock.

7.3 Suction Test

1. Close stopcock if using instrument with stopcock.
2. Press the suction START/STOP key (LED lights up).
The device generates a negative pressure inside of the suction system. The vacuum pump is deactivated once a negative pressure of max. -60 kpa is reached.

CAUTION!

The suction system has a leak if the vacuum pump does not stop. Check the suction system for leaks and repeat function test. Do not use device, if vacuum pump is not operating or not operating correctly.

3. Open the instrument. Suction starts. The vacuum pump is automatically re-activated if the min. negative pressure is exceeded.
4. Press the suction START/STOP key again (LED goes out).

The function test is complete. The device is tested and ready for surgery.

WARNING!

Device defect

Do not use this device if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage of the power plug and power cable.

8 Using the Device during Surgery

WARNING!
Function test
The function test must be performed prior to each device use.

CAUTION!
Equipment should be positioned such that power cable can be easily disconnected.

Before surgery

1. Tube set and instrument are filled with irrigation fluid. The suction system is connected. The stopcocks are closed.
2. Press the suction START/STOP key (LED lights up).
The suction system is pre-evacuated.
Only then can the suction procedure be started with the instrument.

During Surgery

1. Press the irrigation START/STOP key (LED lights up).
2. Insert instrument into abdomen.
Start the irrigation process:
Operate the trumpet valve or open the stopcock.
3. Starting the suction process:
Operate the trumpet valve or open the stopcock.

After Surgery

NOTE!
The manufacturer recommends emptying the tube set using the pump when employing saccharine irrigation fluids (HF use). Remove fluid bags and instrument first.

1. Press the irrigation START/STOP key and the suction START/STOP key (LED goes out).
2. Use the power switch to turn device off.
3. Remove the tube sets.

NOTE!
Comply with hygiene rules and regulations when disposing of the tube set, collected fluid, and the suction container.

9 Care and Maintenance

Special care is necessary when servicing, maintaining, and storing the device and its accessories to maintain the functionality of the device and its accessories.

9.1 Cleaning the Device

- 1. Use the **ON/OFF** switch to turn the device off.
- 2. Remove the power cable.
- 3. Wipe the surface of the device with a soft cloth moistened with the surface disinfectant based on alcohol (for example Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the device.

CAUTION!
Cleaning the device
Do not sterilize the device.

9.2 Maintenance Intervals

WARNING!
Maintenance and calibration
Do not open the device. The device may not be opened or calibrated by the user. Only authorized service technicians are permitted to repair, calibrate, or modify the device or its equipment.

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. This inspection must be carried out once a year. The tests are described in chapter 10 "Annual Inspection".

Manufacturer's specifications

Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

NOTE!
Service or maintenance work may not be carried out during surgery.

9.3 Maintenance by Authorized Service Technician

An authorized service technician has to inspect and service the device at appropriate intervals to ensure the safety and functionality of the unit. The minimum service interval is two years, depending on frequency and duration of use. If the service interval is not maintained, the manufacturer does not assume any liability for the functional safety of the device.

Two-year maintenance interval

A sticker located on the rear panel of the device will remind you of the latest date for the next service or maintenance check.

Authorized service technicians are only trained and certified by the manufacturer.

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians.

Authorized trained personnel

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

Unauthorized personnel

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer

Liability

Technical documents

er of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

Certification

Ask the service technician for a certificate after he or she has inspected the unit or performed any service tasks. This certificate lists the type and scope of the service as well as the date and name of the servicing company together with the signature of the service technician.

WARNING!
Risk of electrical shock
To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

WARNING!
Modification of the device
This device may not be modified without the permission of the manufacturer.

WARNING!
Modified device
If the device is modified, suitable examinations and tests must be carried out to ensure further safe use of the device.

9.4 Replacing the Fuse

WARNING!
Replacing fuse
Replace the fuse only with a fuse of the same type and rating.

CAUTION!
Before replacing the fuse, check the values of the fuse to be inserted acc. to chapter 11 "Technical Data".

The fuse may be defective and is in need of replacement if:

- displays and LEDs (if available on your equipment) do not light up,
- the device does not function.

Check to make sure

- the main power supply cable is properly connected to the power supply input and to a safety socket,
- the house power supply fuse is functioning.

WARNING!
Unplug the power cable from the device before checking the fuse.

The device does **not** have to be opened to replace the fuse.

1. Switch device off.
2. Disconnect device from power supply by pulling the mains plug from the mains socket.
3. Remove power connection cable from mains socket.
4. The fuse holder is located next to the mains socket.
5. Remove fuse holder as depicted in Fig. 9-1 "Opening the fuse holder".
6. **A** Undo the latch of the fuse holder with a small screwdriver.
7. **B** Remove the fuse holder.
8. **C** Check fuse.
9. Insert a new fuse. Use only the specified type of fuse (see chapter 11 "Technical Data").
10. Insert the fuse holder until it can be heard snapping into place.


Fig. 9-1 Opening the fuse holder

9.5 Care of the Reusable Tube Set (NOT FOR SALE IN USA)

All reusable inflow tubes must be cleaned, disinfected, and sterilized before use. This applies in particular to the first use after delivery since the product is delivered non-sterile. Effective cleaning and disinfection are an essential prerequisite for effective sterilization.

The device PP110 can be used with a reusable tube set listed as an accessory (see chapter 12 "Accessory list"). Please comply with the following information when using a reusable tube set.

9.5.1 General Notes

WARNING!

Reprocessing

Only reprocess after training has been completed. The reusable tube set may only be reprocessed by trained personnel.

Within the scope of your responsibility for the sterility of the product in the application of material compatibility, please note the following:

- Only cleaning/disinfection and sterilization methods sufficiently validated and specifically designed for the tube set are to be used.
- The equipment used (washer-disinfector, sterilizer) must be regularly main-

- tained and checked.
- Adherence to the validated parameters during each cycle must be ensured.

Please also comply with the legal or statutory requirements applicable to your country as well as the hygiene rules of the physician's office or hospital. This applies in particular to the differing requirements regarding effective prion deactivation.

Products for cleaning and disinfection

WARNING!
Cleaning and disinfection
The reusable tube set (NOT FOR SALE IN USA) consists of silicone, polysulfone (PSU), POM (Polyoxymethylene) and stainless steel. For cleaning and disinfection, use pH-neutral or slightly alkaline cleaners and disinfectants approved for the utilized materials. If using unsuitable substances (e.g., Neodisher® Mediklar rinsing agent), the tube set and especially connectors made of PSU may become damaged.

Examples for cleaning agents and disinfectants which have been used for validation by the manufacturer:

Steps		Manufacturer	Product	Concentration [%]
Pretreatment		Dr. Weigert	Neodisher MediZym	0.5
Automated	Cleaning	Dr. Weigert	Neodisher MediClean	0.5
	Neutralization	Dr. Weigert	Neodisher Z	0.1
Manual	Cleaning	Dr. Weigert	Neodisher MediZym	0.5
	Disinfection	ASP	Cidex OPA	-

When selecting the cleaning and disinfecting agents, please ensure that the following ingredients are not included:

- Organic, mineral, and oxidizing acids with pH-value < 6.5
- Alkalis (pH value > 9.5, neutral/enzymatic or slightly alkaline cleaner recommended)
- Organic solvents (alcohols, ethers, ketones, petroleum)
- Oxidizing agents (such as hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons

It is imperative to comply with the concentrations, temperatures, and exposure times, as well as the post-rinsing requirements specified by the manufacturer of the cleaning agent and/or disinfectant.

Make sure that the entire tube set is rinsed and dried. If necessary, rinse or dry the tube set both from the Luer lock side and from the side of the spikes.


Never clean the product with metal brushes or steel wool but with a soft cloth.

Rinse aids and / or acid neutralizers are not permitted.


The product shall not be exposed to temperatures exceeding 142 °C (288 °F)!

Instrument oils or greases must not be used to maintain the product.

Workflow Automated Reprocessing

Workflow Manual Reprocessing

9.5.2 Reprocessing

9.5.2.1 Preparation

WARNING!

Membrane

Remove the membrane before cleaning.

1. Before cleaning: Carefully remove the membrane ② from the pressure chamber ① (see Fig. 9-2 "Pressure chamber membrane and its position"). This is done by pulling the strap of the membrane ⑥ towards the top. Make sure not to damage the membrane. Damaged membranes must be replaced.
2. Before starting sterilization, the membrane must be reattached to the pressure chamber according to the instructions in chapter 9.5.2.6 "Reinstalling the Membrane".

Removing the membrane from the product



Fig. 9-2 Pressure chamber membrane and its position

- ① Pressure chamber
- ② Membrane
- ③ Lip of the membrane
- ④ Groove
- ⑤ Notch
- ⑥ Strap
- ⑦ Spikes with protective caps
- ⑧ Tube clamps

9.5.2.2 Pretreatment

Pretreatment must be carried out both during automated and manual cleaning and disinfection. Please note that the disinfectant used for the pretreatment serves only personal protection against contamination and cannot replace the subsequent disinfection carried out after the cleaning step.

Directly after the application (no later than within 2 hrs), large impurities must be removed from the reusable tube set.

Instruction

1. If necessary, remove the protective caps of the spikes (7) (see Fig. 9-2 "Pressure chamber membrane and its position").
2. Open the tube clamps (8).
3. Rinse the tube set and the membrane for at least 1 min under running water (temperature < 35 °C/95 °F). Rinse all lumina of the tube set three times using a disposable syringe filled with water (minimum volume 100 ml).
4. Place the tube set and membrane in a freshly prepared cleaning solution. Observe relevant requirements as listed in chapter 9.5.1 "General Notes" concerning the selection, concentration, temperature, and exposure time.
5. Rinse all lumina of the tube set three times using a disposable syringe filled with cleaning solution (minimum volume 100 ml).
6. Manually remove all externally visible particles with a clean soft cloth and soft brush that is used only for this purpose and that does not scratch the surface.
7. Rinse the tube set and the membrane for at least 1 min under running water (temperature < 35 °C/95 °F). Rinse all lumina of the tube set three times using a disposable syringe filled with water (minimum volume 100 ml).

If you use a combined cleaning and disinfecting agent for this purpose - e.g., for reasons of occupational safety - please take into account the following:

- They should be aldehyde-free (otherwise fixation of blood contamination)
- They must have a tested efficacy (e.g., VAH/DGHM or FDA/EPA approval/ clearance/registration, or a CE marking)
- They must be suitable for disinfecting the product
- They must not contain any of the components listed under chapter 9.5.1 "General Notes" and are compatible with the product.

9.5.2.3 Cleaning and Disinfection

An automated process (cleaning and disinfecting equipment) should be used if possible. A manual process - even when using an ultrasonic bath - should be used only if an automated process is not available or not possible due to the significantly reduced efficiency and reproducibility.

Suitability of the cleaning device used in the context of automated cleaning and disinfection, the selected program, and the used cleaning or disinfecting agent, as well as a guarantee with respect to the attainable hygiene status and material compatibility are the sole responsibility of the operator!

Automated Cleaning and Disinfection

When selecting a washer-disinfector, make sure of the following:

- The washer-disinfector has a tested effectiveness (e.g., DGHM (German Society for Hygiene and Microbiology and Compliance with EN ISO 15883-1) - or FDA approval/ clearance/registration; CE marking according to medical device directive or regulation).

If possible, a tested thermal disinfection program is used (see Vario TD Program of Miele & Cie. KG, type: G 7836 CD with mobile injection unit E450 for MIS instruments; A0 value (scale for the inactivation of microorganisms of a disinfection process; > 3000 or - for older devices - at least 5 min at 90 °C/ 194 °F) (chemical disinfection may leave residues of disinfectant on the prod-

ucts).

The used program is suitable for the product

The used program contains sufficient rinsing cycles.

Only sterile or low-germ (max. 10 germs/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g., purified water/high purity water) is used for rinsing.

The air used for drying is filtered (oil-free, low in particles and germs).

The washer-disinfector is regularly maintained and checked.

When choosing the cleaning agent system, make sure of the following:

- It is generally suitable for the cleaning of instruments made of metal and plastics.
- If no thermal disinfection is used, additionally a suitable disinfectant with tested effectiveness (e.g., VAH (Association for Applied Hygiene) / DGHM - or FDA/ EPA (environmental protection agency) approval/clearance/registration or CE-marking) is used and that this is compatible with the cleaning agent used.
- The chemicals used are compatible with the products.

Instruction

1. Place tube set and membrane into the cleaning and disinfecting equipment. Make sure the tube set and membrane do not touch one another. Use a tubing connector to connect all lumina of the tubing set to the rinsing connectors of the washer-disinfector.
2. Start the verified thermal disinfection program (see Vario TD Program from Miele; AO value > 3000 or - for older devices - at least 5 min at 90 °C/194 °F).
3. Remove the tube set from the washer-disinfector after the end of the program.
4. If necessary, check and dry the tube set and the membrane (see chapter 9.5.2.4 "Visual Inspection" and chapter 9.5.2.5 "Subsequent Drying").

Manual Cleaning and Disinfection

Manual cleaning and disinfection should only be used if an automated process is not available.

When choosing the cleaning and disinfecting agent to be used, make sure of the following:

- It is generally suitable for the cleaning or disinfection of instruments made of metals and plastics.
- A disinfectant with tested effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/Registration or CE marking) is used which is compatible with the cleaning agent used
- The chemicals used are compatible with the products.

Combined cleaning/disinfecting agents should not be used if possible. Combined cleaning/disinfecting agents can be used only in cases of very low contamination (no visible impurities).

It is imperative to comply with the concentrations, temperatures, and exposure times, as well as the post-rinsing requirements specified by the manufacturer of the cleaning agent and disinfectant. Use only freshly prepared solutions, only sterile or low-germ (max. 10 germs/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g., purified water/high-purity water), or only a soft, clean and lint-free cloth and/or filtered air for drying.

Instruction

Manual cleaning

1. Place the tube set and the membrane in the cleaning bath for the specified exposure time of the cleaning agent so that they are sufficiently covered and do not touch each other.
2. Rinse all lumina of the products at least five times at the beginning and at the

Manual disinfection

- end of the exposure time using a disposable syringe filled with the cleaning agent (minimum volume 100 ml). Clean the exterior surfaces with a clean, soft cloth and a soft brush that is only used for this purpose.
3. Remove the tube set and the membrane from the cleaning bath and rinse them at least three times thoroughly (at least 1 min) with water.
4. Rinse all lumina of the tube set at least five times using a disposable syringe filled with water (minimum volume 100 ml).
5. Inspect the tube set and the membrane (see chapter 9.5.2.4 "Visual Inspection").

Instruction

1. Place the cleaned and inspected tube set and the membrane in the disinfection bath for the specified exposure time so that they are sufficiently covered and do not touch each other.
2. Rinse all lumina of the tube set at least five times at the beginning and at the end of the exposure time using a disposable syringe filled with the disinfection solution (minimum volume 100 ml).
3. Remove the tube set and the membrane from the disinfecting bath and rinse at least five times thoroughly (at least 1 min) with water.
4. Rinse all lumina of the tube set at least five times using a disposable syringe filled with water (minimum volume 100 ml).
5. Inspect and dry the product if necessary (see chapter 9.5.2.4 "Visual Inspection" and chapter 9.5.2.5 "Subsequent Drying").

9.5.2.4 Visual Inspection

After cleaning or cleaning/disinfecting, inspect the tube set and the membrane for corrosion, damaged surfaces, splinters, dirt and discoloration, and inspect for damaged products (see chapter Reusability for numerical restrictions on reuse). If the chrome-plated components of the Luer lock adapter or the spike are damaged (scratches, detachments), the product must be replaced immediately (no further use permitted).

A product that is still dirty must be cleaned and disinfected again.

9.5.2.5 Subsequent Drying

After cleaning, the tube set and membrane must be dried before packing. For this purpose, the tube set and membrane must be placed in an appropriate tray and drained or dried with an air pistol and filtered compressed air; if necessary, the individual components must be dried with a soft, clean and lint-free cloth. Alternatively, the components can be treated in an air dryer (10 minutes at 100 °C/ 212 °F).

Check the tube set and the membrane after drying as described in chapter 9.5.2.4 "Visual Inspection".

9.5.2.6 Reinstalling the Membrane

Reinstalling membrane before sterilization (see Fig. Pressure chamber and its position): Place the lip of the membrane (3) into the ring groove (4) of the pressure chamber. The strap (6) must be positioned in the notch (5) provided for this purpose. Press the lip of the membrane into the groove (4). A properly inserted membrane is flush with the pressure chamber and exhibits no wrinkles.

9.5.2.7 Packaging

Pack the reassembled, dry product (see chapter 9.5.2.6 "Reinstalling the Membrane") in disposable sterilization packaging (single or double packaging) which meets the following requirements (material/process):

- Conform to EN ISO 11607-1 and -2 /ANSI AAMI ISO 11607-1 and -2.
- Suitable for steam sterilization (temperature resistance up to at least 142 °C

(288 °F) with sufficient steam permeability).

- Sufficient protection of the product or the sterilization packaging against mechanical damage (e.g., from the spikes).

9.5.2.8 Sterilization

Only a clean, dry and disinfected as well as correctly assembled product (see chapter 9.5.2.5 "Subsequent Drying", chapter 9.5.2.6 "Reinstalling the Membrane" and chapter 9.5.2.7 "Packaging") may be sterilized. Only the sterilization procedures listed below are to be used for sterilization. Other sterilization methods are not permitted.

Steam Sterilization

- Fractionated vacuum process with at least three vacuum steps (with sufficient product drying >20 minutes).
- Steam sterilizer according to EN 13060 and EN 285. ANSI AAMI ST79.
- Validated according to EN ISO 17665-1 (valid commissioning (installation/operational qualification/IQ/OQ, and product-specific performance evaluation (PQ)).
- Maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665-1).
- Sterilization time (exposure time at sterilization temperature) between 5 and 18 min (for prion deactivation) at 132 °C (270 °F)/134 °C (273 °F).

Use of the less effective gravitation procedure is allowed only in case of non-availability of the fractionated vacuum procedure; requires significantly longer sterilization times and must be validated specifically for the products, devices, and parameters in sole responsibility of the user.

The actual required drying time depends on parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer condition, etc.) and must therefore be determined by the user. However, drying times should never be less than 20 minutes.

CAUTION!

Drying time

The indicated drying time depends on several variables, including the following: Altitude, humidity, type of packaging, preconditioning, size of chamber, mass of load, and placement in chamber. Users must verify that the drying time set in their autoclave results in dry surgical equipment when using the method of saturated steam sterilization described here.

The flash sterilization process is generally not permitted.

Other sterilization methods

Do not use hot air sterilization, radiation sterilization, formaldehyde or plasma sterilization

Individual automated processing methods must be independently validated by the operator.

Ethylene oxide sterilization is possible but not validated by the manufacturer

9.5.2.9 Storage

After sterilization, the products in the sterilization packaging must be stored dry and dust-free.

9.5.2.10 Reusability

The products can be reused up to 20 times, provided they are undamaged and uncontaminated, and used with due care; any further use or the use of damaged and/or contaminated products is the responsibility of the user.

If disregarded, any liability is excluded.

CAUTION!**Maximal cleaning / reprocessing cycles**

Do not clean / reprocess the reusable tube sets for more than 20 cycles as indicated on the label.

WARNING!**Signs of damage**

Please check the reusable tube (NOT FOR SALE IN USA) for signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation

CAUTION!**Check of the reprocessed tube set**

After each cleaning / reprocessing, the cleaned / reprocessed tube set must be checked for possible damages caused by the process before it is stored / used properly and professionally for the following use!

WARNING!**Contamination**

Do not use the tube set if signs or indications of a contamination are detected.

10 Annual Inspection

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. These inspections have to be carried out on an annual basis. Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

Manufacturer's specification

The tests described in this chapter are designed specifically for trained personnel or a hospital technician. The operation of the device as well as its functionality and serviceability are easily checked. Each test conducted must be documented with date and signature in chapter 15.1 "Test Log".

Inspection tests

Each test conducted must be documented with date and signature on the test log (chapter 15 "Test Log").

Measured Values and Tolerances

The following measuring tools and resources were used by the manufacturer to determine the measurements and tolerance levels:

Tube set	
Beaker/meas. cup	with scaling up to 2.0 l
2 fluid bags	
Suction tubes	Vacuum tube, drain tube
Suction container	3.0 l
Instrument	Trumpet valve or stopcock
Stop-watch	

WARNING!


An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

10.1 Electrical Safety Test

1. Perform a visual inspection. Make sure
 - the fuse corresponds with the specifications indicated by the manufacturer,
 - labels and stickers on device are legible,
 - the mechanical condition of the device allows for its safe use,
 - the device is clean to ensure proper and safe functionality.
2. Carry out the measurements for the device leakage current and the protective conductor resistance as per IEC 62353 in the current version or according to the applicable national standard (possible measuring points: metal edges around the non-heating device plug for mains connection. See (17) Fig. 5-2 "Device rear").

10.2 Irrigation Test

Preparation

1. Switch device on using the power switch.
2. Insert the tube set.
3. Connect the irrigation tube with both fluid bags. The fluid bags have to be suspended 1.5 m above the device.
4. Start the irrigation process:
Press the irrigation START/STOP key (LED lights up).
5. Open the flushing valve to fill the tube set completely with liquid.
6. Close the flushing valve.

Testing

1. Start test
Empty measuring cup.
2. Hold tube into measuring cup, open outlet, and stop at 30 seconds.
3. End test:
Press the irrigation START/STOP key (LED goes out). The measuring cup must contain approx. 1,5 l ($\pm 10\%$) of water.

10.3 Suction Test

1. Connect the suction system including suction container and suction and irrigation instrument.
Close the stopcock.
2. Press the suction START/STOP key (LED on). The device generates a negative pressure inside the suction system.
The vacuum pump stops when a negative pressure of a maximum of -60 kPa (-0.60 bar) has been reached.
The suction system has a leak, if the vacuum pump does not stop after 1 minute. Check the suction system for leaks and repeat the test.
3. Starting the suction:
Open the instrument.
The vacuum pump automatically restarts after a maximum of 30 seconds, depending on the instrument used, if the minimum negative pressure is exceeded.
4. Completing the suction:
Press the suction START/STOP key again (LED off).

Each successfully conducted test must be documented in the test log.

11 Technical Data

en

Type or model designation:	PP110			
Manufacturer information:	W.O.M. WORLD OF MEDICINE GmbH Salzufer 8, 10587 Berlin			
Mains voltage range	100-240 V~			
Supply frequency range	50/60 Hz			
Fuse designation	2 x T 3.15 AH, 250 V			
Power consumption	Upper voltage range	Current (A)	Voltage (V)	Power consumption (VA)
Upper voltage range				
	Normal operation	0.31	240	73
	Peak	0.33	240	79
Lower voltage range				
	Normal operation	0.61	100	58
	Peak	0.65	100	65
Protection class (I, II, III)	I			
Application part type (B, BF, CF)	BF			
Defibrillator protection (yes/no)	No			
Protection class (IP code)	IP41			
Classification (I, IIa, IIb, III) acc. to Appendix IX of the Europ. MDD	IIa			
Conformity with the following standards: (in their current versions)	IEC 60601-1, EN 60601-1, ANSI/AAMI ES 60601-1, CAN/CSA C22.2 No. 60601-1			
Operating conditions	10 to 40 °C/50 to 104 °F 30 to 75 % relative humidity 70 to 106 kPa air pressure 3000 m maximum altitude above sea level for device use			
Possible use with explosive anesthetic gases	This device is not designed for use with flammable anesthetic gases (Class AP) or flammable anesthetic gases with oxygen (Class APG).			
Storage conditions	-20 to +60 °C/-4 to +140 °F 10 to 85 % relative humidity 70 to 106 kPa air pressure			
Transportation conditions	-20 to +60 °C/-4 to +140 °F 10 to 85 % relative humidity 70 to 106 kPa air pressure			
Maximum sound level	< 80 dB(a)			
Maximum negative suction pressure	-60 kPa			
Maximum pump pressure	450 mmHg			
Maximum irrigation performance	2.0 l/min (with T0505-01/T0506-01/T0500-01 and T0502-01) 3.0 l/min (with T0508-01)			
Precision (repeatability)				
	Flow	±10 % (for maximum flow)		
	Pressure	±50 mmHg (for maximum pump pressure)		
	Deficit	n/a		

Technical Data

en

Accuracy		
	Flow	±10 % (for max. flow)
	Pressure	±50 mmHg (for maximum pump pressure)
	Deficit	n/a
Dimensions	Width x Height x Depth	210 mm x 148 mm x 260 mm
Weight		3.7 kg
Interfaces/ports		
IN/OUT signal for components		Service interface (D-Sub connector 9-pin, only for service, not ESD protected.)
Mains connection		IEC-60320-1 C14
RFID Transponder Technology		Transmit/Receive Frequency Range: 13.56 MHz ± 0.424 MHz Transceiver class: Class I RF Output Power: -10.83 dBμA/m at 10 m/ 32.8 ft Type of Antenna: Inductive Loop Antenna Antenna loop area: 0.00032 m² Modulation: Amplitude-shift keying (ASK) Mode of Operation (Simplex / Duplex): Duplex
Essential performance		Providing a irrigation pressure to generate a flushing flow. The maximum irrigation pressure is 450 mmHg +-10 %. Providing a suction pressure to generate fluid suction. The maximum suction pressure is 60 kPa.

12 Accessory list

en

Table 1: Tubings

Item No.	Packing units	Description
T0500-01	10 + 1	Irrigation day patient tube set, disposable/one-time use (NOT FOR SALE IN USA AND CHINA)
T0501-01	10	Irrigation day tube set, disposable/one-time use (NOT FOR SALE IN USA AND CHINA)
T0502-01	10	Suction tube set (1 connection), disposable
T0504-01	10	Vacuum tube set, incl. filter, can be used 30 days
T0505-01	10	Irrigation tube set, disposable
T0506-01	1	Irrigation tube set, reusable (NOT FOR SALE IN USA)

Table 2: Additional accessories

Item No.	Packing units	Description
Z0521-01	10	Pressure chamber membrane
Z0155-01	1	Cannula for fluid bag, reusable (1 sales unit = 1 individual units)
Z0708-01	1	Tube connection luer lock male, 8 mm (1 sales unit = 1 individual units)
Z0101-01	1	Power supply cord EU, 2 m
Z0176-00	1	Power supply cord UK, 2 m
Z0102-01	1	Power cable EU 2.5 m
Z0178-00	1	Power supply cord CN, 1.8 m
Z0560-01	1	Power supply cord BR, 1.8 m

Precautionary measures

13 Electromagnetic Compatibility

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

The device is to be used only for the purposes described in the Instructions for Use and is intended for use in environments in Professional Healthcare Facility Environment. This applies even if individual requirements meet the conditions for deviating electromagnetic environments. During installation and commissioning as well as during operation of the device, the compliance with the notes and instructions for EMC must be strictly observed.

CAUTION!

Accessories

To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device PP110 must be used only with the accessories listed in chapter 12 "Accessory list".

To ensure the basic safety and essential functionality concerning electromagnetic interference over the life of the device, the device must be restarted no later than after 24 hours so that a diagnostic self-test can be performed. In addition, the specified maintenance intervals must be observed (see chapter 9.2 "Maintenance Intervals").

This device complies with the electromagnetic compatibility (EMC) requirements for medical electrical devices as defined by IEC 60601-1-2. The limits used in testing provide a basic level of safety against typical electromagnetic interference likely to occur in professional health care facilities. Nevertheless, it can happen that individual performance features are no longer available or only to a limited extent due to the presence of EM interference.

13.1 Electrical Connections

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

Hospital employees should be informed about and trained in ESD precautionary measures.

ESD (Electrostatic Discharge) precautionary measures

13.2 Guidelines and Manufacturer's Declaration – Electromagnetic Emissions

The PP110 is intended for use in an environment as described below. The user/operator of the PP110 should make sure the device is operated within such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The PP110 uses RF energy solely for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PP110 is suitable for use in all facilities including those in residential areas and those directly connected to a public utility network supplying buildings used for residential purposes as well.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

13.3 Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity

The PP110 is intended for use in an electromagnetic environment as described below. The user/operator of the PP110 must make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	IEC 60601 test level	Compliance levels	Electromagnetic environment guidelines
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/bursts according to IEC 61000-4-4	± 2 kV for power lines Modulation 100 kHz	± 2 kV for power lines Modulation 100 kHz	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Surges according to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brownouts, and fluctuations of the power supply according to IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	
	0 % UT; 250/300 cycles	0 % UT; 250/300 cycles	
Supply frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.
Conducted RF IEC 61000-4-6	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	
Transients RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM by 1 kHz	3 V/m 80 MHz to 2.7 GHz 80 % AM by 1 kHz	

IMMUNITY to proximity fields from RF wireless communications equipment

Test method IEC 61000-4-3

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
5500						
5785						

WARNING!

Portable HF communication devices can have an effect on the performance characteristics of the device PP110. Therefore, such devices must be kept at a minimum distance of 30 cm (independent of any calculation) from the device PP110, its accessories, and the cables.

14 Error and Warning Messages

Error and Warning Messages	Cause	Troubleshooting
Device error		
All LEDs flash rapidly Short warning beeps, 1x or 2x, are continuously repeated	Electronics defective	Turn the device off and back on after approx. 10 seconds. If the error message is displayed again, the device may not be used any longer. Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.
All LEDs flash rapidly Short 4 warning beeps are continuously repeated	Sensor errors: Impermissible deviation or error in the electronics components measuring the pressure. Pressure on sensor while no tube is inserted.	
All LEDs flash rapidly Short 5 warning beeps are continuously repeated	Motor error: Defective motor control. The motor error warning can also occur during the operation of the device due to a jammed roller wheel.	
All LEDs flash rapidly Short 7x warning beeps are continuously repeated	RFID error: Errors in the transponder electronics	
All LEDs flash rapidly Short 10x warning beeps are continuously repeated	Calibration error: The device is not calibrated properly.	
Suction does not start Suction status LED flashes slowly Short 3x warning beeps repeated 3 times	Vacuum pump defective	Turn the device off and back on after approx. 10 seconds. If the error message is displayed again, the device can be used without the integrated vacuum pump. In this case, use external suction. After the surgical procedure, have the device checked by an authorized service technician.
Tube warning		
Irrigation cannot be started Red and green tube status LEDs flash slowly Irrigation and suction status LEDs light up Warning beeps, 3x short beeps	The pump was switched on while a tube was inserted.	Remove tube. Reinsert tube after a successful device test. A tube may not be inserted while starting.
The red tube status LED lights up after inserting the tube Warning beeps, 3x short beeps	The tube has no remaining uses left or the tube is not approved for this pump.	Reinsert the tube. If the tube is not detected again, remove tube and insert a new one into the pump.
The red LED flashes slowly after pressing the start key while a tube is inserted Warning beeps, 3x short beeps Irrigation cannot be started	The tube has no remaining uses left or the tube is not approved for this pump. Start key was pressed.	Reinsert the tube. If the tube is not detected again, remove tube and insert a new one into the pump.
Roller wheel stops Red and green tube status LEDs light up Warning beeps, 3x short beeps	The tube was detected by the transponder but the roller tube is not correctly tensioned over the roller wheel. Then the pump was started.	Remove tube and reinsert. Check the correct fit of the tube on the left side of the roller wheel. The roller wheel must be tensioned across the roller wheel.
Green tube status LED is flashing Warning beep, 1x short	Last cycle for reusable tube After the surgical procedure is finished, the inserted tube cannot be used any longer.	Make sure a new tube is available for the subsequent surgical procedure.

Error and Warning Messages	Cause	Troubleshooting
After pressing the START key: Red and green tube status LEDs flash slowly Warning beep, 1x short	The pump was started without an inserted tube or the tube was not detected by the pump, e.g. defective transponder.	Insert a new tube into pump before starting.
Device information		
Green tube status LED lights up	Tube OK - irrigation can be started	
All LEDs light up briefly and then go out Warning beep, 1x long	Self-test concluded successfully. The pump is ready	
Suction does not automatically shut off when instrument is closed	Leakage in tube system or defective vacuum suction	Check the tube system and suction container for leaks. If there is no leakage, suction is defective and must be switched off. It is possible to continue using the device without integrated vacuum pump. In this case, use external suction. After the surgical procedure, have the device checked by an authorized service technician.

15 Test Log

15.1 Test Log

[illegible]

15.2 Return Form

Please fill out this form when returning the device:

Name of owner:

Sales partner:

Address of person returning unit:

Street:

House number:

ZIP/Postal code:

City:

Country:

IMPORTANT!

Serial number (see identification plate):

Device type:

Description of defect:

Contact

Signature

Date

Index**A**

Authorized service technician 4

C

Care and maintenance 4

Certification 34

Connect patient tube 25

Connect the instrument 20, 22

Contamination 4

Contraindications 5

D

Delivery inspection 11

Disposal 4

E

ESD (Electrostatic Discharge) precautionary measures 48

Exclusion of liability 4

F

Federal Law 4

G

Grounding contact 13

I

Insert new day tube 25

Insert tube set 20, 22

Inspection tests 43

Intended use 5

Invalidating a tube set 18

M

Manufacturer's specification 43

Manufacturer's specifications 33

Measured Values and Tolerances 44

O

Only for U.S. operators 13

P

Potential equalization 13

Precautionary measures 49

R

Remove day patient tube set 24, 24

Remove the tube set 22

Removing the tube set 20

Returning the device 11

RFID technology (transponder technology) 18

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Deze gebruiksaanwijzing bevat auteursrechtelijk beschermde informatie waar copyright op bestaat. Alle rechten voorbehouden. Het is verboden om deze gebruiksaanwijzing zonder uitdrukkelijke, schriftelijke toestemming van W.O.M. WORLD OF MEDICINE GmbH geheel of gedeeltelijk door middel van fotokopieën, microfilm of met andere middelen te vermenigvuldigen of te verspreiden.