




正面

背面

<h3>Troubleshooting</h3> <table> <tr> <th>Trouble</th><th>Solution</th></tr> <tr> <td rowspan="2">Start up and no response</td><td>Please try several more times according to the correct boot-up method</td></tr> <tr> <td>Please use the attached charging base for charging. The charging time is about 1h</td></tr> <tr> <td>Camera does not turn on</td><td>Shut down and restart the ear cleaning rod</td></tr> <tr> <td>APP quits quickly</td><td>Delete the existing APP and re-scan the QR code for installation.</td></tr> <tr> <td>The power indicator of the ear cleaning rod flashing rapidly in red light</td><td>Low power mode, please charge the ear cleaning rod with the provide usb cabel from the original package</td></tr> <tr> <td>WiFi connection is unstable, sometimes it is disconnected in seconds</td><td>Restart ear cleaning rod, reconnect bebird APP</td></tr> <tr> <td>The product failed to display the real-time image after successfully connecting the bebird APP</td><td>Disconnect the phone from WiFi, restart the product and reconnect the WiFi</td></tr> <tr> <td colspan="2">If the above operation still fails to solve the problem, please contact Customer Service for further assistance</td></tr> </table>		Trouble	Solution	Start up and no response	Please try several more times according to the correct boot-up method	Please use the attached charging base for charging. The charging time is about 1h	Camera does not turn on	Shut down and restart the ear cleaning rod	APP quits quickly	Delete the existing APP and re-scan the QR code for installation.	The power indicator of the ear cleaning rod flashing rapidly in red light	Low power mode, please charge the ear cleaning rod with the provide usb cabel from the original package	WiFi connection is unstable, sometimes it is disconnected in seconds	Restart ear cleaning rod, reconnect bebird APP	The product failed to display the real-time image after successfully connecting the bebird APP	Disconnect the phone from WiFi, restart the product and reconnect the WiFi	If the above operation still fails to solve the problem, please contact Customer Service for further assistance																								
Trouble	Solution																																								
Start up and no response	Please try several more times according to the correct boot-up method																																								
	Please use the attached charging base for charging. The charging time is about 1h																																								
Camera does not turn on	Shut down and restart the ear cleaning rod																																								
APP quits quickly	Delete the existing APP and re-scan the QR code for installation.																																								
The power indicator of the ear cleaning rod flashing rapidly in red light	Low power mode, please charge the ear cleaning rod with the provide usb cabel from the original package																																								
WiFi connection is unstable, sometimes it is disconnected in seconds	Restart ear cleaning rod, reconnect bebird APP																																								
The product failed to display the real-time image after successfully connecting the bebird APP	Disconnect the phone from WiFi, restart the product and reconnect the WiFi																																								
If the above operation still fails to solve the problem, please contact Customer Service for further assistance																																									
<div>  </div> <h3>Product parameter</h3> <table> <tr> <td>Brand</td><td>bebird</td></tr> <tr> <td>Product Name</td><td>bebird smart visual ear cleaning rod</td></tr> <tr> <td>Model</td><td>MI3</td></tr> <tr> <td>Rod weight</td><td>24.3g</td></tr> <tr> <td>Network standard</td><td>IEEE 802.11b/g/n</td></tr> <tr> <td>Antenna</td><td>Internal FPC antenna</td></tr> <tr> <td>Working frequency</td><td>2.4Ghz</td></tr> <tr> <td>Image transmission rate</td><td>20fps</td></tr> <tr> <td>Image sensor</td><td>CMOS</td></tr> <tr> <td>Working temperature</td><td>-10 ~ 50°C</td></tr> <tr> <td>Battery</td><td>300mAh lithium battery</td></tr> <tr> <td>Battery life</td><td>About 75 minutes</td></tr> <tr> <td>Charging time</td><td>1h</td></tr> <tr> <td>Input parameters</td><td>5V --- 0.4A</td></tr> <tr> <td>Working voltage</td><td>4.2V</td></tr> <tr> <td>Rated power</td><td>1.5W</td></tr> <tr> <td>Lens diameter</td><td>3.5mm</td></tr> <tr> <td>Optimum focal length</td><td>1.5 ~ 2cm</td></tr> <tr> <td>Pixel</td><td>1000W</td></tr> <tr> <td>Gyroscope</td><td>6 axis</td></tr> </table> <div>  </div>		Brand	bebird	Product Name	bebird smart visual ear cleaning rod	Model	MI3	Rod weight	24.3g	Network standard	IEEE 802.11b/g/n	Antenna	Internal FPC antenna	Working frequency	2.4Ghz	Image transmission rate	20fps	Image sensor	CMOS	Working temperature	-10 ~ 50°C	Battery	300mAh lithium battery	Battery life	About 75 minutes	Charging time	1h	Input parameters	5V --- 0.4A	Working voltage	4.2V	Rated power	1.5W	Lens diameter	3.5mm	Optimum focal length	1.5 ~ 2cm	Pixel	1000W	Gyroscope	6 axis
Brand	bebird																																								
Product Name	bebird smart visual ear cleaning rod																																								
Model	MI3																																								
Rod weight	24.3g																																								
Network standard	IEEE 802.11b/g/n																																								
Antenna	Internal FPC antenna																																								
Working frequency	2.4Ghz																																								
Image transmission rate	20fps																																								
Image sensor	CMOS																																								
Working temperature	-10 ~ 50°C																																								
Battery	300mAh lithium battery																																								
Battery life	About 75 minutes																																								
Charging time	1h																																								
Input parameters	5V --- 0.4A																																								
Working voltage	4.2V																																								
Rated power	1.5W																																								
Lens diameter	3.5mm																																								
Optimum focal length	1.5 ~ 2cm																																								
Pixel	1000W																																								
Gyroscope	6 axis																																								
<h3>Hazardous Material</h3> <table> <tr> <th rowspan="2">Part name</th><th colspan="6">Hazardous Material</th></tr> <tr> <th>(Pb)</th><th>(Hg)</th><th>(Cd)</th><th>(Cr(VI))</th><th>(PBB)</th><th>(PBDE)</th></tr> <tr> <td>Printed Circuit Components</td><td>×</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr> <tr> <td>Case</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr> <tr> <td>Plastic rubber accessory</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr> </table> <p> <small>This form is formulated in accordance with SJ/T 11364</small> <small>○ : represent the hazardous material ratio is below GB/T26572-2011</small> <small>× : represent the hazardous material ratio is over GB/T26572-2011</small> </p>		Part name	Hazardous Material						(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)	Printed Circuit Components	×	○	○	○	○	○	Case	○	○	○	○	○	○	Plastic rubber accessory	○	○	○	○	○	○						
Part name	Hazardous Material																																								
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)																																			
Printed Circuit Components	×	○	○	○	○	○																																			
Case	○	○	○	○	○	○																																			
Plastic rubber accessory	○	○	○	○	○	○																																			
<h3>Warranty Card</h3> <table> <tr> <td colspan="2">User information</td></tr> <tr> <td>User name:</td><td>Telephone number:</td></tr> <tr> <td colspan="2">Purchase date:</td></tr> <tr> <td>Seller:</td><td>Telephone number:</td></tr> </table>		User information		User name:	Telephone number:	Purchase date:		Seller:	Telephone number:																																
User information																																									
User name:	Telephone number:																																								
Purchase date:																																									
Seller:	Telephone number:																																								
<div>  </div> <h3>Visual Ear Cleaner</h3> <p>FCC Warning</p> <p>This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.</p> <p>Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.</p> <p>NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:</p> <ul style="list-style-type: none"> • Reorient or relocate the receiving antenna. • Increase the separation between the equipment and receiver. • Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. • Consult the dealer or an experienced radio/TV technician for help. <p>Radiation Exposure Statement</p> <p>The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction</p>																																									