Therabody

# ThermBack LED

User Manual & Unit Warnings

Table of contents

1. Product Overview	EN	3
2. What’s in the Box	FRCA	16
3. Download the Therabody App		
4. Getting to Know the Device		
5. Using the Device		
6. How the Device Works		
7. Care, Cleaning, and Charging		
8. Troubleshooting Guide		

# ThermBack LED

## 1. Product Overview

Utilizing a powerful, 4-in-1 combination of heat, vibration, near-infrared LED, and far-infrared therapies, this cordless, wearable device delivers deliver fast, long-lasting relief from lower back pain. Designed for maximum therapeutic benefits, it targets specific areas of the lower back by positioning infrared LED lights along the spine, vibration motors on muscles and tissues, and carbon fiber-conducted heat with far-infrared energy across the entire lower back. Convenient and effective, ThermBack LED revolutionizes how you manage lower back pain.

## Intended Use

The ThermBack LED is intended to provide vibration and topical heating through mechanically powered heat and the infrared spectrum of light to temporarily increase circulation and relieve minor discomfort in the lower back.

## 2. What's in the box

**A** ThermBack LED device

**C** Travel pouch

**B** USB-C to USB-C cable

**D** Instruction manual

A power adapter is not included with this device.  
For unlimited use without battery life restriction, use while plugged in with a USB-C power adapter of at least 40W. Ensure it is from a certified manufacturer and has not suffered any structural damage.





Figure 1. What's in the box

### 3. Download the Therabody App

*If you currently have the Therabody app, proceed to step 3.*

1. Download the Therabody app from the App Store or Play Store.
2. Follow the in-app instructions to create a new app account and verify your email.
3. Once logged in, follow the onboarding steps to connect your device to the app. Be sure that the Bluetooth settings are enabled for your phone.
4. Select "ThermBack LED" as your Therabody device when prompted.
5. Press and hold the power button on your device to turn it on. The device must be powered on to locate and pair to the app.
6. Tap "Pair Device" and select "ThermBack LED" to confirm the device and app connection.
7. The Therabody app allows you to easily control and customize treatments from your phone. In the app you can control device settings, learn the best way to use each therapy, and get routine suggestions that show you exactly how to use the device based on your unique needs. Plus, you can save your favorite settings for easy access to the treatments you use most.

Figure 2. QR code to download the app

### 4. Getting To Know the Device

- 1 Power Button
- 2 Infrared LED Light Button
- 3 Heat Button
- 4 Vibration Button
- 5 USB-C Charging Port
- 6 Outer Compression Straps



Figure 3. Corresponding pictures/diagrams of the device

## 5. Using the Device

### Putting on the Device

1. With a sleek, comfortable design, ThermBack LED can be worn under or over your clothes. For maximum LED light therapy benefits, we recommend wearing it directly on bare skin.
2. Detach the outer compression straps and let them hang down. These come folded and attached to themselves.
3. Reach the device around the back of your body and hold one strap in each hand.
4. Wrap the left strap around your abdomen, then attach the right strap on top.
5. Reach back and grab the outer compression straps. Pull them tightly to your desired level of support and attach them in front. These can be attached over the control panel.
6. Once comfortable, make sure that the Therabody logo is centered on your back. The control panel does not need to be centered in the front.



Figure 4. Corresponding picture/diagram of using the device

## Using the Device

EN

1. **Turn on the device.** Firmly press and hold the power button. The LED indicator on the button will turn on when the device is on.
2. **Select a heat + far infrared setting.** Press the heat button to toggle between the three heat settings: low (default), medium, and high. Far infrared light therapy turns on when the heat setting is on. Heat and far infrared treatment will turn off after 20 minutes.
3. **Select a vibration pattern.** Press the vibration button to toggle through the three vibration patterns: low (default), high, and wave. Vibration treatment will turn off after 20 minutes.
4. **Turn near infrared LED light therapy on or off.** Near infrared LED light therapy will turn on when the device is on (default). Press the near infrared LED light button to turn off this therapy. For this therapy to be effective, allow the device to directly contact your skin. Near infrared LED light treatment will turn off after 10 minutes (the recommended dosage for maximum treatment benefits).
5. **Turn off the device.** The device will automatically shut off after the treatment is completed. Press the power button at any point during the treatment to turn off the device. The device will automatically shut off after 10 minutes if not in use. The LED light on the button will turn off when the device is off.
6. **Lock the device.** When not in use or during transport, lock the device to prevent turning it on accidentally. Quickly press the power button three times to initiate the lock function. Each LED light will flash sequentially to indicate that the device is locked. Quickly press the button three times again to unlock. Each LED light will flash sequentially (in the opposite direction) to indicate that the device is unlocked.

1

2-6


Figure 5. Corresponding picture/diagram of using the device

## 6. How the Device Works

MODALITY	INTENDED BENEFIT	SETTINGS		
Heat Therapy	Heat enhances blood flow to reduce soreness and stiffness, improve mobility, and accelerate recovery.	Low* 39 °C 102 °F (Default)	Medium* 42 °C 107 °F	High* 45 °C 113 °F
FIR (Far Infrared Light Therapy)	Traditional heat offers relief and relaxation by warming the skin and superficial muscles. Far infrared light uses long wavelengths of non-visible light (energy) for a warming effect to temporarily increase circulation and ease discomfort.	10-12 J/cm2 of energy generated across the entire device in a 20 minute treatment when heat therapy is on		
Vibration Therapy	Vibration therapy stimulates the muscles to relieve soreness and tension.	Low (Default)	High	Wave
NIR (Near Infrared Therapy)	NIR uses shorter wavelengths of non-visible light (energy) for a warming effect to temporarily increase circulation and relaxation.	ON Science-backed dose generated while NIR in use (Default)	OFF	

\*Actual skin temperature may vary depending on how the device is worn.

To learn more about the science and benefits behind the modalities included in this device, visit [therabody.com/science](https://therabody.com/science)

## 7. Care, Cleaning, and Charging

### Device Maintenance

The following maintenance instructions are important to ensure that your device continues to work as it was designed. Failure to follow these instructions may cause your device to stop working.

#### Cleaning the Device

Your device should be cleaned before its first use. For ongoing care, clean at least once a week. Carefully follow the cleaning instructions below to ensure longevity of the device and its components.

- The device is NOT waterproof or machine washable. Do not allow excess liquids or water to enter the product.
- Visually inspect the device for any obvious signs of debris build-up.
- Spot clean only. To clean, wipe the inner surface of the strap with a dampened microfiber cloth.
- Clean the LEDs with a 70% IPA or Clorox-type wipe.
- After cleaning, allow the device to dry thoroughly before storing or beginning another treatment.
- Do NOT submerge the device in water or clean it under running water. Do not allow the device to come into contact with corrosive solutions, which can damage the device's appearance and function.

#### Charging the Device

1. The device is USB-C enabled and comes with a charging cable. Before charging the device, ensure that the device is powered off.
2. Connect one end of the charging cable to the USB-C charging port on the device. It is located on the bottom of the control panel. Connect the other end of the charging cable to a power adapter.
3. The LED light will indicate the battery status when the device is on and charging. Press the button to check the battery status while the device is off.
  - Flashing orange: Low battery
  - Flashing blue: Half battery
  - Flashing green: High battery
  - Solid green: Full battery

The battery life of the device is up to 180 minutes, depending on the settings used. For unlimited use without battery life restriction, use while plugged in with a USB-C power adapter of at least 40W.

**\*\*Note:** If using a charger from a third party, ensure that it is from a certified manufacturer and has not suffered any structural damage. Do not try to use the device while charging.

## 8. Troubleshooting Guide

SCENARIO	ERROR CODE DISPLAYED	POTENTIAL CAUSES	POTENTIAL SOLUTIONS
The device stops working	The power button LED indicator blinks orange 10 times and then turns off	The device is not charged	Connect the USB-C cable (included in box) to the device and plug in the other end to an adapter (above 40W). Connect to a power outlet. Verify that the device is charging by checking proper connection between the power adapter, the USB-C cable, and the USB-C port on the device. The battery LED indicator will display the battery level. If using the device while charging, battery level will not increase, but it will remain powered on.
	The power button LED indicator blinks orange, blue, and green 5 times and then turns off	The battery pack has overheated during charging or use	Ensure the device is powered off, remove from charger, and let the device rest for 60 minutes. If the issue still persists after verifying all of the above, contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
The LED light therapy pad stops working	The dark orange LED indicator blinks 10 times	The LED light therapy pad has overheated	Ensure the device is powered off, remove from charger, and let the device rest for 60 minutes. If the issue still persists after verifying all of the above, contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
The heat therapy pad stops working	The 3 orange heat intensity indicators blink 10 times	The heat therapy pad has overheated	Ensure the device is powered off, remove from charger, and let the device rest for 60 minutes. If the issue still persists after verifying all of the above, contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
Vibration therapy stops working	The 3 white vibration indicators blink 10 times	The motor has overheated	Ensure the device is powered off, remove from charger, and let the device rest for 60 minutes. If the issue still persists after verifying all of the above, contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
Heat function is not working	The device is on, heat is not working, and the heat LED indicator is off	Heat treatment is turned off	Press the heat button to turn on the treatment.
	Heat does not turn on even after heat button is pressed, heat LED indicator is off	Heat button issue	If the heat treatment does not turn on after verifying the above, the heat connection may be damaged. Contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
	Heat is low	Low heat level may be selected	Increase the heat intensity by pressing the heat button

## 8. Troubleshooting Guide

EN

SCENARIO	ERROR CODE DISPLAYED	POTENTIAL CAUSES	POTENTIAL SOLUTIONS
The vibration function is not working	The device is on, but the vibration is not working	Vibration treatment is turned off	Press the vibration button to turn on the treatment.
	The vibration does not turn on even after the vibration button is pressed	Vibration button issue	If the vibration treatment does not turn on after verifying the above, the vibration connection may be damaged. Contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
Cannot charge the device	N/A	The power adapter cable is not well connected to the adapter socket on the console	Verify that the power adapter cable is well connected to the adapter socket on the console, and the power adapter is connected to the 100-240 Volt power outlet. If this doesn't work, contact Therabody. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
		The power outlet connected is beyond the range of 100-240 Volt	
The device is charging very slowly	Slow charging or incompletely charging	There may be an adapter issue, or the adapter power is too low	Ensure you are using the USB-C to USB-C cable provided, or one from a certified manufacturer.  If the issue still persists, Contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
Mechanical noise	There is a rattling noise when the device is turned on	Internal malfunction	Contact Therabody for further assistance. Visit <a href="https://www.therabody.com/us/en-us/support/support.html">https://www.therabody.com/us/en-us/support/support.html</a> for the available contact methods.
Bluetooth connectivity issue	Cannot connect to the app or the device isn't working properly once connected	The firmware is not updated to the latest version	Update the device firmware and try to connect to Bluetooth again.
		The phone version or system cannot support the connectivity	Update your phone to the latest version. If this doesn't work, try to install the app on another phone device.

DO NOT attempt to troubleshoot by disassembling the device as this will void the warranty and may cause injury.

### Store and transport the device

To store, gently roll or fold the device.

- To roll: Lay the device flat with the internal side facing you. Starting at the left outer strap, gently roll towards the middle. Avoid rolling the LEDs. Continue to roll, completing the right outer strap.
- To fold: Lay the device flat with the internal side facing you. Fold the left and right outer straps over the device. Lightly fold in half, avoiding folding the LEDs.

**Avoid tightly rolling the device. Do not fold over the LEDs, as this may damage the device's internal components.**

To store and transport the device, lay it flat in its travel pouch to not damage the internal components. Store the device in a cool and dry place (Temperature: 0°C/32°F - 35°C/95°F Relative humidity: 10~95% RH). Do not store the device where temperatures may exceed 35°C/95°F, such as in direct sunlight or in a vehicle.

## Warnings and Guidance (Precautions and Contraindications)

### Background

Therabody products are designed to unlock the body's natural ability to achieve health and well-being. Through science and technology, the Therabody portfolio of products allows people to access the therapeutic benefits of different natural phenomena to meet both their needs and preferences. There will be times when it is advisable to modify how devices are used (precautions) or times when it is not appropriate to use certain devices (contraindications). Read the following safety information for the device in its entirety prior to use.

### Important Safety Information

#### General Device Use

This device is intended for use by people in good health. This device is contraindicated against and should not be used by or on anyone with a history of epilepsy, seizures or cardiopathy.

The device is not recommended for anyone with an electronic implanted device (such as a pacemaker), cardiac arrhythmia, tumors, or acute episodes of inflammatory diseases. The device is not recommended for those who have arteriosclerosis, thromboses, or implants in the body region being treated. The device should not be used if you have dark brown or black spots, such as large freckles, birthmarks, moles, or warts, on the area being treated. The device is not recommended if you have eczema, psoriasis, lesions, open wounds, or active infections other than mild to moderate acne, such as cold sores, in the area being treated. Wait for the infected area to heal before using the device. The device should not be used if you have abnormal skin conditions caused by diabetes or other systemic or metabolic diseases. If you have a history of herpes outbreaks in the area of treatment, use of the device is not recommended unless you have consulted with your physician and have received preventive treatment.

Please consult your physician prior to using the device if you are pregnant and/or nursing. **Immediately stop using the device at the first sign of discomfort.** If you have any medical considerations please consult your doctor before using the device.

### Safety, Precautions, and Contraindications

These recommendations are derived from consultation with medical experts and published research regarding precautions and contraindications as of the printing date. For up-to-date information, please visit us online at <https://www.therabody.com/us/en-us/precautions-and-contraindications.html>.

#### Precautions:

Due care is required in these circumstances and device use may need to be modified. Consult with a medical professional if you currently have or suspect you may have any of the following conditions or if you have any questions.

- Recent injury or surgery
- Hypertension (controlled)
- Sensitivity to heat or pressure
- Mild peripheral neuropathy
- Osteopenia
- Osteoporosis
- Diabetes
- Varicose veins
- Bony prominences or regions
- Scoliosis or spinal deformity

- Medications that may alter client sensations
- This device should fit tightly but should NOT restrict circulation in any way
- Do not wear this device while sleeping, including during naps

#### Contraindications:

The following are circumstances where the potential risks may outweigh the benefits. Consult a medical professional before use.

- Skin rash, open wounds, blisters, local tissue inflammation, infections, bruises, or tumors on the treatment area
- Pregnancy/nursing
- Deep vein thrombosis, osteomyelitis
- Herniated disks
- Bone fracture or myositis ossifications
- Hypertension (uncontrolled)
- Abnormal sensations (e.g., numbness)
- Cancer/tumors
- Epilepsy
- Cardiopathy (heart disease)
- Vascular/arterial and circulation diseases, including blood clots, phlebitis, thromboses, and varicose veins
- Acute or severe cardiac, liver, or kidney disease
- Neurologic conditions resulting in loss or altered sensation
- Direct application to the face, throat, or genitalia
- Recent surgery or injury
- Peripheral vascular insufficiency or disease
- Severe peripheral neuropathy or other causes of sensory impairment
- Bleeding disorders
- Connective tissue disorders
- Medications that may thin the blood or alter sensations
- Direct pressure over surgical site or hardware
- Direct pressure over eyes or throat
- Severe scoliosis or spinal deformity
- Pacemaker or ICD
- Implanted devices such as IUDs, metal pins/plates
- History of embolism
- Photo allergy or disorder (ex. Lupus, porphyria)
- Medications that cause light sensitivity
- Suspicious lesions or skin cancer, please visit your physician
- If taking or using any retinol or/ sun sensitive medications or /products or benzoyl peroxide do not use infrared LED light
- Extreme discomfort or pain felt by client
- Do not use in the presence of unexplained pain
- Allergy to device material: Polyurethane, Spandex, Polyester, Nylon, Plastic, Silicone, Artificial PU Leather

## Limited Warranty and Customer Service

For full warranty information, please visit [www.therabody.com/warranty](http://www.therabody.com/warranty).

Customers who are in need of product support should visit <https://www.therabody.com/us/en-us/support/support.html> for the available contact methods.

FDA-Registered  
Limited Warranty Only With Authorized Retailer Purchase

© 2025 Therabody, Inc. All Rights Reserved.  
Patents at [www.therabody.com/patents](http://www.therabody.com/patents)

Manufactured and Distributed by:  
Therabody, Inc.  
1640 S Sepulveda, Suite 300  
Los Angeles, CA 90025, USA

In the Australian market, the Therabody ThermBack LED device warranty is 1 year.  
For Australian Market Only: The purchase of this device is accompanied by a limited warranty. For complete details about the rights afforded to you under this warranty, visit <https://therabody.com/anz/en-au/warranty-support/>.

Imported and Distributed By:  
Therabody Australia Pty Ltd  
Waterman Business Centre  
175 Maroondah Highway  
Suite 16, C307, Level 3  
Ringwood Vic 3134 Australia

Importers:  
Therabody International Limited  
5th Floor, 40 Mespil Road,  
DO4C2N4 Dublin, Ireland

Therabody UK Limited  
Portoken House 155-157  
The Minorities, London,  
United Kingdom, EC3N1LJ

The weight: 1.94lbs, 882g  
The capacity: 3200mAh  
The chemistry: Inner Fabric- 30% PU, 30% Spandex, 20% Polyester, 20% Nylon  
Outer Fabric: 40% PU, 60% Nylon

## UNIT WARNINGS

READ ALL INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL, ON THE CHARGER, AND ON THE DEVICE BEFORE USING OR CHARGING THE DEVICE AND ITS ACCESSORIES.

When using the device, the following precautions should always be adhered to:

1. **USE ONLY AS INSTRUCTED.** Use only Therabody recommended accessories and replacement parts. This device is not intended to diagnose, mitigate, or treat diseases. The device is an advanced mechanical

tool with electric components. If the device and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Failure to use and maintain the device and corresponding attachments in accordance with the instructions outlined in this manual will void the warranty and may result in product damage or physical injury. The device is intended for commercial and home use.

2. **NOT FOR CHILDREN.** The device is not intended for use by young children or persons with reduced physical, sensory, or reasoning capabilities or lack of experience and knowledge, unless given supervision or instruction concerning use of the device by a person responsible for their safety. The device is not a toy. Do not play with, bend, or pull the electrical components. Advise children not to play with the device or accessories. Adult supervision should be provided for those under 18 using the device.

3. **CHARGING.** If your device doesn't turn on or the battery indicator displays a low battery level, please charge before first use. Charge the battery using the USB-C cable (included). When you use the charging cable (included), make sure its USB-C connector is fully inserted into a compatible power adapter before you plug the adapter into a power outlet. You can also charge the battery with third party cables that are compliant with USB-C and with applicable country regulations and international and regional safety standards. Other adapters may not meet applicable safety standards, and charging with such adapters could pose a risk of death or injury. Do not force the cable into place. Fully charge the device at least once every six months to prevent device and battery damage. Unplug the power adapter and charging cable when not in use.

4. **DO NOT OVERCHARGE.** Do not leave the device connected to the power adapter for more than one hour after the battery has been fully charged. The battery includes a protection system to avoid the risk of overcharging. However, overcharging may reduce its life over time.

5. **CHARGING LOCATIONS.** Ensure that the battery status LED is OFF before charging. Insert the USB-C cable into the port (behind the power button) to charge the device. Do not use the device while charging. Avoid skin contact when the USB-C is plugged into a power source as it may cause discomfort or injury. Do not sleep or sit on the USB-C connector. The device should be charged indoors in a well-ventilated, dry location. Do not charge the device outdoors, in a bathroom, or within 10 feet (3.1 meters) of a bathtub or pool. Do not use the device or charger on wet surfaces, and do not expose the charger to moisture, rain, or snow. Do not use the device or its compatible charger in the presence of explosive atmospheres (gaseous fumes, dust, grain, metal powders, or flammable materials). Sparks may be generated, possibly causing a fire.

6. **UNPLUGGING THE CHARGER.** Pull the plug, not the cord, to reduce the risk of damage to the electrical plug and cord. Never carry the charger by its cord. Keep the cord away from heat, oil, and sharp edges. Do not stretch the charger cord or place the cord under strain. Do not handle the charger, terminals, or the device with wet hands. Always unplug this appliance from the electrical outlet immediately after using and before cleaning. Store cable to ensure it is not stepped on, tripped over, or otherwise subjected to damage or stress. Do not use a power adapter or charging cable that has received a sharp blow, been dropped, run over, or damaged in any way. For long-term storage, store with a fully charged battery.

7. **USAGE AND SKIN SENSITIVITIES.** Adjust the device to fit comfortably on your body. Always adjust the device while the device is OFF. Press and hold the Power Button to turn off the device. Make sure the Battery status LED is OFF. Discontinue the use of the device if you feel any discomfort. Some people may experience reactions to certain materials used in wearable items that are in prolonged contact with their skin. This can be due to allergies, environmental factors, extended exposure to irritants like sweat or other causes. You may be more likely to experience irritation from any wearable device if you have allergies or other sensitivities. If you have known sensitivities, take special care when using the device and accessories. If you experience redness, swelling, itchiness, or any other irritation or discomfort on your skin around, or beneath, your device, remove the device and consult your physician before resuming wear. Continued use, even after symptoms subside, may result in renewed or increased irritation.

8. **DO NOT OPERATE UNDER BLANKET AND PILLOW.** Excessive heating can occur and cause fire, electric shock, or injury.

9. **DEVICE CARE AND SERVICE.** Handle the device with care. The device contains electronic components that can be damaged if dropped, burned, punctured, or crushed. Do not disassemble the device or attempt to repair the device yourself. Disassembling the device may damage it, result in loss of water resistance, and may cause injury. If the device is damaged or if malfunctions occur, contact Therabody Customer service. The warranty will be void if the device, batteries, or charger are

disassembled or if any parts have been removed. Do not use if damaged as these may cause injury. Avoid heavy exposure to dust or sand. Do not puncture or damage the device. Puncturing the device can lead to incorrect operation.

10. **CLEANING INSTRUCTIONS.** Visually inspect the device for any obvious signs of debris build-up. Wipe your device with a damp microfiber cloth. After cleaning, allow the device to dry thoroughly before storing or beginning another treatment. A properly cleaned device should have no visible signs of debris or moisture. Note: This product is not waterproof. Only clean with a dampened cloth wipe by wiping the device. Do not submerge the device in water or clean it under running water. Do not allow the device to come into contact with corrosive solutions, which can damage the device's appearance and function.

11. **STORING THE DEVICE AND ACCESSORIES.** Store in a cool, dry place. Only charge the device when the ambient temperature is between 35°C/95°F and 0°C/32°F. Do not store the device or charging cable where temperatures may exceed 35°C/95°F, such as in direct sunlight, in a vehicle, or in a metal building during the summer.

12. **DEVICE DISPOSAL.** This device contains a lithium-ion battery, and care must be taken upon disposal of the device. Before disposal of this device, please review your local laws and requirements surrounding Lithium-Ion Battery disposal. The preferred method of disposal is recycling the whole device.

### Reporting adverse events to FDA

MedWatch is the Food and Drug Administration's (FDA) program for reporting serious reactions, product quality problems, therapeutic inequivalence/failure, and product use errors with human medical products, including drugs, biologic products, medical devices, dietary supplements, infant formula, and cosmetics. If you think you or someone in your family has experienced a serious reaction to a medical product, you are encouraged to take the reporting form to your doctor. Your health care provider can provide clinical information based on your medical record that can help the FDA evaluate your report. However, we understand that for a variety of reasons, you may not wish to have the form filled out by your healthcare provider, or your health care provider may choose not to complete the form. Your health care provider is not required to report to the FDA. In these situations, you may complete the Online Reporting Form yourself.

You will receive an acknowledgement from the FDA when your report is received. Reports are reviewed by FDA staff. You will be personally contacted only if we need additional information.

### Submitting adverse event reports to FDA

Use one of the methods below to submit voluntary adverse event reports to the FDA:

1. Report Online at: [www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home](http://www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home)
2. Consumer Reporting Form FDA 3500B. Follow the instructions on the form to either fax or mail it in for submission. For help filling out the form, see MedWatchLearn. The form is available at: [www.fda.gov/downloads/aboutFDA/reportsmanualsforms/forms/ucm349464.pdf](http://www.fda.gov/downloads/aboutFDA/reportsmanualsforms/forms/ucm349464.pdf)
3. Call FDA at 1-800-FDA-1088 to report by telephone.
4. Reporting Form FDA 3500 commonly used by health professionals. The form is available at: [www.fda.gov/downloads/aboutFDA/reportsmanualsforms/forms/ucm163919.pdf](http://www.fda.gov/downloads/aboutFDA/reportsmanualsforms/forms/ucm163919.pdf)

### GARMENT CARE.

Meanings of symbols

No.	Icon	Description	Location of Label
1		CE Mark	On rating label
2		Please read the user manual before use.	In packaging
3		Separate collection for waste electrical and electronic equipment. Note: For more information about disposal of equipment, its parts and accessories, please contact your local distributor.	On rating label
4	IP22	IP Classification: The first number 2: Protected against access to hazardous parts with a finger, and the jointed test finger of 12 mm $\Phi$ , 80 mm length, shall have adequate clearance from hazardous parts, and protected against solid foreign objects of 12.5 mm $\Phi$ and greater.	On rating label
5		RCM Mark	On rating label and packaging
6		Type BF apply part	On rating label
7		Manufacturer information	On rating label
8		UKCA Mark for UK	On rating label
9		Bluetooth	On packaging
10	IC	IC Canada	On rating label
11	 UDI : (01) 00840295408514 (10) 0000 (21) 00000	Unique Device Identification (UDI)	On rating label and packaging
12	 UDI : (01) 00840295408514 (10) 0000 (21) 00000	Date of Manufacture	On rating label

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## Medical EMC + Safety

IEC 60601-1 and IEC60601-11 Requirement ( Warning )---For Medical device apply for US, CA  
This device is Class III equipment with type BF applied. It complies with Medical Electrical Safety Standards (IEC 60601-1).

This device also complies with Medical EMC Standard (IEC 60601-1-2).

The has been tested and found to comply with the electromagnetic compatibility (EMC) limits for medical devices to IEC 60601-1-2: 2014. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

1\* WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally."

2\* WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

3\* WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME EQUIPMENT, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result."

## GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSION - FOR ALL EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emission		
The <b>ThermBack LED</b> is intended for use in the electromagnetic environment specified below. The customer or the user of <b>ThermBack LED</b> should assure that it is used in such an environment.		
EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The <b>ThermBack LED</b> uses RF energy only for its internal function. There for, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The <b>ThermBack LED</b> suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage function / flicker emissions IEC 61000-3-3	Not applicable	

## GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY - FOR ALL EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity			
The <b>ThermBack LED</b> is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ThermBack LED</b> should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact  ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines  ± 1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode  ±2kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Voltage dips, short interruptions and voltage variations on power supply input lines  IEC 61000-4-11	0 % UT; 0,5 cycle g) At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°  0 % UT; 250/300 cycle 0 % UT; 250/300 cycle	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the <b>ThermBack LED</b> requires continued operation during power mains interruptions, it is recommended that the <b>ThermBack LED</b> be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz magnetic field  IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U <sub>T</sub> is the a.c. mains voltage prior to application of the test level.			

**Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM**

Guidance and manufacturer's declaration – electromagnetic immunity			
The <b>ThermBack LED</b> is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ThermBack LED</b> should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF  IEC 61000-4-6	3 Vrms  150kHz to 80MHz  6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the <b>ThermBack LED</b> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:  $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{12}{V_2} \right] \sqrt{P}$
Radiated RF  IEC 61000-4-3	10 V/m  80MHz to 2.7GHz  385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	10 V/m  80MHz to 2.7GHz  385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz  $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.7 GHz  where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). <sup>b</sup> Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:  
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.			
<p>a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.</p> <p>b Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the <b>ThermBack LED</b> is used exceeds the applicable RF compliance level above, the <b>ThermBack LED</b> should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the <b>ThermBack LED</b>.</p> <p>c Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.</p>			

**Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS**

Recommended separation distances between portable and mobile RF communications equipment and the <b>ThermBack LED</b>				
Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz outside ISM and amateur radio $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	150 kHz to 80 MHz in ISM and amateur radio $d = \left[ \frac{12}{V_2} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{12}{V_2} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.00	0.35	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	35	70

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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#### FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

#### ISED Warning

This device complies with Innovation, Science, and Economic Development Canada licence-exempt RSS standard(s).

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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The device is compliance with RF exposure guidelines, users can obtain Canadian information on RF exposure and compliance.

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# Therabody

Born in Los Angeles, CA.  
Designed for everybody.

@Therabody