29.4 Mounting the transportation lock on the counterweight of the SCARA robotic arm

The units described below must be assembled and adjusted by Evident. If these units are assembled or adjusted by the customer, the operations are not ensured.

Transportation lock (rear)


ATTENTION

The transportation locks must be remounted each time before the unit is transported. Make absolutely sure not to lose the transportation locks and the mounting material.

» Hex screwdrivers (size 3 mm and 5 mm)

1. Switch the VS200 system off using the main power switch and disconnect the system from the power supply. To do so, disconnect the external power supply unit from the power supply.


2. Remove the following elements of the housing. You can remove more if required.
 - » Right housing side cover
3. Mount the transportation lock to the SCARA robot arm first. See [Mounting the transportation lock on the SCARA robotic arm on page 281](#).

4. Screw the locking plate on the back of the loader to the counterweight hand-tight using the middle screw. Then screw the plate to the rear of the housing using the four smaller screws. Tighten the five screws hand-tight with the tool (hex screwdriver size 3 mm / 5 mm).
5. Reconnect the VS200 scanner to the power supply.

29.5 Mounting the transportation lock on the SCARA robotic arm

The units described below must be assembled and adjusted by Evident. If these units are assembled or adjusted by the customer, the operations are not ensured.

Transportation locks (3) and (4)


ATTENTION

The transportation locks must be remounted each time before the unit is transported. Make absolutely sure not to lose the transportation locks and the mounting material.


» Hex screwdriver (size 3 mm)

1. Switch the VS200 system off using the main power switch and disconnect the system from the power supply. To do so, disconnect the external power supply unit from the power supply.


29 Preparing for transportation

2. Manually move the SCARA robot arm to the right. Move the three parts of the arm so that they are stacked vertically. To move the SCARA robot arm, grasp the segments of the arm.
3. Use the 4 hex screws (size 3mm hex screwdriver) to attach the upper part of the transportation lock.

4. Use the 4 hex screws (size 3 mm hex screwdriver) to attach the lower part of the transportation lock.

30 Proper selection of the power supply cord

If no suitable power supply cord has been provided, please select an appropriate power supply cord with a certification mark by referring to the specifications and the table below.

Evident is not responsible for damage caused by the use of uncertified power cords with Olympus devices.

Specifications

| | |
|------------------------|--|
| Voltage Rating | 125V AC (for 100-120V AC area) or 250 V AC (for 220-240V AC area) |
| Rated current | 9.5 A minimum |
| Temperature Rating | 60 °C minimum |
| Length | 3 m maximum |
| Fittings configuration | Grounding type attachment plug cap. Opposite terminates in molded-on IEC configuration appliance coupling. |

30 Proper selection of the power supply cord

Table 1 - Certification marks for power cords

The power cord has to have a certification mark from one of the bodies listed in table 1 or it must use a cable that has been tested by a body that is listed in table 1 or table 2. The plug has to have at least one certification mark shown in table 1. If you are unable to acquire a cord tested by one of the bodies listed in table 1 in your country, please use a cord that has been tested by a comparable body in your country.

| Country | Agency | Certification Mark | Country | Agency | Certification Mark |
|----------------|-------------|---|-------------|----------|---|
| Argentina | IRAM |  | Japan | JET, JQA |  |
| Australia | SAA |  | Canada | CSA |  |
| Belgium | CEBEC |  | Netherlands | KEMA |  |
| Denmark | DEMKO |  | Norway | NEMKO |  |
| Germany | VDE |  | Austria | ÖVE |  |
| Finland | FEI |  | Sweden | SEMKO |  |
| France | UTE |  | Switzerland | SEV |  |
| United Kingdom | ASTA BSI |  | Spain | AEE |  |
| Ireland | NSAI |  | U.S.A. | UL |  |
| Italy | IMQ |  | | | |

Table 2 - HAR Flexible cables

| Certification body | Printed or embossed harmonization mark (on the plug or cable insulation) | |
|--|--|-------|
| Comité Électrotechnique Belge (CEBEC) | CEBEC | <HAR> |
| VDE Verband der Elektrotechnik Elektronik Informationstechnik e.V. | <VDE> | <HAR> |
| Union Technique de l'Électricité (UTE) | USE | <HAR> |
| Istituto Italiano del Marchio di Qualità (IMQ) | IEMMEQU | <HAR> |
| British Approvals Service for Cables (BASEC) | BASEC | <HAR> |
| N.V. KEMA | KEMA-KEUR | <HAR> |
| SEMKO AB Svenska Elektriska Materielkontroll-anstalten | SEMKO | <HAR> |
| Österreichischer Verband für Elektrotechnik (ÖVE) | <ÖVE> | <HAR> |
| Danmarks Elektriske Materielkontrol (DEMKO) | <DEMKO> | <HAR> |
| National Standards Authority of Ireland (NSAI) | <NSAI> | <HAR> |
| Norges Elektriske Materielkontroll (NEMKO) | NEMKO | <HAR> |
| Asociación Electrotécnica Española (AEE) | <UNED> | <HAR> |
| Hellenic Organization for Standardization (ELOT) | ELOT | <HAR> |
| Instituto Português da Qualidade (IPQ) | np | <HAR> |
| Schweizerischer Elektrotechnischer Verein (SEV) | SEV | <HAR> |
| Elektriska Inspektoratet | SETI | <HAR> |

Underwriters Laboratories Inc. (UL)

SV, SVT, SJ oder SJT, 3 X 18AWG

Canadian Standards Association (CSA)

SV, SVT, SJ oder SJT, 3 X 18AWG

31 Declarations of conformity and disposal

31.1 CE conformity (Europe)

This system complies to the requirements of the following European directives:

- » Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU
- » Low Voltage Directive 2014/35/EU
- » EMC Directive 2014/30/EU
- » Machinery Directive 2006/42/EC
- » Radio Equipment Directive (RED) 2014/53/EU

This system complies with the requirements of Directive 2014/30/EU concerning electromagnetic compatibility according to Standard IEC/EN61326-1.

- » Emission: Class A
- » Immunity: Applied to industrial environment requirements. Operation of this equipment in a residential area may cause interference.

31.2 WEEE declaration (Europe)

In accordance with the European directive on Waste of Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste but should be collected separately. Refer to your local authority in the EU for return and/or collection systems available in your country.

31.3 RoHS conformity (Europe)

This product conforms with the European Union directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment 2011/65/EU.

31.4 FCC (USA)

Part 15 Clause 15.21

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC Part 15.19(a) [interference compliance statement], unless the following statement is already provided on the device label:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF Exposure guidance to the user on what separation distance to the body/limbs the device must be operated in order to comply with FCC RF Exposure requirements.

FCC Supplier's Declaration of Conformity
Hereby declares that the product
Product name: Optical Microscope
Model Number: VS200
Confirms to the following specifications:
FCC part 15, Subpart B, Section 15.107 and Section 15.109
Supplementary Information:
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
Responsible Party Name: Evident Scientific Inc.
Address: 48 Woerd Ave Waltham, MA 02453, U.S.A.
Phone Number: 781-419-3900

This device generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

1. Reorient or relocate the receiving antenna.
2. Increase the separation between the equipment and receiver.
3. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
4. Consult the dealer or an experienced radio/TV technician for help.

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31.6 RFID (Canada)

ISED notice:

This device complies with Industry Canada's licence-exempt RSSs. Operation is subject to the following two conditions:

- (1) This device may not cause interference; and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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| | |
|-------------------------------|--|
| Manufacturer and Registrant: | EVIDENT Technology Center Europe GmbH, Germany |
| 기자재명칭(제품명칭) Equipment Name | 미약 전계강도 무선기기 Weak electric field strength wireless device |
| Basic Model Number: | VS200-BU |
| Series Model Number: | VS200-BU-L |
| Registration No.: | R-R-OIS-VS200-BU |
| Basic Model Number: | VS200-LOADER |
| Registration No.: | R-R-OIS-VS200-LO |

31.8 IMDA (Singapore)

Complies with IMDA Standards.

31.9 UKCA (United Kingdom Conformity Assessed)

Complies with the relevant UK legislation.

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Without permission granted by the NCC, any company, enterprise, or user is not allowed to change frequency, enhance transmitting power or alter original characteristic as well as performance to a approved low power radio-frequency devices. The low power radio-frequency devices shall not influence aircraft security and interfere legal communications; If found, the user shall cease operating immediately until no interference is achieved. The said legal communications means radio communications is operated in compliance with the Telecommunications Management Act. The low power radio-frequency devices must be susceptible with the interference from legal communications or ISM radio wave radiated devices.

VS21-BU (VS200-BU)

CCAO21LP047BT9

31 Declarations of conformity and disposal

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32 Support

If you have any questions or problems that you can not resolve with the help of this manual, please contact your local service organization.

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EVIDENT Technology Center Europe GmbH
Wilhelm-Schickard-Straße 3, 48149 Münster, Germany

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Evident Scientific Private Ltd.
201, 202, 203 (2nd Floor), Tower C, Unitech Cyber Park, Sector 39,
Gurgaon-122022, India

Evident Australia Pty Ltd.
Level 4, 97 Waterloo Road, Macquarie Park, NSW 2113, Australia

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Evident Europe GmbH
Evident Scientific, Inc.
48 Woerd Ave, Waltham, MA 02453, U.S.A.
Evident Korea Co. Ltd.
Evident Scientific Singapore Pte Ltd.