



10 SLEEP REPORT

- After leaving the bed for 5-15 minutes, the recorded data will automatically be transferred to the app when connected the device with the app, you can check the recorded data in the app or recall historical data.
- Have to sleep for at least 30 minutes to have a sleep report.
- The sleep reports can be reviewing on the app all the time.



11 ALARM

Press the “” icon on “Alarm” page to set different alarm. The alarm clock can ring on the Phone or the Sleep Monitor.



12 FREQUENT ISSUES

Issue	Reason	Solution
APP: Device is not connected	<ul style="list-style-type: none">Device is not turned on.Bluetooth on the phone is not activated.The phone is beyond device's Bluetooth range.	<ul style="list-style-type: none">Turn on the device.Turn on the Bluetooth on the phone.Make sure the distance between the phone and the device is within 10 meters.Please charge the device.
Can't view the analysis data last night	<ul style="list-style-type: none">Device detection not started.Bluetooth is not connected.	<ul style="list-style-type: none">Turn on the device before fall asleep.Position the Sleepee correctly in the height of your chest between mattress and sheet.Make sure that the battery is fully charged, that the distance between device and your smartphone is less than 10 meters and that Bluetooth and Sleepee are connected.
Sleeping time not correct on the reports	<ul style="list-style-type: none">Device time is not local time.	<ul style="list-style-type: none">When rest time use the Sleepee, need connected the app to synchronize the local time to avoid the time difference leading to sleep report error.

Issue	Reason	Solution
Account cannot be registered	<ul style="list-style-type: none">Phone is not connected to the network.	<ul style="list-style-type: none">Connect your phone on WiFi, GPRS and other data communication networks.
When nobody is lying on the bed, heartbeat and respiratory signals are still being monitored.	<ul style="list-style-type: none">The sensor strap detecting vibrations or signals other than those of the person being monitored.Data shows delay.	<ul style="list-style-type: none">Avoid touching the user's bed.Remove all large vibration sources that are transmitted by the user, such as washing machines, electric fans, etc. Make sure the Sensors Strip is placed directly below the chest when the monitor is lying down.Please waiting for ten seconds.

Declaration

Please keep in mind that this product is not a medical device or a device that prevents sudden death. Users should not rely solely on and ignore the care of the monitored person. For users with diseases, the family should personally check the condition of the test.
The device is not suitable for people with a pacemaker, please remember!

13 NOTICES

- The best use environment: 8-24cm thick mattress;
- This product is not suitable for hard board beds;
- The use of soft cotton as a mattress may cause the sensor belt to capture signals inaccurately;
- Using a cold mat may affect the sensor's capture signal in the summer. It is not recommended to use a mahjong mat or a bamboo mat. Laying ice silk mats does not affect the use of this product.
- Do not bend this product excessively. This product has a built-in sensor. Excessive bending may damage the sensor. Do not use while charging, which may damage the life of the product.
- Do not contact to high temperature equipment, and do not use it with the electric blankets.
- Do not throw this product into water or wash it directly with water, which may cause damage inside the product.

- The sensor strip is suitable for flat and stable surfaces. Do not place the sensor belt on a movable bed (such as a hammock), which may cause the sensor belt can't work well.
- Make sure that no vibration sources are in close proximity to ensure that the sensor strip accurately track the user's sleeping. The vibration caused by electric fans, loud music, etc. may lead to misjudgment of this product.
- The sensor belt will not cause any harm to the user.
- If you do not use the device for a long time, please turn off the device, you can roll up the sensor belt and store it. Do not fold or compress the sensor belt, which may cause damage to the product.



14 WARRANTY STATEMENT

- Welcome to use Sleepee. Please read the instruction manual before using to ensure correct and safe use.
- Please use the original accessories to avoid unpredictable damage. If you do not follow the correct procedure to use or link incompatible accessories, this behavior will automatically invalidate the warranty and may even endanger the safety of you and others. The company does not assume any responsibility for this.
- If you disassemble or tamper with the outer casing of this product, the product will be out of warranty.
- The contents of this manual are provided according to the conditions at the time of manufacture of the product. Unless otherwise stipulated by applicable law, no clear or tacit guarantee of any kind is given to the accuracy, reliability and content of the document.
- The contents of this manual are protected by copyright laws. Without the written permission of the company, any copy, spread, distribution or application in any forms of this manual are strictly prohibited.

15 FCC Statement

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no

guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.


The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

INSTRUCTION MANUAL

Smart Sleeping Monitor

01 FUNCTION OPERATION

Includes: USB charging cable x 1, User manual x 1



02 PARAMETERS

Product Model	Sleepee
Sensor	PVDF
Weight	143 g (0.315 lb)
Battery	900mAh rechargeable lithium-ion battery
Wireless Connection	Bluetooth 4.0
Communication Distance	10m
CPU	nRF52832
Work Time	25 Days
Charge Current	500mA

03 INSTALL APP

Method 1: Scan the QR code to download the app.



Method 2: Search for "Sleepee" in App Store or Google Play Store to download the app.




iOS: Version 7.0 and newer.
Bluetooth 4.0

Android: Version 4.3 and newer.


04 POWER ON / OFF

Power on: Long press the power button until the blue light comes on.


Power off: Long press the power button until the red light comes on.



05 CONNECTING DEVICE




Turn on the device.




Turn on the phone Bluetooth.







Open the app and follow the instructions given by the app to complete the registration then login.



Press the device ID, then connect the sleep monitor via Bluetooth, when you see the connection icon like this  means connection completed.

Please ensure that the Bluetooth of the mobile phone is turned on during the connection process.
Keep the device close to the phone during the connection.


06 USER'S GUIDANCE



07 INDICATOR DESCRIPTION

The Bluetooth is not connected	The blue light flashes automatically for 30 seconds and then goes out. Press the button to turn off/on the indicator light.
Bluetooth connection successful	The blue light stays on for 5 seconds and then goes out.
Low power alarm	The red light flashes for 3 seconds accompanied by a BtB cue sound. Please charge in time.
Charging	The red light is usually on.
Charging completed	The green light is usually on.

08 FUNCTION OVERVIEW



09 REAL TIME SLEEP

- Press the “Start Monitor” button on the “Sleep” page of the app to start monitoring.
- The device also support offline monitoring, just switch on the Sleepee and place it on the correct positioning, then you can fall asleep, it will recording the data of you sleep automatically.
- The offline sleep data was stored for up to 7 days, and then covered in sequence after 7 days.

