

Hubei Huatong Kangda Medical Technology Co., Ltd.

Intelligent temperature and humidity data Collector

Model No.:HK-90

FCC ID: 2AUKL-HK-90

user's manual

一： Intelligent temperature and humidity data Collector APP Downloading

Use your phone to scan the QR code below:



二： Fix the instrument and turn it on

Use tape to secure the collector to the underarm or to the skin surface where the part is to be measured. Long press the back switch of the collector for three seconds

三： connection

Connected to the receiving device for the first time:


- a) Open the “Switch” button and fix the position, and the blue light will light to indicate the power on;
- b) Open an account to accept the device and turn on the Bluetooth to ensure that the device is connected to the network;
- c) Launch the APP software, register your account and log in
- d) Establish a connection between the data collector and the receiving device in the APP software. Connected to a connected receiving device;

Connect to a connected receiving device:


- a) Turn on the receiving device and turn on Bluetooth;
- b) The data collector is powered on and fixed in position
- c) Start the APP software and automatically establish a connection

四: installation battery (replace the battery)


1. Open the battery compartment cover (can be opened with a 1 dollar coin)



2. Open the battery compartment and replace the new button battery (marked "+" facing up)



3. Close the battery compartment cover



FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception,

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC RF Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction