

Welcome to KANG Smart Sleep

Assist to deep sleep with brain EEG technology

Dear users, thank you for choosing wearable sleep aid products. To ensure safety and proper use of the product features, please read this manual before trying.

Product Specifications

- | | | |
|------------------------------|-------------------------------------|----------------|
| The BCI Smart Sleep wearable | Noise canceling wireless headphones | Product Manual |
| USB charging cable | | |

1. Download Kang APP

Go to APP store and search for: KANG SLEEP



2. Register an account


Complete registration



3. Check battery is full

Press the button to turn on the device.

Connect the wireless noise canceling headphones with cell phone through Bluetooth.



Turn on the device, green light indicates full battery; consistently red & green lights indicate battery low.

4. Connect the device to APP

Open KANG Shui APP.



Click 'Connect' and find the device from the device list named "BrainUp xxx".



Put on the device properly following the instruction page.

5. One-click to fall asleep fast



Click the Sleep button to initiate the smart sleep function.

● Sleep button



Long press the end of sleep something sleep function can be withdrawn

Long press end sleep button

6. Sleep report



Wearable EEG device for brain signal monitoring




Brain-friendly Games

Brain games for memory enhancement, reaction improvement and concentration level training.

Mental health examination: anxiety, dementia and 9 other mental health indicators (coming soon)

Brain ability test: brain energy consumption, brain chaos, left and right brain laterality and 16 other brain ability indicators (coming soon)



FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF Exposure Information

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.