

Quick use manual v1.0

Product name: Intelligent thermometer

Product model: X1-19

1. Sweep code



Scan the QR code above, download it to the mobile phone, and install the intelligent thermometer app

2. Connect

Take out the smart thermometer, start the app, and click the search button to establish the connection.



3. Stick it up.

Apply the medical double-sided adhesive to the metal sheet on the bottom of the thermometer, and gently press it on the chest or armpit of the measured human body.

4. Observe.

Open the smart body temperature stick app and observe the body temperature.

5. how to conduct remote monitoring?

1A mobile app sends invitations to B mobile app;

2. B mobile app accepts the invitation. If a mobile app connects to the device and opens the network, the online status will be displayed on the B mobile app. Click to conduct remote monitoring;

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 5mm between the radiator and your body.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.