


# 1、 Main interface



**Menu content:** Temperature statistics, Personnel management, Temperature record, Settings, Log

**Prompt of Normal temperature:** When the detected temperature is less than 37.3 ° C, unregistered personnel will be prompted by a voice "Visitor's temperature is normal", the prompt of registered personnel will include their name.

**Prompt of Abnormal temperature:** When the detected temperature is greater than or equal to 37.3 ° C, unregistered personnel will be prompted by a voice "Visitor's

temperature is too high", and the abnormal body temperature and personnel information will be stored in the temperature statistics.

**Range of temperature display:** In any mode, the temperature display range of the screener is not narrower than 28.0 °C ~ 42.0 °C


**Tips:** Normal temperature voice prompts only in single player mode

### Abnormal Temperature Detection:



## 2、Function Description


① Temperature statistics: This function can query the temperature statistics and abnormal proportion of the day, query the temperature detection results by date, and export statistical reports by date, etc.




② **Personnel management** : Can register personnel, query and edit personnel information

**1、 Enter the person's name in the user name search box to query the temperature record**

**2、 There are two ways to add people by clicking the "+" in the lower right corner: single person taking pictures (pictures) and batch import**



### ③ Temperature record: Select time range to view historical temperature records



## 3、 Settings and Log

- **Recognition Mode:** Can be adjusted to single or multiplayer mode.
- **Device information:** In the device information, you can see the current app version, thermal imaging version, IP address, MAC address, etc. You can view the device identification by identifying the QR code, and you can update the current app version.
- **Online authorization:** The device is officially authorized before leaving the factory, and the button displays as "Authorized"
- **Advance :** The equipment is calibrated before leaving the factory. If the temperature deviation is too large due to environmental factors, please calibrate

under the guidance of a technician

- **Precaution:** View the usage environment of the temperature measuring device, installation height, identification distance, and common problems
- **Log:** If the system fails, please upload a log file so that technicians can analyze the cause and solve the problem.

## 4、 Equipment Photo



## Warning:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE: This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter

## RF Exposure Statement

To maintain compliance with FCC's RF Exposure guidelines, This equipment should be installed and operated with minimum distance of 20cm the radiator your body. This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter