



*This product is certified to conform to the Class IIa requirements of the 93/42/EEC Medical Devices Directive.*

# **Install & Service Manual Prudent**

## **Flat-panel Digital X-ray Detector**

**Model: Prudent1717/1417/1212**

**ver. 1.0.8**

---

Be sure to read and understand this manual thoroughly before using the product, and keep this manual in an easily accessible location for quick reference when required.

## INTRODUCTION

PRUDENT is a high-resolution digital X-ray imaging device commonly referred to as a flat panel detect. A built-in a-Si TFT flat panel type digital sensor receives X-ray and converts to digital image. X-ray photons are converted to digital output signals. The digital signals are then read out by TFTs. The image data file is saved at the computer for display, which can be linked with PACS and printed by DICOM printer through acquisition software. It is specifically designed to deliver x-ray imaging for medical, veterinary and industrial applications.

## ATTENTION

1. This manual guides the PRUDENT user to perform all installation and set-up procedures. Be sure that the user reads this manual thoroughly.
2. This includes the instructions of NETWORK card installation.
3. The use of calibration data and the method of creating calibration data are demonstrated.
4. Guidelines for the manual mode when connecting the generator directly to the PRUDENT is included.
5. No Modifying warning statement  
Never disassemble or modify the product as it may result in fire or electric shock. Also, since the instrument incorporates parts that may cause electric shocks and other hazardous parts, touching them may cause death or serious injury
6. FCC 15C Compliance statement  
Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
7. RF Exposure SAR statement  
RF exposure compliance  
The available scientific evidence does not show that any health problems are associated with using low power wireless devices. There is no proof, however, that these low power wireless devices are absolutely safe. Low power wireless devices emit low levels radio frequency energy (RF) in the microwave range while being used. Whereas high levels of RF can produce health effects (by heating tissue), exposure of low-level RF that does not produce heating effects causes no known adverse health effects. Many studies of low-level RF exposures have not found any biological effects. Some studies have suggested that some biological effects might occur, but such findings have not been confirmed by additional research. This device has been tested and found to comply with FCC/IC radiation exposure limits set forth for an uncontrolled equipment and meets the FCC radio frequency (RF) Exposure Guidelines in Supplement

## REVISION HISTORY

PRUDENT Manual Revision History

Revision	Revisions	Revised Date	Reviser
1.0.0	Initial drafting	2015.03.20	MoonKi Choi NamYeol Lee JiHong Jung
1.0.1	Add related content: Introduction, Warning, Micromagnetic ect.	2015.04.09	NamYeol Lee
1.0.2	Modified layout: 6.2 Environment, 6.4 PICB	2015.04.14	NamYeol Lee
1.0.3	Add related content: FCC and SAR test	2015.05.06	JiHong Jung
1.0.4	Add related content: Revision Hestory, DOC No.	2015.08.21	JiHong Jung
1.0.5	Revision by change of Company name	2018.03.02	JongMyung Shin
1.0.6	Modified some contents and typos : Page 1,2,9,11	2018.06.11	JongMyung Shin
1.0.7	A model(PIXX2430) with a miniaturized image sensor was added.	2020.06.30	JongMyung Shin
1.0.8	A model(Prudent) with a miniaturized image sensor was added.	2020.07.06	JongMyung Shin

## **Contents**

INTRODUCTION.....	2
ATTENTION.....	2
REVISION HISTORY .....	3
<b>1. Safety Information.....</b>	<b>5</b>
<b>2. Detector Installation .....</b>	<b>11</b>
<b>2.1. Wired Detector .....</b>	<b>11</b>
<b>2.1.1. Network Setting .....</b>	<b>11</b>
<b>2.1.2. Initial Connection Setting .....</b>	<b>12</b>
<b>2.1.3. Multi Connection Setting .....</b>	<b>13</b>
<b>2.2. Wireless Detector.....</b>	<b>14</b>
<b>2.2.1. The settings between Bridge (Repeater) and PC.....</b>	<b>14</b>
<b>2.2.2. Network adapter setting-Network Adaptor (Windows 7,8,10).....</b>	<b>18</b>
<b>2.2.3. Advanced Network Set-up .....</b>	<b>19</b>
<b>Appendix – Guidelines for Pediatric Subjects .....</b>	<b>23</b>
<b>Appendix - Check Vaccine programs in your Workstation PC.....</b>	<b>26</b>
<b>Appendix - Virtual Memory / DEP.....</b>	<b>27</b>
<b>Virtual Memory Windows 7 .....</b>	<b>27</b>
<b>DEP Setup in Windows 7 .....</b>	<b>29</b>
<b>DEP Setup in Windows 8.....</b>	<b>31</b>
<b>Appendix - How to Disable Driver Signature Verification .....</b>	<b>34</b>

## 1. Safety Information

### 1.1. Safety Information

PRUDENT passed all legal safety requirements as an electromedical equipment to guarantee the user's safety. Do not ignore following cautions while handling the products, and read thoroughly this page before use.



#### CAUTION!

- ◆ Do not use the instrument, if a malfunction has occurred, until the problems are solved by qualified personnel.
  - ◆ Do not install the instrument in a location with the conditions listed below.
    - Otherwise, it may result in failure or malfunction, fall or cause fire or injury.
    - Close to facilities where water is used.
    - Where it will be exposed to direct sunlight.
    - Close to heat source such as a heater.
    - Prone to vibration.
    - Insecure place.
    - Dusty environment.
    - Saline or sulfurous environment.
    - High temperature or humidity.
  - ◆ Do not use the instrument unless designated.
  - ◆ Do not touch any screws fixed in the instrument. Otherwise, loosened screws will result in the deterioration of image quality or the damaged instrument.
- |  |
|--|
| <ul style="list-style-type: none"><li>◆ Only authorized engineers from PIXXGEN are qualified for installation. Be sure to follow the instructions in this manual. Any inquiries related to the maintenance should be in touch with PIXXGEN Service Team at <a href="mailto:tech@pixxgen.com">tech@pixxgen.com</a>.</li><li>◆ Approach us if the instrument did not respond as shown in the instructions.</li></ul> |
| <ul style="list-style-type: none"><li>◆ Disclaimer<ul style="list-style-type: none"><li>- Manufacturer is not liable to accidents or breakdown caused by the use of Detector by legally unqualified personnel.</li><li>- The manufacturer is not liable to any accidents or technical problems caused by modification or repair of the device by agents or unspecified engineer.</li></ul></li></ul>               |

**CAUTION!****FCC Statement**

## Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

For product available in the USA/Canada market, only channel 1~11 can be operated. Selection of other channels is not possible.

This device is going to be operated in 5.15~5.25GHz frequency range, it is restricted in indoor environment only.

**CAUTION!****IMPORTANT NOTE:****FCC Radiation Exposure Statement:**

This device meets FCC limits for exposure to radio waves. It is designed not to exceed the limits for exposure to radio waves (radio frequency electromagnetic field) adopted by the

Federal Communications Commission.

This equipment should be installed and operated with minimum distance 0mm between the radiator & your body.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

**WARNING!**



















- ◆ To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth
- ◆ Do not touch signal input, signal output or other connectors, and the patient simultaneously.
- ◆ External equipment intended for connection to signal input, signal output or other connectors, shall comply with relevant IEC Standard (e.g., IEC60950 for IT equipment and IEC60601-1 series for medical electrical equipment). In addition, all such combination-system shall comply with the standard IEC60601-1 and/or IEC60601-1-1 harmonized national standard or the combination. If, in doubt, contact qualified technician or your local representative.
- ◆ No modification of this equipment is allowed
- ◆ Do not modify this equipment without authorization of the manufacturer.
- ◆ In the event of any serious accident with respect to the device, it shall be reported to the manufacturer and the relevant authority.
- ◆ Remove the battery if the ME EQUIPMENT is no likely to be used for some time.

**ELECTROMAGNETIC!**

- ◆ This equipment has been tested and found to comply with the limits for medical devices in EN 60601-1-2:2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
  - Reorient or relocate the receiving device.
  - Increase the separation between the equipment.
  - Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
  - Consult the manufacturer or field service technician for help.



## 1.2 Marks

	Attention, see instruction manual.
	The CE mark shows that the instrument obtained EU countries' requirements. CE number can be followed by the CE mark.
	Manufacturer's name and address.
	EU representative's name, address and contact details.
	Consult instructions for use
	Fragile, handle with care
	This symbol shall be accompanied by the manufacturer's serial number.
	Temperature limit.
	Symbol for the marking of electrical and electro-network equipment that must be recycled.
	Do not disassemble or open the instrument.
	Do not place the instrument near where liquid is present.
	Use a dry cloth only when cleaning the instrument.
	This way up
	Stacking Limited by number
	Keep away from rain
	Alternating Current
	Direct Current
	Device is switched on in order to bring it into the Stand-by condition

### 1.3 Condition for installation

- Do not install this equipment in any of the locations listed below.
  - Where the power supply is unstable.
  - Where temperature and humidity are high.
  - Where the room is without air-conditioner or ventilation.
  - Where it is exposed to direct sunlight.
- It is highly recommended to separate X-ray power in the distributing board of the electricity of the building.
- It is highly recommended to take internet connection nearby workstation pc in the room
- Check the distance between workstation PC in the operation room and the area which detector will be placed.
- The wall outlet or the circuit breaker shall be installed near the equipment and shall be easily accessible when problems occur.
- Turn off the detector after the operation.
  - For wired Detector: Press the stand-by button and check the power display (Green ->Red)
  - For wireless Detector: Press the Stand-by button and remove the battery pack. (Lift up the battery pack and remove it.)

### 1.4 Environment

- Be sure to use and store this equipment under the conditions described below.

	Temperature	Humidity
Storage	-20 to 70°C	10 to 95 % RH(Non-condensing)
Operation	10 to 35°C	20 to 75% RH(Non-condensing)

- Do not expose this equipment to high temperature and humidity since it makes detector malfunctions.

### 1.5 Etc.

- Type of protection against electric shock
  - Class I and/or Internally powered equipment.
- Degree of protection against electric shock
  - Not classified - No applied parts
- Classification according to the degree of protection against ingress of water
  - IPX0, ordinary equipment
- This equipment is not suitable for use in the presence of flammable anesthetics or oxygen
- Mode of operation:
  - continuous operation


## 2. Detector Installation

### 2.1. Wired Detector


#### 2.1.1. Network Setting

The settings between Prudent and Workstation


(1) Open the network and sharing center.



(2) Open the Network Adaptor properties



- (3) Open TCP/IPv4 and set the static IP as  
192.168.100.101~150. 192.168.100.111 is the wired network IP  
factory default of single detector.




## 2.1.2. Initial Connection Setting

Check **Show Network Window** message before connecting. Click **Connect** and then **IP and PORT**

window will appear.

### Select network window

This is to set up connecting ports between PIXX and the workstation computer.



**IP address:** Set the IP address for detector.

**Port number:** Set the PORT number for detector

### 2.1.3. Multi Connection Setting

A setting method for using more than two detectors.

192.168.100.111 is the wired network IP. (Factory default of single detector)


	Single	Dual	Triple
Detector IP	192.168.100.111	192.168.101.112	192.168.102.113
Port	7000	7000	7000
Network Adaptor IP	192.168.100.101	192.168.101.102	192.168.102.103

For changing Detector IP, it needs see  **Appendix Wired, Bridge, Direct, and Wi-Fi Router (DHCP server) setting for Multiple detector.**

## 2.2. Wireless Detector

### 2.2.1. The settings between Bridge (Repeater) and PC

(1) Open a web browser (e.g., Internet Explorer, Chrome, Firefox, or Safari) and enter <http://tplinklogin.net/> and then login to the system as an administrator (ID: admin / PW: admin).



#### CAUTION!

The recommended product for Bridge of Wi-Fi Router is TP Link Archer C7. In this instruction, it has been instructed based on Archer C7 product. The rest of Wi-Fi Router are able to be used. However, if you like to use it, additional test is required.

#### Tip!


When you set up Bridge at <http://tplinklogin.net/> and if it is not able to enter to set up page, disable all Network Adopter at Network Adopter Setting except the Network Adopter that is connected Bridge.

After completed the setting for Bridge, enable all the Network Adopter at Network Adopter Setting that you have disabled previously.


#### Tip!

When setting the Bridge (Repeater), set IP address as automatically in the TCP/IPv4

- (2) Go to the **Advanced > Network > LAN**, and change the IP Address to **192.168.120.2** in order to prevent IP conflicts and then click the **"Save"** button.



- (3) Go to the **Advanced > Wireless Settings**, and change the **"Wireless Network Name"** to **"PIXXAP"**. And select **"Enable WDS Bridging"** And then click the **"Survey"** button.






(4) Select "**PIXXAP120**" SSID from AP list, click the "**Choose**" button.

#### Survey

AP Number: 3

 Refresh

ID	MAC Address	SSID	Signal	Channel	Security	Operation
1	64-E5-99-63-CD-10	PILAIM5G		149	PSK	<a href="#">Choose</a>
2	C4-12-F5-6B-3A-56	PIAL501		36	PSK	<a href="#">Choose</a>
3	C0-25-E9-18-F2-58	PIXXAP120		36	PSK	<a href="#">Choose</a>

(5) Input the password "**1234567890**" and then click the "**Save**" button and reboot.

#### 5GHz WDS

WDS Bridging:

☒ Enable WDS Bridging

SSID (to be bridged):

PIXXAP120

[Survey](#)

MAC Address (to be bridged):

C0-25-E9-18-F2-58

Example: 00-1D-0F-11-22-33

WDS Mode:

Auto

Security:

☐ No Security

☒ WPA-PSK/WPA2-PSK

☐ WEP


Password:

1234567890


[Save](#)




(6) Open the network and sharing center.



(7) Open the Network Adaptor properties




(8) Open TCP/IPv4 and set the static IP as 192.168.120.3~254.



## 2.2.2. Network adapter setting-Network Adaptor (Windows 7,8,10)

(1) Installation

- Insert driver CD
- Run install file following each OS version



### 2.2.3. Advanced Network Set-up


This is to secure data streaming from disconnection.

#### CAUTION!


For the user using own network adapter, it is essential to update the driver of network adapter and adjust the advanced network setting.  
Realtek network adapter in particular is sensitive to set-up.  
This set-up is to protect the network from freezing or slowing down.  
Improper set-up may result in data loss or damaged images or intermittent connectivity.

(1) To check the information of current network adapter, proceed to **Control Panel > Network and Sharing Center > Change adapter settings**.


The network adapter connected to Prudent system appears as '**Unidentified Network**'.



- (2) Run **Unidentified Network** and proceed to **Properties > Composition**.



- (3) Proceed to **Driver** tab and check the version of current network adapter.




- (4) The old version of network adapter driver that is already installed on the computer may influence the speed of network system.


- (5) Update the network adapter driver. The information of network adapter driver is available on the manufacturers' web sites below.

<b>Tip!</b>	<a href="http://www.intel.com">http://www.intel.com</a> Intel <a href="http://www.realtek.com.tw/">http://www.realtek.com.tw/</a> Realtek <a href="http://www.broadcom.com/">http://www.broadcom.com/</a> Broad.com.
-------------	--

- (6) Uncheck all lists in the **Power Management** option.




- (7) Disable 'Wake on~' and 'WOL' options in the **Advanced** option.



(8) Change Receive Buffers value to **Maximum** value.

(Receive Buffers setting option of network adapter can be found in **Properties of Advanced/Performance Option.**)



## Appendix – Guidelines for Pediatric Subjects

### 1) Positioning the pediatric patient

Pediatric patients are not as likely as adults to understand the need to remain still during the procedure. Therefore it makes sense to provide aids to maintaining stable positioning. It is strongly recommended the use of immobilizing devices such as bean bags and restraint systems (foam wedges, adhesive tapes, etc.) to avoid the need of repeating exposures due to the movement of the pediatric patients. Whenever possible use techniques based on the lowest exposure times.

### 2) Shielding

We recommend you provide extra shielding of radiosensitive organs or tissues such as eyes, gonads and thyroid glands. Applying a correct collimation will help to protect the patient against excessive radiation as well. Please review the following scientific literature regarding pediatric radiosensitivity: GROSSMAN, Herman. "Radiation Protection in Diagnostic Radiography of Children". Pediatric Radiology, Vol. 51, (No.1): 141-144, January, 1973: <http://pediatrics.aappublications.org/cgi/reprint/51/1/141>.

### 3) Technique factors

You should take steps to reduce technique factors to the lowest possible levels consistent with good image acquisition. For example if your adult abdomen settings are: 70-85KVp, 200-400mA, 15-80mAs, consider starting at 65-75KVp, 100-160mA, 2.5-10mAs for a pediatric patient. Whenever possible use high KVp techniques and large SID (Source Image Distance).

Table 1 below can be used to estimate technique factors for various body builds. As the patient size increases, the KVp generally increases. Also depicted in Table 1 are the corresponding values of x-ray beam cross-sectional area and the estimated patient thickness in terms of water equivalence.

Table 1

Age	Head	Chest	Abdomen	Extremity (Forearm)
Newborn	67 KVp/2.0 mAs (110cm <sup>2</sup> /9.0cm)	60 KVp/2.0 mAs (140cm <sup>2</sup> /8.0cm)	66 KVp/2.0 mAs (200cm <sup>2</sup> /10cm)	N/A
1-yr-old	72 KVp/2.0 mAs (160cm <sup>2</sup> /12cm)	66 KVp/2.0 mAs (250cm <sup>2</sup> /9.0cm)	70 KVp/4.0 mAs (300cm <sup>2</sup> /13cm)	56 KVp/5.0 mAs (35cm <sup>2</sup> /1.8cm)
5-yr-old	75 KVp/2.0 mAs (210cm <sup>2</sup> /14cm)	70 KVp/2.0 mAs (430cm <sup>2</sup> /10cm)	72 KVp/5.0 mAs (540cm <sup>2</sup> /15cm)	60 KVp/5.0 mAs (84cm <sup>2</sup> /3.3cm)
10-yr-old	77 KVp/2.0 mAs (240cm <sup>2</sup> /15cm)	74 KVp/3.0 mAs (670cm <sup>2</sup> /13cm)	75 KVp/6.0 mAs (820cm <sup>2</sup> /17cm)	62 KVp/6 mAs (140cm <sup>2</sup> /5.0cm)
15-yr-old	79 KVp/2.0 mAs (270cm <sup>2</sup> /16cm)	78 KVp/4.0 mAs (780cm <sup>2</sup> /14cm)	78 KVp/7.0 mAs (900cm <sup>2</sup> /20cm)	65 KVp/6.0 mAs (200cm <sup>2</sup> /6.2cm)
Adult	75 KVp/15 mAs (320cm <sup>2</sup> /20cm)	120 KVp/2.0 mAs (1300cm <sup>2</sup> /15cm)	75 KVp/15 mAs (1200cm <sup>2</sup> /22cm)	65 KVp/8.0 mAs (200cm <sup>2</sup> /7.9cm)

## 4) Dosimetry

Table 2 summarizes the key dosimetry parameters for the four types of radiographic examination for patients ranging from newborn to the adult. In each cell, the first value is the entrance skin kerma(free-in-air) in  $\mu\text{Gy}$ . The second term gives the energy imparted to the patient, expressed in  $\mu\text{J}$ . In parentheses on the second line are the corresponding values of patient effective dose in  $\mu\text{Sv}$ .

Table 2

Age	Head 100 $\mu\text{Gy}/78.2\mu\text{J}$ (10 $\mu\text{Sv}$ )	Newborn Chest 77 $\mu\text{Gy}/66\mu\text{J}$ (19 $\mu\text{Sv}$ )	Abdomen 100 $\mu\text{Gy}/140\mu\text{J}$ (62 $\mu\text{Sv}$ )	Extremity (Forearm) NA
1-yr-old	120 $\mu\text{Gy}/165\mu\text{J}$ (7.3 $\mu\text{Sv}$ )	96 $\mu\text{Gy}/160\mu\text{J}$ (16 $\mu\text{Sv}$ )	230 $\mu\text{Gy}/580\mu\text{J}$ (90 $\mu\text{Sv}$ )	130 $\mu\text{Gy}/9.5\mu\text{J}$ (0.21 $\mu\text{Sv}$ )
5-yr-old	140 $\mu\text{Gy}/260\mu\text{J}$ (5.9 $\mu\text{Sv}$ )	110 $\mu\text{Gy}/340\mu\text{J}$ (18 $\mu\text{Sv}$ )	320 $\mu\text{Gy}/1500\mu\text{J}$ (120 $\mu\text{Sv}$ )	160 $\mu\text{Gy}/44\mu\text{J}$ (0.5 $\mu\text{Sv}$ )
10-yr-old	150 $\mu\text{Gy}/320\mu\text{J}$ (4.3 $\mu\text{Sv}$ )	190 $\mu\text{Gy}/1100\mu\text{J}$ (33 $\mu\text{Sv}$ )	420 $\mu\text{Gy}/3300\mu\text{J}$ (160 $\mu\text{Sv}$ )	200 $\mu\text{Gy}/130\mu\text{J}$ (0.87 $\mu\text{Sv}$ )
15-yr-old	150 $\mu\text{Gy}/400\mu\text{J}$ (3.1 $\mu\text{Sv}$ )	280 $\mu\text{Gy}/2100\mu\text{J}$ (36 $\mu\text{Sv}$ )	550 $\mu\text{Gy}/5100\mu\text{J}$ (140 $\mu\text{Sv}$ )	220 $\mu\text{Gy}/240\mu\text{J}$ (0.92 $\mu\text{Sv}$ )
Adult	1100 $\mu\text{Gy}/3200\mu\text{J}$ (19 $\mu\text{Sv}$ )	150 $\mu\text{Gy}/2500\mu\text{J}$ (34 $\mu\text{Sv}$ )	1100 $\mu\text{Gy}/13000\mu\text{J}$ (290 $\mu\text{Sv}$ )	300 $\mu\text{Gy}/360\mu\text{J}$ (1.1 $\mu\text{Sv}$ )

## 5) Other References for Pediatric Dosimetry

1. Size measurements are based on approximate mean values (averaged across males and females) from: McDowell, M.A., C.D. Fryar, C.L. Ogden, and K.M. Flegal. 2008. Anthropomorphic Reference Data for Children and Adults, United States, 2003-2006. National Health Statistics Reports, 10, 1-48. Available for download at: <http://www.cdc.gov/nchs/data/nhsr/hsr010.pdf>. The weight given for the neonate subgroup is lower than the average to ensure that a broad range of sizes is adequately covered.

2. These suggested subgroups fall within the age groups identified in the guidance entitled "premarket Assessment of Pediatric Medical Devices"

(<http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/UCM089742.PDF>): neonate(birth-1 month), infant(1 month-2 yrs.) child(2-12 yrs), adolescent(12-21 yrs). For design and evaluation of radiological devices, patient size(i.e. height, weight, thickness) is a better indicator.

3. Based on reports of the lifetime incidence of cancer vs. age of exposure data, these pediatric subgroups defined by the Agency cover the region where the largest age dependence is expected for cancer risk. Risk decreases much less steeply as a function of age for individuals over 21 years old. [See NAS National Research Council Committee to Assess Health Risks from Exposure to Low



Levels of Ionizing Radiation.2006. Health risks from exposure to low levels of ionizing radiation: BEIR VII phase 2. Washington, D.C.: National Academy of Sciences, National Academies Press.]

4. The following reference gives current data for anteroposterior and transverse body diameter for pediatric patients ranging in age from 0.5 to 20 years: Kleinman, P.L.,K.J. Strauss, D. Zurakowski, K.S. Buckley, and G.A. Taylor. 2010. Patient size measured as a function of age at a tertiary care children's hospital. American Journal of Roentgenology, 194, 1611-1619.

5. The following reference used cylindrical phantoms with diameters of 8,16,24, and 32cm to represent a neonate, 5year old, 12 year old, and adult patient respectively: Siegel, M. J., et al. 2004. Radiation dose and image quality in pediatric CT: effect of technical factors and phantom size and shape. Radiology, 233(2), 515-522.

#### 6. Summary


- ✓ Image only when there is a clear medical benefit.
- ✓ Image only the indicated area.
- ✓ Use the lowest amount of radiation for adequate imaging based on size of the child(reducing tube output- KVp and mAs)
- ✓ Try to use always short exposure times, large SID values and immobilizing devices.
- ✓ Avoid multiple scans and use alternative diagnostic studies(such a ultrasound or MRI) when possible.

## Appendix - Check Vaccine programs in your Workstation PC.

Some vaccine programs might block detector signal since they recognize our software as a virus or malware.




We strongly recommend removing vaccine programs except MS Essential.


When using MS Essential, please set up "Excluded files and locations" as shown below.




## Appendix - Virtual Memory / DEP


### Virtual Memory Windows 7










### DEP Setup in Windows 7






## DEP Setup in Windows 8


- This function may not work exactly.
- There are CMOS disturbs this function to use, at that case you need to set DEP release in OS.

(1). System

Control Panel → System and Security → System




(2). Click the Advanced system Setting and then the System Properties window will be shown. And then move the Advanced Tab and click the Settings button of Performance






(3). When the Performance Options window show, move the Data Execution Prevention Tab and then select 'Turn on DEP for all programs and services except those I select'




(4). Click the 'Add...' button and select two files as like below

ZView.exe




WBMaker.exe





(5). Confirm the 'ZView' and 'WBMaker' was included in List and then click 'Apply' button




## Appendix - How to Disable Driver Signature Verification


Press the Windows key + C keyboard combination to bring up the Charms Bar, and then click on the Settings Charm.




We need to head into the Modern Control Panel, so go ahead and click on the Change PC settings link.




When the Control Panel opens, switch over to the "Update & recovery" section.




Then click on the Recovery option on the left hand side.




Once selected, you will see an advanced startup section appear on the right hand side. You will need to click on the "Restart now" button.




Once your Computer has rebooted you will need to choose the Troubleshoot option.




Then head into Advanced options.




Then Startup Settings.



Since we are modifying boot time configuration settings, you will need to restart your Computer one last time.




Finally, you will be given a list of startup settings that you can change. The one we are looking for is "Disable driver signature enforcement". To choose the setting, you will need to press the F7 key.



That's all there is to it. Your PC will then reboot and you will be able to install unsigned drivers without any error message.



PIXXGEN Corporation  
5F, SMARTBAY, 123, Beolmal-ro, Dongan-gu  
Anyang-si Gyeonggi-do, 431-804 Korea  
Tel: +82-70-4846-8888 FAX: +82-2-6455-2905  
Web: <http://www.pixxgen.com>



CMC Medical Devices & Drugs S.L.  
C/Horacio Lengo N° 18, CP 29006, Málaga, Spain.