



Pocket Portable
Nebulizer
User Guide

Model fn2000m

✓TM Fill

✓TM Click

✓TM Breathe

**USER: READ THIS GUIDE BEFORE OPERATING THIS DEVICE.
SAVE THIS GUIDE FOR FUTURE REFERENCE.**





SYMBOLS

The symbols that appear on the box and the Flyp™ nebulizer are described below:

SYMBOL	MEANING
	Caution
	Consult instructions for use
	Keep dry
	Do not place in dishwasher

2

	Transient storage temperature limits -10°C to 45°C (14°F to 113°F)
	Federal law (USA) restricts this device to sale by or on the order of a physician or licensed practitioner.
	Type BF applied part
	FCC Approved
	Protected against insertion of fingers and will not become damaged or unsafe during a test in which it is exposed to vertically dripping water when held at an angle.
	Bluetooth enabled technology

3



This page left blank intentionally.



Important information
to prevent damage to your
Flyp™ nebulizer

Important safety information
regarding hazards that may cause
personal injury

Before using your Flyp nebulizer for the first time, be aware of all warnings and safety information. Only use accessories approved by the manufacturer and referenced within this manual. If you do not fully understand all the warnings, safety precautions, and operating instructions, contact the Flyp customer service team at 844.FLYPNEB (359.7632) for technical support. ALWAYS KEEP THIS USER GUIDE HANDY.



TABLE OF CONTENTS

Section 1: Introduction	
Intended Use	8-9
Section 2: Safety Guidelines	
Safety Guidelines	10-15
Section 3: Product Description	
Features and Benefits.....	16-17
Section 4: Flyp™ at a Glance	
Names and Functions of Parts	18-20
Section 5: Using Flyp	
About Flyp's Battery	21-23
Filling Medication Reservoir with Prescription.....	24-25
Turning Flyp ON and OFF	26-27
Inhaling Prescribed Medication.....	28

Section 6: Cleaning and Disinfecting	
Cleaning Parts.....	29-33
Disinfecting Parts.....	34-35
Section 7: Troubleshooting	
Troubleshooting.....	36-37
Section 8: Support	
Learning More.....	38
Section 9: Information	
Technical Data.....	39
Specifications	39-40
Aerosol Performance	41-42
Electromagnetic Compatibility.....	43-52
Disposal and Recycling	53
Warranty	53
Other information: FCC and bluetooth	54-57
Section 10: Replacement Parts and Accessories.....	58



Section 1: INTRODUCTION

INTENDED USE

Indications for use:

The Flyp™ nebulizer, for use by adolescent and adult patients, is intended to aerosolize healthcare provider-prescribed solutions for inhalation that are approved for use with a general-purpose nebulizer. Flyp is intended for use at home or a medical facility, such as a hospital or doctor's office.

Intended user:

Flyp nebulizer, for use by adolescent and adult patients, should only be used by a patient under the supervision of a qualified medical expert, such as a doctor, nurse, or respiratory therapist. The user, or their actively assisting caregiver, should be capable of understanding all of the User Guide's contents. Flyp is intended for use by a single user. Flyp is not a life-saving device. Patients who are in severe respiratory distress, who are unconscious, or who

are not breathing spontaneously should not use this device.

In an emergency, call 911 immediately for medical assistance.

Caution:

Federal law (USA) restricts this device to sale by or on the order of a physician or licensed practitioner.

Precautions:

All warnings and cautions described in the User Guide should be observed.

Service life:

The Flyp main unit has an expected service life of 36 months.

***Note:** The Medication Reservoir (HypersonIQ™ Cartridge) should be replaced as needed, or every 3 months.

***Note:** Charge only with Flyp provided charger.

DO NOT Alter or modify this device.

***Note:** The mouthpiece of this device is removable.



Section 2: SAFETY GUIDELINES

SAFETY GUIDELINES

Read this section to learn how to use Flyp™ safely and correctly and to prevent risks and injuries to you and others.

Keep this User Guide handy for future reference.

WARNING

Failure to follow these instructions could result in serious injury or damage to the device or other property. Read all the safety

information below before using the device.

Medical device:

Flyp is a medical device, available only by prescription and for prescribed medications.

Be sure to follow a qualified medical expert's instructions.

Direct exhaled medication away from others.

If you are using Flyp to treat a serious condition, a back-up device is recommended.

However, in a life-threatening situation, **NEVER** rely on a nebulizer. Call 911 immediately for emergency medical assistance.

Use by others:

For single patient use only.

If others use it, infection may spread.

Medical use:

Flyp is intended solely to deliver prescribed medication to treat a respiratory condition. Using Flyp for any other purpose is dangerous, and neither the distributor, manufacturer, nor

their affiliates can be held liable for any damage or injury caused by improper use or misuse.

Supervision required:

Adult supervision is required when the device is used by children and the infirm. Small pieces may present a potential choking hazard. If a piece is accidentally swallowed always consult a doctor immediately. You should also be aware that the USB Wall Charger presents a potential choking hazard.



Cleaning:

Be sure to clean all components before using them for the first time and after each subsequent use (refer to pages 29-32).

The Main Unit should be cleaned every day, but it should **NEVER** be placed in water or washed in a dishwashing machine. **ALWAYS** disconnect and unplug the USB Wall Charger before cleaning the Main Unit.

Do not drop medication on the Main Unit or into its USB port. If you drop medication on either area, immediately wipe it off with clean dry cloth.

Battery:

The battery is not replaceable. Don't attempt to replace the battery yourself, as you may damage it. Overheating and injury could result. The lithium-ion battery must be recycled or disposed of separately from household waste.

NEVER incinerate the battery. Incineration may cause the battery to rupture. If an ignition source exists, then fire and even an explosion could result. **NEVER** immerse the battery in water, as this may cause the battery to rupture.

Electronic device:

Flyp complies with all applicable electromagnetic compatibility (EMC) standards. You should, however, avoid operating it near other electronic devices.

If you are not going to use the unit for a long period of time, disconnect the USB Wall Charger.

Charging:

Flyp contains an internal, lithium-ion rechargeable battery that cannot be removed.

Charge Flyp only with the USB Wall Charger provided. **NEVER** use the device when it is charging.

All batteries deteriorate over time if they are not used or charged. Do not store Flyp for long periods of time without charging it periodically. If Flyp is not used or charged for a long period of time, the battery may create a hazardous condition.

Proper handling:

Flyp contains sensitive components. Do not drop, crush, puncture, bend, heat, incinerate, or apply strong shock to the device or its parts.

• **NEVER** unscrew or open the Main Unit.



- Do not disassemble, repair, or modify Flyp. You may injure yourself or render the device ineffective.
- Keep the device away from heated surfaces and extreme heat and cold. Do not leave Flyp in a car if it will be subjected to significant heat or cold.
- Do not expose the device to direct sunlight for extended periods of time. Keep the device away from children, pets, pests, lint and dust. You may render the device ineffective.

- Flyp is intended for use at home or a medical facility, such as a hospital or doctor's office.

If you're concerned about dirtying Flyp™, you can carry it in the Draw-String Bag.

Using cords and ports:

NEVER force the USB Wall Charger into the USB Port. Check for obstructions in the USB Port. If the USB Wall Charger and USB Port don't

CAUTION

join with reasonable ease, they probably don't match. Make sure that you are using the USB Wall Charger provided and that you have positioned it correctly.

Disposing of Flyp properly:

For information about the proper disposal of Flyp, including other important regulatory compliance information, see "Disposal and Recycling" on page 53.



Section 3: PRODUCT DESCRIPTION

FEATURES AND BENEFITS

Read this section to learn about Flyp's features and benefits.

Small and convenient:

Flyp™ fits in your pocket. There is no separate control unit, compressor, mask, hose or cup. And Flyp's convenient Draw-String Bag makes it easy to carry with you.

Flyp is rechargeable with a USB Wall Charger, like a cell phone. There are no disposable batteries.

Simple:

Flyp is designed for easy operation and delivers all medications approved for use with general-purpose nebulizers.

Quiet and Virtually silent:

Barely a whisper thanks to HypersoniQ™ Technology. The internal disk vibrates 2000x faster than a hummingbird's wings.

Fast and efficient:

The Flyp nebulizer delivers a standard, 3ml dose of medication in about 7 minutes. So you can get back to doing what love.

Cleaning function:

To active this mode, press and hold the power button for 3 seconds. Flashing blue indicator lights appear. Flyp will produce a pulsing aerosol. This special pulse action helps to loosen dried particles that may reside on the disk. Operate the unit for 1-2 minutes or until the reservoir is empty.




Section 4: FLYP™ AT A GLANCE

NAMES AND FUNCTIONS OF PARTS


Read this section to learn the names and functions of the parts.

If any items are missing, contact the store where you purchased Flyp. You may also contact Flyp nebulizer customer service at 1-844-FLYPNEB




NEBULIZER DISK

ON/OFF
BUTTON
& INDICATOR
LIGHT






USB WALL CHARGER



Note: the Draw-String Bag is not intended to protect Flyp™. It is provided for convenience only.

DRAW-STRING BAG



Section 5: Using Flyp™

ABOUT Flyp's BATTERY

Read this section to learn how to use Flyp correctly.

About Flyp's battery:

Flyp is charged through a USB port, just as many cellular phones and portable electronic media devices are charged.

Flyp has an internal battery that you cannot replace.

Charging Flyp's battery:

You can charge Flyp's battery by connecting Flyp™ to a wall outlet using the USB Wall Charger provided. Only use the charger provided.


NEVER use the device while it is charging.

***Note:** Flyp is designed to shut off automatically after 10 minutes to conserve battery life.



To charge the battery with a wall outlet:

1. Open the USB Cover to reveal the USB Port.
2. Gently insert the USB Wall Charger's Cord into the USB Port.
3. Insert the USB Wall Charger into a wall outlet.




Connection
via wall outlet

INDICATOR LIGHT	= STATUS
No Light	Power off
Solid BLUE light	Power on/Treatment Mode
Flashing BLUE light	Cleaning Mode
Flashing Yellow & Red light	HypersoniQ™ Cartridge not properly attached
Flashing Yellow light (plugged-in)	Battery is charging
Solid Yellow light (plugged-in)	Battery is fully charged
Flashing Yellow light	Low Battery
Solid Blue light & flashing Yellow light	Power on, Low battery
Flashing Blue and Yellow light	Cleaning, Low Battery
Flashing Red light	Battery requires charging



FILLING MEDICATION RESERVOIR WITH MEDICATION

1. There is no need to remove the Medication Reservoir to fill it with medication.




STEP 1

24

2. Insert the medication ampule into the Medication Reservoir.

Squeeze to empty contents and fill reservoir.



STEP 2

3. Close the Reservoir. Be careful not to touch any part of the Stopper that may come in contact with your medication. Then close the Reservoir Cover.

Note: the Medication Reservoir's maximum capacity is 6mL.



STEP 3

25





TURNING FLYP ON & OFF

Gently raise the Mouthpiece using your index finger, revealing the On-Off Button and On-Off Indicator Light.

To turn Flyp™ on, press the On-Off Button. The blue On-Off Indicator Light will light. Visually confirm that an aerosol mist is flowing from the Mouthpiece's end.

To turn Flyp off, press the On-Off Button again. The blue On-Off Indicator Light will no longer be lit.


Flyp will automatically turn off after 10 minutes.





INHALING PRESCRIBED MEDICATION

Place the Mouthpiece between your lips. Inhale and breathe in a calm manner at a normal rate.



Section 6: CLEANING & DISINFECTING

CLEANING PARTS

Read this section to learn how to clean Flyp™.

***Note:** Failure to clean after every use can impact nebulization time and device performance.


ALWAYS disconnect and unplug the USB Wall Charger before cleaning or disinfecting

If, after proper cleaning and charging, Flyp takes longer than 10 minutes to deliver either 1 unit-dose of medication or 3 mL of normal saline (0.9%), the HypersonIQ™ cartridge should be cleaned or replaced. Consult "Troubleshooting" on pages 36-37 to see if corrective steps can be taken.



Cleaning the Mouthpiece and Medication Reservoir:

Remove the Mouthpiece by pulling on it gently. To remove the Medication Reservoir, take hold on each side and lift up.



Prepare a solution of warm water and mild liquid dishwashing soap. Wash the Mouthpiece in the soap solution. Rinse the Mouthpiece and Medication Reservoir under running tap water. These are the only parts that can be submerged in water. Shake off excess water and dry on clean, dry towel.

NEVER rinse the Main Unit under running water.

Cleaning the Main Unit's exterior:

Prior to cleaning the Main Unit's exterior, disconnect and unplug the USB Wall Charger, make sure the USB Cover is closed tightly, and remove the Mouthpiece. Prepare

a solution of warm water and mild dishwashing soap. Then clean the Main Unit by wiping it with dry clean cloth lightly moistened by the soap solution. Use a dry clean cloth moistened with warm tap water to remove any soap residue. To dry the Main Unit, wipe it with a dry clean cloth. Make sure it is completely dry before reinserting the Mouthpiece and storing it.

NEVER place the Main Unit or USB Wall Charger in water or hold them under running water. The lithium-ion battery inside the Main Unit may be dangerous when wet.



CLEANING *HypersoniQ™* Cartridge

1. Open the stopper on the HypersoniQ cartridge (Medication Reservoir).
2. Fill the reservoir half-way with ~3ml distilled water. Close the stopper.
3. Gently raise the Mouthpiece with your index finger.
4. **Press and hold the power button for 3 seconds.** Flashing blue indicator light will appear. Flyp will produce a pulsing aerosol.
5. Allow the aerosol mist to flow from the raised Mouthpiece for 1-2 minutes or until empty.

6. Power Flyp off by pressing the power button.
7. Remove HypersoniQ cartridge and shake out excess water.
8. Let the main unit, cartridge, and mouthpiece dry on a clean, dry towel.
9. Allow the Medication Reservoir interior to dry fully prior to using again.

It is important to clean the HypersoniQ cartridge after each use in order to prevent the buildup of residual medication. If these steps are not followed, the disk may become blocked – impacting performance. Please see troubleshooting on page 36 in case of a suspected blockage.

DEEP CLEANING

Residual medication can accumulate over time and cause buildup on the disk which may slow aerosol production. To prevent this, a deep cleaning with distilled white vinegar (5%) or Flyp cleaning solution on a weekly (or more frequently as-needed) basis is recommended.


1. Check the HypersoniQ cartridge to make sure it is empty.
2. Open the stopper on the HypersoniQ cartridge (Medication Reservoir).
3. Fill the Medication Reservoir approximately half-way with ~3ml distilled white vinegar (or Flyp cleaning solution). Close the stopper.
4. Gently raise the Mouthpiece with your index finger.

5. **Press and hold the power button for 3 seconds.** Flashing blue indicator light will appear. Flyp will produce a pulsing aerosol.
6. Allow the aerosol mist to flow from the raised Mouthpiece for 1-2 minutes or until empty.
7. Power Flyp off by pressing the power button.
8. Remove HypersoniQ cartridge and Mouthpiece. Rinse both under running water, ensuring the stopper of the HypersoniQ cartridge is in the open position.
9. Let the cartridge and mouthpiece dry on a clean, dry towel.
10. Allow the interior to dry fully prior to using again.



OPTIONAL DISINFECTION

Read this section to learn how to disinfect Flyp. Although disinfection is optional, it is recommended that you do it at the end of every day.



DISINFECTING HYPERSONIQ CARTRIDGE (MEDICATION RESERVOIR) & MOUTHPIECE:

Step 1. Remove the cartridge and mouthpiece from Flyp. Remember to open stopper on cartridge.

Step 2. Soak the mouthpiece and cartridge in 70% ethyl alcohol for 10 minutes. Make sure both are completely submerged in the alcohol.

Step 3. Remove mouthpiece and cartridge from the alcohol bath and rinse under running water.

Step 4. Let dry on a clean, dry towel.

DISINFECTING MAIN UNIT'S EXTERIOR:

Step 1. Wipe exterior with clean, dry towel moistened with ethyl alcohol.

Step 2. Allow Flyp to sit for 5 minutes.

Step 3. Moisten with a clean dry towel with water and wipe exterior to remove the ethyl alcohol.

Step 4. Let dry on a clean, dry towel.



Section 7: TROUBLESHOOTING

TROUBLESHOOTING

Most problems with the device can be solved by following this advice.

The following table describes possible troubles, their causes, and what corrective action to take.

If these corrective actions do not return Flyp™ to full working order, read the next section, "Support."

TROUBLE	POSSIBLE CAUSE	CORRECTIVE ACTION
The nebulization rate is extremely low.	The battery is low.	Charge Flyp. Refer to pages 22-23.
	An air bubble has formed in the medication's pathway.	Tap the back of the Main Unit with your index finger until normal flow resumes.
	The Nebulizer Disk has become clogged with medication.	See deep cleaning instructions. Refer to page 33.
The On-Off Indicator Light is on, but there is no mist.	The Medication Reservoir is not correctly inserted.	Remove the Medication Reservoir and install it correctly.
	The Medication Reservoir is not filled.	Fill with medication. Refer to pages 24-25.
	Check that Flyp has been properly cleaned.	See deep cleaning instructions. Refer to page 33.



Section 8: SUPPORT

You can find more information about using Flyp™ on our website. To learn about service and support, and to view tutorials, go to: www.flypnebulier.com


NOTES:

Section 9: INFORMATION

TECHNICAL DATA

Technical Specifications:

Length: 119 mm
Width: 54 mm
Depth: 27 mm
Weight: 102 g
(with battery)



SPECIFICATIONS

Indications for use:

The Flyp™ nebulizer, for use by adolescent and adult patients, is intended to aerosolize healthcare provider-prescribed solutions for inhalation that are approved for use with a general-purpose nebulizer. Flyp is intended for use at home or a medical facility, such as a hospital or doctor's office.

**Power supply:**

Model 1VP0500-1000U

Rated 100-240V, 50/60hz, 0.5A

Patient population:

Adolescent to adult

Method of operation:

Piezoelectric/Ultrasonic

Power source:

Lithium-ion battery (rechargeable)

Weight:

Flyp: 102 g (with battery)

Dimensions:

27 x 119 x 54 mm (Main Unit)

Sound level:

<35 dBA at 1 meter

Operating conditions:Temperature range: 5°C to +40°C
(+41°F to +104°F)

Humidity: up to 95% RH

Storage conditions:Temperature range: -10°C to
+45°C (14°F to +113°F)

Humidity: up to 93% RH

Aerosol Performance Compared to 510K Testing

	Albuterol sulfate 1.0mg/ml @2.5mL	Ipratropium bromide 0.2mg/mL @2.5ml	Cromolyn sodium 1.0mg/ml @2mL
MMAD (μ)	2.31 +/- 0.10	2.17 +/- 0.05	2.26 +/- 0.11
GSD (geometric standard deviation)	1.60 +/- 0.07	1.48 +/- 0.03	1.54 +/- 0.06
Total delivered dose (μ)*	1038.8 +/- 54.6	230.85 +/- 19.25	958.45 +/- 61.45
Total respirable dose (μ)	921.6 +/- 70.3	209.4 +/- 14.9	857.7 +/- 84.4
Respirable fraction (% mass 0.4-4.7 μ)	87.8%	90.3%	88.1%

* Sum of Coarse (74.7 μ), Respirable dose (0.4 – 4.7 μ) + Fine (<0.4 μ).

All data shown at a 95% confidence level.

Medications were tested 3 times each on 3 devices for a total of 27 sample data points.



Aerosol Performance Compared to 510K Testing			
	Albuterol sulfate	Ipratropium bromide	Cromolyn sodium
Coarse particle dose (μ)	69.2 +/- 36.1	7.4 +/- 3.8	85.7 +/- 37.2
Respirable particle dose (μ)	921.6 +/- 74.2	209.4 +/- 14.9	857.7 +/- 84.4
Fine particle dose (μ)	48.0 +/- 20.9	14.0 +/- 5.7	15.0 +/- 7.4
Coarse fraction ($>4.7\mu$)	6.7%	3.1%	8.9%
Respirable fraction ($<4.7\mu$)	87.8%	90.3%	88.1%
Fine fraction ($<1.0\mu$)	4.7%	6%	1.6%

All data shown at a 95% confidence level. Medications were tested 3 times each on 3 devices for a total of 27 sample data points.

ELECTROMAGNETIC COMPATIBILITY

With the increased number of electronic devices, such as computers and cellular phones, medical devices in use may be susceptible to electromagnetic interference. This interference may result in incorrect operation of the medical device and create a potentially unsafe situation. Do not use cellular phones and other electronic devices that generate electromagnetic fields near the medical device.

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this User Guide.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The use of accessories, transducers and/or cables other than those specified, with the exception of those sold by the manufacturer as replacement



parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

EMC Tables. The following tables are provided in accordance with IEC 60601-1-2: 2014, ed. 4.0.

Guidance and Manufacturer's Declaration - Emissions

The Flyp™ fn2000m is intended for use in the electromagnetic environment specified below. The customer or user of the Flyp fn2000m should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
Radiated Emissions: CISPR 11	Class B	
Conducted Emissions: CISPR 11	Class B	



EMC Tables. The following tables are provided in accordance with IEC 60601-1-2: 2014, ed. 4.0.

Guidance and Manufacturer's Declaration - Emissions

The Flyp™ fn2000m is intended for use in the electromagnetic environment specified below. The customer or user of the Flyp fn2000m should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
Harmonics: IEC61000-3-2	Not Applicable	
Flicker: IEC61000-3-3	Not Applicable	The Flyp fn2000m is suitable for use in all establishments, including domestic, and those directly connected the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacturer's Declaration – Immunity

The Flyp fn2000m is intended for use in the electromagnetic environment specified below. The customer or user of the Flyp fn2000m should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD: IEC61000-4-2	±8kV Contact ±2kV, ±4kV, ±8kV, ±15kV Air	±8kV Contact ±2kV, ±4kV, ±8kV, ±15kV Air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be 30%.
EFT: IEC61000-4-4	+/- 0.5 KV L1 and L2 +/- 1 KV L1 and L2	+/- 0.5 KV L1 and L2 +/- 1 KV L1 and L2	Mains power quality should be that of a typical commercial or hospital environment.
Surges: IEC61000-4-5	+/- 0.5 KV L1 and L2 +/- 1 KV L1 and L2	+/- 0.5 KV L1 and L2 +/- 1 KV L1 and L2	Mains power quality should be that of a typical commercial or hospital environment.



Guidance and Manufacturer's Declaration – Immunity

The Flyp™ fn2000m is intended for use in the electromagnetic environment specified below. The customer or user of the Flyp fn2000m should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Voltage dips and interruptions: IEC61000-4-11	0 % UT: 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT: 1 cycle 70% UT: 25 cycles for 50 Hz Single phase: at 0° 0 % UT; 250/300 cycle for 50 Hz and 60 Hz	0 % UT: 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT: 1 cycle 70% UT: 25 cycles for 50 Hz Single phase: at 0° 0 % UT; 250/300 cycle for 50 Hz and 60 Hz	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Flyp fn2000m requires continued charging during mains power interruptions, it is recommended that the Flyp fn2000m be powered from an uninterruptible power supply or battery.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Power Frequency 50/60Hz Magnetic Field: IEC61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.



Guidance and Manufacturer's Declaration – Immunity

The Flyp™ fn2000m is intended for use in the electromagnetic environment specified below. The customer or user of the Flyp fn2000m should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level
Conducted RF: IEC61000-4-4	3 Vrms 0.15-80 MHz 6 Vrms in ISM and Amateur radio bands	3 Vrms 0.15-80 MHz 6 Vrms in ISM and Amateur radio bands
Radiated RF: IEC61000-4-3	10 V/m 80 MHz - 2.7 GHz 1 kHz @ 80% AM	10 V/m 80 MHz - 2.7 GHz 1 kHz @ 80% AM

Electromagnetic Environment – Guidance

Portable and mobile communications equipment should be separated from the Flyp fn2000m by no less than the distances calculated/listed below:

$$D = (3.5/V1)(\text{Sqrt } P)$$

150kHz to 80 MHz

where P is the max. power in watts and D is the recommended separation distance in meters.

$$D = (3.5/E1)(\text{Sqrt } P)$$

80 to 800 MHz

Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1).

$$D = (7/E1)(\text{Sqrt } P)$$

800 MHz to 2.5 GHz

Interference may occur in the vicinity of equipment containing a transmitter.



Recommended Separation Distances for the Flyp fn2000m

The Flyp™ fn2000m is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the Flyp fn2000m can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the Flyp fn2000m as recommended below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $D=(3.5/V1)(\text{Sqrt } P)$	Separation (m) 80 to 800MHz $D=(3.5/E1)(\text{Sqrt } P)$	Separation (m) 800MHz to 2.5GHz $D=(7/E1)(\text{Sqrt } P)$
0.01	0.116667	0.116667	0.233333
0.1	0.368932	0.368932	0.737865
1	1.166667	1.166667	2.333333
10	3.689324	3.689324	7.378648

DISPOSAL & RECYCLING

CAUTION

Your Flyp must be disposed of properly, according to local laws and regulations. Because this product contains a battery, the product must be disposed of separately from household waste. When your Flyp reaches its end of life, contact your local authorities to learn about recycling options.

WARRANTY

The Flyp main unit is warrantied for a period of thirty-six (36) months from the date of purchase against defects in manufacturing. All warranties are based on typical usage and following care instructions.

We will, at our option, repair or replace without charge your device. Repair or replacement is our only responsibility, and your only remedy, under the warranty.

HypersoniQ™ Cartridge is warrantied for a period of three (3) months from the date of purchase against defects in manufacturing.



OTHER INFORMATION

Federal Communication Commission (FCC) Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and

can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: Do

not use cellular phones and other electronic devices that generate electromagnetic fields near the medical device.

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

- Consult the dealer or an experienced radio/TV technician for help.



This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

BLUETOOTH SPECIFICATIONS	FLYP NEBULIZER
Bluetooth Compliance	Version 5.0 Single mode low energy
Operating Frequency	2.360 Ghz to 2.500 Ghz
Max Output Power	TX: +4 Dbm Max
Operating Range	200 meter (maximum)
Network Topology	Star
Operation	BLE peripheral device running GATT Server
Antenna Type	Integrated
Modulation Technology	FHSS (Frequency Hopping Spread Spectrum)
Modulation Type	GFSK
Data Rate	1 Mbit/second
Data Latency	400 ms



BLUETOOTH CON'T	FLYP NEBULIZER
Data Integrity	Adaptive Frequency hopping, 24-bit CRC. All packets have a 24-bit CRC value calculated for them by the transmitter and appended to the packet. The receiver recalculates the CRC and compares the calculated value with the value appended to the packet. If they are not the same, the receiver will not acknowledge that packet, resulting in a retransmission of that packet.
Format	Sends clock time every minute. Send treatment notification and time stamp on button press. Sends data logs when requested from mobile App.
Quality of Service	This device uses Bluetooth Smart technology for wireless communications, which allows for reliable communications in electrically noisy environments. If the connection is lost, the device starts to advertise immediately and the App can reconnect.
Bluetooth Supported Profiles	GATT based custom profile for Flyp
Authentication & Encryption	Supported
Encryption Key Size	128-bit AES HW

Section 10: REPLACEMENT PARTS AND ACCESSORIES

PART	REORDER NUMBER
Mouthpiece	MP-2000
HypersonIQ Replacement Head	hs2000m
USB Wall Charger	DC-1000
Mask Adapter	MA-2000
Draw-String Bag	AIR-4000
User Guide	AIR-5000
Cleaning Kit	CB-2000



85 Willow Street,
New Haven, CT, 06511
(844) 359-7632
www.flypnebulizer.com

v2.13

