


Smart Visual Ear Cleaner

Model:YCCEQ

UserManual

Software


1 Scan code to download app




Process

2 Remove the lid, start automatically


3 Set WIFI




4 Connect YCEQ-XXXX




5 Return to the APP interface, Click the icon



6 If you get this prompt, do not switch WIFI



7 You can record as screen or take a picture in the user interface



Troubleshooting

Malfunction	Solution
Unclear	Wipe the lens
Cannot be turned on	Fully charged
Flashback	Clear memory
light isn't working	Fully charged
Failed to display images after successful connection	Disconnect WIFI, restart the machine, and reconnect

If the above operations do not solve the problem, please contact the online customer service or dealer, we can provide further solutions

Product parameter

Product:	Smart visual ear cleaner
Model:	YCCEQ
Network standard:	IEEE802.11B/g/n
Operating frequency:	2.4GHz
Image transfer rate:	30fps
Operating temperature:	-10°C-55°C
Battery:	380mAh
Service time:	90min
Charging time:	40min
Input current:	DC5V300mAh
Lens diameter:	3.5mm
Optimum focal length:	1.5-2cm
Pixel:	UHD camera
Gyroscope:	three-axis

Indicator light

Machine State	Unconnect	Connect	Charge	Full battery
Light State	Green flash	Green On	Red On	Red Off

Pay attention

1. Not suitable for children under 3 years old.

2. Wipe the camera with a cotton swab if the image is not clear.

3. Please observe around before use, do not use in the presence of people running, to avoid collision.

4. When cleaning the lens, please wipe it carefully with a special alcohol swab.

5. Do not use this product alone by children to avoid accidental injury.


6. Do not use the endoscope in liquid to avoid water damage.

7. Avoid direct exposure of products to sunlight, especially accessories, to avoid softening.

8. During use, the temperature of the product will rise slightly to 35°C. Please feel free to use it.

9. To ensure the best use experience, please keep your child in a sitting and standing position.

10. Replace the ear spoon if it is loose or has been used more than 30 times for cleanliness and safety.



Hazardous substance

In order to reduce environmental pollution and protect the earth, it is specially declared that this product complies with the limited amount of substances in the electronic and electrical products of the People's Republic of China standard GB/T26572-2011 in terms of the content of the following substances

The name and content of harmful substances in the product

Component	Poisonous or harmful substance or element					
	Pb	Hg	Cd	[Cr(VI)]	PBBs	PGDES
Plastic parts	○	○	○	○	○	○
Electronic component	○	○	○	○	○	○
Printed circuit element	○	○	○	○	○	○
Metal assembly	○	○	○	○	○	○
Powercord	○	○	○	○	○	○

CE Indicates that the content of the substance in the part is below the GB/T26572 standard

X Indicates that the content of the substance in the part is over the GB/T26572 standard

Earchon®

ENT and minimally invasive surgical endoscopic camera system solutions provider

Earchon®

Create national medical brand to benefit the health of people around the world

Shenzhen Earchon Bio-technique Co., Ltd

Add.: 1702, Building A, No.9 Furong Road, Tan Tou Community, Song Gang Street, Bao'an District, Shenzhen City, Guangdong Province, China

Earchon®

Smart Visual Ear Cleaner

Model: YCCEQ

Features

Water /dust /fog proof lens

Automatic induction

UHD camera pixels

Status indicator light

Larger display range

Fast charging via Type-C port

Silicon earspoon

Better safety

Smart Visual Ear Cleaner

Model: YCCEQ

Image transfer rate: 30fp

Lens diameter: 3.5mm

Optimum focal length: 1.5-2cm

Pixel: Ultra HD camera

Gyroscope: three-axis

Input current: DC 5V 300mAh

Battery capacity: 380mAh

Service time: 90min

Operating temperature: -10~55°C

CE FC FDA RoHS

Earchon®

Smart Visual Ear Cleaner

Model: YCCEQ

WIFI direct connection

Ear scoop head storage

Material safety

Comfortable to the ear

HD display

High-end chip

Multi-axis gyroscope

Ultra-long life battery

Scan QR code to download APP



不干胶正面

不干胶背面

FCC Statement

1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference.
(2) This device must accept any interference received, including interference that may cause undesired operation.
2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE:
This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.
This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
Reorient or relocate the receiving antenna.
Increase the separation between the equipment and receiver.
Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
Consult the dealer or an experienced radio/TV technician for help.

SAR Information Statement

Your Smart Visual Ear Cleaner is a radio transmitter and receiver. It is designed and manufactured not to exceed the emission limits for exposure to radiofrequency (RF) energy set by the Federal Communications Commission of the U.S. Government. These limits are part of comprehensive guidelines and establish permitted levels of RF energy for the general population. The guidelines are based on standards that were developed by independent scientific organizations through periodic and thorough evaluation of scientific studies. The standards include a substantial safety margin designed to assure the safety of all persons, regardless of age and health. The exposure standard for Smart Visual Ear Cleaner employs a unit of measurement known as the Specific Absorption Rate, or SAR. The SAR limit set by the FCC is 4.0W/kg. * Tests for SAR are conducted with the phone transmitting at its highest certified power level in all tested frequency bands. Although the SAR is determined at the highest certified power level, the actual SAR level of the phone while operating can be well below the maximum value. This is because the phone is designed to operate at multiple power levels so as to use only the power required to reach the network. In general, the closer you are to a wireless base station antenna, the lower the power output. Before a phone model is available for sale to the public, it must be tested and certified to the FCC that it does not exceed the limit established by the government adopted requirement for safe exposure. The tests are performed in positions and locations (e.g., at the ear and worn on the body) as required by the FCC for each model. The highest SAR value for this model phone when tested for use at the ear is 0.260W/Kg and extremity, as described in this user guide, is 0.070W/Kg.

The FCC has granted an Equipment Authorization for this model phone with all reported SAR levels evaluated as in compliance with the FCC RFexposure guidelines. SAR information on this model phone is on file with the FCC and can be found under the Display Grant section of <http://www.fcc.gov/oet/fccid> after searching on FCC ID:2BB6VYCCEQ. Additional information on Specific Absorption Rates (SAR) can be found on the Cellular Telecommunications Industry Association (CTIA) web-site at <http://www.wow-com.com>. * In the United States and Canada, the SAR limit for mobile phones used by the public is 1.6 watts/kg (W/kg) averaged over one gram of tissue. The standard incorporates a substantial margin of safety to give additional protection for the public and to account for any variations in measurements.