Biotage® Liquid Level Sensor for Biotage® Selekt Family

User Manual

Biotage® Selekt Liquid Level Sensor for Biotage® Selekt Family

User Manual

CONTENTS

1	Introduction	
1	Biotage® Liquid Level Sensor	
1	Solvent Notifications and Stop Limits Based on Distance	
1	Safety Stop Limit for Waste	
1	Included in the Kit	
2	Technical Specifications	
3	Safety	
3	Intended Use	
3	Labels Used on the Biotage® Liquid Level Sensor Kit	
3	MIC Giteki Marking on Biotage® Liquid Level Sensor	
3	FCC and ISED for Biotage® Liquid Level Sensor	
4	Safety Requirements	
4	Safety in Other Languages	
5	Connect the Liquid Level Sensors to the System	
5	Connect the Hub to the System	
5	Mount the Liquid Level Sensors on the Reservoirs	
6	Cable and Tube Routing	
7	Connect the Liquid Levels Sensors to the Hub	
7	Connections	
7	Use the Universal Cap Adapter	
8	Enable the Liquid Level Sensors	
9	Administrate Reservoirs	
9	Reservoir Parameters	
9	Calibrate a Reservoir	
10	Transfer Reservoirs	
11	Assign Solvents and Reservoir Types	
12	Status of the Liquid Level Sensors	
13	Solvent and Waste Estimations and Notifications	
13	Run Estimation	
13	Solvent and Waste Notifications	
13	Safety Stop Limit for Waste	
14	Refill a Solvent and Empty Waste	
14	Refill a Solvent Reservoir	
14	Empty the Waste Reservoir	
14	Replace a Solvent or Waste Reservoir	
15	Maintenance and Troubleshooting	
15	Maintenance	
15	Troubleshooting	
16	General Information	
16	Consumables and Accessories	
16	Biotage® Selekt Family User Documentation	
16	Manufacturer	
16	Contact Us	
16	Intellectual Property	
17	Safety Translations	
17	安全について	
18	Sécurité	
19	Seguridad	
20	Sicherheitshinweise	
21	Sicurezza	
22	Sikkerhed	
23	Säkerhet	
24	Turvallisuus	
25	Veiligheid	

Introduction

Biotage® Liquid Level Sensor

Using radar technology, Biotage® Liquid Level Sensor measures the distance between the sensor (mounted on the reservoir cap) and the liquid level, enabling the system to present the actual amount of liquid available in the reservoir in real time.

With the liquid level sensor hub installed on your Selekt or Selekt Enkel system, and its functionality enabled in the system settings, the sensors can be used to monitor the liquid levels in the solvent reservoirs and/or the waste reservoir connected to your system (see Figure 1).

For the reservoirs to be used with the liquid level sensors, they must be defined and calibrated in the **Reservoir Administration** view. The software comes with two preconfigured reservoirs supplied by Biotage. Custom reservoirs can be defined and calibrated on one system and transferred to other Selekt and Selekt Enkel systems using a USB memory device or network share (if enabled).

Reservoirs are assigned to the solvent inlets and waste outlet in the **Solvent Setup** view; see page 11. The status of the connected liquid level sensors is displayed both in the **Solvent Setup** view and the **About** view (see page 12).

Figure 1. Selekt system with Biotage® Liquid Level Sensor used on all solvent reservoirs and the waste reservoir. All the liquid level sensors are connected to Biotage® Liquid Level Sensor Hub.

Solvent Notifications and Stop Limits Based on Distance

The solvent notifications and stop limits are based on the distance between the solvent surface and the bottom of the reservoir. This optimizes how much solvent that can be used without the risk of aspirating air. For more information, see page 13.

Safety Stop Limit for Waste

To ensure that the waste reservoir cannot be overfilled, the pump operation in progress is immediately paused if the distance between the sensor and the waste liquid becomes smaller than allowed. This is useful if e.g. the operator by mistake is using a waste reservoir that is smaller than the one that has been assigned to the solvent inlet in the software.

Included in the Kit

Part Number	Description	Qty
418758	Biotage® Liquid Level Sensor Hub	1
418574	CAN Power Cable	1
418536	Biotage® Liquid Level Sensor	5/3*
418777	Tube Holder	5/3*
418569	Waste Tube Insert	4
418839	Markings Kit	1
418859	Cable Support Clip, Solvent	20
418890	Cable Support Clip, Waste	5

* Five are delivered with the Biotage® Selekt Liquid Level Sensor Kit (P/N 418776) and three with the Biotage® Selekt Enkel Liquid Level Sensor Kit (P/N 418843).

Note that the required caps are not included in the kit and have to be ordered separately, see “Consumables and Accessories” on page 16.

We recommend that the boxes and packing materials are kept by the customer incase the kit needs to be shipped in the future.

Note: All the system photographs and software images in this manual show the Selekt system.

Technical Specifications

Fume Hood/ Ventilation System	A laboratory fume hood with the capacity to handle leakage of solvents and provide suitable airflow to prevent the build-up of flammable vapors.
Location	Indoor use only. No wet locations. Biotage® Liquid Level Sensor must be placed on a stable, and flat level surface.
Ambient Temperature	Operating: 15°C to 32°C (59°F to 90°F). Storage and transportation: -25°C to 60°C (-13°F to 140°F).
Humidity	20% to 90% RH (non-condensing).
Altitude	Operation: up to 2000 m (6562 ft).
Pollution Degree	The Biotage® Liquid Level Sensor kit shall be installed and used in a level 2 environment in accordance with EN 60664-1.
Electrical Supply	24 V
Max Power Consumed	500 mA
Materials	PEEK (tube holder), PPS (sensor house, cap, and universal cap adapter), PP (hub), TPE partly coated with PO (sensor cable)
Weight	Hub: 133.5 g (0.29 lbs) Sensor: 64.2 g (0.14 lbs)
Dimensions	Hub: 157 x 85 x 32 mm (6.2" x 3.3" x 1.3") Sensor: 63 x 55 x 27 mm (2.5" x 2.2" x 1.1")
IP Classification	IP54 for the sensor
RF Type (Sensor)	Transceiver
Operating Frequency Range	57 to 64 GHz
Radar Peak Output Power	18 dBm
Sensor Cable Length	1.8 m (70.9")

Safety

All operations on Biotage® Selekt and Biotage® Selekt Enkel systems must be performed according to the user documentation supplied with the system and its accessories; see “Biotage® Selekt Family User Documentation” on page 16.

Intended Use

Biotage® Liquid Level Sensor from Biotage is intended solely to be used with Biotage® Selekt and Biotage® Selekt Enkel systems to detect solvent and waste levels.

Labels Used on the Biotage® Liquid Level Sensor Kit

All symbols used on the labels on the Biotage® Liquid Level Sensor kit can be referred to in the Biotage® Selekt Family Installation and Safety document (P/N 417511) except the following symbol on Biotage® Liquid Level Sensor:

This device is approved by the Ministry of Internal Affairs and Communications in Japan; see “MIC Giteki Marking on Biotage® Liquid Level Sensor” below.

MIC Giteki Marking on Biotage® Liquid Level Sensor

Japanese Radio Law and Japanese Telecommunications Business Law Compliance.

This device is granted pursuant to the Japanese Radio Law (電波法) and the Japanese Telecommunications Business Law (電気通信事業法).

This device should not be modified (otherwise the granted designation number will become invalid).

FCC and ISED for Biotage® Liquid Level Sensor

This device contains licence-exempt transmitter(s)/ receiver(s) that comply with Innovation, Science and Economic Development (ISED) Canada’s licence-exempt RSS(s) and complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation of the device.

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Changes or modifications made to this equipment not expressly approved by Biotage may void the FCC authorization to operate this equipment.

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- » Reorient or relocate the receiving antenna.
- » Increase the separation between the equipment and receiver.
- » Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- » Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC/ISED radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 20 cm between the radiator and your body.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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Safety Requirements

- » Never connect or disconnect the liquid level sensor hub when the system is turned on.
- » Never use a liquid level sensor assembly that is damaged. Inspect before use.
- » Never lift a reservoir by its liquid level sensor assembly.
- » Always ensure that the solvent inlet tubes, waste outlet tube, and the sensor cables are routed in a safe manner to reduce the risk that the solvent and waste reservoirs tip over. It is recommended that the optional secondary solvent containment is used.
- » Ensure that the power cord, cables, and tubing connected to the system cannot come in contact with chemicals.
- » If spillage occurs on the liquid level sensor assembly, immediately dry it off.
- » Always ensure that the current waste level is correct in the **Solvent Setup** view after disabling **Biotage® Liquid Level Sensor** in the **System Settings** view.
- » Read and follow the safety requirements in the Biotage® Selekt Family Installation and Safety document (P/N 417511) that is supplied with the system.

Safety in Other Languages

See “Safety Translations” on page 17.

Connect the Liquid Level Sensors to the System

Prepare the System

Warning

» Read and follow the safety requirements on page 4.

Connect the Hub to the System

1. Ensure that the system is turned off.
2. Slide the hub into position between two rails on the right side of the system; see A in Figure 2. It is recommended that it is placed into the bottom rail.
3. Connect the CAN power cable supplied with the liquid level sensor kit between the **CAN IN** port on the hub and the **AUX** port at the rear of the system; see B in Figure 2.

Note: If the **AUX** port is plugged, remove the plug.

Figure 2. Connecting the liquid level sensor hub to the system.

Mount the Liquid Level Sensors on the Reservoirs

Note: Only reservoirs with the GL45 or 38-430 neck thread are currently supported. An universal cap adapter that may be used for other reservoirs can be ordered separately; see “Use the Universal Cap Adapter” on page 7.

1. Attach a GL45 cap (P/N 418528SP) or a 38-430 cap (P/N 418529SP) to the reservoir; see Figure 3.

Figure 3. The GL45 and 38-430 caps that can be used with the liquid level sensor kit.

2. Attach the solvent inlet/waste outlet tube to the tube holder; see Figure 4. Two tube inserts are required to attach the waste tube.


Figure 4. A: Two tube inserts are required for the waste outlet tube. **B:** The waste outlet tube attached to the tube holder.

3. Snap the tube holder onto the cap as shown in Figure 5.

Figure 5. Attach the tube holder to the cap.

4. Adjust the position of the tube by gently sliding it up or down the tube holder groove:
 - » **Solvent reservoir:** Ensure that the solvent inlet filter is at the bottom of the reservoir and that the filter and tube will not block the sensor's line of sight to the bottom of the reservoir; see Figure 6.
 - » **Waste reservoir:** Ensure that the end of the waste outlet tube is positioned at the bottom of the reservoir and that the tube will not block the sensor's line of sight to the bottom of the reservoir; see Figure 6.


Figure 6. A: Ensure that the solvent inlet filter, or the waste outlet tube, is at the bottom of the reservoir. **B:** Ensure that nothing is blocking the sensor's line of sight to the bottom of the reservoir.

5. Snap the sensor onto the tube holder as shown in Figure 7.


Figure 7. Snap the sensor onto the tube holder as shown in Figure 7.

6. Close the sensor by flipping it down. Ensure that it clicks into its closed position; see Figure 8.


Figure 8. **A:** The locking mechanism on the sensor. **B:** The sensor in its closed position.

7. Repeat step 1 to 6 for all reservoirs.
8. Ensure that the reservoirs are placed on a flat and stable surface. If using the optional secondary solvent containment, place the metal plate that can be ordered separately (P/N 418565SP) at the bottom of the containment (see Figure 9).

Figure 9. If using the optional secondary solvent containment, place the metal plate that can be ordered separately at the bottom of the containment.

Cable and Tube Routing

Link each sensor cable to its solvent/waste tube using five clips and three markers:

1. Link the sensor cable to its solvent/waste tube using five clips with a distance of approximately 15 cm (5.9"); see C1 to C5 in Figure 10. Place the first clip (C1) after the heat shrink tube; see Figure 11.
- Note:** Twenty (20) solvent clips, and five (5) waste clips are included in the support kit (P/N 418956SP). See the difference in Figure 12.
2. Link the sensor cable to both ends of the solvent/waste tube using two markers with the correct label (S1, S2, S3, S4, or W); see M1 to M2 in Figure 10. Place the first marker (M1) on the heat shrink tube; see Figure 11.
3. Attach a marker with the correct label to the end of the sensor cable; see M3 in Figure 10.

Note: Two markers per solvent/waste tube are included in the support kit (P/N 418956SP). The third marker is supplied with the system.


Figure 10. The cable and tube routing using clips and markers with S1-S5 and W.

Figure 11. Place the first marker (M1) on the heat shrink tube and the first clip (C1) after the heat shrink tube. (The photograph is showing a close-up of M1 and C1 in Figure 10.)


Figure 12. The solvent and waste clips.

Connect the Liquid Levels Sensors to the Hub

Connect each liquid level sensor to the hub port that corresponds to the reservoir it is mounted on; see Figure 13 and Table 1.

Note: The sensor setup has to be confirmed in the software when the **Biotage® Liquid Level Sensor** option is enabled (see “Enable the Liquid Level Sensors” on page 8) and each time a sensor is connected or disconnected to/from the hub. You can at anytime see the sensor setup in the **About** view; see Figure 27 on page 12.

Figure 13. The liquid level sensors connected to the hub.

Connections

Note: The **S3** and **S4** ports on the hub are not supported on Selekt Enkel, i.e. sensors connected to these ports will not be used.

System Port	Hub Port
S1	S1
S2	S2
S3	S3
S4	S4
Waste	W
AUX	CAN IN

Table 1. Connections on the system and the liquid level sensor hub.

Use the Universal Cap Adapter

There is an optional universal cap adapter for reservoirs that do not have the GL45 or 38-430 cap that can be ordered separately

(P/N 418571SP). See the requirements for using the adapter in Figure 14.


Figure 14. The universal cap adapter can be used with reservoirs that have a cap with 1) an inner diameter of 42 mm or larger and 2) a maximum top thickness of 6 mm.

1. Drill a hole of 38 mm in diameter in the middle of the cap; see Figure 14.
2. Insert the universal cap adapter through the hole; see A and B in Figure 15.
3. Attach the 38-430 cap supplied by Biotage to the universal cap adapter; see C in Figure 15.
4. Tighten the 38-430 cap and attach the whole cap assembly to the reservoir.

Figure 15. **A:** The universal cap adapter. **B:** The cap with the 38 mm hole on top of the adapter. **C:** The 38-430 cap attached to the adapter and the modified cap.

Figure 16. The liquid level sensor mounted on a reservoir that requires the universal cap adapter. Do not mount the sensor on a cap that is in a tilted position.

Enable the Liquid Level Sensors

Prepare the System

When the liquid level sensor hub has been connected to the system (see page 5), its functionality has to be enabled in the system settings (see Figure 17).

Note: Only users with system owner privilege can change the system settings.

1. Turn on the system using the power switch located on the left side of the system.
2. Press **Menu** and then **System Settings**.
3. Enable **Solvent Monitoring** and/or **Waste Monitoring**.
4. Enable the **Biotage® Liquid Level Sensor** option.
5. An information dialog opens, showing the sensor setup (see Figure 17). Confirm that the setup is according to your expectations. Otherwise, check the connections at the hub on the right side of the system (see Figure 13 and Table 1 on page 7).

Note: The sensor setup has to be confirmed each time a sensor is connected or disconnected to/from the hub. You can also at anytime see the sensor setup in the **About** view; see Figure 27 on page 12.


Figure 17. The System Settings view when the Biotage® Liquid Level Sensor option has been enabled. An information dialog opens, showing which sensors that are connected to the system.

Administrate Reservoirs

Prepare the System

The software comes with two preconfigured reservoirs supplied by Biotage that cannot be edited or deleted. To add, edit, calibrate, transfer, and delete user-defined reservoirs, press **Menu**, **Data Administration** and then **Reservoirs** (see Figure 18).

Monitoring using liquid level sensors is only possible when using 250 mL or larger glass or plastic reservoirs that have an inner height of 15 to 90 cm (5.9" to 35.4") and have been calibrated on a Selekt or Selekt Enkel system.

Only reservoirs with the GL45 or 38-430 neck thread are currently supported. An universal cap adapter supplied with the liquid level sensors that can be ordered separately may be used for other reservoirs; see page 7.

Note: Only users with system owner privilege can administrate reservoirs.

- » **Description:** A description of the reservoir, e.g. material, manufacturer, and part number. (Optional.)
- » **Calibrated:** The date and time when the reservoir was calibrated. The result of the calibration can be viewed by pressing ⓘ. To calibrate/recalibrate the reservoir, see the instructions below. (Optional.)

Note: Only calibrated reservoirs can be used with the liquid level sensors.


Figure 18. The Reservoir Administration view.

Reservoir Parameters

The following reservoir parameters are available (see Figure 19):

- » **Visible:** If enabled, the reservoir is available for use.
- » **Name:** The name of the reservoir.
- » **Capacity:** The reservoir volume in mL.
- » **Icon:** The icon of the reservoir; see Figure 20. The icon can help you to identify the correct reservoir when assigning reservoirs to the solvent inlets and the waste outlet (see Figure 24 on page 11). (Optional.)


Figure 19. The Edit view.


Figure 20. The available reservoir icons.

Calibrate a Reservoir

During the calibration, water will be dispensed into the reservoir through the waste channel while the reservoir is monitored by the waste sensor (port **W** on the hub). The Reservoir Calibration Wizard will guide you through the whole process.

Note: The last solvents used on the system must be miscible with water.

1. Press **Menu, Data Administration** and then **Reservoirs**.
2. Edit or create the reservoir that you want to calibrate.
3. Press **Calibrate**. Read and follow the instructions that appear on the screen; see Figure 21.

Reservoir Calibration Wizard (Page 1 of 3)

Assign water to one of the solvent inlets and prime with 60 mL of water:

1. Assign water to one of the solvent inlets. The capacity of the reservoir to be calibrated is 1000 mL.
Note: The last solvents used on the system must be miscible with water.
2. Select the solvent inlet with water assigned to it.
Solvent Reservoir: S3: Water **Current Level:** 5000 mL **Required:** ≈ 2366 mL
3. Prime the system with 60 mL of water. (Press **Prime** once.)

Reservoir Calibration Wizard (Page 2 of 3)

Calibrate the "1000 mL Clear Glass GL45":

1. Move the liquid level sensor connected to the W port on the sensor hub and the waste outlet tube to the empty reservoir to be calibrated.
2. Close the liquid level sensor and start the calibration.

Estimated remaining time (hh:mm):

Figure 21. The first and second page in the Reservoir Calibration Wizard.

Transfer Reservoirs

Reservoirs added on one system can be transferred to other Selekt and Selekt Enkel systems using a USB memory device or network share (if enabled).

Note: When exporting reservoirs from the system, they are saved at biotage\selekt\common\reservoirs\ (on USB) or \common\reservoirs\ (in network share folder). Ensure to use the same location when importing reservoirs to the system.

Note: If remote control is enabled on your system, it is not possible to export and import reservoir types using a USB memory device on a remote client.

1. Press **Menu, Data Administration** and then **Reservoirs**.
2. Press **Transfer**. The **Transfer Reservoirs** view opens; see Figure 22.
3. If applicable, select whether you want to transfer reservoirs from/to a USB memory device or the network share folder from the **Available on** drop-down list.

4. Select the reservoir(s) to be transferred.
To only list calibrated reservoirs, enable the **Only Calibrated** option.

To transfer all listed reservoirs, enable the **Select All** option.

Note: If a reservoir with the same name and identical settings exists on the system and the USB/network share, it is grayed out and cannot be selected.

5. Transfer the selected reservoir(s) by pressing **Selected**.
6. If a reservoir with the same name but different settings exists on the destination location, the **Reservoir Already Exists** dialog opens (see Figure 23). The settings for both reservoirs are listed. Select whether to replace the reservoir in the destination (**Replace**), keep both (**Create Copy**), or keep the reservoir in the destination (**Skip**).

Figure 22. The Transfer Reservoirs view. In this example, the selected reservoirs will be copied from the system to the connected USB memory device when **Selected** is pressed. The reservoir "2000 mL Amber Glass" is grayed out (and cannot be selected) as its settings are identical on the system and the connected USB memory device.

Figure 23. The Reservoir Already Exists dialog appears when a reservoir with the same name but with different settings exists in the destination location. In this example, the reservoir "2000 mL Clear Glass 38-430" on the system has a newer calibration than the one on the connected USB memory device.

Assign Solvents and Reservoir Types

Prepare the System

Warning

» If using the optional secondary solvent containment on the top of the system, it must be empty when using columns larger than 400 grams.

Note: For the liquid level sensors to work, they have to be used with reservoirs that have been calibrated. See “Administrate Reservoirs” on page 9.

1. Ensure that each liquid level sensor is connected to the hub port that corresponds to the reservoir it is mounted on; see see Figure 13 and Table 1 on page 7. For instructions on how to mount the liquid level sensors, see page 5.
2. In the software, press **Menu** and then **Solvent Setup**.
3. Assign a solvent to each solvent inlet using the **Solvent** drop-down lists in the **Solvent Setup** view.
To add solvents to the list, press **Solvent Administration**. For more information, see the user manual supplied with the system.
4. Assign a reservoir to each solvent inlet and the waste outlet from the **Reservoir** drop-down lists; see Figure 24. Reservoirs that are calibrated are marked with .
To add reservoirs to the list, press **Reservoir Administration**. For more information, see page 9.

The **Custom** option is available when you want to use a reservoir without adding it to the system; you can enter its capacity manually in the **Solvent Setup** view (see solvent inlet S₄ in Figure 24).

Note: If you select the **Custom** option or a reservoir that is not calibrated, the monitoring will be based on the pumped volume. The current liquid level must then be entered each time you refill a solvent or empty a waste reservoir (see Figure 25).

5. Prime the solvent inlets that have been assigned a new solvent.


Figure 24. The Reservoir drop-down list in the Solvent Setup view. Reservoirs that are calibrated are marked with .


Figure 25. The Solvent Setup view when solvent and waste monitoring using liquid level sensors is enabled. In this example, solvent inlet S₃ and S₄ have not been assigned a calibrated reservoir and the monitoring of the reservoirs will be based on the pumped volume. See the explanation of the sensor icons on page 12.

Status of the Liquid Level Sensors

Prepare the System

The status of the connected liquid level sensors is displayed both in the **Solvent Setup** view (see Figure 26) and in the **About** view (see Figure 27).

Icon	Status
No icon	No liquid level sensor is connected to the system or monitoring using liquid level sensors is not enabled. Note: If solvent/waste monitoring is enabled, it is based on the pumped volume. The current liquid level must be entered each time you refill/empty the reservoir.
	The connected liquid level sensor is open (flipped up) and/or the assigned reservoir is not calibrated (see reservoir S3 and S4 in Figure 25). The monitoring is based on the pumped volume. The current liquid level must be entered each time you refill/empty the reservoir.
	The connected liquid level sensor is used.
	The connected liquid level sensor failed to detect the liquid level. Check the following: » Ensure that the reservoir's liquid level is stable. » Ensure that the reservoir is not tilted. » Ensure that the tubing is properly inserted into the reservoir and that there is nothing disturbing the sensor; see steps 2 through 4 on page 5. » Ensure that there is no liquid on the bottom of the sensor house that can disturb the sensor readings. Dry of any liquid using a clean, lint-free cloth. » Ensure that the sensor is properly closed; see Figure 8 on page 6. » Rotate the sensor (and tube holder) into another position, e.g. rotate it 90 degrees, to see if that solves the problem. The monitoring is based on the pumped volume until the problem has been solved. If the problem has not been solved after using 20% of the set reservoir capacity, the run will be paused and the user notified. Resolve the problem or run without the sensor by flipping it up. Note: To be able to flip up the sensor, push on the tube holder as shown in Figure 29 on page 14.
	Sensor failure. The sensor must be disconnect from the liquid level sensor hub. (The user is notified.)


Figure 26. The status of the liquid level sensors shown in the Solvent Setup view.


Figure 27. The status of the liquid level sensors shown in the About view. Here is also the hub's serial number listed.

Solvent and Waste Estimations and Notifications

Setup and Run a Purification

Run Estimation

You can at any time view the estimated solvent consumption and the required waste capacity:

- » When setting up the run by pressing **Run Estimation** below the gradient table.
- » During the run by pressing **Run Estimation** in the **Solvent Setup** view; see Figure 28.

If there is not enough solvent or waste capacity for the run when the run is started, the system will notify the user.

Note: The available volume displayed in the **Run Estimation** dialog is equal to the current level (in the **Solvent Setup** view) minus the reservoir's stop volume; see Figure 28.


Figure 28. To view the estimated solvent consumption and the required waste capacity for the rest of the run, press Run Estimation in the Solvent Setup view. The estimations are updated in real time.

Solvent and Waste Notifications

When using the liquid level sensors for solvent and waste monitoring:

Notification	Description
 S1, S2	The solvent surface is 60 mm from the bottom of the reservoir. In this example the notification is for solvent inlets S1 and S2. When the solvent surface is 30 mm from the bottom of the reservoir, the run is automatically paused and the user will be prompted to refill the solvent(s).
 OFF	The solvent monitoring has been turned off for the rest of the run by the user. The monitoring is automatically turned on again when the run has been completed.
 W	The waste level is above 85% of the set waste reservoir capacity. When the waste level is at 95% of the set capacity, the run is automatically paused and the user will be prompted to empty the waste reservoir. When the distance between the sensor and the waste liquid becomes smaller than allowed, the run is automatically paused and the user will be prompted to empty the waste reservoir. See "Safety Stop Limit for Waste" below.

Safety Stop Limit for Waste

To ensure that the waste reservoir cannot be overfilled, the pump operation in progress is immediately paused if the distance between the sensor and the waste liquid becomes smaller than allowed. This is useful if e.g. the operator by mistake is using a waste reservoir that is smaller than the one that has been assigned to the solvent inlet in the software.

Refill a Solvent and Empty Waste

Setup and Run a Purification

Refill a Solvent Reservoir

Refill a solvent reservoir by flipping up the sensor; see Figure 29 and Figure 30. This can be done at any time without pausing. When the sensor is closed again, the amount of liquid available in the reservoir is automatically updated in the software.

Figure 29. Push on the tube holder to be able to flip up (open) the sensor.

Figure 30. Refill a solvent reservoir at any time using a funnel.

Empty the Waste Reservoir

Always ensure that the system is idle or paused when emptying the waste reservoir. Empty the waste reservoir by unscrewing the whole liquid level sensor assembly; see Figure 31.

The run is automatically paused when the waste level is 95% of the set capacity.

When the liquid level sensor assembly is attached to the reservoir again (with the sensor closed), the amount of liquid available in the reservoir is automatically updated in the software.

Replace a Solvent or Waste Reservoir

Always ensure that the system is idle or paused when replacing a solvent or waste reservoir. Either move the whole liquid level sensor assembly to the new reservoir by unscrewing the cap (see Figure 31) or detach the tube holder from the cap (see Figure 32) and attach the assembly to the new reservoir using a GL45 cap (P/N 418528SP) or a 38-430 cap (P/N 418529SP) supplied by Biotage.

Note: When replacing the reservoir with another type of reservoir, ensure to reassign the reservoir type in the **Solvent Setup** view (see page 11).

Figure 31. Remove the whole liquid level sensor assembly by unscrewing the cap. The cap rotates freely, independent of the attached tub holder and sensor.

Figure 32. Detach the tube holder from the cap.

When the liquid level sensor assembly is attached to a reservoir again (with the sensor closed), the amount of liquid available in the reservoir is automatically updated in the software.

Maintenance and Troubleshooting

Maintenance

Warning

- » If spillage occurs on the liquid level sensor assembly, immediately dry it off using a clean, non-abrasive cloth.

Clean the Tube Holder

The tube holder must be cleaned when spillage occurs or before moving it to a reservoir with a different solvent. Clean the tube holder using ethanol or similar.

Clean the Sensor House and Cable

Clean the sensor house using a cloth dampened with ethanol or similar.

If necessary, the sensor cable can be cleaned with a cloth lightly dampened with water or ethanol.

After cleaning, wipe dry with a clean, non-abrasive cloth.

Troubleshooting

If a liquid level sensor fails to detect the liquid level (🔔), check the following:

- » Ensure that the reservoir's liquid level is stable.
- » Ensure that the reservoir is not tilted.
- » Ensure that the tubing is properly inserted into the reservoir and that there is nothing disturbing the sensor; see steps 2 through 4 on page 5.
- » Ensure that there is no liquid on the bottom of the sensor house that can disturb the sensor readings. Dry of any liquid using a clean, lint-free cloth.
- » Ensure that the sensor is properly closed; see Figure 8 on page 6.
- » Rotate the sensor (and tube holder) into another position, e.g. rotate it 90 degrees, to see if that solves the problem.

General Information

Consumables and Accessories

Only genuine Biotage consumables and accessories must be used in the system. To order consumables and accessories, see contact information on the back of this document or visit our website www.biotage.com.

Biotage® Liquid Level Sensor Accessories and Spare Parts

Part No.	Description	Qty
418776	Biotage® Selekt Liquid Level Sensor Kit	1
418843	Biotage® Selekt Enkel Liquid Level Sensor Kit	1
418758SP	Biotage® Liquid Level Sensor Hub	1
418536SP	Biotage® Liquid Level Sensor	1
418777SP	Tube Holder	5
418523SP	Reservoir 5L, GL45	5
418524SP	Reservoir 10L, GL45	1
418528SP	GL45 Cap	5
418529SP	38-430 Cap	5
418571SP	Universal Cap Adapter (38-430)	5
418841SP	Funnel	1
418956SP	Cable and Tube Support Kit	1
418565SP	Secondary Solvent Containment Plate	1

For a complete list, please visit our website www.biotage.com.

Biotage® Selekt Family User Documentation

The following user documentation is supplied with the Biotage® Selekt and Biotage® Selekt Enkel systems and/or can be downloaded at www.biotage.com.

- » Biotage® Selekt Family Installation and Safety, P/N 417511
- » Biotage® Selekt User Manual, P/N 416183, or Biotage® Selekt Enkel User Manual, P/N 417510
- » Biotage® Selekt Family Safety Translations, P/N 417512
- » Biotage® Selekt Monitor API User Manual, P/N 418955 (only at www.biotage.com)
- » Biotage® Liquid Level Sensor for Biotage® Selekt Family User Manual, P/N 419144 (this document)

The user manuals supplied with the systems are also available in Chinese (P/N 416183ZH/417510ZH), Japanese (P/N 416183JP/417510JP) and German (P/N 416183DE/417510DE) at www.biotage.com.

Manufacturer

Biotage Sweden AB

Contact Us

Biotage Sweden AB

Box 8
SE-751 03 Uppsala
SWEDEN

Visiting address: Vimpelgatan 5

Phone: +46 18 56 59 00

Fax: +46 18 59 19 22

E-mail: info@biotage.com

Website: www.biotage.com

Please contact your local Biotage representative. See contact information on the back of this document or visit our website www.biotage.com.

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Your Complete Partner for Effective Chemistry

Biotage is a worldwide supplier of instruments and accessories designed to facilitate the work of scientists in life sciences. With our deep knowledge of the industry, academic contacts and in-house R&D teams, we can deliver the best solutions to your challenges. We take great pride in our flexibility and ability to meet our customer's individual needs. With strong foundations in analytical, organic, process, and biomolecule chemistry, we can offer the widest range of solutions available on the market.

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: + 46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: + 82 31 706 8500
Fax: + 82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

INDIA

Tel: +91 11 45653772
india@biotage.com

Distributors in other regions are listed on www.biotage.com

Part Number: 419144-A

© 2023 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.