

13” Bedside Tablet” User Manual

The Bedside Tablet is a Telemedicine, Smart ward Product. The utility model relates to the field of multimedia equipment technology, with state-of-the-art technology with all-in-one machine, including a Sturdy chassis , terminal control, and motherboard. It is used to transmit and display information, and multimedia can be used to publish videos and video on demand broadcasting.

Product clinical application

The patient's room is equipped with a corresponding bedside tablet primarily used for visual calls, remote family visits, and interaction. In cases where patients are in intensive care or have contagious diseases, and it's not convenient for family members to visit in person, they can conduct visits through the platform.

Through the Tele-conference platform , patients and visitors can achieve their visiting goals without direct physical contact, eliminating concerns about cross-infection between the two parties. Visitors do not need to enter the intensive care unit, which eliminates the inconvenience of changing into special visitation clothing and shoes upon entering the room. The system imposes no special requirements on visitors, making it accessible to people of all ages and physical conditions, whether they are elderly or children, and whether their health is good or poor.

This platform combines video and voice , and all a patient's family members need to do is scan the QR code provided during admission to conduct remote video visits with the patient. This system significantly meets the psychological needs of both patients and their families, alleviates family members' anxiety, and, at the same time, reduces the workload of healthcare personnel and minimizes visitation management costs



Product installation

Description of power supply network (The power adapter and POE are forbidden to be plugged into the device interface at the same time. The abnormal power connection will cause device damage immediately. Choose one of the following connection model)

Model 1: The first type of power adapter is powered by DC12V 2A can be directly plugged into the host DC interface, and the network input RJ45 interface;

Model 2: The second type of POE power supply, the front end needs to use the POE separator, the POE switch 48V input, and 12V output to the host through the POE separator;

FCC STATEMENT

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with minimum distance 20cm between the radiator & your body.

FCC NOTICE

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference; and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC CAUTION

Caution the user that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.