

故障排除

故障	解决方法
开机没反应	请使用配套的充电线进行充电，充电时长45分钟
摄像头不亮灯	采耳棒关机重启
APP闪退	卸载现有APP，重新扫描二维码进行安装
电源指示灯快闪	请在剩余电量用完前，使用配套的充电线为采耳棒充电
产品成功连接APP后未能正常显示实时影像	将手机与采耳棒WIFI的连接断开，重启采耳棒，重新连接WIFI
WIFI连接不稳定，有时候几秒钟就断开了	重启采耳棒，重新连接APP
若按照以上操作仍未能解决问题，请与O`haitai app客服取得联系，我们会进一步为您提供切实可行的解决方案	

产品参数

品牌	O`haitai
产品名称	智能可视采耳棒
产品型号	K1
主机重量	60g
网络标准	IEEE 802.11b/g
天线	内置天线
工作频率	2.4GHz
图像传输速率	20fps
图像传感器	CMOS
工作环境温度	-10~45°C
电池	200mAh锂电池
电池续航时间	约60分钟
充电时长	约45分钟
输入电流	DC 5V0.5A
镜头直径	3.5mm
最远焦距	15~20mm
景深	10~50mm
清晰度	1000万
陀螺仪	6轴



注意事项

1. 清洁镜头时，建议使用专业酒精棉签小心擦拭。
2. 使用前请观察周围，请勿在有人跑动的环境下使用产品，以免撞击。
3. 本产品不适用于3周岁以下儿童使用。
4. 禁止儿童单独使用本产品，以免不慎造成伤害。
5. 请勿将主机机身置于液体中使用，以免浸水造成损坏。
6. 产品内置聚合物电池，如不经常使用，也须每月充电一次。
7. 请使用5V的充电器给设备充电，在充电状态下，禁止使用此设备。
8. 避免产品暴露在阳光下直晒，以免软化。
9. 使用时，产品温度会轻微上升(最高不超过35°C)，请放心使用。

有害物质

部件名称	有害物质					
	铅(Pb)	汞(Hg)	镉(Cd)	六价铬(Cr(VI))	多溴联苯(PBB)	多溴联苯醚(PBDE)
印刷电路组件	×	○	○	○	○	○
外壳	○	○	○	○	○	○
塑胶零件	○	○	○	○	○	○

本表格根据SJ/T 11365的规定编制。
O`haitai公司声明该部件所有均质材料中所含的有害物质均符合GB/T26573-2012规定的限量要求。以下X表示该有害物质至少在该部件的某一批材料中含量超过GB/T26573-2012规定的限量要求。

售后保修

用户信息	
用户名：	联系电话：
收件地址：	
销售商家：	联系电话：

O`haitai

O`haitai

高品质生活，采耳享受

开启 O`haitai app

1 扫描二维码下载APP



2 打开APP
进入WIFI连接操作



3 长按电源键
开启采耳棒



O`haitai

High quality life, ear picking enjoyment

open O`haitai app

1
Enable scanning QR code download app



2
Open APP
Enter WIFI connection operation



3
Long press the power button
Open the ear picking rod



Press and hold for 1 second to turn on the device.
Blue light at the charging port

连接产品wifi


- 1 手机连接产品WIFI
选择“HBY-*****”
- 2 返回操作界面
点击“开始”
- 3 进入采耳操作界面



更换耳勺硅胶套



产品配件



产品尺寸



Troubleshooting

fault	resolvent
No response when turned on	Please use the matching charging cable for charging. Charging duration 45 minutes
Camera does not light up	Ear picking rod shutdown and restart
APP Flash Back	Uninstall the existing app and rescan the QR code to install
Power indicator light flashing fast	Please use the matching device before using up the remaining power. The charging line is for charging the ear rod.
The product failed to sync successfully after connecting to the APP. Display real-time images	Disconnect the phone from the ear pick stick WiFi. Restart the ear pick stick and reconnect to WiFi
WiFi connection is unstable, sometimes for a few seconds. It's disconnected	Restart the ear picking stick and reconnect to the APP
If the problem cannot be resolved by following the above steps, please contact O`haitai app customer service. We will further provide you with practical and feasible solutions	

Product parameters

brand	O`haitai
Product Name	Smart visual ear pick
Product model	K1
Host weight	60g
Network standards	IEEE 802.11b/g/n
antenna	Internal antenna
working frequency	2.4GHz
Image transmission rate	20fps
image sensor	CMOS
Working environment temperature	-10~45°C



Note

1. When cleaning the lens, it is recommended to use a professional alcohol cotton swab to wipe it carefully.
2. Before use, please observe the surroundings and do not use the product in an environment where people are running, to avoid impact.
3. This product is not suitable for use by children under the age of 3.
4. Children are prohibited from using this product alone to avoid accidental injury.
5. Do not use the host body in liquid to avoid damage caused by immersion.
6. The product is equipped with a built-in polymer battery, which must be charged once a month if not frequently used.
7. Please use a 5V charger to charge the device. It is prohibited to use this device while it is in a charging state.
8. Avoid exposing the product to direct sunlight to avoid softening.
9. When using, the product temperature may slightly increase (up to 35°C), please rest assured of use.

After sales warranty

user information	
User Name:	Contact number:
Delivery address:	
Sales merchant:	Contact number:

Hazardous substances

parts name	Hazardous substances					
	lead (Pb)	mercury (Hg)	cadmium (Cd)	hexavalent chromium (Cr(VI))	polychlorinated biphenyls (PBB)	polychlorinated dibenzodioxin (PCDD)
printed circuit assembly	×	○	○	○	○	○
housing	○	○	○	○	○	○
plastic Rubber accessories	○	○	○	○	○	○

This table is prepared in accordance with the revision of SJ/T 11365. O indicates that the content of the harmful substance in all homogeneous materials of the component is below the limit requirements specified in GB/T26573-2012. X indicates that the content of the hazardous substance in at least one homogeneous material of the component exceeds the limit requirements specified in GB/T26573-2012.

O`haitai

O`haitai

高品质生活，采耳享受

开启 O`haitai app

1 扫描二维码下载APP



2 打开APP
进入WIFI连接操作



3 长按电源键
开启采耳棒




Press and hold for 1 second to turn on the device.
Blue light at the charging port

连接产品wifi

- 1 手机连接产品WIFI
选择“HBY-*****”
- 2 返回操作界面
点击“开始”
- 3 进入采耳操作界面



产品配件



Product size



This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction