



## **BioPot3 Device**

### **Instructions For Use**

**Sensomedical Labs**

**LTD**



Manufacturer:

Sensomedical Labs LTD.

Nazareth Industrial Park, Mount Precipice

Nazareth 1612102, Israel

+972 (0)4 6800668

## Contents

1. BIOPOT3 Kit .....	5
1.1 BIOPOT VERSIONS, PRODUCT FAMILY DESCRIPTION .....	5
1.2 KIT .....	5
1.3 BIOPOT: .....	6
1.3 Charging .....	6
1.4 Powering .....	7
1.5 Reference and ground switch .....	7
1.6 BIOPOT CAP .....	8
1.7 BIOPOT Name Description .....	8
1.8 BIOPOT LEDs mode .....	9
2. Android Application .....	9
2.1 Features or functional requirements .....	9
2.2 Installation .....	9
2.3 Launch application & set permissions .....	11
2.4 Patient, Session & user details .....	12
2.5 Connect to BIOPOT .....	17
2.6 Main screen .....	19
2.6.1. BIOPOT information: .....	20
2.6.2. Tool bar: .....	20
2.6.3. EEG/EMG channels: .....	22
2.6.4. Accelerometer channels .....	23
2.6.5. Recording files .....	23
2.7 Channels Zoom & Montage .....	23
2.8 Impedance measurement .....	26
2.9 Additional settings, filters and device readout .....	26
2.10 Recording data .....	28
2.10.1. Recording data types .....	30
2.12.1. Connection to BIOPOT was lost .....	35
3. BLE Protocol .....	36
4.1 Characteristics 1 - [Read & Write] .....	36
4.2 Characteristics 2 - [Read & Write] .....	38
4.3 Characteristics 3 - [Read & Write] .....	38
4.4 Characteristics 4 [Notify] .....	39
4.5 Characteristics 5 [Read & Write] .....	40

4.6 Characteristics 6 [Notify] .....	41
4. Requirements .....	42
5.1. Temperature and Humidity: .....	42
5.2. The BIOPOT3 V2 shall withstand during operation the following conditions: Temp: 0°C – 40°C and Humidity: 20 -90 % .....	42
5.3. The BIOPOT3 V2 shall withstand during storage the following conditions: Temp: -10°C – 50°C and Humidity: 20 -90 % .....	42
5.4. Atmospheric Pressure (kPa) .....	42
5.5. ....	42
5.6. Materials in contact with the patient or operator: .....	42
All substances in contact with the skin are typical medical grade plastics and other biocompatible materials which are commonly used in adhesive electrode systems. ....	42
5.7. Cleaning and disinfection: .....	42
5.8. Waste product .....	42
6. Instructions for use .....	43
6.1. The device may be used by a researcher/professional in a specialized lab for data acquisition. It can also be used by a home user together with a mobile device or PC for data acquisition. ....	43
6.2. The device is used by a researcher/professional in a specialized clinical lab for data acquisition. It can also be used by a home user together with a mobile device or PC for data acquisition. ....	43
6.3. The primary user may need to collect biopotential data, bioimpedance data, accelerometer data, or some combination of all. Data collection requirements may scale up over time (i.e. start basic with increased functionality with project progress) .....	43
6.4. The primary user may be working on their own software, or they may require just a basic application for data collection only with no software development required. BOTH potential requirements must be covered: easy software for data collection without development, SDK for customized software development. ....	43
6.5. The primary user likely requires flexibility in the sensor part of the project, and may need to use a basic off-the-shelf option, or may be working on their own sensors. Multiple options may be required in a single customer/project. ....	43
6.6. The secondary user is a student that requires an easy platform for simple data collection for basic research purposes. This does not require programming. ....	43
6.7. The secondary user requires basic software and hardware interfacing and may be encouraged to do their own software development. SDK with examples would be useful, compatibility with MatLab or Labview, and or Python required. ....	43
6.8. In order to start using the device, the operator need to attach the electrode sticker to the device (as described in section 1.6). To attach the electrode to the patient, see "Electrode patch user manual". After the preparation you can start recording data (as described in section 2.10). ....	43
6.9. Hazards related to use of the medical device .....	43
6.10. General: .....	44

6.11. Technical Description: .....	45
7. General information .....	45
7.1. Manufacturer: Senso Medical LTD .....	45
7.2. Costumer: Senso Medical Labs Ltd. ....	45
7.3. Costumer's Address: .....	45
7.4. Equipment Model No.: BP-30-08-00061 .....	45
7.5. Equipment Serial No.: 61 .....	45
7.6. For labels explanation look into Biopot labels file. ....	45
8. List of accessories, detachable parts, and materials .....	45
9. Declaration – electromagnetic emissions .....	46
10. Declaration – electromagnetic immunity .....	46
11. Recommended separation distances between portable and mobile RF communications equipment and the BioPot SYSTEM.....	48
12. Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment	48

# 1. BIOPOT3 Kit


## 1.1 BIOPOT VERSIONS, PRODUCT FAMILY DESCRIPTION

### 1.2 KIT

KIT includes BIOPOT3 device, docking charging station, and USB-C charging cable:




### 1.3 BIOPOT:



#### 1.3 Charging


You need to enter the docking station as picture bellow connected to USB cable:

When power cable is connected the yellow LED should turn on and docking station green LED on:



#### 1.4 Powering

- Short press on power button turns BIOPOT ON, BIOPOT green LED should start blinking.
- Long press on power button turns BIOPOT OFF (close to 1.5 seconds).




#### 1.5 Reference and ground switch

- If switch to the right side that means reference and ground are shorted.
- If switch to the left side that means reference and ground are NOT shorted.



## 1.6 BIOPOT CAP



## 1.7 BIOPOT Name Description

Index	Description	Values
2-9	“SML BIO” [8 chars]	Unchangeable
10-11	Channels number [2 chars]	08/16/19
12	Accelerometer [1 char]	X- without accelerometer A-with accelerometer
13	Bio-impedance [1 char]	X- without bioimpedance I-with bioimpedance
14-19	Serial number [4 chars]	Serial number. Value from 0 to 65535

Example: SML BIO 08AX00005: This means, that the EEG device is having 8channels with accelerometer and no Bio-impedance function with serial number “00005”



## 1.8 BIOPOT LEDs mode

		Green	Blue	Red	Time On [mSec]	Time Off [mSec]
Advertising	Not connected to Bluetooth.	Flashing	Off	Off	100	1500
Idle	while BIOPOT3 connected to Bluetooth and there is no data acquisition.	Off	Flashing	Off	100	100
Data acquisition	while BIOPOT3 acquiring data (also connected to Bluetooth).	Off	Flashing	Off	100	500
impedance	In impedance measurement	Flashing	Flashing	Off	100	500
Low Battery	Low battery while Idle or Acquiring data	Off	Off	Flashing	Same as data acquisition or Idle	

## 2. Android Application

### 2.1 Features or functional requirements

Android device:



- Android device OS version 6.1 and above.
- Android device Bluetooth version 4.2 and above.

### 2.2 Installation

In order to install android application, follow these steps:


- Download APK file that you get from SensoMedical Labs.
- If you already have BIOPOT3 application on your device, you need to uninstall it.
- Move APK file to your phone or tablet.

- Location APK file and press on it




press “install” and follow installation process till installation done.



- After installation done, BIOPOT3 application should appear on your home screen:



## 2.3 Launch application & set permissions

Press on BIOPOT3 icon to open application, on first application launch you need to allow android permission in order:





NOTE: make sure that location is Enabled on your android device.

## 2.4 Patient, Session & user details

Note: this part can be skipped, press “Skip all” and default parameters will be set ( as shown below).

This is Patient Details screen:

11:33


Patient   Session   User

1   2   3

### Patient Details

Patient Code (scan or enter)

0000



Place patient barcode inside the blue frame

Skip all

NEXT

III   O   <

you can Scan patient barcode or enter patient code:

Note: Next button will be disabled if patient code not filled.

11:33


PatientSessionUser

123

Patient Details

Patient Code (scan or enter)

0000



Place patient barcode inside the blue frame

Skip all

NEXT

11:33


PatientSessionUser

123

Patient Details

Patient Code (scan or enter)

051111407592

 Scan Again

Skip all

NEXT

Press “Next” to set Session details:

1:42 100%

Patient Session User

1 2 3

### Session Details

Saving File Format  
**EDF/BDF Format**

File Maximum Size  
**100 MB**

Saving Target  
**Internal Memory**

Folder Name  
**Biopot data files**

Folder will be saved under root in saving target

File Name  
**BioPot**

Date ☒

Time ☒

Example  
BioPot\_20180624\_22.30.41\_0001

[Skip all](#)

Back Next

Here you can set BIOPOT recording file location and files name.

Based on example above, recording application files will be on:  
phone internal storage >> Biopot data files.

Files name start with BioPot followed by time, date and number of file same as shown in Example.

NOTE: these parameters decide the location of stored files and files names.

Press “Next” to set User Details:

The screenshot shows a mobile application interface with a blue header bar. At the top left, the time is 11:34. The header bar contains three tabs: 'Patient' with a checkmark, 'Session' with a checkmark, and 'User' with a blue circle containing the number 3. Below the header, the title 'User Details' is centered. The form contains four fields: 'Name' with the value 'Dr, Anonymous', 'ID' with the value '0000', 'Role' with the value 'Doc', and 'Email' with the value 'No Data'. At the bottom of the form, there is a link labeled 'Skip all'. Below the form, there are two buttons: 'Back' and 'Done'. The 'Done' button is highlighted in blue. At the very bottom, there is a navigation bar with three icons: a square, a circle, and a triangle.

Field	Value
Name	Dr, Anonymous
ID	0000
Role	Doc
Email	No Data

Here you can set user details: Name, ID, Role and Email.

This information will be shown on recording files.


NOTE: this information will be saved on recording text files only.

Press “Done” to “Connect to BIOPOT”.




## 2.5 Connect to BIOPOT

Application locating BIOPOT devices:



If application couldn't find BIOPOT devices, you shall see this screen:




You need to press “Try again” to start location BIOPOT devices again, till you find your BIOPOT device:



If you have more than one device, select the BIOPOT needed to connect and press “Start” Button.

## 2.6 Main screen


If application succeeds to connect to BIOPOT you should see channels main screen:



Main screen divided to several sections:

#### 2.6.1. BIOPOT information:

BIOPOT information located on the top left on main screen:





Here you can find the BIOPOT name, EEG/EMG channels number, accelerometer status and Bio-impedance status.


If you press on  icon, application navigate to connection screen.


#### 2.6.2. Tool bar:


BIOPOT tool bar location on to right on main screen:




Omega icon  navigate to impedance measurement screen.

Flash icon  navigates to samples per second's graph:



Folder icon  navigates to saved recording files.

Setting icon  navigate to:



Here you can set the number of channels you need to see on main screen on each channels types: EEG/EMG, Bio-impedance or Accelerometer.

## Additional Settings

Additional setting label  
setting screen.

navigates to additional

## Patient/Session/User Setup


This label  
& user details screen that set before.

navigate to Patient, Session

### 2.6.3. EEG/EMG channels:


In EEG/EMG section you should see BIOPOT channels, depends on BIOPOT channels number.



Using , you can change charts amplitude values.



Using , you can change charts period value.

Pressing  will reset amplitude to default value 300 micro volt / 10 seconds.

#### 2.6.4. Accelerometer channels

On Accelerometer section you should see BIOPOT accelerometer X, Y, and Z accelerometer axis's if the accelerometer enabled in BIOPOT device.

#### 2.6.5. Recording files




Recording button start data recording.

### 2.7 Channels Zoom & Montage


By pressing on channel chart, you will be navigating to chosen “Zoom channel” channel screen:




You can switch between channels using “Next” and “Prev” buttons and channels types:



All channels referenced to same reference/ground, in order to montage channels you need to change reference to wanted channel on the bottom of the zoom channel screen:



After choosing wanted montage channel press “Apply” button to save



changes:



Montage Channel

View this channel over: Ch. 1 APPLY



When choosing “montage channel”, channel label will change based to chosen montage channel:

From

EEG/EMG Bio Impedance Accelerometer

Zoom Channel: < Prev 5 Next >


to

Zoom Channel: < Prev 5(1) Next >

that's mean that 5 is montaged by channel 1.

2.8 Impedance measurement


From main screen, omega icon  navigate to impedance screen:





Based on impedance colors and “impedance” you can tell the impedance value/interval.

Based on “channels mapping” you can recognize channel number of impedance measurement.

2.9 Additional settings, filters and device readout

To enter “Additional setting” press on  from main screen ,then

 Additional Settings on the bottom of screen appears:



Here you can set/change the BIOPOT parameters.


Hardware low pass filter, hardware high pass filter and sampling rate parameters for A2D chip.


Software low pass filter – this filter process after the A2D HW LPF and HPF.

Notch filter – this filter process for removing noise on application channel and not applied on recording data.

Note: recording data filtered by hardware and software filter, notch filter not changed saved data.



After changing additional setting, press “Apply” button to set changes or press back button on the left top “back” button to go to main screen.


“Device Readout” next to “Additional Setting” tab  navigate to “Device Readout screen:





Here you can see BIOPOT battery percentage and BIOPOT temperature

## 2.10 Recording data


To start recording data, press on  from main screen, recording button should change and timer for recording period should appears next to it .

If you need to stop recording, press on  button.

You can see recording files from application by pressing Folder icon  on main screen:



Extract saved recording files from android device based on “session details” (by default, files saved on: internal Storage >Biopot data files) to:



## 2.10.1. Recording data types


### 2.10.1.1. EDF/BDF files

Based on BIOPOT version, data recording files type changes, BIOPOTs 16bits in EDF and BIOPOTs 24bits in BDF format.

EDF and BDF files can be opened on any EEDF and BDF readers like EEGLAB and BDF

### 2.10.1.2. Text files

On the header of the text file you can find user and patient and session details followed by blocks of data:



#### 2.10.1.2.1. Data format

There 7 types of characteristics 4 buffer data, which depends on the number of channels, number of samples per channels, accelerometer and bio impedance, explained in the table below:

Type	BioPot	Data size [bits]	Channel numbers	Samples per channel	Accelerometer	Bio-impedance	Buffer bytes description	Buffer length
1	8	16	8	15	-	-	4 [TS] + 8*15*2	244
2	8	16	8	14	Activated	-	4 [TS] + 8*14*2 + 12[acc]	240
3	8	16	4	28	Activated	Activated on 4 channels	4 [TS] + 4*28*2 + 4[bio] + 12[acc]	244
4	16	16	16	7	-	-	4 [TS] + 16*7*2	228
5	16	16	16	7	Activated	-	4 [TS] + 16*7*2+ 12[acc]	240
6	16	16	12	9	Activated	Activated on 4 channels	4 [TS] + 12*9*2+ 4[bio] +	236

							12[acc]	
7	19	16	19	6	-	-	4 [TS] + 19*6*2	232
8	19	16	19	6	Activated	-	4 [TS] + 19*6*2+ 12[acc]	244
9	19	16	15	7	Activated	Activated on 4 channels	4 [TS] + 15*7*2+ 12[acc]+ 4[bio]	230
10	8	24	8	10	-	-	4 [TS] + 8*10*3	244
11	8	24	8	9	Activated	-	4 [TS] + 8*9*3+ 12[acc]	232

Example of type 5 from table above:

	Start index	End index	index			Low/high
Timestamps	0	3	0	TS0		low
			1	TS1		
			2	TS2		
			3	TS3		high
Sample 1 of 12 channels	4	27	4	S1	CH1	low
			5			high
			6		CH2	low
			7			high
			8		CH3	low
			9			high
			10		CH4	low
			11			high
			12		CH5	low
			13			high
			14		CH6	low
			15			high
			..		...	...
			Sample 2 of 12 channels		28	51
....	...	...	...			
Sample 9 of 12 channels	196	219	S9			
Bio-Impedance	220	223	220-221	S1	...	...
			222-223	S2		
Accelerometer	224	225	224	X1	Point 1	low
			225			high
	226	227	226	Y1		low
			227			high
	228	229	228	Z1		low
			229			high





			243			high
--	--	--	-----	--	--	------

Example of type 11 from table above:

	Start index	End index	index		Low/high
Timestamps	0	3	0	TS0	low
			1	TS1	
			2	TS2	
			3	TS3	high
Sample 1 of 9 channels	4	27	4	S1	low
			5		
			6		high
			7		low
			8		
			9		high
			10		low
			11		
			12		high
			13		low
			14		
			15		high
			16		low
			17		
			18		high
			19		low
			20		
			21		high
			22		low
			23		
			24		high
			25		low
			26		
			27		high
Sample 2 of 9 channels	28	51	Same	S2	...
....				...	
Sample 9 of 9 channels	196	219		S9	
Accelerometer	220	221	220	X1	low
			221		high
	222	223	222	Y1	low
			223		high
	224	225	224	Z1	low
			225		high
	226	227	230	X2	low
			231		high
	228	229	232	Y2	low
			233		high
	230	231	234	Z2	low
			235		high


General format:

Data start with 4 bytes of timestamps followed by A2d data sample by sample for all channels based on number of samples followed by accelerometer data if accelerometer is enable.


NOTES:

- all data saved in little endian format (low to high).
- Sample size is 2 bytes (BIOPOT 16bits) OR 3 bytes (BIOPOT 24bits).
- Number of samples calculated based on BIOPOT sample size and if accelerometer enabled/disabled.
- See BLE Protocol characteristic 1.

### 2.12.1. Connection to BIOPOT was lost



This message appears when BIOPOT disconnected, application stuck in “locating connection” screen in order to locate same BIOPOT and automatically connect it.



In order to exit this screen, press the back button on the left/right bottom device bar.

### 3. BLE Protocol

BLE data protocol:

#### 4.1 Characteristics 1 - [Read & Write]

Byte Index		
0-4	Registers - Reserved	Reserved
5	Channels number	8

6	Accelerometer mode/status	0 – Accelerometer off 1 – 2g dynamically selectable 2 – 4g dynamically selectable 3 – 8g dynamically selectable 4 – 16g dynamically selectable
7	Bio impedance mode/status	0- off 1- on
8	External memory mode/status	0- off 1- on
9	Number of samples per channels	This value changed based on bytes 5,6 and 7, which equal to:
	$\#samples_{channel} = \frac{244 - 4[ts] - 4 * hasBioImp - 12 * hasAccelerometer}{\#channels * SampleDataSize_{InBytes}}$	
10	BIOPOT data size	0x18 – 24 bits Else -16 bits
11	Sync Enable Not implemented	1 – Sync Enabled (in the test device when this byte is written as 1 a sync simulation will start after 3 seconds) 0 – Sync disabled
12	Serial number	BIOPOT3 serial number
13		
14		
15		

#### 4.2 Characteristics 2 - [Read & Write]

Index		
0	BIOPOT data acquisition	0 – Turn data acquisition off 1 - Turn data acquisition on 2- Read external Memory 3- Erase External Memory

#### 4.3 Characteristics 3 - [Read & Write]

Index		
0	Impedance measurement	0.off 1.On 2.Offset Calibration 3.Impedance Calibration
1	Reserved: Bioimpedance	0 – off 1 – on
2	Reserved	Reserved
3-22	Impedance values	20 impedance values in 100 $\Omega$ resolution. Byte for each channel.

#### 4.4 Characteristics 4 [Notify]

There 7 types of characteristics 4 buffer data, which depends on the number of channels, number of samples per channels, accelerometer and bio impedance, explained in the table below:

Type	BioPot	Data size [bits]	Channel numbers	Samples per channel	Accelerometer	Bio-impedance	Buffer bytes description	Buffer length
1	8	16	8	15	-	-	4 [TS] + 8*15*2	244
2	8	16	8	14	Activated	-	4 [TS] + 8*14*2 + 12[acc]	240
3	8	16	4	28	Activated	Activated on 4 channels	4 [TS] + 4*28*2 + 4[bio] + 12[acc]	244
4	16	16	16	7	-	-	4 [TS] + 16*7*2	228
5	16	16	16	7	Activated	-	4 [TS] + 16*7*2+ 12[acc]	240
6	16	16	12	9	Activated	Activated on 4 channels	4 [TS] + 12*9*2+ 4[bio] + 12[acc]	236
7	19	16	19	6	-	-	4 [TS] + 19*6*2	232
8	19	16	19	6	Activated	-	4 [TS] + 19*6*2+ 12[acc]	244
9	19	16	15	7	Activated	Activated on 4 channels	4 [TS] + 15*7*2+ 12[acc]+ 4[bio]	230
10	8	24	8	10	-	-	4 [TS] + 8*10*3	244
11	8	24	8	9	Activated	-	4 [TS] + 8*9*3+ 12[acc]	232

#### 4.5 Characteristics 5 [Read & Write]

Index																																																						
0	16bits: Channels off/on. 24bits Not used	Bit [32,31, 30.....25]																																																				
1		Bit [24,23, 22.....17]																																																				
2		Bit [16,15, 14.....9]																																																				
3		Bit [8,7,6,5, 4.....1]																																																				
4	Not used																																																					
5	24bits: SW <b>High</b> pass filter																																																					
		<table><tr><td><b>Index.</b></td><td><b>Cutoff [Hz]</b></td><td><b>Index.</b></td><td><b>Cutoff [Hz]</b></td></tr><tr><td><b>0</b></td><td>No HPF</td><td><b>12</b></td><td>15</td></tr><tr><td><b>1</b></td><td>0.1</td><td><b>13</b></td><td>20</td></tr><tr><td><b>2</b></td><td>0.25</td><td><b>14</b></td><td>25</td></tr><tr><td><b>3</b></td><td>0.3</td><td><b>15</b></td><td>30</td></tr><tr><td><b>4</b></td><td>0.5</td><td><b>16</b></td><td>50</td></tr><tr><td><b>5</b></td><td>1</td><td><b>17</b></td><td>75</td></tr><tr><td><b>6</b></td><td>2</td><td><b>18</b></td><td>100</td></tr><tr><td><b>7</b></td><td>2.5</td><td><b>19</b></td><td>150</td></tr><tr><td><b>8</b></td><td>3</td><td><b>20</b></td><td>200</td></tr><tr><td><b>9</b></td><td>5</td><td><b>21</b></td><td>250</td></tr><tr><td><b>10</b></td><td>7.5</td><td><b>22</b></td><td>300</td></tr><tr><td><b>11</b></td><td>10</td><td><b>23</b></td><td>500</td></tr></table>	<b>Index.</b>	<b>Cutoff [Hz]</b>	<b>Index.</b>	<b>Cutoff [Hz]</b>	<b>0</b>	No HPF	<b>12</b>	15	<b>1</b>	0.1	<b>13</b>	20	<b>2</b>	0.25	<b>14</b>	25	<b>3</b>	0.3	<b>15</b>	30	<b>4</b>	0.5	<b>16</b>	50	<b>5</b>	1	<b>17</b>	75	<b>6</b>	2	<b>18</b>	100	<b>7</b>	2.5	<b>19</b>	150	<b>8</b>	3	<b>20</b>	200	<b>9</b>	5	<b>21</b>	250	<b>10</b>	7.5	<b>22</b>	300	<b>11</b>	10	<b>23</b>	500
		<b>Index.</b>	<b>Cutoff [Hz]</b>	<b>Index.</b>	<b>Cutoff [Hz]</b>																																																	
		<b>0</b>	No HPF	<b>12</b>	15																																																	
		<b>1</b>	0.1	<b>13</b>	20																																																	
		<b>2</b>	0.25	<b>14</b>	25																																																	
		<b>3</b>	0.3	<b>15</b>	30																																																	
		<b>4</b>	0.5	<b>16</b>	50																																																	
		<b>5</b>	1	<b>17</b>	75																																																	
		<b>6</b>	2	<b>18</b>	100																																																	
		<b>7</b>	2.5	<b>19</b>	150																																																	
		<b>8</b>	3	<b>20</b>	200																																																	
		<b>9</b>	5	<b>21</b>	250																																																	
		<b>10</b>	7.5	<b>22</b>	300																																																	
<b>11</b>	10	<b>23</b>	500																																																			
6	sampling rate	500[default]/1000/2000																																																				
7																																																						
8	Not used																																																					
9																																																						
10	SW Low pass filter																																																					
		<table><tr><td><b>Index.</b></td><td><b>Cutoff</b></td><td><b>Index.</b></td><td><b>Cutoff</b></td></tr><tr><td><b>0</b></td><td>No LPF</td><td></td><td></td></tr><tr><td><b>1</b></td><td>5hz</td><td><b>11</b></td><td>55hz</td></tr><tr><td><b>2</b></td><td>10hz</td><td><b>12</b></td><td>60hz</td></tr><tr><td><b>3</b></td><td>15hz</td><td><b>13</b></td><td>65hz</td></tr><tr><td><b>4</b></td><td>20hz [Default]</td><td><b>14</b></td><td>70hz</td></tr><tr><td><b>5</b></td><td>25hz</td><td><b>15</b></td><td>75hz</td></tr><tr><td><b>6</b></td><td>30hz</td><td><b>16</b></td><td>80hz</td></tr><tr><td><b>7</b></td><td>35hz</td><td><b>17</b></td><td>85hz</td></tr><tr><td><b>8</b></td><td>40hz</td><td><b>18</b></td><td>90hz</td></tr><tr><td><b>9</b></td><td>45hz</td><td><b>19</b></td><td>95hz</td></tr><tr><td><b>10</b></td><td>50hz</td><td><b>20</b></td><td>100hz</td></tr></table>	<b>Index.</b>	<b>Cutoff</b>	<b>Index.</b>	<b>Cutoff</b>	<b>0</b>	No LPF			<b>1</b>	5hz	<b>11</b>	55hz	<b>2</b>	10hz	<b>12</b>	60hz	<b>3</b>	15hz	<b>13</b>	65hz	<b>4</b>	20hz [Default]	<b>14</b>	70hz	<b>5</b>	25hz	<b>15</b>	75hz	<b>6</b>	30hz	<b>16</b>	80hz	<b>7</b>	35hz	<b>17</b>	85hz	<b>8</b>	40hz	<b>18</b>	90hz	<b>9</b>	45hz	<b>19</b>	95hz	<b>10</b>	50hz	<b>20</b>	100hz				
		<b>Index.</b>	<b>Cutoff</b>	<b>Index.</b>	<b>Cutoff</b>																																																	
		<b>0</b>	No LPF																																																			
		<b>1</b>	5hz	<b>11</b>	55hz																																																	
		<b>2</b>	10hz	<b>12</b>	60hz																																																	
		<b>3</b>	15hz	<b>13</b>	65hz																																																	
		<b>4</b>	20hz [Default]	<b>14</b>	70hz																																																	
		<b>5</b>	25hz	<b>15</b>	75hz																																																	
		<b>6</b>	30hz	<b>16</b>	80hz																																																	
		<b>7</b>	35hz	<b>17</b>	85hz																																																	
		<b>8</b>	40hz	<b>18</b>	90hz																																																	
		<b>9</b>	45hz	<b>19</b>	95hz																																																	
<b>10</b>	50hz	<b>20</b>	100hz																																																			



--	--	--

#### 4.6 Characteristics 6 [Notify]

0	Reserved	Reserved
1	Reserved	Reserved
2	Reserved	Reserved
3	Read Percentage	0-100
4	Read Finished	0/1
5	Sector Number being erased	0-2048
6		
7	Read Page Number	0-32768
8		
9	Write Page Number	0-32768
10		
11	Reserved	Reserved
12	Battery level value	0-100
13	Temperature value	Current Device temperature
14	Battery charging	0. Battery is charging 1. Battery not charging

## 4. Requirements

- 5.1. Temperature and Humidity:
- 5.2. The BIOPOT3 V2 shall withstand during operation the following conditions: Temp: 0°C – 40°C and Humidity: 20 -90 %
- 5.3. The BIOPOT3 V2 shall withstand during storage the following conditions: Temp: -10°C – 50°C and Humidity: 20 -90 %
- 5.4. Atmospheric Pressure (kPa)
- 5.5.
- 5.6. Materials in contact with the patient or operator:  
All substances in contact with the skin are typical medical grade plastics and other biocompatible materials which are commonly used in adhesive electrode systems.
- 5.7. Cleaning and disinfection:  
The reusable electronic component shall routinely be cleaned and disinfected immediately after they have been used, the cleaning procedure is according to the standard OR cleaning procedure.  
The electrode sticker is for one time use so no need for cleaning.
- 5.8. Waste product
  - 5.8.1. The electrode sticker is for single use, after using it you can throw it to a regular trash can.
  - 5.8.2. The electrical Biopot device and the charger must be taken to an approved household hazardous waste site or battery-recycling center.

## 6. Instructions for use

- 6.1. The device may be used by a researcher/professional in a specialized lab for data acquisition. It can also be used by a home user together with a mobile device or PC for data acquisition.
- 6.2. The device is used by a researcher/professional in a specialized clinical lab for data acquisition. It can also be used by a home user together with a mobile device or PC for data acquisition.
- 6.3. The primary user may need to collect biopotential data, bioimpedance data, accelerometer data, or some combination of all. Data collection requirements may scale up over time (i.e. start basic with increased functionality with project progress)
- 6.4. The primary user may be working on their own software, or they may require just a basic application for data collection only with no software development required. BOTH potential requirements must be covered: easy software for data collection without development, SDK for customized software development.
- 6.5. The primary user likely requires flexibility in the sensor part of the project, and may need to use a basic off-the-shelf option, or may be working on their own sensors. Multiple options may be required in a single customer/project.
- 6.6. The secondary user is a student that requires an easy platform for simple data collection for basic research purposes. This does not require programming.
- 6.7. The secondary user requires basic software and hardware interfacing and may be encouraged to do their own software development. SDK with examples would be useful, compatibility with MatLab or Labview, and or Python required.
- 6.8. In order to start using the device, the operator need to attach the electrode sticker to the device (as described in section 1.6). To attach the electrode to the patient, see "Electrode patch user manual". After the preparation you can start recording data (as described in section 2.10).
- 6.9. Hazards related to use of the medical device
  - 6.9.1. Inadequate operating instructions, such as:
    - 6.9.1.1. Inadequate specification of accessories to be used with the medical device may lead to device malfunction.
    - 6.9.1.2. Inadequate specification of pre use checks may cause inability to perform data recording resulting with inconvenience to the patient and staff or physical harm to patient resulting from heating up of the electrodes
    - 6.9.1.3. Inadequate specifications of service and maintenance may lead to device

- 6.9.2. Use by unskilled/untrained personnel may lead to wrong usage and wrong results that may pose a risk to the patient's life or Harm to the device.
- 6.9.3. Inadequate warning of hazards likely with re-use of single use medical devices. Like the electrode sticker.
- 6.9.4. Incorrect measurement and other metrological aspects may cause delay in treatment or unnecessary treatment due to hardware failure.
- 6.9.5. When using the device be aware of Sharp edges or points that may cause to operator or patient injury.
- 6.9.6. Replacement of the components is forbidden.
- 6.9.7. To prevent electrical shocks or burns that may happen as a result of mains or other high voltage from other device finds a current return path through BPT V1 SFC's test should be done.
- 6.9.8. For electromagnetic interference risks look into RMF.
- 6.9.9. Power supplied from an external DC power source, have hazardous situation. HW damage, may result in patient injury (burns)
- 6.9.9.1. Non continuous Connection may lead to Device misuse, interfering of other devices, Device not in range and RF noise
- 6.9.10. Lithium batteries:
  - 6.9.10.1. are classified as Class 9 hazardous material UN3480, a power source with high energy density and dangerous materials in a closed metal case!
  - 6.9.10.2. Lithium batteries must be handled only by qualified and trained personnel.
  - 6.9.10.3. Don't replace the lithium battery.
  - 6.9.10.4. Instruction for fire extinguishing:
    - Extinguish with water, if possible, cover battery completely in water
    - Extinguishing with water will produce fluoride, phosphate, fluoride-oxide and carbon-oxide.
    - Alternatively extinguish with a CO2 extinguisher
- 6.10. General:
  - 6.10.1. Do not stack and locate the system close to other equipment
  - 6.10.2. Other cables and accessories may negatively affect EMC performance
  - 6.10.3. Portable RF communications equipment including antennas, can affect medical electrical equipment. Any equipment should not be used closer than 30 cm (12 inches) to any part of the BioPot system, including cables specified by manufacturer"
  - 6.10.4. For BioPot system that are classified as class A according to CISPR 11, the instructions for use include the following note: "The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment."

#### 6.11. Technical Description:

- 6.11.1. The technical description describes precautions to be taken to prevent adverse events to the patient and Operator due to electromagnetic disturbances.
- 6.11.2. Compliance for each emissions and immunity standard or test specified by this collateral standard, e.g. emissions class and group and immunity test level
- 6.11.3. Any deviations from this collateral standard and allowances used
- 6.11.4. All necessary instructions for maintaining basic safety and essential performance with regard to electromagnetic disturbances for the expected service life.

## 7. General information

7.1. Manufacturer: Senso Medical LTD

7.2. Costumer: Senso Medical Labs Ltd.

7.3. Costumer's Address:

P.O. Box:	2653
Industrial Park Building, Mount	St.\ Industrial Zone
precipice.	
Nazareth Galilee	, City
Postal Code:	1612102
Israel	
Fax No.:	046800668
Tel No.:	046800668

7.4. Equipment Model No.: BP-30-08-00061

7.5. Equipment Serial No.: 61

7.6. For labels explanation look into Biopot labels file.

## 8. List of accessories, detachable parts, and materials


1. Biopot device
2. Electrode sticker
3. Charger (Docking + USB to Type C cable)
4. Box includes all the above parts

## 9. Declaration – electromagnetic emissions

Declaration – electromagnetic emissions		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group1 Class A	The BioPot SYSTEM uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Class A	The BioPot SYSTEM is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the BioPot SYSTEM or shielding the location.
Voltage Fluctuations And Flicker IEC 61000-3-3:2013	Complies	

## 10. Declaration – electromagnetic immunity

Declaration – electromagnetic immunity			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	8 kV contact 2, 4, 8, 15kV air	8 kV contact 2, 4, 8, 15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	2 kV for power supply lines N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s)  2 kV line(s) to earth  2 kV Signal input/output) to earth	1 kV line(s) to line(s) 2 kV line(s) to earth N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5cycle at 0°, 45°, 90°, 135°,180°, 225°, 270° and 315° 0% UT; 1cycle and 70% UT; 25/30 cycles Single phase at 0° 0% UT; 250/300 cycle	0% UT; 0.5cycle at 0°, 45°, 90°, 135°,180°, 225°, 270° and 315° 0% UT; 1cycle and 70% UT; 25/30 cycles Single phase at	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BioPot SYSTEM requires continued operation during power mains interruptions, it is recommended that the BioPot SYSTEM be powered from an uninterruptible power supply or a battery.

		0° 0% UT; 250/300 cycle	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 (A/m)	30 (A/m)	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE</b> UT is the a.c. mains voltage prior to application of the test level.			
Conducted RF IEC 61000-4-6	3V, 6V	3Vrms, 6V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BioPot SYSTEM, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{12}{V_2} \right] \sqrt{P}$ $d = \left[ \frac{12}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{23}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. D Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3V/m	3V/m	
	3V from 0.15 to 80MHz; 6V from 0.15 to 80MHz and 80% AM at 1kHz          3V/m from 80MHz to 2.7GHz	3V from 0.15 to 80MHz; 6V from 0.15 to 80MHz and 80% AM at 1kHz          3V/m from 80MHz to 2.7GHz	

## 11. Recommended separation distances between portable and mobile RF communications equipment and the BioPot SYSTEM

Recommended separation distances between portable and mobile RF communications equipment and the BioPot SYSTEM				
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz outside ISM bands $d = [\frac{3,5}{V_1}] \sqrt{P}$	150 kHz to 80 MHz in ISM bands $d = [\frac{12}{V_2}] \sqrt{P}$	80 MHz to 800 MHz $d = [\frac{12}{E_1}] \sqrt{P}$	800 MHz to 2,5 GHz $d = [\frac{23}{E_1}] \sqrt{P}$
0.01	0.12	0.2	0.4	1
0.1	0.37	0.64	1.3	2.6
1	1.17	2	4	8
10	3.7	6.4	13	26
100	11.7	20	40	80

## 12. Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment							
Test frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)	Compliance level (V/m)
385	380 – 390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1.8	0.3	27	27
450	430 – 470	GMRS 460, FRS 460	FM <sup>c)</sup> ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
710 745 780	704 – 787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9	9
810 870 930	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0.3	28	28



1720	1 700	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulati on <sup>b)</sup> 217 Hz	2	0.3	28	28
1845	–						
1970	1 990						
2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulati on <sup>b)</sup> 217 Hz	2	0.3	28	28
5240	5 100	WLAN 802.11 a/n	Pulse modulati on <sup>b)</sup> 217 Hz	0.2	0.3	9	9
5500	–						
5785	5 800						

## FCC Warning

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses And can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement.

The device can be used in portable exposure condition without restriction.