

LED-B7S Equipment

Installation Manual

Version: V1.0

Version Notes:

Version	Revision Notes	Revised By	Date
V1.0	First Draft	Guo Shuai	2024 .01. 30

Table of contents

1. Overview.....	3
2. Specifications & Features.....	3
2.1 Product Specifications Sheet.....	3
2.2 Main Uses and Scope of Use.....	5
2.3 Product View	5
2.4 Product Dimensions.....	5
2.5 Product Interface Description.....	6
3. Description Of Power Supply and Network Supply Mode.....	8
4. Installation Instructions.....	8
5. Cleaning Instructions.....	9

1 Overview

LED-B7S Bedside Terminal is a product designed for telemedicine and smart ward call systems.

As a networked LCD all-in-one device, this product includes a casing, terminal control mainboard, and mounting board. It is used for nurse calling, information transmission and display, and also supports multimedia broadcasting and video-on-demand services.

2 Specifications and features

Product features: appearance, stability, high definition, practicality, maintenance and scalability, etc.

2.1 Product Specifications Sheet

1	Basic parameters	Material: plastic structure /tempered glass	
2		Installation: wall mounting	
3		color: White	
4	Motherboard parameters	CPU: A55, 4 cores	
5		GPU: G52	
6		Memory : 1G	
7		External storage: 8G	
8		System: Android 12	
9	Screen parameters	Size: 7 inches	
10		Screen ratio: 16: 9	
11		Resolution: 1024 × 600	
12		Brightness: 22 0cd/m²	
13		Uniformity: >80%	
14		Viewing angle: U/D/R/L (CR>10): 8 5 /8 5 /8 5 /8 5	

15		Screen life: 20 000 hrs	
16	I/O interface	Power interface: 5pin 1.25, with lock, DC 12-24V	
17		Network interface: 4pin 1.25, with lock, 100M LAN	
18		Handle interface: Fusion TYPE-C interface	
19		Reset button: Hidden machine reset button	
22	Supported formats	Audio format: MP3/WMA/AAC	
23		HD video format: RMVB/AVI/MPG/MKV/VOB/MP4	
24		Image format: JPEG/BMP/PNG	
25	power	Rated power: $\leq 5\text{ W}$	
26		Standby power: $\leq 1\text{W}$	
27		Working voltage: DC12V/2A	
28	Peripheral parameters	size : 207.3 * 106.6 * 15.4mm	
29		Net weight of the whole machine: 0. 28 KG	
30		Packing instructions: bubble bag + carton	
31	Overall performance	Microphone: 1	
32		Sound output: 8 ohms 1 watt	
34		WIFI/Bluetooth , NFC card reader , camera: 200W	Optional
35		Power supply method : non-standard POE (12-24V) DC power supply 12V Standard POE (48-54V)	Optional
36	touchscreen	Capacitive 10-point touch screen, using industrial-grade touch chip	
37	other	Anti-static: contact $\geq 4\text{KV}$, non-contact $\geq 8\text{KV}$	
38		Heat dissipation method: aluminum profile heat sink passive heat dissipation	


2.2 Main Purposes and Scope of Use

Telemedicine, smart ward

2.3 Product view



2.4 Product Dimensions (mm)




2.5 Product interface description




Front view

- 1 TYPE-C call button interface
- 2 Hidden reset button
- 3 Light sensitive area
- 4 Camera
- 5 Call/End button
- 6 Locate button
- 7 Function button
- 8 Microphone
- 9 Call button
- 10 Nightlight button
- 11 Nightlight
- 12 Microphone





Rear View



- 1 Network data input interface, 4pin 1.25 pitch connector
- 2 Power input interface (12V-24V), 5pin 1.25 pitch connector

3. Description of power supply and network methods


- ◆ POE power supply, the POE switch outputs 48V to the POE splitter, and the POE splitter reduces the voltage to 12V and outputs the powered product;
- ◆ Power over Ethernet - Ethernet power supply technology refers to the transmission of data signals for some IP- based terminal equipment without any changes to the existing Ethernet Cat.5 cabling and infrastructure. Technology that provides DC power supply to such equipment. POE technology can ensure the security of existing structured cabling while ensuring the normal operation of existing networks and minimizing costs.
- ◆ Power over Ethernet (PoE) is a power distribution technology that transmits power and data over twisted pairs to any device connected to an Ethernet network;




Note: When setting up a static IP network, DNS cannot be 0 or blank. An IP address must be filled in, otherwise the save will fail.

4. Installation Instructions

- ◆ Hanging board size chart (mm)
- ◆ Drill tail screws for locking plate (PA3*12mm)*4pcs



- ◆ Anti-theft screw installation diagram
- ◆ Anti-theft screws (PM2.5*18mm)*1pcs



5. Product cleaning instructions



Cleaning products Use a soft lint-free cloth or cotton lint-free cleaning tool. Do not squeeze the display screen hard when cleaning. When wiping, it is recommended to wipe from one side of the display screen to the other until it is completely wiped clean.

Caution: Do not use rough cloths or abrasive paper for wiping. Do not use cleaners containing chemical components, and never spray liquid directly into the product's ventilation holes to avoid damaging the product through electrical short-circuiting!

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a

particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC Radiation Exposure Statement

This equipment must be installed and operated in accordance with provided instructions and the antenna(s) used for this transmitter must be installed to provide a separation distance of at least 20 cm from all persons and must not be co-located or operating in conjunction with any other antenna or transmitter. End-users and installers must be provided with antenna installation instructions and transmitter operating conditions for satisfying RF exposure compliance.