50mm

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Intelligent visual ear sampler  
(Product Multi Color General Manual)

### APP download method

Please use the scan function on your mobile browser to scan below QR code download and installation

Software QR code

desktop icon

#### Cooler

[IOS app store](#) (苹果商店)

[Android Google](#) (安卓国外)

[Android china](#) (安卓国内)

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Select after scanning the code  
[Android Google](#)  
(安卓国外) option  
After entering the  
interface click to select  
a third-party browser to  
download

Select after scanning the code  
[IOS app store](#) option  
open [App store](#)  
app store.

You can also directly  
[App store](#) search for  
[Cooler](#) app store.

#### Cooler

[IOS app store](#) (苹果商店)

[Android Google](#) (安卓国外)

[Android china](#) (安卓国内)

**Kind reminder:** Scan and download  
code on Android phones Reminder for  
safety, please leave it on your mobile  
desktop updates - find again Pure  
mode - click close to download.

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### Connection process

**Step 1.** First ensure the battery level when receiving the product Adequate, alternating blue and red lights flashing is sufficient insufficient battery, please charge in a timely manner

**Step 3.** Open the middle box after opening the APP interface After the following X displays letters and numbers, indicates successful connection, then click start. The key will display the image

**Step 2.** Open the phone settings - WIFI function and connect it Cooler xxxx is a hot topic  
(Note: This WIFI is machine transmitted and cannot be used for internet access)

**Step 4.** If you encounter this prompt, please do not switch to WIFI Click (Use) to proceed. If you make a mistake in not using the option, You can go back to step 2 to reconnect the hotspot

### Product parameters

Product Name:	Intelligent visual ear sampler
Model:	S7
Weight:	~15 grams
Transmission method:	CMOS
Battery capacity:	~180mAh
Usage duration:	~50 minutes
Input current:	DC5V-300mA
Lens pixels:	Ultra high definition camera
Connect hotspots:	Cooler XXXX
Network standards:	IEEE802.11b/g/h
Best distance:	1.2-1.6CM
Working frequency:	2.4GHz
Working temperature:	-5~40 degrees Celsius

### points for attention

- When the image is unclear, please wipe the lens with a cotton swab
- When cleaning the lens, please use a professional alcohol swab and gently wipe it
- Before use, please observe the surroundings and use in a safe and undisturbed environment, Avoid collisions
- This product is not suitable for children under the age of 5
- Children are prohibited from using this product alone. It must be operated under adult guidance
- Please keep this product out of reach of children
- Do not use the endoscope in liquid to avoid damage caused by immersion in water
- Please put the silicone case provided on the ear spoon for use
- This product is only used as a personal care item and cannot be used as medical equipment
- Please place the product in a cool and dry place, away from sources of fire and corrosive materials
- The product comes with a built-in battery, which needs to be charged once a month if not frequently used
- The temperature of the product will gradually increase slightly (up to a maximum of 38 °C) during use It has no impact on the human body, please use with confidence

### Indicator light status

- Press and hold the power button for two seconds, and the blue light will flash
- The WIFI successfully transmitted by the mobile software connected to the machine is displayed with a constant blue light
- When the battery is low, the red and blue lights flash alternately
- Flashing red light during charging (charging when the machine is turned off)
- The fully charged blue light remains on continuously

### Troubleshooting

#### 1. Screen lag

- The data transmission distance is too far, and the product is not compatible with the phone during use Need to maintain optimal distance
- In the case of insufficient battery, it can also cause lagging, Please charge in a timely manner before use

#### 2.WIFI signal cannot be transmitted/intermittent

- The lower the battery, the weaker the WIFI signal. Please fully charge before use
- Beyond the WIFI signal transmission distance, the product and phone must maintain Best distance
- There are too many other interfering signals in the surrounding area, which can affect the transmission quality, Please try to use it in a relatively pure environment
- Ensure that the product has successfully connected to the phone signal

#### 3.No image display

- Confirm if the phone has successfully connected to the product hotspot
- After confirming the connection, if there is still no display, please close the app Reconnect and open the app

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference, and  
(2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception,

which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:  
-- Reorient or relocate the receiving antenna.  
-- Increase the separation between the equipment and receiver.  
-- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.  
-- Consult the dealer or an experienced radio/TV technician for help. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction

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50mm宽\*110mm高

110mm