



# User Guide


# Introduction

Thank you for choosing Keepr! Keepr uses a combination of this breath alcohol measuring device and a paired mobile app interface in order to capture, track, and share your breath samples with your Contacts.

Contacts include any trusted accountability partners who are either helping you make the best, safest decisions whenever you consume alcohol, or who are helping you maintain a particular level or frequency of sobriety.

**Note:** Keepr™ is not intended to treat, cure, prevent, mitigate or diagnose a disease.

Model number: HHB-P000 Rev 1

 **keepr** 11035 Aurora Avenue,  
Des Moines, IA 50322

# Table of Contents

- 1** Introduction
- 2** First Time Setup
  - 2.1** Perform Initial Keepr Device & App Startup
  - 2.2** Add Contacts for Sample Sharing
  - 2.3** Set Up Sharing Conditions
- 3** Taking Your Sample
  - 3.1** Taking Your Sample
- 4** Contacts
  - 4.1** Add a Contact
  - 4.2** Request a Sample
  - 4.3** Set Up Sharing Conditions
  - 4.4** View a Contact's Current Location
  - 4.5** Demote a Contact
  - 4.6** Remove a Contact
- 5** History
  - 5.1** View Sample History
  - 5.2** Share Your History
- 6** Schedules
  - 6.1** Create a Schedule
- 7** Additional Procedures
  - 7.1** Unpair Bluetooth Connection
  - 7.2** Use Device Without App
- 8** Symbols Glossary
- 9** Maintenance and Care

# Table of Contents

- 9** Maintenance and Care
  - 9.1** Cleaning
  - 9.2** Storage
  - 9.3** Recalibration
- 10** Specifications
- 11** Warnings
- 12** Warranty and Disclaimers
- 13** Troubleshooting

# First Time Start Up

We recommend completing the following procedures before you start taking and sharing samples using Keepr.

# First Time Start Up

## Perform Initial Keepr Device and App Startup

To use your Keepr device with your mobile phone for the first time, complete the following steps.

1. Use the USB-C power cord to connect your Keepr device into a compatible power block (refer to the Specifications section for more details).

**Result:** The device LED screen displays the Charging  icon.

2. Install and open the Keepr app on your phone and complete the following substeps.
  - A. Select the option to log in.
  - B. Follow the steps to create a new account.
  - C. Ensure you are logged in to your account.

# First Time Start Up

## Perform Initial Keepr Device and App Startup

---

**D. Optional.** When prompted, for best accountability results, grant the app access to your camera and/or location when the app is in use.

**Note:** This is not required in order to use the Keepr app and device.

**3.** To connect the device to your phone app, complete the following substeps.

**A.** In the phone app, go to the **My Device** tile and start the device pairing process.

**B.** Press and hold the button for 5 seconds to enter **Pairing Mode**.

**Result:** The bottom light flashes blue for either 60 seconds, or until it successfully pairs with the app.


**C.** At the same time, using the Keepr app on your cell phone, activate device pairing.

# First Time Start Up

## Perform Initial Keepr Device and App Startup

---

**Result:** One of the following results occurs:

- If a Bluetooth connection is established, the screen displays the Bluetooth Connected  icon. In addition, the name of the device displays on the list of available devices in the app.
- if a Bluetooth connection is not established in 60 seconds, the bottom light stops flashing blue. The app on your phone also displays a notification that the pairing was not successful.
- Repeat the previous substeps as many times as necessary to achieve connection.

**Result:** You are now ready to take Breath Alcohol Content samples using Keepr, whenever you want.



# First Time Start Up

## Add Contacts For Sample Sharing

---

For your Keepr to be the most effective accountability tool it can be, it is crucial to set up sharing with your accountability partners, called Contacts in the app.

As part of the First Time Setup for your Keepr app, you may add one or more Contacts using the steps below. After you add them, Contacts can log into the app and become one of your Authorized Users. Your Authorized Users are your primary accountability partners, and therefore have additional capabilities. For more detailed information about Contacts and Authorized Users, refer to the Contacts tile section later in this guide.

To add a Contact in the Keepr app, complete the following steps.

1. On the home screen, navigate to the **Contacts Tile**.
2. On the **Contacts Tile**, select the option to manage requests.

**Result:** The Requests screen displays.

3. Select the option to invite a new Contact.
4. At the prompt, enter the email address of the person you want to add as a Contact.

# First Time Start Up

## Add & Promote Contacts For Sample Sharing

**Result:** The invited email address displays under Outgoing Requests on the Requests screen, and the intended recipient receives an emailed invitation to download the app.

- 5.** The other user must log in to the app, navigate to their Requests screen and accept your invitation under their Incoming Requests.

**Result:** The invited user displays in your list of Contacts.

- 6. Optional.** If you want the other user to become your Authorized User, allowing them to better act as your accountability partner, the other user must complete the following substeps.

- A.** The other user must log into the app and navigate to the **Contacts Tile**, select the desired Contact.

- B.** The other user must select your name from their list of Contacts.

- C.** On the subsequent screen, they must request Authorized User privileges.

# First Time Start Up

## Add & Promote Contacts For Sample Sharing

---

**Results:** The following results occur:

- On your **Contacts Tile**, the other user now displays as an Authorized User, denoted by the text “Authorized” appearing below their portrait.
- On the other user’s Contacts tile, you now display a Monitored User, denoted by the text “Monitored” appearing below your portrait.

You are now connected to one or more Contacts and can freely share samples with them.

# First Time Start Up

## Set Up Sharing Conditions

---

For best results when using Keepr, we recommend opening each of your Contacts' sharing conditions and ensuring that your information is being shared to a degree that you are comfortable with. Your Authorized Users are always notified whenever you update your sharing conditions with them, so be sure to communicate openly about what Keepr information you intend to share with them and whether you ever need to change what you are sharing.

To set up what information you share with each of your Contacts, complete the following steps.

- 1.** On the home screen, navigate to the **Contacts Tile**.
- 2.** On the **Contacts Tile**, select the desired Contact.
- 3.** On the subsequent screen, to begin adjusting sharing conditions, select **Edit** and then toggle any of the following conditions on or off.
  - **Breath Alcohol Content.** If enabled, the Contact can see the Breath Alcohol Content value of your shared sample results. Breath Alcohol Content value of your shared sample results.
  - **Sample Location.** If enabled, the Contact can see the GPS location the app captured for each of your shared breath samples.

# First Time Start Up

## Set Up Sharing Conditions

---

- **Sample Photo.** If enabled, the Contact can see the photo the app captured for each of your shared breath samples.
- **Live Location.** If enabled, the Contact can see your current GPS location whenever they select your name in their list of Contacts.
- **Sample Requesting.** If enabled, the Contact can request that you take a breath sample.
- **Scheduling.** Only Authorized Users can enable this sharing conditions. They can only enable or disable this condition when they first become your Authorized User. If enabled, your Authorized User can create, modify and view your sample-taking schedule.

**4.** Save your changes.

# Taking Your Sample

Take regular samples using your Keepr device. All samples you take are stored in the Keepr app for later viewing.

# Taking Your Sample

## Take a sample using the device

To take a sample using the device, complete the following steps.


**Prerequisites:** Ensure that your Keepr device and your phone app are paired and have sufficient battery life. For pairing instructions, refer to the First Time Setup section.


To avoid contaminating your sample, wait 20 minutes after eating, drinking, or smoking before completing these steps.

1. Open the Keepr app on your cell phone and log in to your account.

**Note:** If you want this sample to be retained in the app and to count towards any existing breath sample requests, you must use the app alongside the Keepr device.

2. To wake up the Keepr device, open the lid.

**Result:** The screen activates, displaying the **Getting Ready**  icon.

3. Wait for the screen to stop displaying the **Getting Ready**  icon and to instead display the **sample now** message.

**Note:** If the **sample now** message does not display, ensure that the device lid is still open.

# Taking Your Sample

## Take a sample using the device

---

**4.** In the Keepr app on your cell phone, complete the following substeps.

**A.** On the main screen, select the large button at the bottom of the screen.

**Note:** If any open sample requests are assigned to you, a green dot appears next to this button.

**Result:** The sample delivery screen displays, showing your camera view, your current sample requests from Contacts, and an option to switch between a private and shared sample.

**B.** Indicate whether your sample is **Shared** or **Private**.

**Note:** Private samples are for your own personal use, and are therefore not visible to any of your Contacts. You can still view and optionally share private samples in the History menu and its exported email reports.



# Taking Your Sample

## Take a sample using the device

---

- C.** Hold your phone so that your face is displayed in the camera frame.
- 5.** While still holding your phone in place and ensuring you aren't covering the exhaust port at the base of the device, take a deep breath, bring the mouthpiece at the top of the Keepr device to your lips, and perform a sustained blow until you hear the device make a click.

**Result:** The app indicates that the device has collected enough volume of your breath in order to calculate your current Breath Alcohol Content.

- 6.** Wait for the device to calculate your Breath Alcohol Content.

# Taking Your Sample

## Take a sample using the device

---

**Results:** The following results occur:

- The device displays your Breath Alcohol Content, up to 2 decimal points (such as **0.00** ).
- The phone app displays your Breath Alcohol Content, current GPS location, your mid-test photo capture, and the current date and time. This information is now shared with Contacts according to your established Sharing Conditions.
- The phone app indicates whether your sample satisfies any of your pending requests with Contacts.

After 10 seconds of displaying your reading, the Keepr device returns to the Getting Ready **...** status.

# Contacts

Interact with your accountability partners in  
Keepr.

# Contacts

For your Keepr to be the most effective accountability tool it can be, it is crucial to set up sharing with your accountability partners, called Contacts in the app. In this tile, you can access a large variety of features related to your Contacts, including: adjusting what information you share with each Contact; requesting samples from them; adding Contacts; becoming an Authorized User for another Contact; removing them; viewing their current location data (if they share it with you); and viewing various ingoing and outgoing requests.

There are 2 levels of Contacts you can have in your Contacts list:

## **Basic Contact:**

Displays in your Contacts list as a Contact with no additional icon. Note the following regarding Basic Contacts:

- Whenever another Keepr user accepts your email invitation, or whenever you demote an Authorized User, they appear as your Basic Contact.
- Your Basic Contacts can potentially request breath samples from you (and get notified if you miss the 20 minute request window), view your sample history, and view your location, but only if you have set up their Sharing Conditions accordingly.
- Only Basic Contacts can be removed from your Contacts list.

# Contacts

## **Authorized User:**

- Denoted with the word “Authorized” in your Contacts list.  
Your Authorized Users have the following unique features in the app:
- Authorized Users can create, modify, and enforce a schedule for your samples, if you allow it. For more information, refer to the Schedules Tile section of this guide.
- Authorized Users receive notifications in the following additional circumstances:
  - When you miss a scheduled sample.
  - When you update their Sharing conditions, such as turning off photo or location sharing.
  - When you demote them to a Basic Contact.

When a Contact first becomes your Authorized User, they have the option to gain a Sample Requesting and Scheduling for your account. Note that you cannot turn the Scheduling condition on or off after the moment your Contact becomes your Authorized User

- Authorized Users must first be demoted back to a Basic Contact before you can delete them from your Contacts list.
- In addition, if you act as an Authorized User for any users, those users appear denoted with the text “Monitored” in your Contacts list.

To complete any action related to Contacts in the app, refer to the following procedures.

# Contacts

## Add a Contact

---

To add a Contact in the Keepr app, complete the following steps.

1. On the home screen, navigate to the **Contacts Tile**.
2. On the Contacts tile, select the option to manage requests.

**Result:** The Requests screen opens.

3. Select the option to invite a new Contact.
4. At the prompt, enter the email address of the person you want to add as a Contact.

**Result:** The invited email address displays under Outgoing Requests on the Requests screen, and the intended recipient receives an emailed invitation to download the app.

5. The other user must log in to the app, navigate to their Requests screen and accept your invitation under their Incoming Requests.

**Result:** In your **Contacts Tile**, the invited user displays in your list of Contacts.

You are now connected to one or more Contacts and can freely share samples with them.

# Contacts

## Request Authorized User Privileges

---

To become an Authorized User, granting you more capabilities you can utilize as an accountability partner for another user, complete the following steps.

- 1.** On the home screen, navigate to the **Contacts Tile**.
- 2.** On the **Contacts Tile**, select the desired Contact.
- 3.** On the subsequent screen, request Authorized User privileges.
- 4. Optional.** If you want the ability to schedule tests for your Monitored User, enable the Scheduling option..

**Results:** The following results occur:

- In your **Contacts Tile**, the other user now displays as a Monitored User, denoted by the text as “Monitored” appearing below their portrait.
- On the other user’s **Contacts Tile**, you now display as an Authorized User, denoted by the text “Authorized” appearing below your portrait.

# Contacts

## Request a Sample

To request a sample from one of your Contacts, complete the following steps.

**Prerequisite:** Your Contact must first allow sample requests from you. Refer to the Set Up Sharing Conditions procedure.

1. On the home screen, navigate to the Contacts tile.
2. On the Contacts screen, on the row of the desired Contact, request a breath sample.

**Results:** The following results occur:

- On your screen, the Contact's row lists what time the request expires, with an option to cancel the sample request. It also displays in your list of outgoing requests.
- The Contact receives a notification that you requested a sample.
- After the Contact takes their next shared sample, you will receive a notification. You can view their test results in your My Device tile and History tile, if the Contact is sharing their Breath Alcohol Content ratings with you.



# Contacts

## Set Up Sharing Conditions

---

To set up what information you share with each of your Contacts, complete the following steps.

- 1.** On the home screen, navigate to the **Contacts Tile**.
- 2.** On the **Contacts Tile**, select the desired Contact.
- 3.** On the subsequent screen, to begin adjusting sharing conditions, select **Edit** and then enable or disable the checkmark for any of the following conditions.
  - **Breath Alcohol Content.** If enabled, the Contact can see the Breath Alcohol Content value of your shared sample results.
  - **Sample Location.** If enabled, the Contact can see the GPS location the app captured for each of your shared breath samples.
  - **Sample Photo.** If enabled, the Contact can see the photo the app captured for each of your shared breath samples.
  - **Live Location.** If enabled, the Contact can see your current GPS location whenever they select your name in their list of Contacts.

# Contacts

## Set Up Sharing Conditions

---

- **Sample Requesting.** If enabled, the Contact can request that you take a breath sample.
- **Scheduling.** Only Authorized Users can enable this sharing condition. They can only enable or disable this condition when they first become your Authorized User. If enabled, your Authorized User can create, modify, and view your sample-taking schedule.

**4.** Save your changes.

# Contacts

## View a Contact's Current Location

---

To view the current GPS location of one of your Contacts, complete the following steps.

**Prerequisite:** Your Contact must first share their current location with you. Refer to the Set Up Sharing Conditions procedure.

- 1.** On the home screen, navigate to the Contacts tile.
- 2.** On the Contacts tile, select the desired Contact.
- 3.** On the subsequent screen, view their GPS location.

# Contacts

## Demote a Contact

---

If you want to demote an Authorized User, or if you want to remove a Contact altogether, you must first demote them to the Basic Contact level. If you want to demote an Authorized User, complete the following steps.

1. On the home screen, navigate to the **Contacts** tile.
2. On the **Contacts** tile, select the desired Authorized User.
3. On the subsequent screen, demote the Authorized User's privileges.
4. When prompted, accept the demotion.

**Result:** The icon denoting the Contact's higher level will be removed. The Contact will be notified of this change in the app.

# Contacts

## Remove a Contact

---

If you no longer want a particular user to be one of your Contacts, complete the following steps. Note that you can only remove Basic Contacts.

1. On the home screen, navigate to the **Contacts Tile**.
2. On the **Contacts Tile**, select the desired Contact.
3. On the subsequent screen, remove the Contact.
4. When prompted, accept the removal.

**Result:** You and the other user are no longer Contacts. Any previously created schedules you or that user made for each other are now deleted.

# History

View and share any of your recorded breath samples

# History

On the History tile, you can review your own breath sample history, as well as any of your Contacts' sample histories, if they elect to share it with you. You can review various data points associated with a user's breath sample history, including:

- Total samples taken during the date range.
- Total sample requests missed during the date range.
- Total samples with a Breath Alcohol Content rating above 0.00 during the date range.
- Your highest recorded Breath Alcohol Content ratings per day during the date range, displayed in the bar graph.
- All recorded Breath Alcohol Content ratings during the date range, displayed in a list under the bar graph.

In addition, this tile allows you to generate and email a PDF report based on your own history. This report changes to match whatever date range you select in the tile. The report includes the following information:

- Date generated
- Keepr device ID
- Original Keepr activation date
- Your name
- Your email

# History

- Date range included in the report
- Total samples taken during the date range
- Total missed samples during the date range
- Number of samples with a Breath Alcohol Content value over 0.00 during the date range (in list format)
- Specific dates and times of any missed samples during the date range (in list format)
- Specific dates, times, and values of any samples with a Breath Alcohol Content value over 0.00 during the date range (in list format)

Whether just viewing your history or sharing your history report, you have the option to include private samples or not.

To view your or another user's sample history, or to share your history, refer to the following procedures.



# History

## View Breath Sample History

---

To optionally view your or another user's sample history, complete the following steps.

**Prerequisite:** If viewing another user's sample history, they must first grant you the Breath Alcohol Content sharing condition.

- 1.** On the home screen, navigate to the **History Tile**.
- 2.** At the top of the **History Tile**, to view any Contact's history (including your own), select the desired user.
- 3.** Specify whether to view results by days, weeks, or months at a time.
- 4.** Select the specific dates you want to view.
- 5.** Specify whether to view either Private samples or Shared samples.

**Note:** You can only view Private samples if viewing your own history.

# History

## View Breath Sample History

---

**Results:** Whenever you modify any of the previously mentioned settings, the data values displayed on screen all update accordingly:

- Whenever you modify any of the previously mentioned settings, the data values displayed on screen all update accordingly.
- Total samples taken during the date range.
- Total sample requests missed during the date range.
- Total samples with a Breath Alcohol Content rating above 0.00 during the date range.
- Your highest recorded Breath Alcohol Content ratings per day during the date range, displayed in the bar graph.
- All recorded Breath Alcohol Content ratings during the date range, displayed in a list under the bar graph.

# History

## Share Your History

---

To share your history with any Contact, Authorized User, or other email address, complete the following steps.

**Note:** You can only share your own history.

- 1.** On the home screen, navigate to the **History Tile**.
- 2.** On the **History Tile**, select the option to share your report.
- 3.** Specify the desired date range.
- 4. Optional.** To include private samples (if any) in the report, enable the Private Samples option.

**Note:** If private samples exist in your history and you do not enable the Private Samples option, your shared report will indicate that you did not share your private samples.

- 5.** In the sharing screen, select any Contacts you want to share this history with.

# History

## View Sample History & Share Results on App

---

**6. Optional.** If you want to share the results with non-Contacts, for each desired email address, select New Email and enter the corresponding email addresses in the text field.

**7.** After all desired Contacts and non-Contacts have been selected, select the button at the bottom of the screen to generate and share your report.

**Result:** Everyone you shared the history with will receive an emailed PDF showing the sample history.

# Schedules

View and share any of your recorded breath samples.

# Schedules

If you act as an Authorized User for another user, you can create a schedule requiring them to take breath samples within specified time ranges. You must have scheduling permissions before you can create the other user's schedules.

When you use a time range (instead of a single, specific time), the app will randomly notify your Monitored User to take a breath sample at some point during that range. For example, if you schedule a test with a time range of 4:00 PM to 5:00 PM on Tuesdays, then the Monitored User will get a notification sometime in that range, such as 4:24 PM that Tuesday. After they receive the notification, they will have the standard amount of time to complete a breath test before it counts as a missed sample for them.

When scheduling multiple test requests for the same day, keep in mind that you cannot send or schedule more than one request every 20 minutes. For example, every Monday, your Monitored User receives two scheduled test requests: one sometime between 1:00 PM and 3:00 PM, and another at 9:00 PM. However, in this example scenario, the app would not allow you to add an additional test request for Mondays at 9:15 PM, because that wouldn't be more than 20 minutes after the existing 9:00 PM test request.

To view and modify your sample schedules, or to create a sample schedule for yourself or your Monitored User, refer to the following procedures.

# Create Schedules

## Create Schedules

To create a schedule for yourself or for your Monitored User (for whom you act as Authorized User), complete the following steps.

**Prerequisite:** If creating a schedule for your Monitored User, you must have previously enabled the Scheduling sharing condition when you first became their Authorized User.

- 1.** On the home screen, navigate to the **Schedules Tile**.
- 2.** To view existing schedules, complete the following substeps.
  - A.** In the **Schedules Tile**, select the desired user from the list.

**Result:** All scheduled tests for that user display.

**B. Optional.** If you want to disable or enable previously existing tests that you created for this user, in the list of scheduled tests, toggle the desired tests to Off or On.

- 3.** For every breath sample that you want to schedule for this user, complete the following substeps.
  - A.** Select the option to schedule a new test.

**Result:** The Add New Recurring Test screen displays.

# Schedules

## Create Schedules

---

**B.** In the list of weekdays, select each day on which you want the user to receive a sample request.

**C.** For **Time**, complete one of the following actions.

- To schedule a request for a specific time on the selected weekdays, enter a single value in the first text field, leaving the second text field blank.
- To schedule a request to occur sometime during a specified time range, enter the desired start time and end time for the time range.

**Note:** Keepr will send the user one single request at a randomly selected time within the specified range and weekday(s).

**D.** To schedule the test, select the button at the bottom of the screen.

As a result, the Monitored User is notified once per test you scheduled for them. The new tests are only visible to your Monitored User and to you, their Authorized User.

From now on, until you disable this Schedule or are demoted from being their Authorized User, the user will be prompted to test at the given time. If this test is ever ignored, it will count as a missed sample for them.



# Additional Procedures


The following procedures are not typically used in day-to-day Keepr operations.

# Additional Procedures

## Unpair Bluetooth Connection


If you ever need to associate your device with an alternate phone, or to simply disassociate it from your current cell phone, you can easily unpair your device. To unpair your device's Bluetooth connection, complete the following steps, starting from a powered down state.

1. Ensure your device has sufficient power and press the button on the device.

**Result:** The LED screen activates, displaying the Getting Ready  icon.

2. In the phone app, navigate to the Settings menu.
3. Find the device in the list of paired devices and select the option to forget the device.

**Results:** The following results occur:

- If the device was successfully disconnected, Bluetooth Connected icon disappears.
- The app indicates that the device is no longer connected.
- The device returns to the Getting Ready  status.

# Additional Procedures

## Use Device without App

---

You can use your Keepr device without being paired with the cell phone app. To do so, complete the Take a Sample procedure but ignore any steps taking place in the app.

**Note:** The device cannot store, log, or share any samples you take when being used in this state. Therefore, for a superior accountability experience, we recommend only using Keepr when connected to the app service.

# Symbols Glossary



## **Charging**

Indicates that the device is currently charging using the USB-C cable.



## **Critically Low Battery**

Indicates that the device has no power remaining.



## **Critically Low Battery**

Indicates that the device has no power remaining.



## **Bluetooth Connected**

Indicates that the device is currently connected to a cell phone using the Keepr app.



## **Getting Ready/Calculating**

Indicates that you should wait to act until the device has finished preparing. This icon displays both during the device's startup process and also while the device calculates your Breath Alcohol Content level.

# Symbols Glossary

sample now

## **Sample Now**

Indicates that the device is ready to collect your breath sample. For instructions, refer to the **Take a Sample** procedure above.

0.00

## **Results**

Indicates the Breath Alcohol Content calculated from your sample. This sample displays up to 2 decimal points, and is automatically logged and shared with your Contacts according to your configured permissions.

E01

## **Error**

Indicates that the device has encountered an error. Refer to the **Troubleshooting** section. If the error starts with the letter “U”, then this was a user error rather than a device error.



# Maintenance and Care

# Maintenance and Care

## Cleaning

---

Do not submerge the Keepr device in water or put in dishwasher. Only clean the device using a damp cloth with water.

## Storage

---

It is crucial to store and operate the Keepr device only in specific conditions. The device is unlike most other consumer electronics, containing a variety of sensitive components used to measure your Breath Alcohol Content. Device abuse and improper storage may result in deteriorated measurement instruments, which may provide inaccurate sample result values.

The instruments should only be stored in the following climate.

- 60 to 80 degrees Fahrenheit
- 30 to 60 percent relative, non-condensing humidity

## Recalibration

---

To ensure the most accurate Breath Alcohol Content results, your device should be recalibrated every 6 months. To calibrate it, you will need to contact CST at (999) 999-9999 and follow their instructions to mail it in for calibration.

# Specifications

<b>Power Source</b>	USB-C  Output: 5V/0.5 – 3A  IEC/UL 62368-1, or IEC/UL 62368-3, or IEC/UL 60950-1 certified
<b>Operating Range</b>	5 ° - 40 °C
<b>BrAC Range</b>	0 - 0.400 g/210 L
<b>BrAC Accuracy</b>	+/-0.005 g/210L at 0.100 g/210L at 25°C.



# Warnings

**Warning:** Do not insert anything other than the USB-C cable provided with the product into the USB-C port.

**Warning:** Inspect the charging cable prior to use. Do not use the charging cable if it has become frayed or otherwise damaged.

**Warning:** Do not inhale through the device.

**Warning:** In the event of lithium ion battery leakage, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash with copious amounts of water and seek medical advice.

**Warning:** Do not immerse the device in any fluid.

**Warning:** The mouthpiece presents a choking hazard for children. Keep the mouthpiece out of the reach of children.

**Warning:** Detachable mouthpiece contains a magnet. Ingesting magnets can cause serious injury or death. Seek immediate medical attention if magnets are swallowed or inhaled.

**Warning:** Users with allergies to plastics should not use the device.

# Warranty and Disclaimers

**Note:** When disposing of this device, do so in accordance with local laws and guidelines.

**Note:** This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:  
**(1)** This device may not cause harmful interference, and **(2)** this device must accept any interference received, including interference that may cause undesired operation.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

# Warranty and Disclaimers

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**Note:** The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc and any use of such marks by Keepr™ is under license. Other trademarks and trade names are those of their respective owners

**Note:** When disposing of this device, do so in accordance with local laws and guidelines.

**Note:** The FCC requires the user be notified that any modifications or changes made to this device that are not expressly approved by Keepr™ may void the user's authority to operate it.

**Note:** Keepr™ is not intended to treat, cure, prevent, mitigate or diagnose a disease.

# Troubleshooting

To troubleshoot errors and issues, refer to the following recommended procedures.

Error or Issue	Recommendation
EXX	'E' followed by any two-digit number, indicating an error that cannot be fixed by the user. Contact Customer Service.
U01	Insufficient Breath – This error displays when you have not exhaled hard enough or for long enough into the Keepr for it to take a sample of your breath.
U02	Temperature Out of Range – This error is displayed when the device detects temperatures above or below its recommended operating range.

