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
60mm

60mm

Smart Ring USER MANUAL


Directory

English	1~7
German	8~13
French	14~19
Italian	20~25
Spanish	26~31
Japanese	32~37


-1-

Component Description

Transparent Resin, Charging port, Oxygen/Heart Rate Sensor, Test Status Indicator

Packing List


Smart Ring, Charging cable, User manual

-2-

Power On


Connect the charging cable to power on the smart ring

-3-

Connection APP

1. Search and download SmartHealth APP or scan the QR code on the right to select the appropriate APP version to download.
2. Open the app and follow the instructions to match the ring.


-4-

Lighting status

The red light flashes when charging and turns off when fully charged

The green light flashes during the heart rate test and the red light flashes during the blood oxygen test

-5-

Wearing style

Wear the ring on your own suitable finger, choose a comfortable size and wear it on the index finger, middle finger, and ring finger to better monitor the data

Wear the sensor on the inside of your finger
Heart rate and blood oxygen measurements are more accurate

-6-

Note

Before using and operating the device, please read and follow the following precautions to ensure optimal device performance and avoid dangerous or illegal situations.

Warning!

Please keep the charging cable and device in a dry environment when charging. Do not touch the charging cable with wet hands and avoid it being immersed in liquids such as water.

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Do not disassemble or modify the battery, insert foreign objects, or immerse it in water or other liquids to avoid causing battery leakage, overheating, fire, or explosion.

Do not drop, squeeze, or puncture the battery. Avoid the battery causing the big external pressure, which can cause internal short circuits and overheating issues.

FC CE RoHS

-37-

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This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure statement

The device has been evaluatec to meel general RF exposure requirement. The device can be used in portable exposure condition without restriction.