




User Manual

Farmer's Hands FHL2 Manual



Product Overview



- Farmer's Hands makes precise farming possible by collecting data on and monitoring calves' suckling, rumination, cough frequency, rest time, eating, physical activity, and more.
- Analyzes delicate movement signals using a sensor on the calf's neck.
- Comfortable-to-attach and reusable sensor provides early disease detection.



Suckling



Rumination



Cough



Resting



Eating

Component

FHL Sensor(Front)



Button

LED


FHL Sensor(Back)



Label

Label

Parts



tags

Sensor

Buckle

Weight

How to use

1. Turn on the power by pressing and holding the button for 3 seconds.
2. Make sure the red LED is blinking.
3. Press the button to turn on the LED and attach the sensor to the neckband.

Precautions

- Do not disassemble or apply excessive shock to the product.
- Do not allow children to chew on the product.
- Do not drop from a high place.
- Do not immerse in water.
- As this is a device that uses BLE wireless communication, do not wrap the sensor with an electromagnetic shielding material such as aluminum foil.
- The product works normally with only one battery installed, but the battery life will be halved.

Product Specification



Model Name	FHL2
Size	130mm x 67mm x 33mm (LxWxH)
Weight	216g
Battery	Li-Ion 18650 / 3.7V 7000mAh (3500mAh*2EA)
Battery use	6 Months
Network	BLE 4.2
Sensor	Accelerometer + Gyroscope (6 axis)
Memory	128Mb Flash Memory

FCC Compliance Statement

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

CAUTION : Any Changes or modifications not expressly approved by the manufacturer could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- **Reorient or relocate the receiving antenna.**
- **Increase the separation between the equipment and receiver.**
- **Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.**
- **Consult the dealer or an experienced radio/TV technician for help.**

Thank You

