
Digital Microscope

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ⓘ When the battery is installed for the first time, the microscope must be plugged into the power cord and energized for 2 seconds before the battery protection function is activated.


1.Introduction to the structure and functions of microscopes	03
2.Installing and turning on and off the microscope	05
3.Image observation usage process	06
4.Microscope function operation instructions	07
5.WIFI function operation instructions	09
6.Microscope connects to computer	10
7.Microscope connects to TV	11
8.Microscope parameters	11
9.Precautions	12
10.Troubleshooting	13

This manual includes the usage instructions of the digital microscope. Please read this manual carefully before use and keep it with the instrument. The product color matching in the manual is for reference only, please refer to the actual object.

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
(Figure 1)

① Install the battery

Open the battery cover, install the battery (18650 battery), and then close the battery cover (Figure 1).

(Figure 2)

② Power active

Each time after installing the battery, you must first use the power cable to power up the microscope and activate the battery management circuit (1-2 seconds is enough) (Figure 2).

(Figure 3)

③ Turn on/off

Turn on: Short press the power button on the back of the microscope until the indicator light on.

Turn off: Press and hold the power button and wait 2-3 seconds until the microscope turns off (Figure 3).

(Figure 4)

④ Charge

Battery indication: normal use; charging state (Figure 4).

(Figure 5)

① Turn on

Press the switch button on the back until the indicator light on (Figure 5).

(Figure 6)

② The use of glass slide

Open the specimen clip of the mobile platform, place the glass slide and fix it (Figure 6).

(Figure 7)

③ Use of stage

Tighten the XY axis knob of the platform by hand and align the specimen to the center of the objective lens. When viewing the specimen, twist the XY axis knob by hand to move the specimen and view other parts of the specimen (Figure 7).

(Figure 8)

④ Use of objective lens and focusing

- Turn the objective lens to the center (vertical direction) (Figure 8)
- According to the type of specimen, choose the appropriate lens to view it (generally view it from a low magnification lens first) (Figure 8)
- After selecting the appropriate objective lens, slowly rotate the focus knob upward or downward until the image is clear (Figure 9)

(Figure 9)

(Figure 10)

① Light source adjustment

Press the light source knob to switch the incident and transmitted light, and rotate it left or right to adjust the brightness of the light source (Figure 10).

(Figure 11)

② Button operation instructions

Press the WiFi icon button once, the display will show the operation instructions of the button, press it again to exit the operation instructions page (Figure 11).

(Figure 12)

③ Digital amplification adjustment

Turn the WiFi icon knob, adjust image digital enlargement or reduction (Figure 12).

(Figure 13)

④ Photograph

Press the photo button to switch to photo mode, and press it again to take photos (Figure 13).

(Figure 14)

⑤ Video

Press the video button to switch to video mode, press it again to start recording (at this time the display will light up with a flashing red dot), press it again to stop recording (Figure 14).

(Figure 15)

⑥ Album

Press the album button and the display will enter the album list to view photos or videos. Rotate the WiFi icon knob to rotate the desired image file, press it again to play the image file; press the album button again to exit the album (Figure 15).

(Figure 16)

⑦ Menu button


Press the menu button to enter the first level menu page. You can rotate the WiFi icon knob to select the parameters and language selection that need to be set. After selection, press the WiFi icon button to turn it into the OK button (Figure 16).

(Figure 17)

⑧ Measure

- Press the menu button to enter the first-level menu page, select the measuring ruler option by rotating the WiFi icon knob, select a different ruler according to different objective lenses (such as: now 4X objective lens, select the 4X ruler), and then press the knob to confirm, and finally press the menu key to exit the menu page; (Figure 17).

(Figure 18)

- Measurement and conversion method : (as shown in the figure, the diameter of the circular hole is measured by 40 small cells from the scale, $40 \times 17.9 (4X \text{ objective}) = 716UM (\text{micron})$, the conversion millimeter is equal to $0.716MMO$) (Figure 18).

- Calculation formula:

The number of small cells measured * base size (micron)

4X objective base 17.9; 10X objective base 7.37;

40X objective base 1.86

(Figure 19)

① Turn on WIFI

Press and hold the WIFI icon button (2 seconds) to enter WIFI mode. At this time, the microscope emits a WIFI signal and displays a password ID (TN300_488XXXXXXX) Default password (12345678) (Figure 19).

(Figure 20)

② Mobile APP connects to WIFI


Click the Connect WIFI button on the APP to enter the WIFI settings page, select the WIFI corresponding to the microscope WIFI signal name (TN300_488XXXXXXX) and connect. Return to the APP main interface and wait for the "Start Observation" text to be displayed (different mobile different response speeds) (Figure 20).

(Figure 21)

③ Start observation

After clicking to start observation, you will enter the real-time image preview page, where you can control the microscope to take photos, record videos, photos enlargement or reduction, and adjust colors (Figure 21).

(Figure 22)

④ APP enter album

- Click the device album icon to enter the device album; (Figure 22).
- Click the local album icon to enter the local album; (Figure 22).


Note: The device album refers to the album in the memory card of the microscope, and the local album is the album in the mobile APP.

⑤ Microscope controls APP

The button control of the microscope can also control the APP to take photos, record videos, enter the album, adjust color, and enter the menu (Figure 23).

(Figure 23)


(Figure 24)

⑥ Other instructions

While the APP takes photos and videos and saves them in the photo album of the mobile phone APP, the memory card of the microscope will also simultaneously save the images and videos data (Figure 24).

6. Microscope connects to computer >>>

EN

(Figure 25)

① Install computer software


Install StrangeView imaging software on the computer (Figure 25).

(Figure 26)

② Connect data cable

Plug one end of the microscope Type-C data cable into the computer USB and the other end into the microscope Type C interface. At this time, the microscope will display 3 options (memory, camera, charging). Rotate the WIFI icon knob, select the camera, and then press the WIFI knob to confirm. (Figure 26).

(Figure 27)

③ Computer viewing images

After connecting to the computer, it will take 1-3 seconds to automatically install the driver, then open the image software MatataXplore, and you can view the images in real time on the computer and edit the images (Figure 27).

(Figure 28)

① Video data cable connects to TV

Plug the video data cable into the HDTV jack of the microscope, and the other end into the TV port, then select the input signal source corresponding to the TV, and the TV can watch (Figure 28).

Product name	Digital microscope
Product color	White
Sensor	Optical Format 1/2.8"
	Pixel Size 2.9 μm x 2.9 μm
	SENSOR resolution 1920(H) x 1080(V)
Pixel specifications	200w
Objective lens	achromatic 4X, 10X, 40X
Screen size	7-inch high-definition IPS screen
Display resolution	1024*600

① Operation

- Components should not be exposed to direct sunlight and should be placed in a dry, clean environment away from high temperatures and severe vibrations.
- The microscope is a precision instrument and should be handled with care to avoid shock and collision during transportation.
- In order not to affect the clarity of the image, avoid dirt or fingerprints from leaving on the lens surface.
- Do not turn the left and right zoom handwheels in opposite directions at the same time, otherwise a malfunction may occur.

② Maintenance

- All lenses should be kept clean. If there is fine dust, you can use a blower to blow it away or gently wipe it with cotton gauze; if there are oil stains and fingerprints, you can use cotton gauze dipped in a small amount of a 3:7 ethanol-ether mixture to gently wipe it away.
- Do not use organic solvents to wipe the remaining surfaces of the microscope, especially the surfaces of plastic products. Neutral detergent should be used for cleaning.
- You should not disassemble or assemble the microscope by yourself to avoid affecting the performance of the microscope.
- When the microscope is not in use, it should be covered with a dust cover to prevent it from coming into contact with dust, and stored in a moisture-proof place to avoid rust or mold.
- To maintain the performance of the microscope, regular inspections are recommended (contact your nearest dealer for details).

① Unable to turn on after installing battery

When the battery is installed for the first time or is installed again after being removed, it needs to be activated by plugging into the microscope power cable interface to start the battery power supply.

② After the microscope is turned on, it automatically turns off

Generally, the battery power is low. You can see the battery power indicator icon in the upper right corner of the display screen, or there is a flashing indicator light on the back of the microscope, which indicating that the battery power is low (Figure 4).

③ Display image has black corners

Generally, the objective lens converter of the microscope is not turned in place. There will be obvious positioning resistance after the objective lens is turned in place. The objective lens must be perpendicular to the stage.

④ Microscope cannot see image

Generally, the specimen is not aligned with the center of the lens. First, turn the 40X objective lens downward, then adjust the XY axis knob to align the specimen with the center of the lens, and slowly adjust the focus up and down until the image is clear.

⑤ The produced specimen viewing image is dark and unclear

This is because the specimen is thick and the light from the bottom light source cannot penetrate the specimen.

⑥ Unable to see plant leaves or object surfaces

This type of specimen is impermeable to light. You need to press the light source knob to switch the incident light to view the specimen.

⑦ Why is there no image when connecting to HDTV TV


When connecting an HDMI signal cable to a TV or other display device with HDMI, pay attention to setting the input signal source of the display device.

⑧ No image when the microscope is connected to the computer via USB


The microscope is connected to the computer via USB. A menu will appear on the microscope's display. You need to select the connection mode and then press the WIFI icon knob.

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FCC Caution:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.