
# ScopeAround

## Ear, Nose & Throat Smart Otoscope

**Model: SA200**

# Contents

Device information	1
Safety information	2
Symbol Description	3-4
Introduction	5
Features	5
Specifications	6-7
What's in the Box?	8
Product Overview	9
Operational Interface	10-11
Steps and warnings for use	12-15
Basic Operations	16-20
Transfer Files	20
Export image/video data	21
Battery charging	22
Safety Warnings	23
Troubleshooting	24
Regulatory Information	25

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# Device information

## 1.1 General description

The Ear, Nose & Throat Smart Otoscope manufactured by Dongguan HeMS Electronic Technology Co., Ltd. (hereinafter referred to as “Dongguan HeMS”) is an active, invasive and non-sterile device. The Ear, Nose & Throat (ENT) Smart Otoscope is a battery-powered and digital hand-held device with WiFi function, used in conjunction with mobile apps, supporting iOS and Android platforms. The high-definition camera is equipped with 6 LED lights, which can present photos or videos on the mobile APP.

## 1.2 Intended purpose

The Ear, Nose & Throat Smart Otoscope is intended for the examination of the external ear canal and tympanic membrane of adult patients. It can also be used for general illumination of oral and nasal cavity and image capture.

## 1.3 Intended user

Medical staff and lay persons.

## 1.4 Intended patient population

Adult patients.

## 1.5 Contraindications

No contraindication.

## 1.6 Working principle

The Ear, Nose & Throat Smart Otoscope is equipped with a high-resolution camera at the tip. When inserted into the ear canal, the camera captures images of the ear canal, eardrum, and any abnormalities present. It has built-in LED lights that illuminate the ear canal, providing a clear view for the camera to capture detailed images. The captured images are transmitted in real-time to a display screen. The display screen is a smartphone, through a wired or wireless connection. These can then be reviewed through the App.

## Safety information

**Warning:** Fire and explosion hazard. Do not operate the device in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide; in oxygen-enriched environments; or in any other potentially explosive environment.

**Warning:** Patient injury risk. Use only approved compatible accessories according to the manufacturer's directions for use and only connect items that are specified as part of the device system. Failure to do so can affect patient and operator safety, compromise product performance and accuracy, and void the product warranty.

**Warning:** Patient injury risk. Do not shine otoscope in patient eyes as this can cause eye injury.

**Warning:** Patient injury risk. Do not use otoscope without attaching a speculum as this can cause ear injury or cross-contamination.

**Warning:** Cross-contamination risk. The speculum is reusable. Before and after reuse, cleaning and disinfection are necessary. Otherwise, it may spread contamination from one patient to another patient.

**Warning:** Do not use the otoscope in an MR environment.

**Warning:** Do not use this otoscope to examine the eye.

**Caution:** It is important to note that while the otoscope can be a valuable tool for examining the ears, it should be used with caution and under appropriate circumstances. In cases of uncertainty or if there are any concerns, it is always recommended to consult a healthcare professional for proper evaluation and management.

**Caution:** Make sure to see a healthcare provider if the patient has the following symptoms: discharge or bleeding from the ears, dizziness, fever, ringing in the ears and severe pain in the ears.

**Caution:** The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


## Symbol Description

Graphical symbols	Description
<div data-bbox="130 275 270 337"> <div>EC</div> <div>REP</div> </div>	Authorized representative in the European Union
	Manufacturer
	UKCA marking
	The operation guide must be read
	CE mark: indicates that the device complies with the EU 2017/745
	Medical device
	Unique device identifier
	Date of manufacture

	Serial number
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice
	Type BF applied part
	The device complies with FCC Part 15
	Recycling

## Introduction

The Scopearound SA200 is an easy-to-use hand-held video otoscope for capturing videos and images of the ear or body. It utilises interchangeable attachments making it ideal for multiple applications.

## Features

- 1.85-inch high-definition screen display with 360\*360 resolution.
- Detachable lens or accessories for easy assembly and storage.
- Supports photo and video capture with touch of a button.
- Supports TYPE-C charging and data transfer.
- Supports app operation with Wi-Fi, comes with Wi-Fi live video streaming capability.

## Specifications

Monitor	Screen Type	1.85–inch IPS LCD Screen
	Screen Resolution	360*360
	Image Resolution	720*720
	Video Resolution	720*720
	Screen Material	Tempered Glass
	Port	Type–C
	Lightness adjustment	bright, medium bright, slightly bright, off
	Magnification	1.0~1.5x

Camera	Diameter of Lens	4.3mm
	Camera Megapixel	1MegaPixel
	Field of View	78°
	Focal Length	1–3cm
	Camera Lights	6 LED lights
	Surface and edges	Free of pores, cracks, and remainders of tooling agents
	Water resistance	The inside of otoscope can withstand a positive pressure of 22kPa for more than 1 minute without air leakage


# Specifications

Others	Working Voltage	3.7V
	Charging Voltage	DC 5V
	Storage Capacity	Built-in 8GB Memory
	Battery Capacity	1380mAh Lithium Polymer Battery
	Video format	AVI

Network	Image and video display	APP, WIFI
	Network standard	IEEE802.11b/g/n
	Antenna	Built-in IPEX antenna
	Operating frequency	2.4GHZ
	Data rate	802.11b/g/n


Working condition	Temperature: 0~45°C (32~113°F)
Transportation and storage condition	Temperature: -20°C~50°C
	Humidity: 15%~95%
	Atmospheric pressure: 50kPa~106kPa
Service life	5 years

## What's in the Box?

1. SA200 Digital Otoscope
2. Storage
3. USB Cable
4. User Manual & Safety Warning Manual
5. Specula Attachment
6. Spatula Attachment
7. Skin Examination Attachment

# Product Overview

## Control Button Descriptions


	Press to scroll upwards in the menus Zoom 1.0~1.5x
	Press to scroll downwards in the menus Left/Right Ear Switch
<b>ok</b>	Confirm Operation Play and Pause in Photo/Video mode Long Press: Enter/Exit Album
	Short Press: Enter Device Settings/Delete Image or Video in Album Long Press: Turn On/Off Wi-Fi
	– Home Page Status Short Press: Take a Photo Long Press: Start/End a Video Recording – Non-Home Page Status Short Press: Back to Home Page
	Short Press: Adjust the Camera Lights <b>By default, the lowest brightness level is used when turning on the device. When looking at your ears, be sure to use the lowest brightness level to avoid burning your ears.</b> Long Press: Power On/Off

# Operational Interface

①	Up
②	Down
③	Date
④	Left/Right Ear Mode
⑤	OK Button
⑥	Button Description
⑦	Settings Button
⑧	Time
⑨	Battery Status Display

# Operational Interface

Settings	Description
Playback	View captured photos and videos
Language	Supports English, German, French, Italian, Spanish, Russian, Japanese, Chinese
Wi-Fi	Turn On/Off the Wi-Fi
Date Stamp	Turn On/Off the Date Stamp
Date/Time	Setting the Date and Time
Format	Formatting the files stored in the TF card (Note: Formatting will delete all files on the TF card, please proceed with caution.)
Auto Power Off	Close — 5mins — 10mins
Factory Reset	Restore all device settings to factory defaults with one click
Version	View device model and version information

## Steps and warnings for use

### Examination of the ear canal, oral cavity, and skin

1.1 First, disinfect the accessories to prevent cross contamination.

1.2 When inspecting the ear canal, use ear the specula accessory (for protection, to prevent excessive penetration and ear injury)

1.3 When examining the oral cavity, use the tongue depressor accessory (pressing down on the base of the tongue can provide a better view of the throat).

1.4 When checking the skin, use the skin accessories.

## Installation/disassembly of otoscopes and accessories

### 2.1 Installation/removal of the ear accessories

When installing, simply align the earbuds with the earpieces and insert them into place. When disassembling, hold the base of the probe with two fingers, and pull out the base of the ear specula accessory with the other hand (**Note: try to hold the base as much as possible to prevent the silicone sleeve on the head from pulling off**).

Align the otoscope and insert it

There is a limited position on the earring accessory, just insert it into place

When disassembling, use two fingers to protect the otoscope, and use the other hand to hold the base of the accessory and pull it out

## 2.2 Installation/disassembly methods for tongue depressor accessories

During installation, simply align the tongue depressor with the earpiece and insert it into place. When disassembling, hold the base of the probe with two fingers and pull out the base of the tongue depressor accessory with the other hand (**Note: try to hold the base as much as possible**).

Align the otoscope and insert it

There is a limited position on the earring accessory, just insert it into place

When disassembling, use two fingers to protect the otoscope, and use the other hand to hold the base of the accessory and pull it out

### 2.3 Installation/disassembly methods for skin accessories

When installing, simply align the skin accessory with the earpiece and insert it into place. When disassembling, hold the base of the probe with two fingers and pull out the base of the skin accessory with the other hand.

Align the otoscope and insert it

There is a limited position on the earring accessory, just insert it into place


When disassembling, use two fingers to protect the otoscope, and use the other hand to hold the base of the accessory and pull it out


## Basic Operations

1. Install and connect the otoscope lens or the accessory you need to use.

2. Long press the power button to turn on the device for further use.

# Basic Operations

## Connecting to the Mobile App


1. Download the app named Smart Otoscope, you can download it from Apple Store/Google Play or scan the QR code below.


2. When installing on an iPhone for the first time, when opening the app, a pop-up window pops up allowing the use of Bluetooth, local network devices, and location information (please allow or else it will not work).

3. When installing an Android phone for the first time, when opening the app, a pop-up window will appear allowing the use of location information (please allow or else it will not work), allowing access to photos (please allow or else it will not work).

4. Turn on the device's Wi-Fi (when turned on, press and hold the settings button to turn on Wi-Fi), and the Wi-Fi icon will flash on the right side of the main interface (as shown in the figure below).

## Connection method 1 – Bluetooth network connection

When the Bluetooth of the mobile phone is turned on, after opening the device's Wi-Fi, the APP will automatically obtain the latest device's Wi-Fi, and there will be a connection pop-up window on the main interface,

1. Click on the Connect icon on the pop-up window, select Add Device Wi-Fi, and the phone will automatically connect to the device Wi-Fi.


2. After successful connection, the Connected green icon will be displayed. Click View to enter the image interface.


## Connection Method 2– Manual Connection

If the connection method 1– Bluetooth network configuration is not successful, please use manual connection method or if the Bluetooth of the phone is not turned on and the device's Wi-Fi is opened, there will be no pop-up window on the main interface of the app, and you also need to manually connect the device's Wi-Fi.

1. Click on Connect your phone to Wi-Fi on the app interface, connect to the device's Wi-Fi (Wi-Fi name: Smart otoscope\_xxxxxx), and enter the Wi-Fi password: 55556666


2. After successful connection, click enter to enter the image interface.

**Note:** After connecting to the app, you can only operate it on your phone.

## Transfer Files

After connecting the device to the computer via a USB data cable, the following interface will appear on the right. There are two options (U-Disk Mode and Camera); use the up and down keys to select, and press the OK key to enter.

Number	Name	Description
1	U-Disk Mode	Selecting this option will cause the device to appear as a USB drive on the computer, allowing for file transfer and copying.
2	Camera	Selecting this option will allow the device to display images on the computer (image capture software such as Amcap is generally used). The procedure is as follows: After selecting this mode, open the Amcap software on the computer to view the images.

To transfer the images or videos to a PC, connect the otoscope to your PC with the USB cable.


## Export image/video data

### 4.1 Device file export method:

The device's photo/video files are saved on the device's memory card. When the device is turned on, it can be connected to the computer through a data cable to copy the photos and video files to the computer.

### 4.2 Data Export Method for Mobile APP (Taking Apple as an Example)

The first step is to click on the photo album on the app's image interface. The second step is to select the file you want to download. The third step is to click on the forward icon in the lower right corner. The fourth step is to select save image from the pop-up prompt to save the file to the Apple phone's photo album.

## Battery charging

Use a conventional 5V/1A or 5V/2A USB interface power supply to charge the device. The charging connection diagram is as follows

Charging line

Power on charging status

Shutdown charging status

## Describe the cleaning and disinfection methods for otoscopes and accessories separately

### Cleaning and disinfection methods for otoscopes and accessories

1. Clean and disinfect the accessories using alcohol wipes from the packaging accessories.

2. Or use 75% medical alcohol for cleaning and disinfection (customers can purchase it themselves)

## Safety Warnings

 **The device's LED lights support three levels of brightness adjustment. When turned on, the default brightness is the lowest level. When looking at the ears, be sure to use the default lowest level of brightness to avoid burning the ears. If the brightness is adjusted back to the lowest level from the highest level, wait for the lens temperature to cool down before inserting it into the ears.**

- Keep away from damp and exposed environments.
- Avoid contact with sharp objects.
- Do not use alcohol, detergents, or other organic cleaning agents to clean the lens. It is recommended to use a soft cloth dampened with water.
- Do not pierce the product shell with nails or other sharp objects. Do not hammer or step on the product.
- Avoid using the product near heat sources or under extreme heat conditions.
- Do not disassemble the product in any way.

### EMC Compliance

- This device complies with EMC standard IEC 60601-1-2. It does not normally affect nearby equipment and devices. It is not normally affected by nearby equipment and devices. It is not safe to operate the device in the presence of high-frequency surgical equipment.
- The use of the device adjacent to or stacked with other equipment or medical electrical systems should be avoided because it could result in improper operation. If such use is necessary, the device and other equipment should be observed to verify that they are operating normally.

# Troubleshooting

## Unable to Power On

- Check if the device has power. Please charge your device.
- Press the reset button to force a reset.

## System Crash, Unable to Operate

- Press the reset button on the device to force a reset.

## No Response When Charging

- Ensure the charging port is free of dirt, foreign objects, moisture, etc.
- Try unplugging and plugging back in the charging cable.
- Try using a new charging cable.

## No Response When Connected to the Computer

- The device needs to be turned on to enter USB mode when connected to a computer.
- Make sure the data cable interface is free of dirt, foreign objects, moisture, etc.
- Try unplugging and plugging back in the data cable (including the computer's USB side).
- Try restarting the device or computer (sometimes the computer's USB may not be recognized).
- Try using a new data cable.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

# Regulatory Information

## FCC Regulatory Compliance

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## RF Exposure Compliance

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End user must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

The portable device is designed to meet the requirements for exposure to radio waves established by the Federal Communications Commission (USA). These requirements set a SAR limit of 1.6 W/kg averaged over one gram of tissue. The highest SAR value reported under this standard during product certification for use when properly worn on the body.

**Please follow the following guidelines.  
Otherwise, it may cause damage to the  
environment or equipment.**

- When removing the device from the environment with the highest storage temperature for use, it should be left in the working environment for at least 30 minutes to cool.
- When removing the device from the lowest storage temperature environment for use, place it in the working environment for at least 30 minutes for preheating.

**This product is suitable for patients to operate  
and use on their own.**

Name: MedPath GmbH

Address: Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

SRN: DE-AR-000000087 Tel.: +49 (0) 89 8130 6837

Email: info@medpath.pro

Name: Dongguan HEMS Electronic Technology Co., Ltd.

Address: Room 301, Building 2, No.3, Jiaoshe Xinxing Third Road,  
Dongkeng Town, Dongguan City, Guangdong Province, China, 523000

SRN: CN-MF-000042670 Tel.: +86-769-81013289 Email: don@hemselec.com

UK Responsible Person

Name: MedPath Limited

Address: 27 Old Gloucester Street, London WC1N 3AX, United Kingdom

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© ScopeAround SA200 D0918V1.3
12-month Limited Warranty