


MEDIEYE-W

User Manual

1. Device Configuration




- ① Therapy Section: This is the part that outputs the therapy light.
- ② Button Section: consists of a power button on the top and a mode button and mode indicator LED.
- ③ Indicating Section: Displays the current device status with LEDs.
- ④ Charging Section: Charges using a USB-C type charging port.

2. How to Use



- ① Press and hold the power button (top button) for 2 seconds to turn on the power
(When powered on, you will hear a sound, and the front top LED will light up)
- ② Press the button at the mode button(bottom button) to change the mode.
(There are 4 modes in total, each with 1, 2, 3 and 5 minutes of treatment. Each mode change will sound)
- ③ Once the mode is set, put it on and press the mode button for 2 seconds to start the treatment.
- ④ When the treatment is complete, press and hold the power button for 2 seconds to power off.

3. Note



- Features are not available while charging for safety reasons.
- While charging, you can check the charging status via the mode LED next to the button
- If it does not light up, it is a battery issue and should be charged before use.
- It is normal for the therapy LED to appear to be flashing.

MEDIAEYE-W

User Manual

4. Caution

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications (including the antennas) to this device that are not expressly approved by the manufacturer may void the user's authority to operate the equipment.

Note : This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Modifications not expressly approved by the manufacturer could void your authority to operate the equipment under FCC rules.

FCC RF Radiation Exposure Statement:

- This equipment complies with FCC RF Radiation exposure limits set forth for an uncontrolled environment. This device and its antenna must not be co-located or operating in conjunction with any other antenna or transmitter.