


Smart ring instruction manual

I. Product Introduction

This smart ring integrates high-precision pulse wave acquisition sensor, gravity sensor, temperature sensor, which can accurately measure physiological indicators such as heart rate, blood oxygen saturation, heart rate variability index, emotional stress, temperature, etc., and it can also monitor the amount of exercise, the number of steps, the workout record, the quality of sleep and so on. The rich capillaries in the belly of the finger facilitate the collection of physiological parameters, thus the smart ring has a higher degree of accuracy compared to other wearable devices.

Smart ring structure diagram:



2. Wearing method.

This product has eight sizes:

U.S. No.	European No.	Japan No.	Diameter(Inner)	circumference(Inner)
#5	#50	#9	15.7mm	49.3mm
#6	#52	#12	16.5mm	51.8mm
#7	# 55	#15	17.3mm	54.3mm
#8	#57	#17	18.1mm	56.9mm
#9	#60	#19	18.9mm	59.4mm
#10	#62	#23	19.8mm	62.2mm
#11	#65	#25	20.6mm	64.7mm
#12	#68	#28	21.4mm	67.2mm

☺ Use a model that fits your finger size and wear it with the sensor window facing your finger.

☺ It is easy to measure the circumference of the finger by wrapping a thin, non-elastic string around the finger.

3. Parameters

Waterproof rating: IP68

Battery capacity: 15-28mAh (the larger the ring size, the larger the battery capacity)

Battery life: 7-15 days (the longer the set monitoring cycle, the longer the battery life)

Bluetooth version: BLE5.0

Charging voltage: DC5V

Memory capacity: 512M NAND Flash

Execution standard: GB9706.1 & YY0505

4. Preservation Methods

- 1、Don't place the ring in the charging compartment when it is not used for a long time;
- 2、The charging compartment is not waterproof, it should be placed in a cool and dry place;
- 4、Try to avoid contact with acid and alkaline substances

5. Precautions

- 1, the set collection period should not be too short, frequent collection will lead to shorten the duration.
- 2、The charging compartment is not waterproof and should be placed in a cool and dry place;
- 3, try to avoid contact with acid, alkaline substances

6. Precautions

This device complies with part 15 of the FCC Rules. Operation is subject

to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception,

which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction

7. App



ChipletRing(iOS)



ChipletRing(Android)



GooglePlay



应用宝