

Heart Rate Monitor

Owners Manual





HORSEPAL

www.horsepal.com



HORSEPAL

Heart Rate Monitor



Using the Heart Rate Monitor and Handle

1. Thoroughly wet the horse in the oval area indicated in the image below. For best results, apply a small amount of ECG gel to the electrode pads.
2. Turn on the HRM sensor and place the Horsepal HRM handle against the horse (in the oval area shown in Image 2), applying slight pressure.
3. The heartbeat detection icon will flash if a good signal is detected. Once a good signal is received, the timer will start counting. The heart rate should be displayed within 8-10 seconds.

Troubleshooting

If the heartbeat detection icon is not flashing, adjust the handle slightly while keeping the electrodes within the oval area in Image 2. Try again and repeat until a connection is established.

If the horse has thick hair, clipping the hair in the oval area will significantly improve the HRM handle's ability to read the heart rate.❏



How to connect your Heart Rate Monitor to the Horsepal APP

1. Download the Horsepal 2.0 App by scanning the QR code below.
2. Open the app and follow the setup instructions. You will need to enter details for a horse and a rider before you can begin training.
3. Once setup is complete and horses and riders have been added, you can start training. Press the Start Training button in the app (red button at the bottom).
4. Ensure your HRM sensor is turned on, then navigate to the Sensor tab. Find your sensor in the list and select it.


Congratulations! You have successfully connected your Horsepal HRM sensor.

Note: Once your sensor is connected, it will automatically reconnect next time.

Download Horsepal 2.0

Access the full Horsepal experience
with our Application.

An ever-expanding universe of powerful
tools for Equine management.





For Instructions and Videos on
how to setup your device
Scan the below QR code



SCAN HERE

HORSEPAL

www.horsepal.com



Tested To Comply with FCC Standards

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment"